

Read manual thoroughly before attempting to install or operate lift.

This manual has been provided to assist you with lift installation and operation.

For further assistance please contact your authorized Harmar Mobility dealer or Harmar's Technical Service Department.

Tel: 800-833-0478 Fax: 866-234-0478

Email: tech@harmar.com

Dealer:

Serial Number:

TABLE of CONTENTS



PREPARATION

PREPARATION	1
	TOOLS REQUIREDp 2
	UN-PACKING THE LIFTp 3
	WIRING THE VEHICLEp 4
	TRIAL FIT PROCEDUREp 5
INSTALLATIO	V
	INSTALLATIONp 6
	ADJUSTMENTSp 9
	DOCKING DEVICESp 10
OPERATION	
	OPERATIONp 11
	SAFETY & MAINTENANCEp 13
TECHNICAL	

TECHNICAL

TROUBLE-SHOOTINGp 14 EXPLODED VIEWS / PARTS BREAKOUTp 15-30 THE FOLLOWING IS A SUGGESTION OF BASIC TOOLS RECOMMENDED TO HAVE ON HAND DURING INSTALLATION

AS WELL AS A GUIDE TO PORTRAY COMPLEXITY OF THIS TYPE OF INSTALLATION

THIS IS, BY NO MEANS, MEANT TO REPRESENT THE REQUIREMENTS AS A WHOLE, AS EACH INSTALLATION WILL

VARY TO SOME DEGREE. HARMAR LIFTS ARE DESIGNED TO INSTALL WITH AS LITTLE ASSEMBLY AS POSSIBLE.





AS ALWAYS, IF YOU HAVE ANY QUESTIONS, CONCERNS, OR COMMENTS, PLEASE FEEL FREE TO CONTACT HARMAR'S TECHNICAL SERVICE DEPARTMENT @ 1-800-833-0478 or tech@harmar.com


UNPACKING THE LIFT

IF ANY PARTS ARE MISSING OR IF ANY DAMAGE IS NOTED, IMMEDIATELY CONTACT THE DISTRIBUTOR FROM WHICH YOU PURCHASED YOUR LIFT. DO NOT ATTEMPT TO INSTALL OR USE THE LIFT WITH ANY MISSING OR DAMAGED PARTS

BE SURE TO CHECK THE CONTENTS OF THE BOX AGAINST THE PACKAGE CHECK LIST TO BE SURE THAT ALL THE PARTS ARE INCLUDED WITH THE LIFT THAT ARE REQUIRED FOR PROPER INSTALLATION









IMPROPER WIRING IS THE #1 CAUSE OF PROBLEMS IN THE OPERATION OF A VEHICLE LIFT.

FOLLOW THE WIRING INSTRUCTIONS CAREFULLY

LOCATED IN THE HARDWARE PACK IS THE VEHICLE WIRING HARNESS. THE HARNESS IS MANUFACTURED TO, AND COMPLIES WITH, THE SAE J1128 REQUIREMENTS, THE WIRE HARNESS IS APPROXIMATELY 23 ft LONG AND WILL ACCOMMODATE MOST VEHICLES.

YOUR WIRING HARNESS HAS BEEN LEFT UN-ASSEMBLED FOR EASE OF INSTALLATION. THE END CONNECTOR HAS BEEN INCLUDED SEPARATELY FROM THE HARNESS TO ALLOW THE INSTALLER TO THE RUN THE WIRE UNDER AND THROUGH THE VEHICLE WITH THE SMALLEST CLEARANCE POSSIBLE. FOLLOW THESE INSTRUCTIONS TO ASSURE PROPER INSTALLATION. DO NOT CONNECT RED WIRE UNTIL THE VERY END.

UNWIND THE HARNESS AND LAY IT FLAT.
ONE END OF THE HARNESS HAS 2 COVERED PINS.
THIS IS THE LIFT END OF THE HARNESS AND
GOES TO THE REAR CARGO AREA WHERE THE LIFT
WILL MOUNT.

- 1 BEGIN ROUTING THE WIRING HARNESS AT THE VEHICLE BATTERY.
 ATTACH THE BLACK WIRE TO THE NEGATIVE TERMINAL ON THE BATTERY.
 DO NOT ATTACH THE RED WIRE UNTIL THE END.
- 2 RUN THE WIRING HARNESS UNDER OR WHEN POSSIBLE THROUGH THE VEHICLE, BACK TO THE REAR CARGO AREA / TRUNK, ALWAYS LOCATE THE WIRING WHERE IT CAN NOT BE SNAGGED BY ROAD DEBRIS AND AWAY FROM THE GAS TANK. THROUGH THE VEHICLE IS BEST, GAINING ENTRY INTO THE VEHICLE THROUGH THE FIREWALL HOWEVER IF YOU WISH TO RUN THE HARNESS UNDER THE VEHICLE, A HOLE WILL MOST LIKELY NEED TO BE DRILLED TO GET THE WIRE INTO THE REAR CARGO AREA / TRUNK, USED SUPPLIED GROMMET IN HOLE TO PROTECT THE WIRE. IF THE HARNESS IS TOO LONG FOR THE VEHICLE, COIL THE EXCESS WIRE AND SECURE IT TO THE VEHICLE FRAME WITH SUPPLIED TIE WRAPS.


DO NOT CUT OR SHORTEN THE HARNESS.

3 ONCE THE HARNESS IS RUN INTO THE VEHICLE, REMOVE PIN'S PROTECTIVE TUBING.



A=RED

4 INSPECT THE PIN'S RETAINING FLANGES.
THE MAY HAVE BECOME DEFORMED
WHILE RUNNING THEM THROUGH THE VEHICLE.
THESE ARE CRITICAL TO SECURE THE PINS INSIDE
THE END CONNECTOR, ADJUST AS NEEDED.



MOLDED INTO

(5)

PLASTIC


5 FLIP OPEN END CONNECTOR'S HINGED RETAINER.
INSERT PINS AS SHOWN: B=BLACK,

VERIFY THE WIRES CAN NOT BE PULLED OUT BY MODERATELY PULLING ON WIRE.


BE SURE WIRE SEALS ARE INSIDE CONNECTOR, AND CLOSE THE END CONNECTOR'S HINGED RETAINER.

PLUG CONNECTOR INTO LIFT

6 LAST: ATTACH RED WIRE TO THE POSITIVE TERMINAL ON THE BATTERY



WHEN THE INSTALLATION
REQUIRES THE WIRING
HARNESS BE RUN ON THE
UNDERSIDE OF THE VEHICLE,
ROUTE THE HARNESS AWAY
FROM THE EXHAUST SYSTEM,
BRAKE LINES, FUEL LINES < GAS
TANK, PINCH POINTS, AND SHARP
EDGES. LOCATE THE WIRING
HARNESS WHERE IT CAN NOT
BE SNAGGED BY ROAD DEBRIS




PROBING FOR 12 VOLTS MAY INDICATE
A CONNECTION, BUT NOT NECESSARILY
A SUPFICIENT CONNECTION.
THE LIFT'S MOTOR CAN DRAW UP TO
30 AMPS AT SOME POINTS, REQUIRING
ALL OF THE AVAILABLE WIRE TO FLOW
PROPER CURRENT. POOR CONNECTIONS
ARE THE ALP PROBLEM ASSOCIATED WITH
A SLOW; INTERMITTENT, WARM.
POORLY PERFORMING MOTOR.
NOT ONLY WILL THE MOTOR PERFORM POORLY,
BUT WILL PREMATURELY DETERIORATE.

IF A SINGLE STRAND OF A MULTISTRAND WIRE IS MAKING CONTACT,

127 WILL APPEAR ON THE METER, OR TEST LIGHT,
BUT WILL NOT ALLOW THE MOTIOR TO OPERATE
IT IS ALWAYS BEST TO TEST BOTH


CURRENT AND VOLTAGE, OR RUN THE MOTION
WITH KNOWN GOOD SHOP BATTERY OR POWER

SOURCE WHEN TESTING FOR POWER.



TRIAL FIT PROCEDURE






THE AL-425 / AL-225 / AL-215 WILL REQUIRE 3 HOLES BE DRILLED IN THE VEHICLE, TO BOLT THE LIFT'S BASE TO THE FLOOR. THE POSITIONING OF THE BASE IS VERY IMPORTANT TO THE SUCCESSFUL INSTALLATION AND OPERATION OF THIS LIFT.


PRIOR TO ANY HOLES BEING DRILLED, IT IS HIGHLY RECOMMENDED THE LIFT BE RUN THROUGH ITS RANGE OF MOTION WHERE THE BASE IS THOUGHT TO BE MOUNTED.







FOLD THE CHAIR OR SCOOTER'S SEAT BACK (or REMOVE SEAT) PREFERABLY (2) WITH DOCKING DEVICE INSTALLED TO ROUGHLY OBSERVE THE CHAIR'S TRANSPORTABLE DIMENSIONS. MEASURE LENGTH, WIDTH, AND HEIGHT








SIMULATE EACH OF THE CHAIR'S DIMENSIONS (LENGTH, WIDTH, and HEIGHT) WHILE RUNNING THE LIFT THROUGH ITS FULL RANGE OF MOTION IN AND OUT OF THE VEHICLE'S CARGO AREA





ONCE IT HAS BEEN VERIFIED EACH OF THE HELD MEASUREMENTS DOES NOT INTERFERE WITH THE VEHICLE'S CARGO AREA / OPENING. MARK THE 3 LOCATIONS IN WHICH TO DRILL

ONCE CONFIDENCE HAS BEEN GAINED ON THE POSITION OF THE LIFT, AND THE "TRIAL FIT" PROCEDURE HAS BEEN COMPLETED, THREE 3/8" HOLES WILL NEED TO BE DRILLED IN THE MARKED LOCATIONS

CAUTION! AVOID THE VEHICLE WIRING, FUEL LINES, FUEL TANKS SPARE TIRES, ETC. WHEN DRILLING

INSTALLATION TIP:


CHOOSE A SAFE LOCATION IN WHICH TO DRILL A SINGLE PILOT HOLE FROM ABOVE IN ONE OF THE 3 MARKED LOCATIONS.


LOOK UNDERNEATH THE VEHICLE AND MEASURE OUT THE SPECIFIED DISTANCE FROM THE PILOT HOLE.

IF ALL IS CLEAR IN THAT AREA, CONTINUE TO DRILL PILOT HOLES IN REMAINING LOCATIONS.


INSPECT ALL THREE PILOT HOLES FROM UNDER THE VEHICLE. IF ALL CLEAR, DRILL 3/8" HOLES IN PILOT LOCATIONS

ONCE HOLES ARE DRILLED, FASTEN BASE TO VEHICLE FLOOR USING SUPPLIED HARDWARE IN THE HARDWARE PACK. BE SURE TO TIGHTEN ALL FASTENERS.





ALL THREAD KIT MAY BE REQUIRED





SOME VEHICLES HAVE STORAGE COMPARTMENTS IN THE REAR CARGO AREA, POTENTIALLY EFFECTING LIFT PLACEMENT. MOST OF THESE CAN BE OVERCOME WITH THE AID OF AN ALL THREAD KIT. SHOWN ABOVE IS THE TYPICAL APPLICATION FOR OVERCOMING THESE STORAGE COMPARTMENTS OR WELLS.

SPARE TIRE STORAGE:
SOME VEHICLES STORE
THE SPARE TIRE UNDER
THE FLOOR OF THE VEHICLE.
MAKE SURE THAT THE BASE
CAN BE BOLTED THROUGH
A STEEL FLOOR OR FRAME.


STORAGE


WELL






BA-01 INSTALLATION







ONCE BASE IS FASTENED TO THE VEHICLE, RAISE PLASTIC COVER, AND ATTACH LIFT TO BASE.

BE SURE TO COMPRESS ALL SPLIT LOCK WASHERS. (TIGHTEN ALL FASTENERS)







AL-425





ATTACH LIFTING ARM AS SHOWN.












The Docking Device is the interface between the lift and the Chair / Scooter.

Most chairs / scooters can be lifted either by the center seat post or the 4 post seat frame.



DD-21 Spreader Bar & Bridles




DD-19 C-Arm

Although most chairs / scooters can be lifted by the above Docking Devices, Many specific devices exist making the application even easier.

Please consult the instruction provided with your Docking Device, or contact our technical service department to find out exactly how this device is right for you.


















Park power chair or scooter for pick-up along side (parallel to) the vehicle.

ALWAYS MAKE SURE THE VEHICLE'S PARKING BRAKE IS FIRMLY SET

Open the hatch of the vehicle Fold the seat-back of the chair or scooter.

KEEP HANDS AND FEET FROM UNDER THE CHAIR OR SCOOTER AS IT IS BEING LIFTED






Press the "OUT" button on the hand controller, to rotate the arm outside the vehicle.

(If using AL-215, *manually* rotate lifting arm outside of vehicle)






Once outside, press the "DOWN" button to lower the strap.


Attach the docking device supplied with your lift.
Attach the strap hook to the docking device.

INSPECT LIFTING STRAP PRIOR TO EACH USE

Press the "UP" button to begin raising your power chair or scooter.

ENSURE THE LIFTING STRAP IS SECURE AND THAT IT POINTS STRAIGHT DOWN. FAILURE TO DO SO COULD RESULT IN THE SCOOTER OR CHAIR SWINGING, CAUSING DAMAGE TO OPERATOR, CHAIR OR VEHICLE.





Once the chair has been raised enough to clear the bumper / threshold, rotate the chair or scooter into the vehicle by pressing the "IN" button.

(If using AL-215, *manually* rotate lifting arm inside of vehicle)

Once the chair or scooter is fully inside the cargo area, Press "Down" until the chair or scooter is firmly placed on the floor

Leave some tension on the strap during transport to eliminate movement of the chair, or docking device.

Unloading is the reverse of loading.



SAFETY:

CAUTION: DO NOT OPERATE THIS LIFT UNTIL YOUR DEALER

HAS SATISFACTORILY INSTRUCTED YOU IN THE PROPER OPERATION OF THE LIFT.

YOUR HARMAR INSIDE LIFT HAS BEEN ENGINEERED AND DESIGNED FOR YEARS OF TROUBLE FREE USE. ALTHOUGH, WITH EVERYDAY USE, SOME PARTS MAY BECOME LOOSE OR WORN. IMPORTANT! CHECK REGULARLY FOR ANY WORN, LOOSE, OR DAMAGED PARTS OF YOUR LIFT. IF ANYTHING IS OBSERVED, <u>DO NOT USE THE LIFT!</u> CONTACT YOUR DEALER OR INSTALLER OF THE LIFT FOR REPAIRS TO BE MADE. FAILURE TO ACT MAY CAUSE SEVERE INJURY

YOUR HARMAR INSIDE LIFT SHOULD ONLY BE USED FOR THE LOADING AND UNLOADING OF SCOOTERS AND POWER WHEELCHAIRS FOR WHICH IT IS DESIGNED. IF YOUR PARTICULAR APPLICATION (VEHICLE & SCOOTER/CHAIR) CHANGES (ex. buy a new car / upgrade to a new power chair) PLEASE CONSULT YOUR INSTALLER OR DEALER, AS THE ORIGINAL EQUIPMENT MAY NOT CONFIGURE THE SAME. DO NOT ADD TO OR MODIFY ANY PART OF THE LIFT SYSTEM WITHOUT FIRST CONSULTING THE MANUFACTURER OF THE LIFT, ANY MODIFICATIONS MAY VOID ANY WARRANTIES AS WELL AS THE STRUCTURAL INTEGRITY OF THE LIFT.

ALWAYS CHECK THE LIFTING STRAP BEFORE EACH USE FOR DAMAGE OR WEAR




IF THERE IS ANY SIGN OF TEARING / LOOSE THREADS / DAMAGE / WEAR **DO NOT USE THE LIFT.** SEVERE DAMAGE TO CHAIR / SCOOTER / VEHICLE AND PERSONAL INJURY MAY OCCUR. PLEASE CONTACT YOUR DEALER OR INSTALLER, AND HAVE THE STRAP REPLACED IMMEDIATELY.

CAUTION: WHEN ATTACHING THE STRAP HOOK TO YOUR CHAIR / SCOOTER, INSURE THE LIFTING STRAP IS SECURE AND TAUGHT AND THAT IT POINTS STRAIGHT DOWN BEFORE LIFTING, FAILURE TO DO SO COULD RESULT IN THE SCOOTER OR CHAIR SWINGING TOWARD THE OPERATOR OR THE VEHICLE.

CAUTION: WHEN USING THE LIFT, KEEP YOUR HANDS AND FEET FROM UNDER THE CHAIR / SCOOTER AS IT IS BEING LOADED OR UNLOADED.

CAUTION: STOP LIFTING BEFORE THE STRAP HOOK CONTACTS THE STRAP ROLLER MOUNTED ON THE LIFTING ARM. CONTACT WILL RESULT IN DAMAGE TO THE LIFT. A PROPERLY ADJUSTED AND OPERATED LIFT WILL AVOID THIS.

CAUTION: INSURE THE CHAIR / SCOOTER IS FIRMLY SITTING ON THE FLOOR OF THE VEHICLE AND NOT SUSPENDED BY THE LIFT DURING TRANSPORTATION. IF THE SCOOTER / CHAIR IS NOT FIRMLY ON THE FLOOR OF THE VEHICLE, IT MAY MOVE DURING TRANSPORTATIONCAUSING DAMAGE TO THE LIFT, TO THE VEHICLE, TO THE SCOOTER / CHAIR, OR TO THE PASSENGERS.



MAINTENANCE:

YOUR HARMAR INSIDE LIFT HAS BEEN DESIGNED AND ENGINEERED TO BE AS TROUBLE FREE AS POSSIBLE TO THE OWNER. BUT, AS WITH ANY MECHANICAL DEVICE, REGULAR CARE SHOULD BE GIVEN WHILE OWNING AN USING THIS DEVICE.

MAINTENANCE IS NO EXCEPTION.

WE RECOMMEND THAT DEALERS SCHEDULE A PREVENTATIVE MAINTENANCE INSPECTION AT LEAST ONCE A YEAR ON MOTORS, LIFT FRAME, WIRING HARNESSES, AND ALL MOVING PARTS OF THE LIFT.

CHECK FOR PAINT CHIPS AND TOUCH UP AND BARE METAL WITH A GOOD GLOSS BLACK ENAMEL OR LACQUER TO INHIBIT RUST. THIS MAY BE NECESSARY MORE FREQUENTLY WHEN SUBJECTED TO SALT AIR OR ROAD SALT.

ALWAYS CHECK The LIFTING STRAP BEFORE EACH USE. IF ANY WEAR OR DAMAGE IN NOTICED **DO NOT USE THE LIFT!** CONTACT THE DEALER OR INSTALLER FOR REPAIR, FAILURE TO DO SO WILL RESULT SEVERE INJURY OR DAMAGE.

HAND CONTROL UNITS SHOULD BE PROPERLY STORED INSIDE THE VEHICLE WHENEVER POSSIBLE. IF LIFT IS INSTALLED IN THE BACK OF A PICK-UP TRUCK, WE STRONGLY RECOMMEND UN-PLUGGING THE CONTROL IN BETWEEN USES, AND STORING IN THE CAB. EXPOSURE TO WEATHER IS OK, BUT RAIN FLOODING THE TRUCK BED, AND THE HAND CONTROLLER BEING SUBMERGED IN WATER, COULD CAUSE THE CONTROLLER TO FAIL

WARRANTY:

THIS LIFT COVERED BY A THREE YEAR FACTORY WARRANTY ON MATERIALS AND WORKMANSHIP. CONTACT YOUR LOCAL DEALER FOR MORE INFORMATION OR GO TO WWW.HARMARMOBILITY.COM

PROBLEM: LIFT WILL NOT OPERATE,


OPERATES SLOWLY, or OPERATES INTERMITTENTLY.

REASON: BAD ELECTRICAL CONNECTIONS / CIRCUIT BREAKER PROCEDURE:

CHECK / CLEAN <u>ALL</u> CONNECTIONS THAT MIGHT BE LOOSE OR DIRTY THIS IS THE #1 CAUSE OF A POORLY PERFORMING LIFT.
THE MOTOR ON THE LIFT REQUIRES <u>QUALITY</u> ELECTRICAL CONNECTIONS

TO OPERATE AT FULL CAPACITY. ANY BREAK IN THE WIRING WILL SLOW DOWN THE MOTOR, OR CAUSE IT TO OPERATE INTERMITTENTLY. WHICH IN TURN WILL PREMATURELY DETERIORATE THE MOTOR.

THE LIFT'S VEHICLE POWER HARNESS RUNNING THROUGH OR UNDER THE VEHICLE IS SUBJECT TO ROAD DEBRIS, AND RUBBING AGAINST THE VEHICLE'S FRAME, CAUSING A SHORT. INSPECT THE FULL LENGTH OF THE HARNESS TO BE SURE THE JACKET OF INSULATION IS NOT TORN.



PROBLEM: LIFTING STRAP OPERATES IN REVERSE
REASON: BUTTON WAS HELD TOO LONG AND UNWOUND LIFTING STRAP
LIKE A YO-YO, THE LIFTING STRAP WOUND IN REVERSE
PROCEDURE: PRESS WHICH EVER BUTTON EXTENDS STRAP, AND
CONTINUE TO HOLD UNTIL STRAP RETURNS TO PROPER DIRECTION


HARNESS.

PROBLEM: SCOOTER OR CHAIR TIPS FRONT OR BACK WHILE BEING LIFTED REASON: CHAIR OF SCOOTER IS BEING LIFTED OFF CENTER PROCEDURE: ADJUST FOR CENTER OF GRAVITY ON DOCKING DEVICE. CHECK TO SEE IF PROPER DOCKING DEVICE IS BEING USED


PROBLEM: LIFTING ARM IS HARD TO ROTATE
REASON: DRY / DIRTY BUSHINGS
PROCEDURE: CLEAN AND LUBRICATE GREY BUSHINGS IN BASE WITH
WHITE/LITHIUM GREASE OR EQUIV.

BE CAREFUL WHEN TROUBLESHOOTING WITH A TEST LIGHT, OR VOLTAGE METER. THEY MAY GIVE FALSE INDICATIONS

THIS SECTION IS RESERVED FOR THE INSTALLER / DEALER HARMAR DOES NOT RECOMMEND THESE PROCEDURES BE ATTEMPTED BY ANYONE WITHOUT PROPER KNOWLEDGE OF AUTOMOTIVE ELECTRICAL CIRCUITRY



THESE INSTRUMENTS WILL DETECT VOLTAGE, BUT MAY NOT INDICATE A TEAR OR POOR CONNECTION IN THE WIRING




BE CAREFUL WHEN TROUBLESHOOTING WITH A TEST LIGHT, OR VOLTAGE METER. THEY MAY GIVE FALSE INDICATION AS SHOWN ABOVE. THEY WILL ONLY INDICATE THERE IS <u>SOME</u> WIRE CONNECTION--NOT NECESSARILY--ENOUGH.

HARMAR STRONGLY RECOMMENDS TROUBLESHOOTING ELECTRICAL PROBLEMS WITH KNOW GOOD / FULLY CHARGED 12 V AUTOMOTIVE BATTERY OF KNOW GOOD / FULLY CHARGED POWER CHAIR/SCOOTER BATTERY









	ASSEMBLY PARTS					
ITEM	QTY	PART NO.	DESCRIPTION			
1	3	3-8_WASHER	3/8 FLAT WASHER ZP STEEL			
2	2	91102A760	3/8 SPRING LOCK WASHER			
3	2	92985A638	HHCS, 3/8-16 x 3-1/2"L ZP STEEL GRADE 5			
4	2	95462A031	HEX NUT, 3/8-16, ZP, GRADE 5			
5	2	95615A140	NUT, NYLOC, 3/8-16 ZP STEEL GRADE 5			
6	3	ALA28050	INSTALLTION WASHER			
7	1	ALA42501	3/8-16 X 10.00 ALL THREAD			
8	2	H205	LEG, EXTENSION			
9	4	H4205	CAP for 1.5sq x .120w Tube			
10	1	H425100	BASE HOUSING WELD			
11	4	HHCS	1/4-28 x 1/2"			

EXPLODED VIEW - BA-01




ITEM	OLA	PART NO.	DESCRIPTION
1	4	90272A828	PhPnHd, 10-32 x 7/16"L, ZP STEEL
2	1	90840A131	HEX NYLOC NUT, 3/9-16, HALF HEIGHT
3	1	91102A760	3/6 SPRING LOCK WASHER
4	1	91247A635	HHCS, 3/8-16 x 2-3/4"L ZP STEEL
5	1	92885A821	HHCS 3/8-16 X 0.625
6	2	92865A626	HHC8 3/8-16 X 1.25
7	1	92665A638	HHCS, 3/8-16 \times 3-1/2°L FULL THREAD, G5 ZP
8	3	95615A140	NUT, NYLOC, 3/8-16, ZP STEEL GRADE 5
9	2	98025A031	WASHER 3/8-ID
10	1	AL206	ARTICULATING LEG ASSY
11	3	ALA28050	INSTALLTION WASHER
12	ĭ	ALA@0133	J-HOOK
13	1	H209	BA01r1 LEG EXTENSION
14	2	H4205	CAP for 1.5sq x .120w Tube
15	1	H43000	BASE HOUSING WELD
16	2	H43001	VERTIGAL JOINER FOOT
17	1	H43002	VERTIGAL JOINER
18	1	H43008	BAMM - VERTICAL JOINER - 2
19	2	H4311111	ANCHOR, BLIND, Expandable, 5/18" LHBW08 #1EG
20	4	HHCS	1/4-28 x 1/2*

EXPLODED VIEW - POWER BASE




		ASSEMBLY PART	S	ASSEMBLY PARTS			
ПЕМ	QTY	PART NO.	DESCRIPTION	ПЕМ	QTY	PART NO.	DESCRIPTION
1	6	90272A148	MACHINE SCREW # 6-32 X 0.50	17	1	H425200	MOTOR HOUSING WELD
2	4	90631A007	NYLOCK NUT#6-32	18	1	H425210	MOTOR HOUSING COVER WELD
3	3	91102A750	1/4" SPLIT LOCK WASHER	19	1	H425230	OUTER POST TUBE
4	4	91255A244	8C8 10-24 X 5/8	20	1	H425231	MOTOR HOUSING GEAR PIN - A
5	3	91255A553	BHCS 1/4-29 X 0.50	21	1	H425232	MOTOR HOUSING GEAR PIN - B
6	2	91772A086	MS #2-58 X 1.00	22	1	H425240	DOUBLE SPROCKET 19-35 , 8-35
7	2	92311A194	SET SCREW 8-32 X 0.50	23	1	H425241	DOUBLE SPROCKET 19-35 , 8-40
8	4	92949A537	BHCS 1/4-20 X 0.50	24	1	H425248	MOTOR HOUSING-GEAR PLATE
9	4	94639A331	SPACER: ID=#6, C/ID=313	25	2	H425250	LIMIT RING
10	1	98407A134	E-CLIP 3/8 SHAFT OID	26	1	H425253REV1	AL425 PCB PLATE
11	1	98407A138	E-CLIP 1/2 SHAFT OD	27	1	H425255	LIMIT/MOTOR HOUSING SPACER
12	1	ALA25-PCB	PRINTED CIRCUIT BOARD	28	1	H425262	BUSHING 1/2ID,5/80D,3/4L
13	1	ALA20034	FLANGE BUSHING	29	1	H425263	BUSHING 3/8/ID,6/8/OD,1/2L
14	1	ALA20035	CUP BUSHING	30	1	H425264	BUSHING 1/2ID,840D,1/2L
15	1	ALA21410	FLANGE BUSHING BRZ 3/6X1/ZX3/6	31	2	OMRON_SS-5GL2	LIMIT SWITCH
16	1	ALA41011	COL MOTOR	32	1	RHSS	1/4-20 x .26" NYLON Flat Head

EXPLODED VIEW - 425 ARM




			ASSEMBLY PARTS	ASSEMBLY PARTS			
ITEM	QTY	PART NO.	DESCRIPTION	ITEM	QTY	PART NO.	DESCRIPTION
1	1	_BHCS	1/4-20 x 1.75 Button Hand Cap	14	2	ALA21410	FLANGE BUSHING BRZ 3/8X1/2X3/8
2	1	_BHC\$	10-32 x 3/8" Button Head	15	1	ALA67898	LIFTING STRAP
3	8	80067A180	8-32 X 0.25 THREAD CUTTING SCREW	16	1	H21072	DETACHABLE HOOK
4	3	81255A553	BHCS 1/4-28 X 0.50	17	1	H425400	HOISTWELD
5	1	92311A558	SET SCREW: 1/4-28 X 0.26	18	1	H425403	HOIST IPMER ARM
6	1	92785A463	CLEVIS PIN & CLIP: 3/8 X 2.26	19	1	H425410	ALAZS COVER-BOTTOM
7	1	92785A486	CLEVIS PIN & CLIP: 3/6 X 2.50	20	1	H425411	ALA25 COVER-BACK
8	2	92949A710	BHCS 1/2-13 X 0.75	21	1	H425420	TURNBUCKLE ASSEMBLY
9	1	ACORN_NUT	1/4-20 Acom nut	22	1	H425430-REV1	HOIST MOTOR COVER
10	1	ALA21013	LIFTING MOTOR	23	2	H425444	ALA25 STRAP ROLLER
11	1	ALA21020	SPOOL CHAIN, 50 LINKS	24	1	H425445	PIN, ROLLER, STRAP, AL425
12	1	ALA21055	STRAP TAB	25	1	TU1-037-049-1_090	ASTM AS13 TUBBNG, 3/8" dla x .049wall - 1.08"L fin BLACK PC
13	1	ALA214	SPOOL	26	1	WASHER	14" SPLIT LOCK WASHER

EXPLODED VIEW - AL225/215 ARM




	ASSEMBLY PARTS			ASSEMBLY PARTS			
ITEM	QTY	PART NO.	DESCRIPTION	ITEM	QTY	PART NO.	DESCRIPTION
1	1	_BHCS	1/4-20 x 1.75 Button Head Cap	10	1	ALA214	SPOOL
2	6	90067A190	8-32 X 0.25 THREAD CUTTING SCREW	11	2	ALA21410	FLANGE BUSHING BRZ 3/8X1/2X3/8
3	1	91265A263	BHCS 10-32 x 3/8" BLACK OXIDE	12	1	ALA67997	LIFTING STRAP
4	3	912554553	BH/CS 1/4-28 X 0.50	13	1	H21072	DETACHABLE HOOK
5	3	92735A466	CLEVIS PIN & CLIP: 3/8 X 2.50	14	1	H215700	AL215/225/ & AG425 COVER-BOTTOM
6	1	ACORN_NUT	1/4-20 Acom nut	15	1	H225400	AL215/225/AG425 ARM WELDMENT
7	1	ALA05111	ROTATION MOTOR	16	2	H225402	AL215 & AL225 ROLLER
8	1	ALA21020	SPOOL CHAIN, 50 LINKS	17	1	TU1-037-049-1_090	TUBING, 3/8" dia x .049wali - 1.09"L fin BLACK PC
9	1	ALA21055	STRAP TAB	18	1	WASHER	1/4" SPLIT LOCK WASHER


EXPLODED VIEW - AL425




ASSEMBLY PARTS					
TEM	QTY	PART NO.	DESCRIPTION		
1	2	905664033	NYLOCK NUT 1/2-13		
2	1	9064QA131	NYLOCK NUT 3/8-16		
3	4	91102A790	3/8 SPRING LOCK WASHER		
4	1	912364636	HHCS 3/8-18 X 3.50		
8	4	91247A836	HHCS 3/8-16 X 2.76		
6	1	91247A721	HHCS 1/2-13 X 2.25		
7	4	91375A619	SET SCREW 3/8-16 X 0.50		
8	1	92949A724	BHSC 1/2-13 X 3		
9	2	98017A210	WASHER 0.50 ID		
10	2	ALA28060	VELCRO, 1" PC from 25 YD ROLL		
11	1	H425002	MOTOR HOUSING ASSEMBLY		
12	1	H425003	COLLAR & POST ASSEMBLY		
13	1	H425004	HOIST ASSEMBLY		
14	1	H425009	ROUND TUBE PLUG		
15	1	H425435	BASE COVER, ALA25		
18	1	POWER_SWITCH	POWER SWITCH		




EXPLODED VIEW - AL225







		ASSEMBLY P	ARTS
TEM	QTY	PART NO.	DESCRIPTION
1	1	9064QA131	NYLOCK HUT 3/8-16
2	2	90640A133	NYLOCK NUT 1/2-13
3	4	91102A760	3/8 SPRING LOCK WASHER
4	1	912364638	HHCS 3/6-16 X 3.50
5	4	91247A835	HHCS 3/8-16 X 2.75
6	4	91375A&19	SET SCREW 3/8-16 X 0.50
7	2	92948A724	BH&C 1/2-13 X 3
8	2	ALA28060	VELCRO, 1" PC from 25 YD ROLL
9	1	H215810C	ANGLED SQUARE POST
10	1	H225004	AL225 ARM ASSEMBLY
11	1	H425002	MOTOR HOUSING ASSEMBLY
12	1	H425435	BASE COVER, AL425
13	1	POWER_SWITCH	POWER SWITCH


EXPLODED VIEW - AL215





		ASSEMBL	Y PARTS
TEM	QTY	PART NO.	DESCRIPTION
1	1	9064QA131	NYLOCK NUT 3/8-16
2	2	90640A133	NYLOCK NUT 1/2-13
3	4	91102A760	3/8 SPRING LOCK WASHER
4	1	91239A838	HHCS 3/8-16 X 3.50
5	2	913754535	SET SCREW 1/4-20 X .375
6	4	91375 <i>85</i> 75	SET SCREW 5/16-18 X 0.3125
7	4	928654621	HHCS 3/8-16 X 0.625
8	2	92649A724	BHSC 1/2-13 X 3
g	1	ALA20034	FLANGE BUSHING
10	1	ALA20036	CUP BUSHING
11	1	H215210C	BASE HOUSING COVER WELDMENT
12	1	H215215C	AL215 OUTER TUBE WELDMENT
13	1	H215810C	AMGLED SQUARE POST
14	1	H225004	AL225 ARM ASSEMBLY
15	1	RHSS	1/4-20 x .25" NYLON Flat Head





THANK YOU FOR MAKING HARMAR AMERICA'S LEADER IN LIFTS AND RAMPS