## RAMP MODEL SELECTION GUIDE

- 1. Always have a qualified assistant present when using any portable ramp.
- Never exceed a slope greater than 2" on 12" with an occupied chair or scooter. Never exceed a slope greater than 3" on 12" with an unoccupied chair or scooter.
- Always make sure that top of ramp is secured on step or landing before
  using. It may be necessary to anchor top of ramp to landing surface
  with steel security pins provided.
- 4. **Make certain** there is adequate head clearance prior to loading an **occupied** chair or scooter into a vehicle.
- 5. Always follow manufacturer recommendations for chair or scooter.
- 6. Always use your lap belt.
- 7. Slope Recommendations:

## 1:12 SLOPE

- 1. The ADA recommended slope for long (up to 30') home and commercial ramps.
- 2. Works for most strong, unassisted manual chair users, although PVI always recommends assistance.

#### 2:12 **SLOPE**

1. Maximum acceptable grade for portable ramp use by occupied chairs and scooters with a qualified assistant.

## 3:12 SLOPE

1. For loading **unoccupied** chairs and scooters.

#### **ANY SLOPE GREATER**

- 1. PVI does **NOT** recommend.
- 8. Watch safety DVD prior to use.

To determine the length of ramp for your application, first you will need to determine the rise. The rise is the vertical measurement between the ground and where the top of the ramp is going to sit.

Prairie View Industries recommends NOT to exceed a 2:12 slope on most applications. Some ramps may need to meet ADA (1:12 slope) requirements that are specified by your state. If your rise is over 24", please contact your local dealer about our modular ramp system. To find a 2:12 slope, take total amount of rise and divide by 2. If you have a 12" rise, divide by 2 and the resulting number would be the length of ramp which is required in feet or you would need a 6' ramp. A 1:12 slope would require 1 foot of ramp for every inch of rise.



**Note:** Consideration must be taken if the ground is not level and slopes away from the rise where ramp will be placed. Please contact your dealer for evaluation on these applications.


# QUICK REFERENCE SLOPE GUIDE CHART

Easy as 1-2.3

**Step #1** - Measure the rise of your application.

Step #2 - Determine your maximum slope.\*

**Step #3** - Pick the proper ramp length, based on your selections, in the graph below.



\*Refer to Slope Recommendations on page to left or visit our website www.pviramps.com/products/calculator to use our ramp length calculator.

