



HELPFUL HINT:

For ease of use, lift the push handle upwards during the entire recline process.

RECLINE POSITION

Step: 1

Rotate Push Handle downward:

Push-in the **Grey Push Buttons** located on the side of the **Push Handle**, at the same time, rotate the **Push Handle** downward. (fig. 1 & 2)



Fig. 1
Grey Push Button



Fig. 2
Rotate Push Handle Downward

Step: 2

Recline the seat back:

Lift **Push Handle** upward (fig. 3) to relieve the occupant's weight. At the same time, pull upwards on the **Recline Pull Loop** (fig. 4), and recline the seat back to desired position. (fig. 5)



Fig. 3
Lift Push Handle Upwards



Fig. 4
Pull Upwards on Recline Pull Loop



Fig. 5
Recline the Seat Back

UPRIGHT POSITION

To return seat back to upright position repeat steps 1 & 2. Push seat back up to bring upright, instead of reclining down.