
INSTRUCTIONS FOR USE

Sara 3000 / Lite / Steady / Flex Active Slings

Including Arjo Wipeable Standing Clip Sling

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使用说明书 · Bedienungsanleitung · Mode d'emploi · Istruzioni per l'uso · Gebruiksaanwijzing

Contents

Foreword	2	Detach the Clips (3 steps)	9
Customer Support	2	Applying the Sling	10
Definitions in this IFU	2	In Chair/Wheelchair (17 steps)	10
Intended Use	3	Removing the Sling	11
Patient/Resident Assessment	3	In Chair/Wheelchair (11 steps)	11
Expected Service Life	3	Cleaning and Disinfection	12
Safety Instructions	4	Flites - Do NOT wash	12
Safety Practices	4	Wipeable Sling - Disinfecting Agents	12
Serious Incident	4	Machine Washable Slings - Cleaning Instruction	
EN Parts Designation	5	(3 steps)	13
Preparations	6	Machine Washable Slings - Cleaning Agents	13
Before First Use (6 steps)	6	Care and Preventive Maintenance	13
Before Every Use (8 steps)	6	Before and After Every Use	13
After every use	6	When Soiled, Stained and Between Patients	13
Select Sling Size	7	Storage	13
Measure	7	Service and Maintenance	13
Sling Selection	8	Troubleshooting	14
Transfer/Transport	8	Technical Specifications	15
Allowed Combinations	9	Parts and Accessories	15
Safe Working Load (SWL)	9	Label on the Sling	16
Clip Attachment & Detachment	9		
Attach the Clips (6 steps)	9		

WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product. Mandatory to read the Instructions for Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2020.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

The content of this publication may not be copied either whole or in part without the consent of Arjo.

Foreword

Thank you for purchasing Arjo equipment. Please read this *Instructions For Use (IFU)* thoroughly.

Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

Customer Support

If you require further information, please contact your local Arjo representative. The contact information is located at the end of this *IFU*.

Definitions in this IFU

WARNING

Means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

CAUTION

Means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

NOTE

Means: This is important information for the correct use of this system or equipment.

Intended Use

The Active Sling is a product intended for assisted transfer and rehabilitation of patients/residents with limited ability to move.

The Arjo Wipeable Standing Clip Sling is intended to be used together with a lifting device to provide body support during raising to a standing position and short transfers of adult patients/residents. It is intended to be cleaned by means of wiping.

The Active Flites Sling has been designed to support hospital or care facility patients/residents while using a patient/resident transfer lifter, under the supervision of trained nursing staff. Active Flites Slings are to be used for a limited period only, and, by nature of their design, must be treated as a disposable and patient/resident specific product.

The Active Sling, the Arjo Wipeable Standing Clip Sling, and Active Flites Sling should be used together with Arjo lift devices in accordance with the Allowed Combinations specified in the Instructions for Use (IFU).

The Active Sling, the Arjo Wipeable Standing Clip Sling, and Active Flites Slings shall only be used by appropriately trained caregivers with adequate knowledge of the care environment, and in accordance with the instructions outlined in the Instructions for Use.

The Active Sling and the Arjo Wipeable Standing Clip Sling is intended to be used in hospital environments, nursing homes, other healthcare facilities and home care.

The Active Flites Sling is intended to be used in hospital environments, nursing homes and other healthcare facilities.

The Active Flites Slings must not be used to support a patient/resident while bathing or showering. The Active Flites Slings are disposable products and cannot be washed.

The Active Sling, the Arjo Wipeable Standing Clip Sling, and Active Flites Sling should only be used for the purpose specified in this Instructions for Use. Any other use is prohibited.

Patient/Resident Assessment

We recommend that facilities establish regular assessment routines.

For Active Sling and Active Flites Sling caregivers should assess each patient/resident according to the following criteria prior to use:

- Patient/resident sits in wheelchair
- Is able to partially bear weight on at least one leg
- Has some trunk stability
- Dependent on caregiver in most situations
- Physically demanding for caregiver
- Stimulation of remaining abilities is very important

Safe Working Load (SWL) for Active Sling:

- MAA3070M, MAA3071M, MAA3061M, MAA3060M: 180 kg (397 lbs);
- TSS.500-504, TSS.500-503S-1C, TSS.501-502SV, TSS.511: 200 kg (440 lbs);
- MFA3000 (Sara 3000 Flites): 200 kg (440 lbs);

The right type and size of slings should be used after proper assessment of each patient/resident's size, condition and the type of lifting situation.

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

Prior to the transfer with the Arjo Wipeable Standing Clip Sling caregivers should assess each patient/resident according to the following criteria:

- Size
- Mobility Level
- Medical conditions
- Ability to understand and respond to instructions

The Arjo Wipeable Standing Clip Sling is intended for patients/residents:

- who are able to bear weight on at least one leg and have some trunk stability.
- whose weight is no more than 200 kg (440 lb).
- has a length in the range of 140 cm and 200 cm (55 1/8" and 78 3/4").

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

Expected Service Life

The expected service time of Active Sling and Active Flites Sling are the maximum period of useful life.

The expected life of the sling is dependent on the actual use conditions. Therefore, before use, always make sure that the sling does not show any signs of fraying, tearing or other damage and that there is no damage (i.e. cracking, bending, breaking). If any such damage is observed, do not use the sling. If you have any doubts about sling safety, as a precaution and to ensure safety, do not use the sling.

The expected service life of the Active Slings:

Service Time: 2 years;

Shelf Time: 5 years.

The expected service life of the Active Flites Slings:

Service Time: 2 weeks (based on 4 transfers per day);

Shelf Time: 5 years.

The expected service life of Arjo Wipeable Standing Clip Sling is the maximum period of useful life. The expected service life of this sling is two (2) years, subject to preventive maintenance being carried out in accordance with the instructions for care and maintenance found in the Instructions for Use. The shelf life, i.e., the period before the first use of the device is 5 years.

EN

Safety Instructions

WARNING

To avoid injury, make sure that the patient is not left unattended at any time.

WARNING

To avoid falling, make sure that the user weight is lower than the safe working load for all products or accessories being used.

WARNING

To avoid injury, always assess the patient prior to use.

WARNING

To avoid injury, never wash a Flites. A Flites is only intended for “single patient use”.

WARNING

To avoid injury, never let the patient smoke when using the sling. The sling can catch fire.

WARNING

To avoid injury, only store the equipment for a short period of time. If stored longer than stated in the *IFU*, it can cause weakness and breakage in the material.

WARNING

To avoid injury, keep the equipment away from sun/UV-light. Exposure to sun/UV-light can weaken the material.

Safety Practices



At any time, if the patient becomes agitated, stop transferring/transporting and safely lower the patient.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

EN

Parts Designation

EN

1. Sling label/print (located on the outside)
2. Waist part
3. Inner attachment clip
4. Outer attachment clip
5. Attachment strap
6. Support belt (buckle)
7. Buckle

8. Leg flap
9. Leg strap (loop)
10. Waist belt
11. Support belt (hook and loop attachment)
12. Flites identification label (located on the outside)

Example of slings. Not all models are represented.

Preparations

Before First Use (6 steps)

WARNING

To avoid injury, always make sure to inspect the equipment prior to use. If any part is missing or damaged - Do NOT use the sling.

EN

1. Check all parts of the sling, see section “Parts Designation” on page 5. If any part is missing or damaged - Do NOT use the sling.
2. Read this *IFU* thoroughly.
3. Check that the sling is clean.
4. Choose a designated area where the *IFU* should be kept and easily accessible at all times.
5. Make sure to have a rescue plan ready in case of a patient emergency.
6. For questions, contact your local Arjo representative for support.

NOTE

To keep track of the service life of the wipeable sling, document the serial number and the date it was first used.

NOTE

Arjo recommends the name of the patient is written on the sling care label to avoid cross contamination between patients.

Before Every Use (8 steps)

1. Patient assessment must always determine the methods used in everyday care, see “Intended Use” on page 3.

WARNING

To avoid the patient from falling, make sure to select the correct sling size according to the *IFU*.

2. Check the allowed combinations against the sling and lift, see “Allowed Combinations” on page 9.

WARNING

To avoid injury, always make sure to inspect the equipment prior to use. If any part is missing or damaged - Do NOT use the sling.

3. Check all parts of the sling, see “Parts

Designation” on page 5. If any part is missing or damaged - Do NOT use the sling.

4. Check for:
 - fraying
 - loose stitching/delamination
 - cracks
 - tears
 - fabric holes
 - soiled fabric
 - damaged clips/loops/buckles
 - unreadable or damaged label/print

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this *IFU*.

5. Check that the sling is clean. If the sling is not clean see “Cleaning and Disinfection” on page 12.
6. Always check the washing symbol on the *Flites* label. If the *Flites* has been washed, throw the *Flites* away. See “*Flites* - Do NOT wash” on page 12.
7.  If this symbol appears, do NOT use the *Flites*. The *Flites* label indicates it has been washed.
8. Before using the sling, read the lift *IFU* for transfer, transport and bathing.

After every use

For Flites

If the patient is leaving the hospital, throw away the patient specific *Flites*.

For Arjo Wipeable Standing Clip Sling

1. Clean and disinfect the sling according to the instructions in this *IFU*.
2. Check all parts of the sling, see “Parts Designation” on page 5. If any part is missing or damaged, dispose the sling.

Select Sling Size

Measure

The patient's physical disabilities, weight distribution and general physique needs to be taken into consideration when selecting a sling.

For Active Slings and Active Flites Sling (3 steps)

The sling size depends on the patient's body shape, what kind of lift is being used, and on the attachment of the sling. For best result try out the sling on the patient.

Each sling size has its own colour, which is found on the binding of the sling. (See Fig. 1).

1. Place the sling around the patient's waist.
2. Close the support belt and make sure it overlaps with enough margin at approximately 10 cm (4 inches).
3. Make sure the sling attachments reach the lift. If the attachment doesn't reach, change to another sling size.

For Arjo Wipeable Standing Clip Sling (2 steps)

The patient's physical disabilities, weight distribution and general physique need to be taken into consideration when selecting a sling.

1. Measure the patient's waist in seated position.
2. Select a sling size that is closest to that measurement (See table below). In case of doubt choose the smaller one. Each sling size has its own colour which is found on the print on the outside of the sling.


To double-check your size selection, align the flap fold with patient's side (See Fig. 2). If the sling is too small, this flap fold will be too far back (See Fig. 3). If the sling is too large, the flap fold will be further in the front (See Fig. 4).

EN

Fig. 1

Article number	Product description	Waist size in cm (in)			
		S	M	L	XL
TSS.50XW	Arjo Wipeable Standing Clip Sling	64-84 cm (25" - 33")	76-104 cm (30" - 41")	88-124 cm (34½" - 49")	108-160 cm (42½" - 63")

Sling Selection

Transfer/Transport

Sling Type	Article number	Safe Working Load	Product Description	Size
Flites	MFA3000	200 kg (440 lbs)	Flites (Disposable and Patient Specific)	One size
Wipeable	TSS.500-503W	200 kg (440 lbs)	Wipeable sling Clip attachment Fixed support belt with buckle	S, M, L, XL
Machine Washable	TSS.500-503S-1C	200 kg (440 lbs)	Sling split hook and loop attachment 1 clip	S, M, L, XL
	TSS.500-504	200 kg (440 lbs)	Fixed support belt	S, M, L, XL, XXL
	TSS.501-502SV	200 kg (440 lbs)	Sling split hook and loop attachment	M, L
	TSS.501-503A	200 kg (440 lbs)	Sling split hook and loop attachment 1 clip	M, L, XL
	TSS.511	200 kg (440 lbs)	Transfer	One size
	MAA3060M	180 kg (397 lbs)	2 clip/loop transfer & buckles (roll)	One size
	MAA3061M	180 kg (397 lbs)	Transfer	One size
	MAA3070M	180 kg (397 lbs)	3P clip standing and raising aid (roll)	One size
MAA3071M	180 kg (397 lbs)	3P clip standing and raising aid (flat)	One size	

Article numbers with "A-suffix" (e.g. MLA7000A) has a BAA (Buy American Act) compliant country of origin.

Allowed Combinations

WARNING

To avoid injury, always follow the allowed combinations listed in this *IFU*. No other combinations are allowed.

Safe Working Load (SWL)

WARNING

Always follow the lowest SWL in the complete system (lift or sling).

- Only use the earlier listed slings with the *Sara 3000/Lite/Stedy/Flex* standing and raising aid.
- The SWL of the *Sara 3000*, *Sara Flex* is 200 kg (440 lbs), the SWL of the *Sara Lite* is 175 kg (385 lbs) and the SWL of the *Sara Stedy* is 182 kg (400 lbs).
- Always follow the lowest SWL in the complete system. E.g. the *Sara 3000* lift has a SWL of 200 kg (440 lbs) and the MAA3060 has a SWL of 180 kg (397 lbs). This means that the MAA3060 Flites has the lowest SWL. The patient is not allowed to weigh more than the lowest SWL.

	Standing/ Raising Aid	Sara 3000, Sara Flex	Sara Lite	Sara Stedy
	SWL	200 kg (440 lbs)	175 kg (385 lbs)	182 kg (400 lbs)
Active Sling	SWL	Size	Size	Size
MFA3000	200 kg (440 lbs)	One size	One size	One size
TSS.500-503W	200 kg (440 lbs)	S, M, L, XL	–	S, M, L, XL
TSS.500-503S-1C	200 kg (440 lbs)	S, M, L, XL	S, M, L, XL	S, M, L, XL
TSS.500-504	200 kg (440 lbs)	S, M, L, XL, XXL	S, M, L, XL, XXL	S, M, L, XL, XXL
TSS.501-502SV	200 kg (440 lbs)	M,L	M,L	M,L
TSS.501-503A	200 kg (440 lbs)	M, L, XL	M, L, XL	M, L, XL
TSS.511	200 kg (440 lbs)	One size	–	–
MAA3060M	180 kg (396 lbs)	One size	–	–
MAA3061M	180 kg (396 lbs)	One size	–	–
MAA3070M	180 kg (396 lbs)	One size	One size	One size
MAA3071M	180 kg (396 lbs)	One size	One size	One size

EN

Clip Attachment & Detachment


Attach the Clips (6 steps)

1. Select the inner or outer attachment clips.
2. Place the clip on the spreader bar lug. (See Fig. 5)
3. Pull the strap down.
4. Make sure the lug is locked at the top end of the clip. (See Fig. 6)
5. Make sure the strap is not squeezed in between the clip and the spreader bar.
6. Make sure the straps are not twisted.

Detach the Clips (3 steps)

Make sure that the weight of the patient is taken up by the receiving surface before removing the clip.

1. Pull the strap up. (See Fig. 7)
2. Make sure the lug is unlocked at the bottom of the clip.
3. Remove the clip.

Applying the Sling

In Chair/Wheelchair (17 steps)

For transfer/transport see the *Sara 3000/Lite/Stedy/Flex lift IFU*.

1. Place *Sara 3000/Lite/Stedy/Flex* nearby.
2. Apply the brakes on the wheelchair.
3. In case of using *Sara Stedy* lift, make sure you:
 - Fold up the seat pads and let the patient step up onto the *Sara Stedy*.
 - Fold down the seat pads and let the patient sit down.
4. Place the sling around the patient's lower back. If needed, slightly lean the patient forward. **(See Fig. 8 – Stedy) or (See Fig. 9 – Lite/3000/Flex)**
5. For machine washable slings and Flites: Position the bottom of the sling, aligned horizontally, about five centimetres (two inches) above the patient's waistline. For wipeable sling: Position the wipeable sling centred around the patient's waistline and make sure the soft lilac surface is in contact with the patient (label on the outside).
6. Make sure that:
 - patient's arms are outside the sling,
 - support belt is positioned loosely around the patient's body and
 - support belt is not twisted or trapped behind the patient's back.
7. In case of using *Sara 3000, Sara Flex* lift, make sure you place the leg flaps around the legs and these are not twisted. **(See Fig. 10)**
8. Fasten the support belt securely by pressing the buckles or hook and loop attachments together. The support belt should be tight but comfortable for the patient.

EN

9. Position the *Sara 3000/ Sara Lite/Sara Flex* lift in front of the patient. Open up the lift legs, if needed. See appropriate lift *IFU*.
10. Place the patient's feet on the lift foot support. In case of using *Sara 3000/Sara Lite/Sara Flex* lifts, make sure the patient's legs are in contact with the knee support.
11. Apply the brakes on the lift.
12. Attach the sling to the lift. See *"Clip Attachment & Detachment"* on page 9 and **(See Fig. 11)**
If the attachment strap is too loose, choose the inner attachment clip. If the attachment strap is too tight, choose the outer attachment clip.
13. **Transfer sling only:** Attach the leg flaps loops to the central lug on the *Sara 3000, Sara Flex* lift. See appropriate lift *IFU*. **(See Fig. 12)**

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

14. Make sure that all clips/loops are attached securely and the patient is comfortable in the sling. **(See Fig. 6)**
15. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clip/loop.
16. Encourage the patient to assist during the raising exercise by holding on to the support grips with both hands.
17. Continue raising until the patient is in a comfortable standing position. Make sure the patient leans back against the sling at all times. See corresponding lift *IFU*.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

Fig. 8 – Stedy

Label on the outside

Fig. 9 – Lite/3000/Flex

Label on the outside

Fig. 10

Fig. 11



Fig. 12

Removing the Sling

In Chair/Wheelchair (11 steps)

1. Apply the brakes on the wheelchair.
2. Position the patient with the his/her back against the chair/wheelchair. Open up the lift legs on *Sara 3000/Lite/Stedy/Flex*, if needed. (See Fig. 13)

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

3. Lower the patient into the chair/wheelchair. See *Sara 3000/Lite/Stedy/Flex lift IFU*.
4. Make sure that the patient's lower back is positioned all the way in against the chair/wheelchair (A) and that the weight of the patient is taken up by the receiving surface before removing the clips/loops. (See Fig. 14)
5. Apply the brakes on the *Sara 3000/Lite/Stedy/Flex* lift.

6. Detach the sling clips/loops from the *Sara 3000/Lite/Stedy/Flex*. (See Fig. 7)
7. **Transfer Sling only:**
 - Detach the leg straps from the *Sara 3000/Lite/Stedy/Flex*.
 - Open and pull them out from under the patient's legs.
 - Tuck the leg straps backwards along the patient's side.
8. Remove the patient's feet from the foot support.
9. Move the *Sara 3000/Lite/Stedy/Flex* away from the patient.
10. Open the sling support belt.
11. Lean the patient forward and remove the sling from behind the patient. (See Fig. 15)

EN

Fig. 13

Fig. 14


Fig. 15

Cleaning and Disinfection

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this *IFU*.

WARNING

To avoid material damage and injury, clean and disinfect according to this *IFU*:

- No other chemicals are allowed.
- Never clean with chlorine (except for the wipeable sling).
- Chlorine will deteriorate the surface of the material (except for the wipeable sling).

WARNING

To avoid injury, always remove the sling before disinfecting the lift.

WARNING

To avoid injury, never wash a *Flites*. A *Flites* is only intended for “single patient use”.

All *Sara 3000/Lite/Stedy/Flex* slings, except the *Flites*, should be cleaned when they are soiled or stained and between patients.

Flites - Do NOT wash

The *Flites* is for single patient use.

Do NOT clean, wash, disinfect, wipe or sterilize the *Flites*. If the *Flites* has been subjected to any such treatment, it shall be discarded.


The *Flites* is marked with the “Do NOT wash” symbol. (See Fig. 16)

If the *Flites* is washed, the “Do NOT use/lift” symbol will appear. (See Fig. 17)

Fig. 16

Fig. 17

Wipeable Sling - Cleaning and Disinfection Instructions (6 steps)

1. Detach the sling from the lift.
2. If needed, rinse to remove any deposits or visible dirt from the sling.
3. Wipe all parts of the sling with one of the disinfecting agents listed under “Wipeable Sling - Disinfecting Agents”.
4. Allow contact time according to the disinfecting agent manufacturer’s instructions for use.
5. If stated in the disinfecting agent manufacturer’s instructions for use, rinse to remove chemical residue.
6. Air dry at room temperature.

Do NOT

- machine wash
- use gas sterilization
- use any mechanical pressure
- tumble dry
- use autoclave
- dry clean
- steam
- iron

Wipeable Sling - Disinfecting Agents

The sling can be disinfected using:

- Ethanol based solution, up to 70%
- Isopropanol based solution, up to 70%
- Hydrogen peroxide based solution, up to 1.5%
- Quaternary ammonium based solution, up to 0.1%
- Bleach (sodium hypochlorite) based solution, up to 10 000 ppm

NOTE

Wipeable sling may show some signs of discoloration and become slightly sticky after repeated contact with bleach (sodium hypochlorite) at 10 000 ppm.

Machine Washable Slings - Cleaning Instructions (3 steps)

1. Detach the sling from the lift.
2. Close all buckles and hook and loop attachments on the slings and belts prior to washing.
3. Machine wash the sling at 70°C (158°F). In order to obtain sufficient level of disinfection, please wash according to local hygiene legislation. If allowed according to the product label, tumble dry in low temperature, max. 60°C (140°F).

Do NOT

- wash with other items that have rough surfaces or sharp objects
- use any mechanical pressure, pressing or rolling during the washing and drying process
- use bleach
- use gas sterilization
- use autoclave
- dry clean
- steam
- iron

Machine Washable Slings - Cleaning Agents

For all washable slings (NOT the *Flites* or the wipeable sling) use common commercial detergent without optical brightener.

No other chemicals are allowed e.g. chlorine, softener, iodine based disinfectants, bromine and ozone.

Care and Preventive Maintenance

WARNING

To avoid injury to both patient and caregiver, never modify the equipment or use incompatible parts.

Before and After Every Use

Visually check all exposed parts

The caregiver shall inspect the sling before and after every use. The complete sling should be checked for all deviations listed below. If any of these deviations are visible, replace the sling immediately.

- Fraying
- Loose stitching/delamination
- Cracks
- Tears
- Holes
- Discolouration or stains from bleaching
- Sling soiled or stained
- Unreadable or damaged label/print
- For *Flites*: Make sure the “Do NOT wash” symbol on the head support is intact. If damaged or missing and the underlying red “Do NOT use/lift” symbol appears, replace the sling immediately. See “*Label on the Sling*” on page 16.

When Soiled, Stained and Between Patients

Clean/Disinfect, All slings except the *Flites*

The caregiver should make sure that the sling is cleaned according to “*Cleaning and Disinfection*” on page 12, when it is soiled or stained and between patients.

Flites

Washing or disinfection is not allowed. Do not wipe, sterilize or soil the *Flites*. If the *Flites* has been subjected to any treatment, it shall be discarded.

Storage

When not in use, the slings should be stored away from direct sunlight where they are not subject to unnecessary strain, stress or pressure, or to excessive heat or humidity. The slings should be kept away from sharp edges, corrosives or other things that could cause damage on the sling.

Service and Maintenance

It is recommended that qualified personnel inspect the condition of the sling twice a year (every 6 month), according to ISO 10535.

Troubleshooting

EN

Problem	Action
The patient is not positioned correctly in the sling.	<ul style="list-style-type: none"> • Make sure the clips/loops are attached at the right length. • Make sure the straps are not twisted. • Make sure the patient is positioned on the inside of the sling. The outside has a label or print for reference. • Make sure the patient is centered in the sling.
The patient feels discomfort in the leg area when seated in the sling.	<ul style="list-style-type: none"> • Make sure that there are no creases on the leg flaps of the sling. • Make sure the leg flaps are attached at the right length.
It is difficult to apply the clips.	<ul style="list-style-type: none"> • Make sure the clip is not caught between the lug on the lift and the clip attachment. • Select the other attachment clip position (inner or outer).
It is difficult to use the sling in combination with the lift.	Use the sling only in accordance with the combinations described in section “ <i>Allowed Combinations</i> ” on page 9.
The position of the patient in the sling should be more reclined or seated.	Use the clips/loops/lift cords to adjust the patient’s position.
The patient wishes to have more support for the back.	Use the inner clips/loops for attachment or tighten the lift cords to adjust the patient’s position.
The patient feels discomfort in the waist area.	Adjust the support belt, tighten or loosen.
The sling is sliding up under the patient’s arms.	Tighten the support belt for a more comfortable fit.
The patient feels discomfort while using the sling.	<ul style="list-style-type: none"> • Make sure the patient is holding on to the handles on the lift. • Adjust the support belt for a more comfortable fit.
The <i>Flites</i> is soiled, stained or wet.	Dispose the <i>Flites</i> and replace it with a new one.
The “Do NOT wash” label is damaged or missing and the underlying red label “Do NOT use/lift” appears.	Dispose the <i>Flites</i> and replace it with a new one.
The patient cannot reach fully standing position.	Select another sling size or another attachment clip position (inner or outer).

Technical Specifications

General	
Safe working load (SWL) = Maximum total load	See "Allowed Combinations" on page 9
Service life - Recommended period of use	Flites: 2 weeks, based on 4 transfers/day* All other slings: 2 years* * Refer to "Expected Service Life" on page 3.
Shelf life - Maximum storage period of new unpacked product	5 years
Model and Type	See "Sling Selection" on page 8
Operating, Transport and Storage Environment	
Temperature	0°C to +40°C (+32°F to +104°F) Operating and Storage -25°C to +60°C (-13°F to +140°F) Transport
Humidity	Max 15–70% at +20°C (68°F) Operating and Storage Max 10–95% at +20°C (68°F) Transport
End of Life Disposal	
Box	The box consists of cardboard, recyclable to local legislations.
Package	The bag consists of soft plastic, recyclable to local legislations.
Sling	Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.

EN

Parts and Accessories

Sling	Part Description	Article number
TSS.500S-1C	Belt	TSB.100-C
TSS.501S-1C	Belt	TSB.100-C
TSS.502S-1C	Belt	TSB.100-C
TSS.503S-1C	Belt	TSB.200-C
TSS.501SV	Belt	TSB.100-V
TSS.502SV	Belt	TSB.100-V
MAA3070M	Belt and Belt retainer	KKX21440
MAA3071M	Belt and Belt retainer	KKX21440.0
MAA3071M	Fleece	KKX20580

Label on the Sling

EN

Care and Washing Symbols	
	Machine wash 70°C (158°F)
	No bleaching allowed
	No tumble drying allowed
	Tumble drying Does not apply to <i>Flites</i> or wipeable sling
	No Ironing allowed
	Do NOT wash Only applies to the <i>Flites</i> and wipeable sling
	No dry cleaning allowed
	Wipe clean Only applies to the wipeable sling

Certificates/Markings	
	CE marking indicating conformity with European Community harmonised legislation
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

Flites Symbols	
	Only on the <i>Flites</i> . (Disposable and patient specific)
	“Do NOT wash” symbol. Located on the outside of the sling.
	“Do NOT use/lift” symbol. The <i>Flites</i> has been washed. Located on the outside of the sling.

Symbol for Clip	
	Use a clip spreader bar.

Article Number	
REF XXXXXX-X	Article number with -X refers to the sling size.
REF XXXXXXX	If the article number is stated without the size letter at the end, it means the sling is a One Size.

Fibre Content	
PES	Polyester
PE	Polyethylene
PU	Polyurethane
PA	Polyamide
TPU	Thermoplastic Polyurethane

Misc. Symbols	
	Safe Working Load (SWL)
	Patient's name symbol
	Record symbol
	Read the <i>IFU</i> before use
	Article number
	Serial number
	Manufacturing date and year
	Manufacturer name and address
	Size of the sling illustrated by colour and letter. S for small, M for medium, L for large, and XL for extra large. Only applies to the wipeable sling.

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Intentionally left blank

Intentionally left blank

AUSTRALIA
Arjo Australia
Building B, Level 3
11 Talavera Road
Macquarie Park, NSW, 2113, Australia
Phone: 1800 072 040

BELGIQUE / BELGIË
Arjo Belgium nv
Evenbroekveld 16
9420 Erpe-Mere
Belgium
T: +32 (0) 53 60 73 80
F: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.com

BRASIL
Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelice, 329 PB02
Galpão - Lapa
São Paulo - SP - Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA
Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA
Arjo Czech Republic s.r.o.
Na Strži 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK
Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND
Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA
Arjo Ibérica S.L.
Parque Empresarial Rivas Futura, C/Marie
Curie 5
Edificio Alfa Planta 6 oficina 6.1-.62
ES-28521 Rivas Vacia, MADRID
Tel: +34 93 583 11 20
Fax: +34 93 583 11 22
E-mail: info.es@arjo.com

FRANCE
Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG
Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA
Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: italy.promo@arjo.com

MIDDLE EAST
Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building, North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND
Arjo BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND
Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@Arjo.com

NORGE
Arjo Norway AS
Olaf Helsets vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

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Address page - REV 25: 01/2020

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

04.SF.00-INT1

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Arjo Inc.