



- Separates into two pieces for easy carrying
- Rotation of pin either locks or allows ramp to be separated into two pieces
- Folds, carries like a suitcase
- Patented hinge eliminates pinch point
- Closure strap locks ramp panels together
- Durable welded construction
- Safety guide indicates if ramp is on a safe slope
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included



"QuickCam" Separates Ramp With No Extra Tools And No Missing Hardware



Included Safety Guide Indicates If Ramp Is On A Safe Slope

Item Number	Length	Width	Folded Size	Weight	Capacity*
BAR536	5'	36"	32¼" x 17½" x 11½"	50 lbs.	800 lbs.
BAR636	6'	36"	38¼" x 17½" x 11½"	58 lbs.	800 lbs.
BAR736	7'	36"	44¼" x 17½" x 11½"	67 lbs.	800 lbs.
BAR836	8'	36"	50¼" x 17½" x 11½"	76 lbs.	800 lbs.

Not sure which ramp length you need? Refer to pages 18 - 19 to help you.

* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).

