



Harmar Pool Lifts quickly and easily meet ADA accessibility requirements with commercial-grade quality that enhances any aquatic setting. Their smooth, durable finish, whisperquiet operation and comfortable designer seating ensure that mobility-impaired swimmers gain safe access to pools and spas.

Features & Benefits

- Capacities up to 400 lbs.
- Blends with outdoor furniture; custom color options available.
- Rust-resistant stainless steel and aluminum construction
- Smooth, whisper-quiet drive system
- 24V battery pack runs 40+ load cycles on one charge; removes for storage
- Reserve power lifts out of water even if battery goes flat.
- · Easy, water-resistant hand control
- Convenient wall-mounted charger
- Folding arms, foot rest and seat belt
- Mounting hardware (lift-specific)
- Public 'lift notice' sign

P300

In-ground anchor for new construction. 350 lbs. (160 kg) capacity



P350

Surface-bolt mounted for an existing concrete deck or new pad. 350 lbs. (160 kg) capacity



P400

High capacity, compact design, surface-bolt mounted for an existing concrete deck or new pad. 400 lbs. (181 kg) capacity



Portable Pool Lifts

Removable wheels and stand for easy transportation and storage. Surface-bolt mounted for an existing concrete deck or new pad.

Includes 'soft touch' thumb screws for easy removal and reinstallation





P400P 400 lbs. (181 kg)capacity

Accessories

P300C Lift Cover P300B Additional Battery Pack for redundancy and 100% pool lift availability.

















Custom Seating Options

Natural

Black Plus

Salsa

Porcelain



Debonair

Adobe Linen

ADA Compliance Now Required For All Public Pools & Spas

Lemon

Cypress



ADA Compliant	P300	P350	P400	P350P	P400P
Application	New Construction or Paver Mounting	Existing Concrete Deck or New Pad	Heavy Duty; Existing Concrete Deck or New Pad	Portable Existing Concrete Deck or New Pad	Portable Heavy Duty; Existing Concrete Deck or New Pad
Mounting	In-Ground Anchor	Surface Bolts	Surface Bolts	Surface Bolts	Surface Bolts
Capacity	350 lbs. (159 kg)	350 lbs. (159 kg)	400 lbs. (181 kg)	350 lbs. (159 kg)	400 lbs. (181 kg)
Profile	Very Slim	Standard	Extremely Compact	Standard	Extremely Compact
Construction	Rust Resistant, Stainless Steel and Aluminum				
Battery / Charger	24V Battery Pack / Wall-Mounted				
Seating	Comfortable Patio Design, Including Seat Belt	Comfortable Patio Design, Including Seat Belt	Comfortable Patio Design, Including Seat Belt	Comfortable Patio Design, Including Seat Belt	Comfortable Patio Design, Including Seat Belt
Seat Color	'Cognac' Standard; Custom Options				
Frame	Powder Coated 'Luna'				
Warranty	3-Year Limited Warranty; 1-Year Battery				





