

# Respiratory Care

Solutions across the continuum of care









# Empower Best Practice.

Medline offers the depth and breadth of products to meet all of your needs with innovative solutions to enhance patient experience, improve outcomes and increase efficiency – all of which is clinically supported by our Respiratory Specialists and delivered by our unparalleled supply chain.

## Table Of Contents

### Airway Management

Endotracheal Tube Holder, mBrace	5
Endotracheal Tubes, Subglottic	6
IntelliCuff® Pressure Monitor	7
Phoneris® Inner Cannula/ Aeroflex® System	8
Tracheostomy Ties	9

### Bronchial Hygiene

Volume Incentive Spirometers	11
------------------------------	----

### Gas Delivery & Suction Devices

Air & Oxygen Flowmeters	13
Medical Gas Hoses	14
Oxygen Conservers	15
Oxygen Regulators	15
Suction Regulators	16–17
Suction Aspirator	18

### Humidification & Filtration

Bacterial/Viral Filters	20
Bypass HME	21
Heat & Moisture Exchangers	22–23
Filtered Heat & Moisture Exchangers	23
Non-Heated Circuits	24
Transport Ventilator Circuits	24
Pulmonary Function Testing Filters & Accessories	25

### Medication Delivery & Asthma Management

Medication Nebulizers	27–29
MDI Spacers	30–31
Peak Flow Meters	31
Medication Delivery Accessories	32
Masks, Aerosol	33
Nebulizer Compressors	34–35

### Non-Invasive Interfaces

DreamPort™ CPAP Mask Solution	37
CPAP Masks AirGel	38–39
CPAP Masks Silicone	40
CPAP Accessories	41
NIV Veraseal	42
Go-PAP™ Systems	43
O <sub>2</sub> -Max™ Systems	43–44

### O<sub>2</sub> & CO<sub>2</sub> Disposables

Southmedic OxyMask™	46–47
Cannulas, O <sub>2</sub>	48–50
Cannulas, CO <sub>2</sub>	51
Tubing, O <sub>2</sub>	52
Masks, O <sub>2</sub>	53
Masks, Aerosol	54
Masks, Tent & Trach	54
Bubble Humidifiers	55
O <sub>2</sub> Accessories	55

### Pulse Oximeters

Fingertip Pulse Oximeters	57
Handheld Oximeters	58
Pulse Oximeter Probes & Accessories	59

### Resuscitation

Resuscitation Bags, Self Inflating	61–67
Colorimetric CO <sub>2</sub> Detectors	68
Resuscitation Masks & Accessories	69
Mouth-to-Mask Resuscitator	69
Hyperinflation Bags	70

### Suction Catheters & Trach Care

Closed Suction	72–74
Open Suction Catheters	75
Open Suction Catheters, Sleeved	75
Open Suction Kits	76
Open Suction Trays with Solutions	77
Sterile Solutions	77
Trach Care & Cleaning Trays	78
Mucus/Specimen Traps	79

### Top Respiratory Products for Seasonal Flu



# Airway Management

Endotracheal Tube Holder, mBrace	5
Endotracheal Tubes, Subglottic	6
IntelliCuff® Pressure Monitor	7
Phoneris® Inner Cannula/Aeroflex® System	8
Tracheostomy Ties	9



## mBrace

The mBrace ET tube fastener provides the enhancements care providers have been looking for. Designed with patient comfort and caregiver efficiency in mind, the mBrace can help health care providers achieve their skin breakdown initiatives.

- Clear silicone adhesive pads provide complete visibility to the skin and are gentle upon removal
- Silicone pads can be peeled back, adjusted and reapplied, giving the caregiver access to clean sensitive skin without removing the device
- Additional strap connection points allow caregiver flexibility for difficult patients and added security to help prevent lip ulcers
- Sliding ET tube track allows the ET tube to be moved side to side in accordance with hospital protocols
- ET tube strap securely locks the ET tube in place
- Adjustable neck strap provides security
- Available with 4 point head strap for improved security
- Available with and without a bite block

Item No.	Description	Pkg.
DYNJMBRC2	mBrace with 2 Point Headstrap	12/cs
DYNJMBRC2B	mBrace with 2 Point Headstrap and Bite Block	12/cs
DYNJMBRC4	mBrace with 4 Point Headstrap	12/cs
DYNJMBRC4B	mBrace with 4 Point Headstrap and Bite Block	12/cs



DYNJMBRC2,  
2 point headstrap



DYNJMBRC4,  
4 point headstrap





## SubG™ Endotracheal Tube


Our SubG subglottic endotracheal tubes are designed to reduce the risk of nosocomial infections due to aspiration of oral and upper airway secretions.

- Polyurethane cuff helps to prevent channel formation
- High volume, low pressure barrel-shaped cuff for optimal sealing
- Kink-resistant tubing with easy-to-read insertion markings
- Subglottic suction port for removing secretions above the cuff
- Suction line cap keeps contaminants out of lumen

Item No.	Description	Pkg.
DYNJ18860	SubG ET Tube, Size # 6.0	10/bx
DYNJ18865	SubG ET Tube, Size # 6.5	10/bx
DYNJ18870	SubG ET Tube, Size # 7.0	10/bx
DYNJ18875	SubG ET Tube, Size # 7.5	10/bx
DYNJ18880	SubG ET Tube, Size # 8.0	10/bx
DYNJ18885	SubG ET Tube, Size # 8.5	10/bx

### Advantages of the SubG ET Tube

- Subglottic suction line and evacuation port for removal of pooled secretions above the cuff
- Barrel-shaped cuff made of polyurethane to aid in the reduction of microaspirations
- Thermosensitive material, a bullnose tip, and easy to read marking for ease in insertion





## IntelliCuff® Cuff Pressure Controller

Continuous and automatic cuff management system.

- IntelliCuff cuff pressure controller continuously measures and automatically maintains set cuff pressure of ET or tracheal tubes during mechanical ventilation
- Secures airway management in intensive care units, operating rooms and air transport
- Large display keeps important data visible at all times, while the convenient buttons make adjusting and verifying settings easy
- Automatic regulation of cuff pressure helps keep pressure consistently within range, which could help reduce the risk of ventilator associated pneumonia (VAP) and tracheal injuries

Item No.	Description	Pkg.
HCS951001	Pressure Controller	1/ea
HCS282016	Disposable Tubing	10/cs
HCS951011	Car Adapter	1/ea
HCS951012	Power Cord	1/ea
HCS951016	Battery Cover	1/ea
HCS951019	Mount, Click	1/ea
HCS951055	Mount, Rail	1/ea
HCS951057	Multi-Angle Mount	1/ea



Optimal cuff pressure with IntelliCuff



Over-inflated cuff



Under-inflated cuff

For more information or to view an In-Service Video to see how IntelliCuff works, visit [www.Intellicuff.com](http://www.Intellicuff.com)



HCS951001



## Phoneris® Inner Cannula and Aeroflex® Ventilator Circuit Connector

The Phoneris/Aeroflex system uses an innovative design to comfortably and securely connect tracheostomy patients to ventilator circuits.

- Secure – A tight connection helps avoid accidental disconnects
- Gentle – Swivel enables comfortable movement and reduces torque on the stoma
- Easy – Puts an end to improvised, unapproved securing devices


Item No.	Description	Pkg.
HCSTIC4A	Phoneris Inner Cannula; Jackson 4	10/cs
HCSTIC6A	Phoneris Inner Cannula; Jackson 6	10/cs
HCSTIC8A	Phoneris Inner Cannula; Jackson 8	10/cs
HCSTIC10A	Phoneris Inner Cannula; Jackson 10	10/cs
HCSTIC4681A	Phoneris AeroFlex	50/cs
HCS3222	15 x 22 mm Flex Connector	50/cs
HCS3215	15 x 15 mm Flex Connector	50/cs




HCSTIC4681A



HCS3222



**Positive-Click jaws** hold tight to help prevent disconnect. Releases easily with a twist



## Tracheostomy Ties

These soft, comfortable straps have adjustable hook-and-loop closures to keep tracheostomy tubes securely in place, ensuring uninterrupted airflow.

### Premium Tracheostomy Ties

- Soft, hypoallergenic, no stretch cotton ensures patient safety and comfort while reducing skin irritation
- Offered in a wide range of sizes

#### 2-Piece

Item No.	Description	Pkg.
HCS540	Adult, Fits Up to 20" Neck	10/cs
HCS541	Pediatric, Fits Up to 18" Neck	10/cs
HCS542	Neonate/Infant, Fits Up to 9" Neck	10/cs

#### 1-Piece

Item No.	Description	Pkg.
HCS597S	Neonate/Pediatric, Fits 7–9" Neck	10/cs
HCS597M	Small Adult, Fits 9–14" Neck	10/cs
HCS597LA	Large Adult, Fits 13–19" Neck	10/cs

### Standard Medline Tracheostomy Ties

- Available in 2 sizes that fit a broad range of patients

#### 2-Piece

Item No.	Description	Pkg.
ZPP501	Adult, 18¾" Long	20/cs
ZPP501P	Pediatric, 13½" Long	20/cs

#### 1-Piece

Item No.	Description	Pkg.
ZPP301	Adult, 15" Long	20/cs
ZPP301P	Pediatric, 11" Long	20/cs
ZPP401	Anti-Disconnect	20/cs




HCS540



HCS597M



HCS597LA



ZPP301



ZPP501



ZPP301



# Bronchial Hygiene

Volume Incentive Spirometers

11

## Incentive Spirometers

Our Volume Incentive Spirometer uses an intuitive design to help patients correctly perform and monitor their post-surgical breathing exercises.

- One way valve ensures patients inhale rather than exhale through the device
- Easy to use Floating ball indicator and guidelines help patients inhale at the proper rate
- Sliding goal indicator can be easily adjusted allowing patients to monitor their progress
- Flexible tubing and dual sided calibrations allow use in various positions
- Available in 2,500 mL and 5,000 mL capacities



Item No.	Description	Pkg.
HCS93250	Incentive Spirometer 2,500 mL	12/cs
HCS93500	Incentive Spirometer 5,000 mL	12/cs



HCS93250



# Gas Delivery & Suction Devices

Air & Oxygen Flowmeters	13
Medical Gas Hoses	14
Oxygen Conservers	15
Oxygen Regulators	15
Suction Regulators	16–17
Suction Aspirator	18

**5** YEAR  
warranty

## Air & Oxygen Flowmeters

### Features:

- Complete line of flowmeters from neonatal to adult
- Available in standard Thorpe style and Dial style
- MRI options are available
- One of the smallest and lightest brass bodied Thorpe style flowmeters available.
- Color coded inner tube and labeling
- No strip "D" shape knob
- Includes an inlet filter and inlet extension
- Dimensions: 5.6 H x 2.7" W
- Weight 6-oz.
- Calibrated at 50 PSI
- Serialized for ease of identification



Dial/Click Style  
Flowmeter



Dual/Y-block Flowmeter  
Assembly



Thorpe Style



Swivel

### Accessories



Swivels/Christmas Trees  
HCS4700G (Green)  
HCS4700Y (Yellow)



Power Takeoff

### Flowmeter Ordering Information

**RTFMXXXXX(YY)\***

Gas	LPM	Connection Type	Options
O = Oxygen A = Medical Air	2X = 200 CC 1L = 1 LPM 3L = 3.5 LPM 4L = 4 LPM (Dial only) 8L = 8 LPM 15 = 15 LPM 25 = 25 LPM (Dial only) 70 = 70 LPM	AA = 1/8" FNPT DH = DISS Handtight DN = DISS Nut OM = Ohmeda Male CM = Chemetron Male PB = Puritan Bennett	PT = Power Takeoff DI = Dial M = MRI Compatible* S = Swivels

\*Some restrictions apply

For dual Y-block flowmeter configurations:  
Replace RTFM prefix with RFMD



## Medical Gas Hose Assemblies

Medline offers custom-configured medical gas hose assemblies in four easy steps:

- ¼" ID Non-Conductive

1. Select hose type
2. Select the inlet fitting (gas service and type)
3. Select the outlet fitting (gas service and type)
4. Specify the hose length

### RT = Hose Assembly

#### 1. Hose Type

OXY = Oxygen

VAC = Vacuum

N2O = Nitrous Oxide

AIR = Medical Air

N2 = Nitrogen

CO<sub>2</sub> = Carbon Dioxide

WAG = Waste Anesthetic Gas Disposal

#### 2. Inlet Fitting (Connects to the Source)

OM: OHIO Male

DC: DISS Male with Check

CM: Chemetron Male

PB: Puritan-Bennett Male Coupler

OF: OHIO Female Coupler

DN: DISS Nut/Nipple Female

DH: DISS Handtight Female

CF: Chemetron Female Coupler

PF: Puritan-Bennett Female

#### 3. Outlet Fitting (Connects to the Medical Device)

OM: OHIO Male

DC: DISS Male with Check

CM: Chemetron Male

PB: Puritan-Bennett Male Coupler

OF: OHIO Female Coupler

DN: DISS Nut/Nipple Female

DH: DISS Handtight Female

CF: Chemetron Female Coupler

PF: Puritan-Bennett Female

#### 4. Hose Length, in feet (e.g. 2 = 2-ft., 4 = 4-ft., 7 = 7-ft., etc.)


#### EXAMPLE:

**RTOXYDHD5** ¼" ID non-conductive oxygen hose assembly with a DISS handtight Female inlet and DISS Male outlet, 5' in length.

For custom configurations, contact your Medline Representative.

Yellow – RTAIRDHDN3

Green – RTOXYDHD5



## EconO2mizer II Single-Lumen Pneumatic Conservor

### Less noise for better compliance

Among the quietest pneumatic conserving devices available.

- Very quiet with no hissing during pulse delivery
- Aids in compliance with ambulatory patients
- 6 conserving flow settings (1, 2, 2.5, 3, 4 and 5)
- 3 continuous settings (2, 3 and 4)
- Features all-brass construction in high-pressure zones
- Weight: 13-oz.
- Size: 4.8 x 2.4 x 1.5" (12.2 x 6.1 x 3.8 cm)
- 5.7:1 conserving ratio



HCS2CON1L

Item No.	Description	Pkg.
HCS2CON1L	870 connection	ea
HCSTHANDLE	Replacement T Handle, Black	ea

## Oxygen Regulators

For all medical and emergency applications.

- Fitted with brass sleeves around key parts
- Laser-etched knob and body

### Economy

Item No.	Description	Flow Range	Flow Increments	Pkg.
HCS8708M	870 CGA Connection	0–8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	1/ea
HCS8715M	870 CGA Connection	0–15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1/ea
HCS8725M	870 CGA Connection	0–25 LPM	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25	1/ea
HCS5408M	540 CGA Connection	0–8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	1/ea
HCS5415M	540 CGA Connection	0–15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1/ea



HCS8715M

HCS5415M

### Super Mini

Item No.	Description	Flow Range	Flow Increments	Pkg.
HCS8715MA	870 Connection	0–15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1/ea

### Accessories

Item No.	Description	Pkg.
HCS870W	Oxygen Cylinder Wrench, Small (CGA-870)	10/bg
HCS870WC	Oxygen Cylinder Wrench with Security Chain, Small (CGA-870)	10/bg
HCS870WLG	Oxygen Cylinder Wrench, Large (CGA-870)	25/bg
HCS166840B	Brass Yoke Washer Seal for CGA-870 Oxygen Regulator	25/bg



HCS8715MA

For specialty oxygen regulators (EMS, all-brass, etc.) please contact your **Medline Sales Representative** or email [ProductSupportRespiratory@Medline.com](mailto:ProductSupportRespiratory@Medline.com) with your request.



NEW


# The Choice is Clear

With a unique transparent and disposable back, Medline regulators enhance efficiency, reduce costs and improve patient safety.

## Identify Contamination Instantly

Up to 37% of suction regulators in hospitals are contaminated due to fluid breaches, and pathogens can spread back through the suction line within 48 hours.<sup>1</sup> The Medline Suction Regulator has a clear back plate offering a view of the internal components of the regulator, so it's easy to see when it's contaminated.

Suction Regulators: A Potential Vector for Hospital Acquired Pathogens by Keith S. Kaye, MD, MPH; Dror Marchaim, MD; Chester Smialowicz, MD; Lauren Bentley, MSBME



## There's no time for downtime

The clear back plate provides an internal view of the regulator, making it possible to see fluid breaches. Plus, cleaning is simple: the unit is easily taken apart, allowing you to discard the contaminated back plate and replace with a new component in minutes.



### Step 1:

Remove contaminated regulator



### Step 2:

Dispose



### Step 3:

Reassemble with only 1 NEW Component in just minutes

## Clear Back Plate

### INFECTION CONTROL MEASURE

Remove from the outlet and identify contamination immediately

For specialty gas delivery products including blenders, analyzers, and more, please contact your **Medline Sales Representative** or email [ProductSupportRespiratory@Medline.com](mailto:ProductSupportRespiratory@Medline.com) with your request.

## Suction Regulators

### Features:

- One of the smallest and lightest suction regulators available (10-oz.)
- ABS polymer casing – impact resistant and durable
- User-adjustable timing can be accessed from the back of the unit
- Easy to read, color-coded gauges
- 10 Year warranty (Analog)
- Only unit on the market with a disposable, see-through back body
- Disposable component for easy cleaning
- Dimensions: 6 H x 3 W x 4" L
- Weight: Continuous: 10-oz., Continuous/Intermittent: 13-oz.



**10YEAR**  
warranty

### Optional Color Coding



### Continuous Models – White Labels

Continuous  
3 Mode



Pediatric  
2 Mode



Continuous  
High/Surgical



### Continuous/Intermittent Models – Grey Labels

Adult Continuous/  
Intermittent




Pediatric Continuous/  
Intermittent



## Suction Regulator Ordering Information

**RTXXXXXX**



## Vacuum Traps

Item No.	Description
RTTRAPMT	1/8" Male NPT Fitting
RTTRAPDH	Female DISS Handtight
RTTRAPDN	Female DISS Nut/Nipple



Vacuum/Safety Traps

For MRI compatible please add "-M" to the end of the part number (available for adult/analog units only).  
For color coding, please specify in sales notes. Default color is white.

## Suction Aspirator and Disposable Canister Kits

Our easy-to-use aspirators are ultra-quiet.

- Weight: 12.1-lb. (5.49 kg)
- Maximum suction: -21 inHg/520 mmHg
- Rugged construction is ideal for bedside carts
- Disposable 850 cc or 1,200 cc suction canister sets with float lids and 850 cc suction canister sets with filter lids.
- HCS7121, HCS7851, and HCS7852 are kits that include suction tubing, 1/each 6' and 18" length.
- Replacement lids are not available

Item No.	Description	Pkg.
HCS7000	Portable Suction Aspirator	1/ea



HCS7000

### Canister with Float Lid

Item No.	Description	Pkg.
HCS7850	Canister with Float Lid: 850 cc	12/cs
HCS7851	Canister with Float Lid & Tubing: 850 cc	12/cs
HCS7120	Canister with Float Lid: 1,200 cc	12/cs
HCS7121	Canister with Float Lid & Tubing: 1,200 cc	12/cs



HCS7120

### Canister with Filter Lid


Item No.	Description	Pkg.
HCS7852	Canister with Filter Lid & Tubing: 850 cc	50/cs
HCS7853	Canister with Filter Lid: 850 cc	100/cs



HCS7850

### Accessories

Item No.	Description	Pkg.
HCS4053	Hydrophobic Filter Set, Barb Both Ends	3/pk
HCS701311	Replacement Gauge Assembly: Includes Pressure Gauge, Pressure Adjustment Knob & Barb Tubing Adapter	1/ea



Tubing included with  
HCS7851, HCS7852, and HCS7121



# Humidification & Filtration

Bacterial/Viral Filters	20
Bypass HME	21
Heat & Moisture Exchangers	22-23
Filtered Heat & Moisture Exchangers	23
Non-Heated Circuits	24
Transport Ventilator Circuits	24
Pulmonary Function Testing Filters & Accessories	25

## Bacterial/Viral Filters

Medline's filters can be used for both respiratory and anesthesia applications. They have translucent housings to aid in viewing built up secretions or other obstructions.

### Performance Specifications

	DYNJAABV1	DYNJAAHEPA
<b>Description</b>	Bacterial/Viral Filter	HEPA Filter
<b>Tidal Volume (mL)</b>	≥250 mL	200–1,500
<b>BFE (%)</b>	99.99	99.9999
<b>VFE (%)</b>	99.99	99.9999
<b>Volume (mL)</b>	34	52
<b>Weight (g)</b>	15.1	40
<b>Flow Resistance @ 30 LPM</b>	0.71	1.48
<b>Flow Resistance @ 60 LPM</b>	1.02	3.33
<b>CO<sub>2</sub> Port</b>	No	Yes
<b>Filter Type</b>	Electrostatic	HEPA
<b>Pkg.</b>	50/cs	50/cs



DYNJAABV1



DYNJAAHEPA

Used for both  
respiratory  
and anesthesia  
applications

## Bypass Heat and Moisture Exchanger

The AARC recommends that ventilator circuits should not be routinely changed and that passive humidifiers can be used for up to 48 hours.<sup>1</sup> It is widely understood that the opening of a respiratory circuit increases the risk of hospital acquired infection. Delivering aerosolized medication creates a challenge in maintaining a closed circuit.

The Medline Bypass Heat and Moisture Exchanger supports AARC clinical recommendations by eliminating the need to open the circuit during aerosol medication delivery for up to 72 hours.

### Features:

- Easy and intuitive to use
- Clear marking of “HME” and “Aerosol” modes
- Use for up to 72 hours

### Performance Specifications

- Tidal Volume 250–1500 mL
- Moisture Output 32.09 mg H<sub>2</sub>O/L @ 500 mL TV at 72 hours
- Resistance to Flow 0.6 cm H<sub>2</sub>O @ 30 L/min flow
- Dead space 29.2 mL
- Weight 39.3 g
- Duration of Use Up to 72 hours

Item No.	Description	Pkg.
DYNJAABHME	15 mm ID/22 mm OD x 22 mm ID	25/cs

Humidification  
Simplified



1. AARC Evidence-Based Clinical Practice Guidelines, Care of the Ventilator Circuit And Its Relation to Ventilator-Associated Pneumonia. Respiratory Care, September, 2003, Volume 48, No 9.



# Heat and Moisture Exchangers and Filtered Heat and Moisture Exchangers

Medline's broad portfolio of HMEs and HMEFs meet your patient's humidification needs throughout the care continuum.

- Multiple tidal volume ranges and sizes to meet your clinical needs
- Non-filtered, electrostatic-filtered, and HEPA-filtered options available
- Angled options available to reduce the need for additional adapters and dead space
- Translucent housing for viewing of possible obstructions



	DYNJAAHME10	DYNJAAHME12	DYNJAAHME14	DYNJAPHME18	DYNJAPHME20
<b>Description</b>	HME, Small with Port	HME, Small, No Port	HME, Angled	HME, Neonatal	HME, Pediatric with Port
<b>Type</b>	Straight	Straight	N/A	Straight	Straight
<b>Connection</b>	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	15 mm OD x 22 mm ID / 22 mm OD x 15 mm ID	15 M x 15 F	22 M / 15 F x 22 F / 15 M
<b>Vt (mL)</b>	150–1500	120–1500	120–1500	10–15	120–1500
<b>Vd (mL)</b>	40.7	32.1	34.5	2	31.3
<b>Weight (g)</b>	14.7	14.7	16.8	2	13.4
<b>CO<sub>2</sub> Port</b>	Yes	No	Yes	No	Yes
<b>Resistance (cm H<sub>2</sub>O)</b>	0.8 @ 30 LPM	0.6 @ 30 LPM	0.5 @ 30 LPM	1.5 @ 10 LPM	2.1 @ 30 LPM
<b>Moisture Output (Mg H<sub>2</sub>O/L)*</b>	34.1 @ 500 mL	33.7 @ 500 mL	33.2 @ 500 mL	28.5 @ 250 mL	33.6 @ 500 mL
<b>HME Media</b>	ElectroStatic	Foam	Paper	Foam	Electrostatic
<b>Pkg.</b>	50/cs	50/cs	50/cs	50/cs	50/cs

	DYNJAAHME30	DYNJAAHME40
<b>Description</b>	HME, Trach	HME, Trach with Suction and Oxygen ports
<b>Type</b>	N/A	N/A
<b>Connection</b>	15 mm ID	15 mm ID
<b>Vt (mL)</b>	> 25	>50
<b>Vd (mL)</b>	14	13.4
<b>Weight (g)</b>	3.5	7.1
<b>CO<sub>2</sub> Port</b>	No	No
<b>Resistance (cm H<sub>2</sub>O)</b>	0.51 @ 30 LPM	0.42 @ 30 LPM
<b>Moisture Output (Mg H<sub>2</sub>O/L)*</b>	25.9 @ 500 mL	26 @ 500 mL
<b>HME Media</b>	Paper	Paper
<b>Pkg.</b>	50/cs	50/cs



DYNJAAHME30



DYNJAAHME40

Our Humidification and Filtration offering is expanding, for the latest options visit our website or call your Medline Representative



	DYNJAAHME1	DYNJAAHME5	DYNJAAHME9	DYNJAAHME11
<b>Description</b>	Filtered HME, Large with Port	Filtered HME, Large, No Port	Filtered HME, Small with Port	Filtered HME, Small, No Port
<b>Type</b>	Straight	Straight	Straight	Straight
<b>Connection</b>	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M
<b>Vt (mL)</b>	150–1500	150–1500	120–1500	120–1500
<b>Vd (mL)</b>	40.7	40.7	31.3	31.3
<b>Weight (g)</b>	18.5	18.5	15	15
<b>CO<sub>2</sub> Port</b>	Yes	No	Yes	Yes
<b>Resistance (cm H<sub>2</sub>O)</b>	0.8 @ 30 LPM	0.8 @ 30 LPM	2.1 @ 30 LPM	2.1 @ 30 LPM
<b>Moisture Output (Mg H<sub>2</sub>O/L)*</b>	34.1 @ 500 mL	34.1 @ 500 mL	33.6 @ 500 mL	33.6 @ 500 mL
<b>Filter Media</b>	ElectroStatic	ElectroStatic	Electrostatic	Electrostatic
<b>BFE(%)</b>	99.999	99.999	99.999	99.999
<b>VFE(%)</b>	99.999	99.999	99.999	99.999
<b>Pkg.</b>	50/cs	50/cs	50/cs	50/cs




	DYNJAAHME13	DYNJAPHME17	DYNJAAHMEHEPA	DYNJ603011
<b>Description</b>	Filtered HME, Small, Angled with Port	Filtered HME, Pediatric with Port	Filtered HME, HEPA with CO <sub>2</sub> Port	Filtered HME, Flat/Rec, Large, Sample Port
<b>Type</b>	Angled	Straight	Straight	Straight
<b>Connection</b>	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 ID x 22 OD / 15 ID
<b>Vt (mL)</b>	120–1500	90–1500	200–1500	200–1500
<b>Vd (mL)</b>	33.7	20.5	52	62
<b>Weight (g)</b>	17.2	13.6	40	37.2
<b>CO<sub>2</sub> Port</b>	Yes	Yes	Yes	Yes
<b>Resistance (cm H<sub>2</sub>O)</b>	1.9 @ 30 LPM	2.4 @ 30 LPM	1.5 @ 30 LPM	<1 @ 30 LPM
<b>Moisture Output (Mg H<sub>2</sub>O/L)*</b>	33.4 @ 500 mL	34.0 @ 250 mL	30	32.5
<b>Filter Media</b>	Electrostatic	Electrostatic	HEPA	Electrostatic
<b>BFE(%)</b>	99.999	99.999	99.9999	>99.99%
<b>VFE(%)</b>	99.999	99.999	99.9999	>99.99%
<b>Pkg.</b>	50/cs	50/cs	50/cs	50/cs

## Non-Heated Breathing Circuits

Medline's family of non-heated breathing circuits support your mechanically ventilated patients with several configurations, which includes "Y" and parallel "Y" options.

### Adult Non-Heated Circuits

Item No.	Description	Pkg.
HCS1830	Adult Dual Limb Circuit, Wye without Ports, 6-ft.	20/cs
HCS72R7002	Adult Dual Limb Circuit, Uncuffed, Wye without Ports, 6-ft.	20/cs
HCS1831	Adult Dual Limb Circuit, Wye with Temperature Sensing Port, 6-ft.	20/cs
HCS1831P	Adult Dual Limb Circuit, Parallel Wye with Temperature Sensing Port, 6-ft.	20/cs
HCS1571215	Adult Dual Limb Circuit, Uncuffed, Parallel Wye with Temperature Sensing Port, Bacterial Filter, 5-ft.	20/cs
HCS72R7401	Adult Dual Limb Circuit, Uncuffed, Parallel Wye without Ports, HMEF, and Bacterial Filter, 6-ft.	20/cs
DYNJ300	Adult Dual Limb Circuit, Wye without Ports, Uncuffed, 6-ft.	20/cs



### Infant and Pediatric Non-Heated Circuits

Item No.	Description	Pkg.
DYNJ100	Infant Dual Limb Circuit, Uncuffed, Parallel Wye without Ports, 5-ft.	20/cs
DYNJ200	Pediatric Dual Limb Circuit, Uncuffed, Parallel Wye without Ports, 5-ft.	20/cs

### MRI-Conditional Transport Ventilator Circuit

Medline offers ready-to-use, transport ventilator circuits that include expandable tubing and an MRI-conditional PEEP valve

- Tubing is expandable up to 120" in length
- Circuit and PEEP valve are MRI-conditional
- Patients can be scanned with the device in operation under the following conditions: static magnetic field up to 3.0 T; maximum spatial gradient field less than or equal to 30 T/m

Item No.	Description	Pkg.
DYNJAA300	Transport Ventilator Circuit with Bacterial Filter	20/cs
DYNJAA300NF	Transport Ventilator Circuit	20/cs



Medline can customize a ventilator circuit to your needs, contact your Medline Representative



## Pulmonary Function Testing (PFT) Filters and Accessories

Medline offers multiple types of pulmonary function filters and accessories to meet your clinical needs. Our PFT filters have a BFE of 99.98% and VFE of 99.92% and mate with various PFT machines. They are packaged separately or kitted with a nose clip and mouthpiece.



DYNJ756SB



DYNJ756ELB

### PFT Filters

Item No.	Description	Fits the Following Devices	Pkg.
DYNJ756BMN	PFT Filter Kit, Standard Round Patient End (Includes Nose Clip and Mouthpiece)	Koko, Kokolegend, Keystone, Koko-Trek, Cybermedic Full Gas Systems, Collins 4g and 3g Full Gas Systems	100/cs
DYNJ756ELB	PFT Filter, Elliptical Patient End		100/cs
DYNJ756NB	PFT Filter Kit, Standard Round Patient End		100/cs
DYNJ756SB	PFT Filter Kit, Small Round Patient End		100/cs
DYNJ757BMN	PFT Filter Kit, Standard Round Patient End (Includes Nose Clip and Mouthpiece)	Sensormedics, 2200, 2400, 2450, 2800 Vmax, Keystone I & II, Pall, Cybermedics, Moose, Gould 5004, New Vyaire Masterscreen (Launched April 2019)	100/cs
DYNJ757ELB	PFT Filter, Elliptical Patient End		100/cs
DYNJ757NB	PFT Filter Kit, Standard Round Patient End		100/cs
DYNJ757SB	PFT Filter Kit, Small Round Patient End		100/cs
DYNJ758B	PFT Filter Kit, Standard Round Patient End	Collins DC II and HDPFT 4000 Plethysmograph	100/cs
DYNJ758BMN	PFT Filter Kit, Standard Round Patient End (Includes Nose Clip and Mouthpiece)		100/cs
DYNJ758ELB	PFT Filter, Elliptical Patient End		100/cs
DYNJ759B	PFT Filter Kit, Standard Round Patient End		100/cs
DYNJ759BMN	PFT Filter Kit, Standard Round Patient End (Includes Nose Clip and Mouthpiece)	Incheck, Incheck Dial, Multispiro, Creative Biomedics, Vitalograph, Schiller, Microdirect, SDI (Pulmogard)	100/cs
DYNJ759ELB	PFT Filter, Elliptical Patient End		100/cs
DYNJ760B	PFT Filter Kit, Standard Round Patient End		100/cs
DYNJ760BMN	PFT Filter Kit, Standard Round Patient End (Includes Nose Clip and Mouthpiece)		100/cs
DYNJ760ELB	PFT Filter, Elliptical Patient End	Brentwood, Medical Graphics All	100/cs
DYNJ780B	PFT Filter Kit, Standard Round Patient End		100/cs
DYNJ780BMN	PFT Filter Kit, Standard Round Patient End (Includes Nose Clip and Mouthpiece)		100/cs
DYNJ780ELB	PFT Filter, Elliptical Patient End		100/cs

### Nose Clips

Item No.	Description	Pkg.
DYNJ750	Nose Clip	100/bx
DYNJ751	Nose Clip with Marshmallow Pads	100/bx
DYNJ752	Large Nose Clip	50/bx
DYNJ755	Nose Clip and Mouthpiece Kit	100/bx



DYNJ751



DYNJ750



DYNJ756BMN

### Mouthpiece

Item No.	Description	Pkg.
DYNJ754	Mouthpiece, White, Thermoplastic	50/bx

DYNJ754

# Medication Delivery & Asthma Management

Medication Nebulizers	27–29
MDI Spacers	30–31
Peak Flow Meters	31
Medication Delivery Accessories	32
Masks, Aerosol	33
Nebulizer Compressors	34–35

# Medication Nebulizers

Breath-Actuated



## Aeromist BAN

The Aeromist breath-actuated nebulizer delivers the performance you require from a high-efficiency nebulizer, with additional features to improve operational efficiency.

- Consistent drug delivery with minimal waste; demonstrated clinical efficacy in lab testing
- Fast treatment time improves patient satisfaction and saves caregiver time
- Limited drug loss to the environment reduces caregiver exposure to aerosol
- Connects directly with a face mask, no additional pieces needed
- Easily switch between breath-actuated and continuous mode
- Versions with masks intended for use in continuous mode

Item No.	Mouthpiece	Mask Type	Filter	Tubing Length	Tubing Connector	Pkg.
HCS4433BAN	•			7-ft.	Standard	50/cs
HCSU4433BAN	•			7-ft.	Universal	50/cs
HCS4435BAN		Adult		7-ft.	Standard	50/cs
HCSU4435BAN		Adult		7-ft.	Universal	50/cs
HCS4436BAN		Pediatric		7-ft.	Standard	50/cs
HCSU4436BAN		Pediatric		7-ft.	Universal	50/cs
HCS4433FBAN	•		•	7-ft.	Standard	50/cs



HCS4433BAN/  
HCSU4433BAN

HCS4436BAN

HCS4435BAN

HCS4433FBAN



# Medication Nebulizers

## Entrainment

### SideStream® Nebulizer

The SideStream® nebulizer offers treatment times up to twice as fast compared to conventional nebulizers with an improved respirable fraction.

- Diamond Jet provides a high respirable fraction and consistent delivery
- Treatment times in approximately 4 minutes

Item No.	Description	Pkg.
HS800ML	SideStream Nebulizer	50/cs
1025529ML	SideStream Pediatric Aerosol Mask	50/cs
1025531ML	SideStream Adult Aerosol Mask	50/cs




# Medication Nebulizers

## Standard

### Handheld Nebulizer Kits

Medline offers a wide selection of reusable and disposable nebulizer kits for medication delivery.

- Ideal particle size of less than five microns, increasing the amount of medication delivered deep within lungs
- Standard fittings accommodate most clinical applications
- High respirable dose, leading to more efficient treatments and less wasted medication



HCSU4483

Item #	Adapter	Valved Tee	Mouthpiece	Reservoir Tube	Tubing Length	Tubing Connector	Mask Type	Disposable/Reusable*	Pkg.
HCS4485					7-ft.	Standard	Adult	Disposable	50/cs
HCSU4485					7-ft.	Universal	Adult	Disposable	50/cs
HCS4483	22 mm OD/ID		●	●	7-ft.	Standard		Disposable	50/cs
HCS448314	22 mm OD/ID		●	●	14-ft.	Standard		Disposable	50/cs
HCSU4483	22 mm OD/ID		●	●	7-ft.	Universal		Disposable	50/cs
HCSU448314	22 mm OD/ID		●	●	14-ft.	Universal		Disposable	50/cs
HCS4482	22 mm OD/ID		●		7-ft.	Standard		Disposable	50/cs
HCSU4482	22 mm OD/ID		●		7-ft.	Universal		Disposable	50/cs
HCS4486					7-ft.	Standard	Pediatric	Disposable	50/cs
HCSU4486					7-ft.	Universal	Pediatric	Disposable	50/cs
HCS4480*					None	None		Disposable	50/cs
HCS4483VT		22 mm OD/OD		●	7-ft.	Standard		Disposable	50/cs
HCS4481VT		22 mm OD/ID			7-ft.	Standard		Disposable	50/cs
HCS4482R	22 mm OD/ID		●		7-ft.	Standard		Reusable	50/cs
HCS4483R	22 mm OD/ID		●	●	7-ft.	Standard		Reusable	50/cs
HCS4482RR	45° Angled		●		7-ft.	Standard		Reusable	50/cs
HCS4493**	22 mm OD/OD		●	●	7-ft.	Standard		Disposable	50/cs
HCS4492**	22 mm OD/OD		●		7-ft.	Standard		Disposable	50/cs
HCS4495**	22 mm OD/OD				7-ft.	Standard	Adult	Disposable	50/cs

Disposable up to 2 weeks of use/reusable up to 6 months of use.

\*Neb cup only, no additional components

\*\*Side air intake design



Upstream Design



Side Air Intake Design

# MDI Spacers

## OptiChamber® Diamond

The OptiChamber® Diamond valved holding chamber's small and intuitive design enhances medication delivery and compliance for all patients.

- Detachable LiteTouch mask has a comfortable, soft-seal interface, designed to create ideal seal with minimal force
- Anti-static valved holding chamber, smaller than most conventional chambers, intuitive design that enhances medication delivery
- Even under low pediatric pressures and flows, valves open freely
- Stepped mouthpiece helps patients transition from pediatric facemasks to mouthpiece; allows attachment to a 22 mm connector
- Mouthpiece and adapter can be removed from the chamber for easy cleaning



1098215ML



1098213ML

Item No.	Description	Pkg.
1098213ML	OptiChamber Diamond Only, 10 Pack	10/cs
HCS1143035	OptiChamber Diamond Only, 50 Pack	50/cs
1098216ML	OptiChamber Diamond with Large Mask	10/cs
1098215ML	OptiChamber Diamond with Medium Mask	10/cs
1098214ML	OptiChamber Diamond with Small Mask	10/cs

## LiteTouch Masks

The LiteTouch masks are designed to provide greater comfort and a more secure seal for increased treatment effectiveness.

- Mask quickly and comfortably conforms to facial contours, providing a soft, secure fit for increased compliance and treatment effectiveness
- Flexible seal makes treatment more comfortable and less technique dependent
- Ensures a proper fit with as little as one pound of force

Item No.	Description	Pkg.
1083783ML	LiteTouch Mask for MDI Spacer, Large	10/cs
1083784ML	LiteTouch Mask for MDI Spacer, Medium	10/cs
1083785ML	LiteTouch Mask for MDI Spacer, Small	10/cs



## Medline MDI Spacers

Medline spacers are designed for patients who have difficulty with metered dose inhalers (MDIs). With spacers, more medication is deposited in the lungs and less is deposited in the mouth and throat.

- A simple and economical way to improve effectiveness of inhaler
- Easy to use and helps avoid some of the unpleasant taste of inhaled medication
- Anti-static and compatible with any standard mask with a 22 mm fitting
- Masks are made of a soft silicone and fit snugly to ensure optimal drug delivery
- The collapsible model is portable and convenient

Item No.	Description	Pkg.
HCSA100	MDI Spacer Only	10/bx
HCSA200	MDI Spacer Only, Collapsible	10/bx
HCSA112	MDI Spacer with Infant Mask	10/bx
HCSA113	MDI Spacer with Pediatric Mask	10/bx
HCSA114	MDI Spacer with Adult Mask	10/bx



HCSA113

HCSA112



HCSA114

## Peak Flow Meters

Peak flow meters help patients manage their asthma and gauge the effectiveness of their treatment program, which may help them predict their next attack.

### Personal Best

Lightweight, portable peak flow meter helps patients monitor their asthma. Handheld device measures Peak Expiratory Flow Rate so patients can initiate preemptive therapy in order to help prevent serious asthma attacks.

- Universal range allows the Personal Best peak flow meter to be used for all patients, from pediatric to adult
- Integrated case for storage
- Mouthpieces available separately

Item No.	Description	Pkg.
1022973UML	Personal Best Peak Flow Meter, Universal Range	10/cs
HS714200ML	Plastic Mouthpiece with One-Way Valve	200/cs
HS715100ML	Plastic Pediatric Mouthpiece	100/cs
HS757100ML	Oval Disposable Plastic Mouthpiece	100/cs



1022973UML



HCSPF100

### Medline

- Individually calibrated to ensure accuracy
- Durable for years of reliable performance
- The pocket model offers superb convenience and portability
- Available in low range and full range versions—low range models are better for geriatric, pediatric or other patients with diminished flow rates

Item No.	Description	Pkg.
HCSPF100	Portable Peak Flow Meter, Full Range, 60–720 LPM	12/bx
HCSPF200	Standard Peak Flow Meter, Full Range, 60–850 LPM	12/bx
HCSPFPEDI	Mini Wright Pediatric Peak Flow Meter, 30–400 LPM	100/cs
HCSPF220	Disposable Pediatric Mouthpiece, Adapter	20/bx
HCSPF262	Disposable One Way Valve Mouthpiece, Adult	100/cs



HCSPF200



HCSPFPEDIH



HCSPF220



HCSPF262



## Valved Tee Adapters

- For inline medication delivery
- Available in 22 mm OD/OD and 22 mm OD/ID

Item No.	Description	Pkg.
HCSA1008	22 mm OD/OD	30/cs
HCSA1009	22 mm OD/ID	30/cs



HCSA1009

HCSA1008

## Southmedic Inline MDI Adapters

- Single-patient-use adult metered dose inhaler (MDI) adapter ensures safe and complete delivery of prescribed aerosol medication
- Adapter remains inline to help reduce cross-contamination and improve patient comfort by reducing manipulation of circuit
- Arrows clearly mark the direction of flow; built-in tethered plug can be used to seal the breathing circuit when no MDI canister is in place
- Compatible with both plastic and metal canister stems
- Can be used with or without dosage counter

Item No.	Description	Pkg.
SMDIA1000ML	22 mm ID → 22 mm OD / 15 mm ID	25/cs
SMDIA1100ML	22 mm OD / 15 mm ID → 22 mm ID	25/cs
SMDIA200ML	22 mm OD → 22 mm OD / 15 mm ID	25/cs
SMDIA1300ML	15 mm OD → 15 mm ID	25/cs



SMDIA1000MLH

## Aerosol Masks

Extra patient comfort comes from an anatomical form and a soft, light and comfortable vinyl material. All aerosol masks feature an elongated fit for optimal patient fit.


- A secure fit with elastic straps
- Standard 22 mm connector
- Not made with natural rubber latex

### Adult

Item No.	Description	Pkg.
HCS4630B	Aerosol Mask with Straight Connector	50/cs
HCS4630	Aerosol Mask with Swivel Connector	50/cs

### Pediatric

Item No.	Description	Pkg.
HCS4631B	Aerosol Mask with Swivel Connector	50/cs



## Aeromist Compact™ Nebulizer Compressors

Our smallest nebulizer compressor. Choose from the standard white model or color options.

- Kits include a disposable nebulizer medicine cup, tee mouthpiece, 6" corrugate reservoir, 7-ft (2.13 m) of tubing, and 5 replacement filters
- Color choices: electric blue, neon green, neon pink, charcoal, or white
- Durable hard plastic shell
- Warranty: 5 years
- Replacement Kits: HCS4483 (Disposable) and HCS4483R (Reusable)


### Specifications

- Method of nebulization: dry piston compressor
- Electrical requirements: 120 VAC / 60 Hz
- Maximum jet nebulizer capacity: 10 mL
- Sound noise level: ≤55 dBA
- Pressure: 30 PSI, Flow Rate: 8 LPM



HCS70004BL

Item No.	Description	Pkg.
HCS70004G	Neon Green Compressor and Nebulizer Kit	8/cs
HCS70004BL	Electric Blue Compressor and Nebulizer Kit	8/cs
HCS70004P	Neon Pink Compressor and Nebulizer Kit	8/cs
HCS70004C	Charcoal Compressor and Nebulizer Kit	8/cs
HCS70004	White Compressor and Nebulizer Kit	8/cs
HCS70004ADLT	White Compressor and Nebulizer Kit with Adult Mask	8/cs
HCS70004RD	White Compressor with Reusable and Disposable Neb Kit	8/cs
HCS70004DP	Sample Pack, Two of each color	8/cs
HCS70004F	Replacement Filters	5/pk
HCS70004CAP	Replacement Cap for HCS70004	8/pk



Components included with each compressor

## Aeromist Trek Portable Nebulizer Compressor

Better compliance with easy treatments, anytime and anywhere.

- Suitable for travel with a sleek design, extra small size, light weight and battery operation
- Durable hard plastic shell
- Compact and stable, with a small footprint
- Kit includes a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, mouthpiece, rechargeable battery, carry case, AC and DC adapters, adult and pediatric masks, and 5 replacement filters
- Warranty: 3 years for the compressor, 6 months for the battery

Item No.	Description	Pkg.
HCS2GO	Portable Battery-Operated Nebulizer Compressor Kit	1/ea
HCS2GONB	Portable Nebulizer Compressor, Battery not included	1/ea
HCS2GOBATT	Replacement Battery for Portable Nebulizer Compressors (Items HCS2GO and HCS2GONB)	1/ea
HCS2GOF	Replacement filters	5/pk
HCS2GOAC	Replacement AC Adapter	1/ea
HCS2GODC	Replacement Car Charger	1/ea



HCS2GO

## Aeromist Buddies™ Nebulizer Compressors

Kid-friendly nebulizers help put young patients at ease, while still providing safe and effective aerosol treatments.

- Kit includes disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, 6" reservoir tube, T-Mouthpiece, pediatric mask, and 5 replacement filters
- Durable hard plastic shell
- Warranty: 5 years

### Specifications

- Method of nebulization: dry piston compressor
- Electrical Requirements: 120 VAC / 60 Hz
- Maximum nebulizer capacity: 10 mL
- Sound noise level: ≤60 dBA
- Pressure: 30 PSI, Flow Rate: 8 LPM

Item No.	Description	Pkg.
HCSMONKEY	Monkey Compressor and Nebulizer Kit	8/cs
HCSMONKEYCAP	Replacement Filter Cap for HCSMONKEY	8/pk
HCSFISHNEB*	Fish Compressor and Nebulizer Kit	8/cs
HCSPUPPY*	Puppy Compressor and Nebulizer Kit	8/cs
HCSPUPPYCAP	Replacement Filter Cap for HCSPUPPY	8/pk
HCSPUPPYRD*	Puppy Compressor with Disposable and Reusable Nebulizer Kits	8/pk
HCSTURTLE	Turtle Compressor and Nebulizer Kit	8/cs
HCSTURTLECAP	Replacement Filter Cap for HCSTURTLE	8/pk
HCS70004F	Replacement Filters	5/pk

\*Comes with Grab-N-Go backpack for convenient storage.



All compressors include a complementary nebulizer kit.

**Kid  
Friendly  
Nebulizers**





# Non-Invasive Interfaces

DreamPort™ CPAP Mask Solution	37
CPAP Masks AirGel	38–39
CPAP Masks Silicone	40
CPAP Accessories	41
NIV Veraseal	42
Go-PAP™ Systems	43
O <sub>2</sub> -Max™ Systems	44

## DreamPort™ CPAP Mask Solution

The DreamPort™ CPAP Mask Solution from Bleep LLC, is an innovative CPAP device designed to eliminate leaks, and improve patient comfort and compliance.

- It is one-size-fits-all and is customizable to each patient's needs
- Adhesive, hypoallergenic surgical tape on the DreamPorts create an airtight seal
- The DreamPorts pair with the DreamWay that ensures a strapless connection to standard-sized CPAP tubing
- Lightweight and easy to use

NEW



Non-Invasive Interfaces

Item No.	Description	Pkg.
CPAP100383	Bleep DreamPort Sleep Solution Kit	24/cs
CPAP100382	Bleep DreamPorts, 32 Ports per Box	54/cs
CPAP100387	Bleep DreamWay CPAP Interface	24/cs
CPAP100434	Bleep DreamPreps, 100 Wipes	8/cs



CPAP100383  
(includes 1 DreamWay Interface and 1 box of DreamPorts)



CPAP100382



CPAP100387



CPAP100434

## AIRgel™ CPAP Masks

Amazing comfort and customized fit technology. The patented AIRgel replaceable cushioning minimizes discomfort and skin irritation while effectively sealing airflow. It is easily detached and reattached for easy cleaning, saving time and money.

- Proprietary AIRgel cushioning minimizes discomfort and skin irritation while effectively sealing airflow
- The flexible and formable mask provides a custom fit—flex, fit and forget

## Aura Nasal Mask

- The shapeable shell of each mask can be easily molded to fit individual needs
- Flexible spacebar is designed to provide superior stability without touching the skin
- Active headgear connector moves with the patient to help maintain an effective seal

Item No.	Description	Mask Size	Pkg.
CPAP50176	Adjustable Gel Cushion	Small	1/ea
CPAP50177	Adjustable Gel Cushion	Medium	1/ea
CPAP50178	Adjustable Gel Cushion	Large	1/ea



## Mojo® Full-Face Gel Mask

- Full-face design in four sizes
- Adjustable tension wheel minimizes discomfort while ergonomic cushions adapt to unique facial features for a comfortable fit
- Magnetic headgear clips hold firmly and disconnects with one hand easy on-and-off

Item No.	Description	Headgear Style	Mask Size	Pkg.
CPAP50835	Vented	One Size	Small	1/ea
CPAP50836	Vented	One Size	Medium	1/ea
CPAP50837	Vented	One Size	Large	1/ea
CPAP50838	Vented	One Size	XL	1/ea
CPAP55045	Headgear Only	One Size	One Size	1/ea



Headgear included

CPAP50177

## AIRgel™ CPAP Masks

### iQ Blue® Nasal Gel Mask

- One-size nasal mask design fits a wide variety of patients
- Clearest field of vision possible with no forehead pad required—patients can keep their glasses on while reading or watching TV
- Easy on-and-off thanks to quick release, breathable headgear

Item No.	Mask Style	Headgear Style	Mask Size	Pkg.
CPAP50655	Vented	Stable Fit	One Size	1/ea
CPAP55044	Headgear Only	Stable Fit	One Size	1/ea



CPAP50655

### MiniMe® 2 Pediatric Gel Mask

- Nasal mask in four sizes, designed for children ages 2 to 12
- Headgear with 5 adjustment points
- Soft, breathable headgear for all-night comfort
- Available in vented and non-vented designs
- Adjustable cushion with memory wire for a true custom fit

Item No.	Description	Headgear Size	Mask Size	Pkg.
CPAP60216	Vented	X-Small/Small	X-Small	10/cs
CPAP602172	Vented	X-Small/Small	Small	10/cs
CPAP602182	Vented	Medium/Large	Medium	10/cs
CPAP60221	Vented	Medium/Large	Large	10/cs
CPAP60256	Non-Vented	X-Small/Small	X-Small	10/cs
CPAP602572	Non-Vented	X-Small/Small	Small	10/cs
CPAP602582	Non-Vented	Medium/Large	Medium	10/cs
CPAP60261	Non-Vented	Medium/Large	Large	10/cs
CPAP55055	Headgear	X-Small/Small	N/A	10/cs
CPAP55056	Headgear	Medium/Large	N/A	10/cs



CPAP60217



### Ascend® Masks

Patented AIRgel replaceable cushion is easily detached and reattached for easy cleaning and replacement. Leak ports are built-in, so a separate exhalation port is not needed.

- Minimizes pressure points and skin irritation
- Ergonomic cushion technology delivers an effective seal for the most comfortable fit
- Flexible spacebar with a unique bendable design provides superior stability, holding a seal without contacting the skin
- Lightweight, breathable headgear slides off easily with one hand
- 360° rotating swivel for optimal tube positioning and easier mobility



CPAP50815

### CPAP Full-Face Mask

Item No.	Description	Pkg.
CPAP50825	Starter kit with Small, Medium, Large Cushions and Headgear	1/ea
CPAP50813	Full-Face Mask, Small	1/ea
CPAP50814	Full-Face Mask, Medium	1/ea
CPAP50815	Full-Face Mask, Large	1/ea
CPAP50207	Replacement Cushion, Small	1/ea
CPAP50208	Replacement Cushion, Medium	1/ea
CPAP50209	Replacement Cushion, Large	1/ea
CPAP55051	CPAP Headgear Replacement, One Size	1/ea

### CPAP Nasal Mask

Item No.	Description	Pkg.
CPAP50174	Starter Kit with Small, Medium, Large Cushions and Headgear	1/ea
CPAP50187	Nasal Mask, Small	1/ea
CPAP50188	Nasal Mask, Medium	1/ea
CPAP50189	Nasal Mask, Large	1/ea
CPAP50201	Replacement Cushion, Small	1/ea
CPAP50202	Replacement Cushion, Medium	1/ea
CPAP50203	Replacement Cushion, Large	1/ea
CPAP55051	CPAP Headgear Replacement, One Size	1/ea



## Silicone CPAP Masks

A consistent seal without excessive tightness, delivering the prescribed therapy with enough comfort to help improve compliance.

- OmniClip and silicone forehead pad for a custom fit at the bridge of the nose for patients of all head shapes and sizes
- SilentVent elbow flushes out exhaled air quietly, without an additional inline exhalation port
- Headgear attaches with convenient clips, includes a handy wash bag to keep it clean and fresh
- Cushions, forehead pads, elbows, cones and clips are individually replaceable
- Available in four sizes, with and without headgear

Item No.	Size	Description	Pkg.
CPAP9400	Small	Full Face Mask with Head Strap	1/ea
CPAP9401	Medium	• SilentVent Anti-Asphyxia Elbow, 22 mm Male	
CPAP9402	Large	• OmniClip with Silicone Forehead Pad	
CPAP9403	X-Large	• OmniClip with Silicone Forehead Pad	
CPAP9404	Small	Full Face Mask, No Head Strap	1/ea
CPAP9405	Medium	• SilentVent Anti-Asphyxia Elbow, 22 mm Male	
CPAP9406	Large	• OmniClip with Silicone Forehead Pad	
CPAP9407	X-Large	• OmniClip with Silicone Forehead Pad	
CPAP9408	Small	Nasal Mask with Head Strap	1/ea
CPAP9409	Medium	• SilentVent Elbow, 22 mm Male	
CPAP9410	Large	• Silicone Forehead Pad	
CPAP9411	X-Large	• OmniClip with Silicone Forehead Pad	
CPAP9412	Small	Nasal Mask, No Head Gear Strap	1/ea
CPAP9413	Medium	• SilentVent Elbow, 22 mm Male	
CPAP9414	Large	• Silicone Forehead Pad	
CPAP9415	X-Large	• OmniClip with Silicone Forehead Pad	




CPAP9403



CPAP9408



CPAP9400



CPAP9408

## Aloha Nasal Pillow Systems

Medline offers the Aloha nasal pillow system for CPAP therapy patients. This lightweight, nasal pillows system is quiet and offers a stable and comfortable seal.

Item No.	Description	Pkg.
SUHCM014	Includes: One Small, Medium and Large Pillow Headgear	1/ea
SUHCMP014SP	Replacement Nasal Pillows for Aloha, Small	1/ea
SUHCMP014MP	Replacement Nasal Pillows for Aloha, Medium	1/ea
SUHCMP014LP	Replacement Nasal Pillows for Aloha, Large	1/ea
SUHHG014	Replacement Headgear for Aloha	1/ea




SUHCM014

## CPAP Accessories

### Tubing

Item No.	Description	Pkg.
SUH0610	6-ft. Tubing with Flex Connector	10/bx
SUHTUB008	8-ft. Tubing with Flex Connectors	1/ea
SUHTUB010	10-ft. Tubing with Flex Connectors	1/ea
SUHTUB0018	18-ft. Tubing with Flex Connectors	1/ea



SUH0610

### Chin Straps

Item No.	Description	Pkg.
SUHCS001	Standard	1/ea
SUHCS003	Adjustable	1/ea
SUHCS004	Deluxe	1/ea
SUHCS006	Premium	1/ea



SUHCS006

## Citrus II® Cleaning Supplies

Item No.	Description	Pkg.
SUH1001	Citrus II Spray Cleaner, 1.5-oz.	24/cs
SUH1002	Citrus II Spray Cleaner, 8-oz.	12/cs
SUH1003	Citrus II Mask Wipes	62 wipes/ dispenser, 12/cs
SUH1004	Citrus II Cleaner Concentrate, 32-oz.	1/ea
SUH1005	Citrus II CPAP Wipes	12 bx/cs
SUHCAP1003S	CPAP Cleaning Wipes, Unscented	62 wipes/ dispenser, 12/cs



## Veraseal 2 Non-Invasive Ventilation Masks

Veraseal 2 NIV Masks feature the advanced, AIRgel cushion technology, which promotes advanced skin protection, superior performance, and comfort. It is available in multiple sizes and configurations to meet your non-invasive ventilation needs.

Item No.	Description	Pkg.
CPAP50827BX	Vented with AAV, Small	10/bx
CPAP50828BX	Vented with AAV, Medium	10/bx
CPAP50829BX	Vented with AAV, Large	10/bx
CPAP50830BX	Vented with AAV, X-Large	10/bx
CPAP50937BX	Non-Vented with AAV, Small	10/bx
CPAP50938BX	Non-Vented with AAV, Medium	10/bx
CPAP50939BX	Non-Vented with AAV, Large	10/bx
CPAP50940BX	Non-Vented with AAV, X-Large	10/bx
CPAP50933BX	Non-Vented, Small	10/bx
CPAP50934BX	Non-Vented, Medium	10/bx
CPAP50935BX	Non-Vented, Large	10/bx
CPAP50936BX	Non-Vented, X-Large	10/bx



CPAP50938BX



CPAP50934BX



CPAP50938BX



CPAP50828BX



CPAP50827BX



CPAP50939BX

## Go-PAP™ CPAP Systems

The Go-PAP CPAP System offers emergency CPAP with integrated nebulization using the flow from an oxygen flowmeter or cylinder.

- Connects to barbed connector and can generate up to 90 LPM flow
- Includes 3 options for PEEP of 5.0, 7.5 and 10.0 cm H<sub>2</sub>O
- Integrated nebulization without interrupting the delivered flow to the patient



Item No.	Description	Pkg.
CPAP3134602	With BiTrac ED Mask and Standard Headgear, Adult Medium	10/cs
CPAP3134602N	With BiTrac ED Mask, Standard Headgear and Nebulizer Canister, Adult Medium	10/cs
CPAP3134602NC	With BiTrac ED Mask, Standard Headgear, Nebulizer Canister and Neb-Connect, Adult Medium	10/cs
CPAP3134603	With BiTrac ED Mask and Standard Headgear, Adult Large	10/cs
CPAP3134603N	With BiTrac ED Mask, Standard Headgear and Nebulizer Canister, Adult Large	10/cs
CPAP3134603NC	With BiTrac ED Mask, Standard Headgear, Nebulizer Canister and Neb-Connect, Adult Large	10/cs

## O<sub>2</sub>-Max™ CPAP Systems

The O<sub>2</sub>-Max CPAP System offers emergency CPAP with integrated nebulization using the flow from an oxygen flowmeter or cylinder.

- Connects to DISS port and can generate up to 140 LPM flow
- Includes 3 options for PEEP of 5.0, 7.5 and 10.0 cm H<sub>2</sub>O
- Integrated nebulization without interrupting the delivered flow to the patient



### O<sub>2</sub>-MAX CPAP System with 3-Set CPAP (5.0, 7.5 and 10.0 cm H<sub>2</sub>O)

Item No.	Description	Pkg.
CPAP3137555X	With BiTrac ED Mask, Adult Large	10/cs
CPAP3137555XN	With BiTrac ED Mask and Nebulizer Canister, Adult Large	10/cs
CPAP3137555X1	With BiTrac ED Mask and Ohmeda Quik-Connect, Adult Large	10/cs
CPAP3137555N1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Large	10/cs
CPAP3137555X2	With BiTrac ED Mask and Chemtron Quik-Connect, Adult Large	10/cs
CPAP3137556X	With BiTrac ED Mask, Adult Medium	10/cs
CPAP3137556XN	With BiTrac ED Mask and Nebulizer Canister, Adult Medium	10/cs
CPAP3137556X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs
CPAP31375561	With BiTrac ED Mask and Ohmeda Quik-Connect, Adult Medium	10/cs
CPAP3137556X2	With BiTrac ED Mask and Chemtron Quik-Connect, Adult Medium	10/cs

### O<sub>2</sub>-MAX CPAP System with 5-Set CPAP (5.0, 7.5, 10.0, 12.5 and 15 cm H<sub>2</sub>O)

Item No.	Description	Pkg.
CPAP3137554X	With BiTrac ED Mask, Adult Large	10/cs
CPAP3137554X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Large	10/cs
CPAP3137554XN	With BiTrac ED Mask and Nebulizer Canister, Adult Large	10/cs
CPAP3137557X	With BiTrac ED Mask, Adult Medium	10/cs
CPAP3137557X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs
CPAP3137557XN	With BiTrac ED Mask and Nebulizer Canister, Adult Medium	10/cs

O<sub>2</sub>-MaxCPAP Systems continues on the next page



## O<sub>2</sub>-Max™ CPAP Systems (continued)

### O<sub>2</sub>-MAX CPAP System

(with Fixed 10 cm H<sub>2</sub>O CPAP)

Item No.	Description	Pkg.
CPAP3137535X	With BiTrac ED Mask, Adult Large	10/cs
CPAP3137538X	With BiTrac ED Mask, Adult Medium	10/cs
CPAP3137538XN	With BiTrac ED Mask and Nebulizer Canister, Adult Medium	10/cs
CPAP31375381	With BiTrac ED Mask and Ohmeda Quik-Connect, Adult Medium	10/cs
CPAP3137538X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs

### O<sub>2</sub>-MAX Trio CPAP System

(Can be set to 30%, 60% or 90% Fio<sub>2</sub>)

Item No.	Description	Pkg.
CPAP3138007X	3-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Medium	10/cs
CPAP3138017X	3-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Large	10/cs
CPAP3138007X1	3-SET PEEP System with BiTrac ED Mask with Head Strap, Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs
CPAP3138107X	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Medium	10/cs
CPAP3138017X1	3-SET PEEP System with BiTrac ED Mask with Head Strap, Ohmeda Quik-Connect and Nebulizer Canister, Adult Large	10/cs
CPAP3138117X	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Large	10/cs
CPAP3138117XN	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Small Volume Nebulizer, Adult Large	10/cs
CPAP31381171	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Large	10/cs
CPAP3138008X	5-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Medium	10/cs
CPAP3138108X	5-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Medium	10/cs
CPAP3138018X	5-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Large	10/cs
CPAP3138118X	5-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Large	10/cs
CPAP3138003X	10 cm Fixed PEEP System with BiTrac ED Mask with Head Strap, Adult Medium	10/cs
CPAP3138013X	10 cm Fixed PEEP System with BiTrac ED Mask with Head Strap, Adult Large	10/cs

## O<sub>2</sub>-Max Accessories

Item No.	Description	Pkg.
CPAP3137028	Small Mask for Go-PAP and O <sub>2</sub> -Max CPAP Systems	10/cs
CPAP3137029	Medium Mask for Go-PAP and O <sub>2</sub> -Max CPAP Systems	10/cs
CPAP3137030	Large Mask for Go-PAP and O <sub>2</sub> -Max CPAP Systems	10/cs
CPAP313111111	Quik-Connect™ Ohmeda DISS	10/cs
CPAP3138000	O <sub>2</sub> -Max Trio with 6" tubing and DISS	10/cs

# O<sub>2</sub> & CO<sub>2</sub> Disposables

Southmedic OxyMask™	46–47
Cannulas, O <sub>2</sub>	48–49
Cannulas, CO <sub>2</sub>	50–51
Tubing, O <sub>2</sub>	52
Masks, O <sub>2</sub>	53
Masks, Aerosol	54
Masks, Tent & Trach	54
Bubble Humidifiers	55
O <sub>2</sub> Accessories	55

# OxyMask™

Everyone breathes easier with OxyMask

Introducing OxyMask—a unique open-style oxygen mask by Southmedic. It's proprietary design allows one mask to provide a broad range of oxygen concentrations while increasing patient comfort, reducing safety risks for patients, and helping to improve overall clinical and operational efficiencies

## OxyMask™

Open mask design that delivers oxygen flows from 1–15 liters per minute at a range of 24–90% FiO<sub>2</sub>

Item No.	Size	Universal Tubing Length	Pkg.
OM11258ML	Adult	7-ft.	25/bx
OP11258ML	Adult Plus	7-ft.	25/bx
OK11258ML	Pediatric	7-ft.	25/bx
OT10258ML*	Infant	7-ft.	25/bx
OM112514ML	Adult	14-ft.	25/bx
OP112514ML	Adult Plus	14-ft.	25/bx
OAPLUS11258**	Adult Plus OxyArm	7-ft.	25/bx

\*Infant Mask delivers from ¼ L/min- 6+ at 22%-65% FiO<sub>2</sub>)

\*\*OxyArm™ delivers 24–50% FiO<sub>2</sub>



## OxyMask™ Non-PVC

Item No.	Size	Universal Tubing Length	Pkg.
OM70258ML	Adult	7-ft.	25/bx
OP70258ML	Adult Plus	7-ft.	25/bx
OK70258ML	Pediatric	7-ft.	25/bx
OT70258ML*	Infant	7-ft.	25/bx

\*Infant Mask delivers from ¼ L/min- 6+ at 22%-65% FiO<sub>2</sub>)

### Patient Experience



Reduce risk of CO<sub>2</sub> re-breathing

Easier communication

Ability to sip fluids to take oral medications

Fewer interruptions in care

### Caregiver Efficiency



Reduced risk of errors

Only one product to learn

Standardization of care

Time saving

### Operational Efficiency



Reduce SKUs

Reduce O<sub>2</sub> consumption

Save storage space

Reduce waste

## OxyMulti-Mask™

- Open mask design that delivers oxygen flows from 1–15 liters per minute at a range of 23–83% FiO<sub>2</sub>
- Delivers oxygen through standard small bore tube using 22mm adapter or can be used with aerosol nebulizers

Item No.	Description	Universal Tubing Length	Pkg.
OHH10258ML	Adult Mask with 22 mm Adapter	7-ft.	25/bx
OPHH10258ML	Adult Plus Mask with 22 mm Adapter	7-ft.	25/bx
OKHH10258ML*	Pediatric Mask with 22 mm Adapter	7-ft.	25/bx
OAT1025ML	Adult Trach Mask	7-ft.	25/bx
OAT1025CSML	Adult Trach Mask with Comfort Neck Strap	7-ft.	25/bx

\*Pediatric mask delivers 22–93% FiO<sub>2</sub>



OHH10258ML

## OxyMask™ Aerosol

Closed mask with open bottom design and right angle swivel adapter for optimal aerosol delivery, even when patient is laying down or sitting.

Item No.	Description	Pkg.
OMN1025ML	Adult Mask with 22 mm Adapter	25/bx
OTN1025ML	Infant Mask with 22 mm Adapter	25/bx
OKN1025ML	Pediatric Mask with 22 mm Adapter	25/bx




OMN1025ML

## OxyMask™ ETCO<sub>2</sub>

Open design mask that delivers oxygen flows from 1–15 liters per minute at a range of 24–65% FiO<sub>2</sub> with end-tidal CO<sub>2</sub> monitoring.

### With Female Pigtail Adapter

Item No.	Size	Universal Tubing Length	Female Pigtail Adapter Length	Pkg.
OM21258SLMML	Adult	7-ft.	8-in.	25/bx
OP21258ML	Adult Plus	7-ft.	8-in.	25/bx
OK21258ML	Pediatric	7-ft.	8-in.	25/bx
OM212514ML	Adult	14-ft.	15-in.	25/bx



OM21258SLMML

### With Male Gas Sampling Line

Item No.	Size	Universal Tubing Length	Male Gas Sampling Line Length	Pkg.
OM21258SLM	Adult	7-ft.	8-ft.	25/bx
OP21258SLMML	Adult Plus	7-ft.	8-ft.	25/bx
OK21258SLMML	Pediatric	7-ft.	8-ft.	25/bx
OM212514SLM	Adult	14-ft.	15-ft.	25/bx

Visit  
[medline.com/  
go/oxymask](https://www.medline.com/go/oxymask)  
for more info



## SuperSoft Oxygen Cannulas

Even short term use of an oxygen cannula can lead to the development of pressure ulcers on the ears and cheeks. The Medline SuperSoft nasal cannula is made of an advanced material that is soft, flexible and non-tacky. Also features curved, non-flared nasal prongs. You will notice a difference right away:

- Maintain skin integrity – Help reduce pressure and friction on the ears and cheeks
- Promote patient compliance with added comfort – Patients are less apprehensive about wearing a cannula that feels soft and comfortable
- Not made with natural rubber latex

### Adult

Item No.	Tubing Length	Connector	Pkg.
HCS4510	No Tubing (Cannula Harness Only)	Male	50/cs
HCSS4504S	4-ft.	Standard	50/cs
HCSU4504S	4-ft.	Universal	50/cs
HCSS4514S	7-ft.	Standard	50/cs
HCSS4516S	14-ft.	Standard	50/cs
HCSU4516S	14-ft.	Universal	50/cs
HCSS4515S	25-ft.	Standard	25/cs
HCSU4515S	25-ft.	Universal	25/cs

### Pediatric

Item No.	Tubing Length	Connector	Pkg.
HCSS4518S	7-ft.	Standard	50/cs
HCSU4518S	7-ft.	Universal	50/cs

### Infant

## Soft-Touch Oxygen Cannulas

- Curved, non-flared nasal prongs
- Available in a wide variety of sizes and lengths
- Crush-resistant tubing
- Not made with natural rubber latex

### Adult

Item No.	Tubing Length	Connector	Pkg.
HCS4510S	No Tubing (Cannula Harness Only)	Standard	50/cs
HCS4504B	4-ft.	Standard	50/cs
HCSU4504B	4-ft.	Universal	50/cs
HCS4514	7-ft.	Standard	50/cs
HCS4516B	14-ft.	Standard	50/cs
HCSU4516B	14-ft.	Universal	50/cs
HCS4515B	25-ft.	Standard	25/cs
HCSU4515B	25-ft.	Universal	25/cs
HCS02KIT1*	Oxygen Set Up Kit	N/A	10/cs

\*Contains two adult soft-touch cannulas, Christmas tree adapter, empty bubble humidifier, oxygen tubing connector and 25-ft. green oxygen tubing



HCSU4517S

### Crush-Resistant Tubing

This multi-channel interior feature helps to maintain a strong airflow

Item No.	Tubing Length	Connector	Pkg.
HCSS4517S	7-ft.	Standard	50/cs
HCSU4517S	7-ft.	Universal	50/cs

### Neonatal

Item No.	Tubing Length	Connector	Pkg.
HCSS4520S	7-ft.	Standard	50/cs
HCSU4520S	7-ft.	Universal	50/cs

### Preemie

Item No.	Tubing Length	Connector	Pkg.
HCSS4521S	7-ft.	Standard	50/cs
HCSU4521S	7-ft.	Universal	50/cs

HCS4514



### Pediatric

Item No.	Tubing Length	Connector	Pkg.
HCS4518	7-ft.	Standard	50/cs
HCSU4518	7-ft.	Universal	50/cs

### Infant

Item No.	Tubing Length	Connector	Pkg.
HCS4517	7-ft.	Standard	50/cs
HCSU4517	7-ft.	Universal	50/cs

## Standard Oxygen Cannulas

- Crush-resistant tubing with standard and universal oxygen connectors
- Straight non-flared nasal prongs
- Not made with natural rubber latex

### Adult

Item No.	Tubing Length	Connector	Pkg.
HCS4511B	7-ft.	Standard	50/cs
HCSU4511B	7-ft.	Universal	50/cs
HCS4512B	14-ft.	Standard	50/cs
HCSU4512B	14-ft.	Universal	50/cs
HCS4513B	25-ft.	Standard	25/cs
HCSU4513B	25-ft.	Universal	25/cs

## High-Flow SuperSoft Oxygen Cannulas

- Maintain skin integrity—Help reduce pressure and friction on the ears and cheeks.
- Adult liter flow from 1–15 LPM
- For use with high-flow bubble humidifiers
- Not made with natural rubber latex

### Clear Tubing

Item No.	Tubing Length	Connector	Pkg.
HCS4514HF	7-ft.	Standard	25/cs
HCSS4516HF	14-ft.	Standard	50/cs
HCS4515HF	25-ft.	Standard	10/cs
HCS4514HFK*	High Flow Oxygen Kit	Standard	7/cs

### Green Tubing

Item No.	Tubing Length	Connector	Pkg.
HCS4514HFG	7-ft.	Standard	25/cs
HCS4516HFG	14-ft.	Standard	25/cs
HCS4515HFG	25-ft.	Standard	10/cs

\*Kit includes 7-ft. High-Flow Cannula and 350 mL empty hi-flow bubble humidifier

## Soft-Touch Oxygen Cannulas with Foam Ear Covers


- Pre-attached foam ear covers to reduce skin irritation and increase patient comfort
- Curved, non-flared nasal prongs
- Crush-resistant tubing
- Not made with natural rubber latex

### Clear Tubing

Item No.	Tubing Length	Connector	Pkg.
HCS4514FM	7-ft.	Standard	25/cs
HCS4516FM	14-ft.	Standard	25/cs
HCSU4516FM	14-ft.	Universal	25/cs
HCS4515FM	25-ft.	Standard	25/cs
HCSU4515FM	25-ft.	Universal	25/cs



HCS4511B



HCS4514HF



HCS4514HFG



Foam Ear Covers

HCS4514FM

## SuperSoft CO<sub>2</sub> Cannulas

Our comprehensive line of CO<sub>2</sub> nasal sampling cannulas features a "Pull N' Peel" design to resist tangling, while also providing accurate, uninterrupted CO<sub>2</sub> sampling and oxygen delivery. Universal O<sub>2</sub> connectors save you from searching for adaptors.

- Universal O<sub>2</sub> Connector is designed to fit any flowmeter without the need for adaptors
- Male and female luer connectors are available to meet caregiver needs
- Maintain skin integrity—help reduce pressure and friction on the ears and cheeks



CO<sub>2</sub> sampling through the top curved nares

Oxygen delivery through the bottom ports

### CO<sub>2</sub> Nasal Sampling Cannulas with Male Luer

Item No.	Description	Gas Sampling Line Length	UOM
HCS4560S	Adult, 7-ft.	4" Pigtail	25/cs
HCS4562S	Pediatric, 7-ft.	7-ft.	25/cs
HCS4568BS	Adult, 7-ft.	7-ft.	100/cs
HCS4568S	Adult, 7-ft.	7-ft.	25/cs
HCS4571S	Adult, 7-ft.	7-ft.	25/cs
HCS4565S	Adult, 10-ft.	10-ft.	25/cs
HCS4565BS	Adult, 10-ft.	10-ft.	100/cs
HCS4566BS	Adult, 14-ft.	14-ft.	100/cs
HCS4566S	Adult, 14-ft.	14-ft.	25/cs



22 mm O<sub>2</sub> Adaptor

Universal O<sub>2</sub> Connector



\*Standard Connector

### CO<sub>2</sub> Nasal Sampling Cannulas with Female Luer

Item No.	Description	Gas Sampling Line Length	UOM
HCS4561S	Adult, 7-ft.	4" Pigtail	25/cs
HCS4563S	Pediatric, 7-ft.	7-ft.	25/cs
HCS4568FBS	Adult, 7-ft.	7-ft.	100/cs
HCS4568FS	Adult, 7-ft.	7-ft.	25/cs
HCS4572S	Adult, 7-ft.	7-ft.	25/cs
HCS4565FS	Adult, 10-ft.	10-ft.	25/cs



Female Luer Connector



Male Luer Connector

### CO<sub>2</sub> Oral Sampling Cannulas with Male Luer

Item No.	Description	Gas Sampling Line Length	UOM
HCS4573S	Adult, 7-ft.	7-ft.	25/cs

### CO<sub>2</sub> Oral Sampling Cannulas with Female Luer

Item No.	Description	Gas Sampling Line Length	UOM
HCS4574S	Adult, 7-ft.	7-ft.	25/cs



Oral scoop sampling for mouth-breathing patients.



4-port nasal sampling allows for dual nares CO<sub>2</sub> sampling and O<sub>2</sub> delivery.

## MicroFilter CO<sub>2</sub> Cannula

- Compatible with reflective-style capnography monitors
- Proven to measure end-tidal CO<sub>2</sub> within 2 mmHg of CapnoLine® Plus<sup>1</sup>
- SuperSoft harness maintains skin integrity and increases patient comfort
- 4-port design available with or without oral sampling scoop


Accurate CO<sub>2</sub>  
sampling you  
can trust

## Oral/Nasal Sampling

Item No.	Size	Universal Tubing Length	MicroFilter Line	Pkg.
HCS4607MO	Adult	7-ft.	7-ft.	25/cs
HCS4614MO	Adult	14-ft.	14-ft.	25/cs

## Nasal Sampling

Item No.	Size	Universal Tubing Length	MicroFilter Line	Pkg.
HCS4607MN	Adult	7-ft.	7-ft.	25/cs
HCS4614MN	Adult	14-ft.	14-ft.	25/cs
HCS4707MN	Pediatric	7-ft.	7-ft.	25/cs
HCS4714MN	Pediatric	14-ft.	14-ft.	25/cs



A unique filter located inside the gas luer prevents snags and reduces moisture levels to capnography machines.

## CO<sub>2</sub> Masks

### Capnoxygen Mask by Southmedic

- Closed mask allows for oxygen delivery and end-tidal CO<sub>2</sub> monitoring

### With Female Pigtail Adapter

Item No.	Size	Universal Tubing Length	Female Pigtail Adapter Length	Pkg.
CO201ML	Adult	7-ft.	8-in.	25/bx
CO203ML	Adult Plus	7-ft.	8-in.	25/bx
CO204ML	Pediatric	7-ft.	8-in.	25/bx

### With Male Gas Sampling Line

Item No.	Size	Universal Tubing Length	Male Gas Sampling Line Length	Pkg.
CO201SLMML	Adult	7-ft.	8-ft.	25/bx
CO203SLMML	Adult Plus	7-ft.	8-ft.	25/bx



CO201SLMML

<sup>1</sup>Data is on file

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™ \*Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 05/2019-16-PM-0190a(2)-[WF#3215496]



## Oxygen Tubing

- Unique, relaxed, crush-resistant, low-memory material helps prevent kinking
- Selected products in green color to provide greater visual contrast to assist in patient safety
- For standard and high flow applications
- Not made with natural rubber latex

### Clear

Item No.	Tubing Length	Connector	Pkg
HCS4507	7-ft. (2.1 m)	Standard	50/cs
HCSU4507	7-ft. (2.1 m)	Universal	50/cs
HCS4524	14-ft. (4.3 m)	Standard	50/cs
HCSU4524	14-ft. (4.3 m)	Universal	50/cs
HCS4525	25-ft. (7.6 m)	Standard	25/cs
HCSU4525	25-ft. (7.6 m)	Universal	25/cs
HCS4540FT	40-ft. (1.0 m)	Standard	10/cs
HCS455010	50-ft. (15.2 m)	Standard	10/cs
HCSU455010	50-ft. (15.2 m)	Universal	10/cs

### Green

Item No.	Tubing Length	Connector	Pkg
HCS4507G	7-ft. (2.1 m)	Standard	50/cs
HCSU4507G	7-ft. (2.1 m)	Universal	50/cs
HCS4524G	14-ft. (4.3 m)	Standard	50/cs
HCSU4524G	14-ft. (4.3 m)	Universal	50/cs
HCS4525G	25-ft. (7.6 m)	Standard	10/cs
HCSU4525G	25-ft. (7.6 m)	Universal	25/cs
HCS4540FTG	40-ft. (1.0 m)	Standard	25/cs
HCS4550GB	50-ft. (15.2 m)	Standard	10/cs
HCSU4550GB	50-ft. (15.2 m)	Universal	10/cs



## Corrugated Aerosol Tubing

Item No.	Description	Pkg
HSK600	Corrugated, Segmented Tubing, 100-ft. (30.5 m), Opaque with Cuttable Cuff Every 6-in. (15.2 cm)	1/cs
HSK626	Corrugated, Segmented Tubing, 6-ft. (1.8 m), Opaque	50/cs
HCS4541	Corrugated Tubing, 6-in. (15.2 cm), Blue	50/cs



HCS4541



HSK600

## Oxygen Masks

Choose from a large selection of masks for all your oxygen-delivery needs. All oxygen masks feature an elongated fit for optimal patient fit.

- Soft, anatomical form for patient comfort
- A secure fit with elastic straps
- Convenient choice of universal or standard connectors
- Not made with natural rubber latex
- All oxygen masks comes with 7' of 3 channel crush-resistant clear tubing



HCS121E

### Adult

Item No.	Description	Connector	Reservoir Bag	Safety Vent	UOM
HCS4600B	Medium Concentration	Standard			50/cs
HCSU4600B	Medium Concentration	Universal			50/cs
HCS4610B	High-Concentration	Standard	●	●	50/cs
HCSU4610B	High-Concentration	Universal	●	●	50/cs
HCS4640B	Partial Non-Rebreathing	Standard	●	●	50/cs
HCSU4640B	Partial Non-Rebreathing	Universal	●	●	50/cs
HCS120E	Venturi Mask with 6 Jet Dilluters	Standard			50/cs
HCSU120E	Venturi Mask with 6 Jet Dilluters	Universal			50/cs
HCS121E	Venturi Mask with Single Jet Dilluter	Standard			50/cs
HCSU121E	Venturi Mask with Single Jet Dilluter	Universal			50/cs
HCS123E	Venturi Mask with Single Jet Dilluter and Corrugate Tubing	Standard			50/cs
HCS4670B	Total Non-Rebreathing	Standard	●		50/cs
HCSU4670B	Total Non-Rebreathing	Universal	●		50/cs
HCS4680B*	3-in-1 Mask	Standard	●	●	50/cs
HCSU4680B*	3-in-1 Mask	Universal	●	●	50/cs



HCS4610B  
High- Concentration  
(Reservoir Bag and Tubing Not Pictured)



HCSU4640B  
Partial Non-Rebreathing  
(Reservoir Bag and Tubing Not Pictured)



HCS4670B  
Total Non-Rebreathing  
(Reservoir Bag and Tubing Not Pictured)

### Pediatric

Item No.	Description	Connector	Reservoir Bag	Safety Vent	UOM
HCS4601B	Medium Concentration	Standard		●	50/cs
HCSU4601B	Medium Concentration	Universal		●	50/cs
HCS4642B	Partial Non-Rebreathing	Standard	●	●	50/cs
HCSU4642B	Partial Non-Rebreathing	Universal	●	●	50/cs
HCS4675B	Total Non-Rebreathing	Standard	●		50/cs
HCSU4675B	Total Non-Rebreathing	Universal	●		50/cs
HCS4685B*	3-in-1 Mask	Standard	●	●	50/cs
HCSU4685B*	3-in-1 Mask	Universal	●	●	50/cs



Comes with bag and tubing

\*3-in-1 masks come packaged with two 1-way valves, reservoir bag, 7' of oxygen tubing, oxygen connector, and oxygen mask



HCS4630B



HCS4631B

## Aerosol Masks

Extra patient comfort comes from an anatomical form and a soft, light and comfortable vinyl material. All aerosol masks feature an elongated fit for optimal patient fit.

- A secure fit with elastic straps
- Standard 22 mm connector
- Not made with natural rubber latex



HCS4630

### Adult

Item No.	Description	Pkg.
HCS4630B	Aerosol Mask with Straight Connector	50/cs
HCS4630	Aerosol Mask with Swivel Connector	50/cs

### Pediatric

Item No.	Description	Pkg.
HCS4631B	Aerosol Mask with Swivel Connector	50/cs

## Tent and Trach Masks

### Adult

Item No.	Description	Pkg.
HCS4632B	Face Tent with 22 mm Connector	50/cs
HCS245E	Trach Mask	50/cs

### Pediatric

Item No.	Description	Pkg.
HCS4621B	Trach Mask	50/cs
HCS4761	22 mm Aerosol Adapter, Works with HCS4632B and HCS4621B	50/cs



HCS4632B



HCS245E

## Bubble Humidifiers

Choose from a wide variety of sizes and configurations.

### Prefilled

Item No.	Description	Pkg.
HCS00350	350 mL, with Adaptor	20/cs
HCS00550	550 mL, with Adaptor	12/cs
HCS00640	650 mL, with Adaptor	10/cs



HCS00350



HCS00550

### Empty Bubble Humidifier

Item No.	Description	Pkg.
HCS64377	300 mL, 6 PSI	50/cs
VTB6900	350 mL, Hi-Flow 6 PSI	10/cs
HCS02KIT1*	Oxygen Setup Kit	10/cs
HCS4514HFK	High Flow Kit: 1 HCS4514HF, 1 VTB6900	7/cs

\*Contains two adult soft-touch cannulas, christmas tree adapter, empty bubble humidifier, oxygen tubing connector, and 25' green oxygen tubing.



VTB6900



HCS64377

Medline can customize an oxygen set up kit to your needs, contact your Medline representative

## Accessories

Medline's complete selection includes a wide variety of everyday connectors and adaptors.

Item No.	Description	Pkg.
HCS4700G	Green Nut Nipple Adapter	50/cs
HCS4700Y	Yellow Nut Nipple Adapter	50/cs
HCS0440	EZ Wrap Cushions for Cannulas	50/cs
HCS0590	Swivel Connector for Oxygen Tubing	25/cs
HCS4701	Rigid Connector for Oxygen Tubing	200/cs
AGISH749	Y Splitter Connector	5/cs
DYNJ750	Oxygen Nose Clip	100/cs
HCS1642	Pressure Line Adapter, 5–7 mm ID	50/cs
HCS2503	Disposable Water Trap	25/cs
HCS974	Aerosol Drainage Bag with "T" Adapter Hanger	50/cs
HCS974Y	Aerosol Drainage Bag with "Y" Adapter Hanger	50/cs
HCS676	15" Humidifier Connector	50/cs
HCS1790	21" Humidifier Connector	50/cs
HCS8986	Humidifier Connector with Rigid Black Elbow Connector	1/ea
HSK962A	Circuit Connector, 22 OD / 15 ID 22 OD	50/cs
WST219	Aerosol Tee Adapter	50/cs



HCS2503



HCS0590



HCS4700G



HSK962A



HCS1790



HCS974Y



# Pulse Oximeters

Fingertip Pulse Oximeters 57

Handheld Oximeters 58

Pulse Oximeter Probes  
& Accessories 59

## Fingertip Pulse Oximeters

Get quick and precise SpO<sub>2</sub> and pulse readings in seconds with our high-quality fingertip pulse oximeters. Small and lightweight, these clip sensors automatically take readings when a finger is inserted.

- Fit a wide range of finger sizes
- Long battery life and are solidly constructed for worry-free use
- All models include AAA batteries, lanyard and warranty



### pulSTAT™ Finger Pulse Oximeter

Patented multi-directional display offers 6 viewing modes.

- SpO<sub>2</sub>, pulse rate, pulse bar, plethysmograph
- Low power consumption provides 30 hours of use
- Bright, digital OLED with adjustable readout sizes
- Soft padding for patient comfort
- Hinged battery cover
- 2 year warranty

HCSM70D, 1/ea



### High-Impact Finger Pulse Oximeter

Rugged casing withstands accidental drops.

- SpO<sub>2</sub>, pulse rate and pulse bar readings
- Multi-directional display with full-color OLED screen
- Durable silicone boot construction resists cracking and chipping
- Hinged battery cover
- Non Rx
- 2 year warranty

HCSFSATSPORT, 1/ea



### Pediatric Finger Pulse Oximeter

Child-friendly frog design makes kids smile.

- SpO<sub>2</sub>, pulse rate, pulse bar, plethysmograph
- Patented multi-directional display has 6 viewing modes
- Bright, digital OLED with adjustable readout sizes
- Soft padding for patient comfort
- Hinged battery cover
- 2 year warranty

HCSM70P, 1/ea



### Soft-Touch Finger Pulse Oximeter

Soft padding enhances patient comfort

- SpO<sub>2</sub>, pulse rate, perfusion index, pulse bar, and waveforms
- Backlit OLED display with 7 viewing modes
- Adjustable brightness
- 1 year warranty

HCSM70S, 1/ea



### Finger Pulse Oximeter

Basic pulse oximeter with soft finger pad

- SpO<sub>2</sub>, pulse rate
- Bright LCD display
- 1 year warranty

HCSM70E, 1/ea

## Handheld Oximeters

Dependable and easy-to-use pulse oximeters with alarms for spot checks and continuous monitoring.

### Medline MD300

- SpO<sub>2</sub> and pulse rate
- Adjustable display: 2 readout sizes and 10 brightness levels
- Includes 1 adult finger probe
- Automatic power off
- Includes protective silicone cover

Item No.	Description	Pkg
HCSMD300	Handheld Spot Check Oximeter with 30" (70 cm) Adult Finger Probe	1/ea



### Medline MD400


- SpO<sub>2</sub>, pulse rate, pulse bar, waveform
- Adjustable audible and visual alarms
- Two display modes
- Comfortable finger clip probe
- Data transfer to computer via USB cable
- Low battery indicator

Item No.	Description	Pkg
HCSMD400	Handheld Continuous Economy Pulse Oximeter with Adult Finger Probe	1/ea



## Replacement Probes

Item No.	Description	Pkg
HCSMD400A	Adult Fingertip Style, for HCSMD400 and HCSMD300	1/ea
HCSMD400P	Pediatric Fingertip Style, for HCSMD400 and HCSMD300	1/ea
HCSMD400I	Infant Wrap Style, for HCSMD400 and HCSMD300	1/ea



## Pulse Oximeter Probes

Medline has probes for a wide variety of pulse oximeters, including handheld and tabletop models, such as the popular 7-pin Nellcor® and 9-pin BCI®.

- Probe choices for unique requirements, materials and applications
- Disposable probes in neonatal through adult sizes
- Reusable probes in pediatric and adult sizes
- Micro-Foam for easy sensor application and repositioning, Medplast for maximum adhesive for long-term use in one location, or Transpore for transparency with repositioning

### Adult

Item No.	Description	Type	Compatible With	Pkg
PCL0603104	Reusable	Finger Clip	BCI 9-PIN	1/ea
PCL0202301	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0603120	Reusable	Finger Clip	Nellcor 7-PIN	1/ea
PCL1003120	Reusable	Soft Boot	Nellcor 7-PIN	1/ea
PCL0342121	Disposable	Medplast Style	Nellcor 7-PIN	24/cs
PCL0322101	Disposable	Transpore Style	Nellcor 7-PIN	24/cs
BCI3078	Reusable	Ear	BCI 9-PIN	1/ea

### Pediatric

Item No.	Description	Type	Compatible With	Pkg
PCL0623104	Reusable	Finger Clip	BCI 9-PIN	1/ea
PCL1023104	Reusable	Soft Boot	BCI 9-PIN	1/ea
PCL0202302	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0322302	Disposable	Transpore Style	BCI 9-PIN	24/cs
PCL0623120	Reusable	Finger Clip	Nellcor 7-PIN	1/ea
PCL1023120	Reusable	Soft Boot	Nellcor 7-PIN	1/ea
PCL0202122	Disposable	Micro-Foam Style	Nellcor 7-PIN	24/cs
PCL0322122	Disposable	Transpore® Style	Nellcor 7-PIN	24/cs

### Infant

Item No.	Description	Type	Compatible With	Pkg
PCL0202303	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0202123	Disposable	Micro-Foam Style	Nellcor 7-PIN	24/cs

### Neonatal

Item No.	Description	Type	Compatible With	Pkg
PCL0202304	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0342124	Disposable	Medplast Style	Nellcor 7-PIN	24/cs
PCL0202124	Disposable	Micro-Foam Style	Nellcor 7-PIN	24/cs
PCL0703120	Disposable	Multi Site Sensor	Nellcor 7-PIN	1/ea

### Universal

Item No.	Description	Type	Compatible With	Pkg
PCL0703104	Reusable	Multi Site Sensor	BCI 9-PIN	1/ea



7 - Pin Nellcor



9 - Pin BCI



Finger Clip  
PCL0603120



Medplast  
PCL0342121



Micro-Foam  
PCL0202301



Transpore  
PCL0322101



Soft Boot  
PCL1003120



# Resuscitation

Resuscitation Bags, Self Inflating	61–67
Colorimetric CO <sub>2</sub> Detectors	68
Resuscitation Masks & Accessories	69
Mouth-to-Mask Resuscitator	69
Hyperinflation Bags	70

## Manual Resuscitators

Innovations provide clinicians with improved ease of use. They are available in infant, pediatric and adult sizes, and come in drawstring bag color-coded to indicate size.

### A Better Bag

- Delivers up to 99% FiO<sub>2</sub>
- Maneuverable – The body of the bag swivels 360° for multiple clinician hand positions, helping to reduce clinician fatigue
- Translucent – The translucent, lightweight bag body and components allow for visualization of vomitus or occlusions in the system
- Textured – Provides an enhanced, non-slip grip
- Repositionable – The patient port swivels 360°, allowing the clinician freedom to move around the patient without placing excessive torque or stress on the tracheal tube or mask
- Manometer allows for a clear visual reading of the centimeters of water pressure being delivered to the patient

### An Improved Mask

- Firmly Sealing—premium material provides a soft cushion and a firm seal

	<b>Adult (Mask Size #5)</b>	<b>Small Adult (Mask Size #5)</b>	<b>Pediatric (Mask Size #4)</b>	<b>Infant (Mask Size #2)</b>
<b>Bag Volume</b>	1,900 mL	1,500 mL	1,000 mL	300 mL
<b>Stroke Volume</b>	950–960 mL	640–660 mL	600–620 mL	140–170 mL
<b>Dead Space with Mask</b>	10 mL	10 mL	10 mL	10 mL
<b>Inspiratory Resistance</b>	<5 cm H <sub>2</sub> O	<5 cm H <sub>2</sub> O	<5 cm H <sub>2</sub> O	<5 cm H <sub>2</sub> O
<b>Exhalation Resistance</b>	<5 cm H <sub>2</sub> O 50 LPM Flow	<5 cm H <sub>2</sub> O 50 LPM Flow	<5 cm H <sub>2</sub> O 50 LPM Flow	<5 cm H <sub>2</sub> O 50 LPM Flow
<b>Attainable Delivery Pressure</b>	>80 cm H <sub>2</sub> O	>80 cm H <sub>2</sub> O	>80 cm H <sub>2</sub> O	>80 cm H <sub>2</sub> O
<b>Ventilatory Frequency</b>	>60 BPM (No Load)	>80 BPM (No Load)	>60 BPM (No Load)	>90 BPM (No Load)
<b>Body Mass Use Range</b>	>40 kg (88-lbs.)	>40 kg (88-lbs.)	<40 kg (88-lbs.)	<10 kg (22-lbs.)
<b>Minimal Tidal Volume</b>	600 mL	600 mL	500 mL	150 mL



## Resuscitation Bags

### Adult Resuscitation Bags

Item No.	Description	Mask Size
CPRM1116	Bag Reservoir, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116B	Bag Reservoir, Pre-Inflated Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116C	Bag Reservoir, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116C2	Bag Reservoir, CO <sub>2</sub> Detector, Medium and Large Adult Masks, and 7-ft. Oxygen Tubing	5, 6
CPRM1116F	Bag Reservoir, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116FC	Bag Reservoir, Filter, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116FM	Bag Reservoir, Manometer, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116FMT	Bag Reservoir, Manometer, Filter, Medium Adult Mask, and 14-ft. Oxygen Tubing	5
CPRM1116FP	Bag Reservoir, PEEP Valve, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116FPC	Bag Reservoir, PEEP Valve, Filter, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116FPM	Bag Reservoir, Manometer, PEEP Valve, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116L	Bag Reservoir, Large Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116M	Bag Reservoir, Manometer, Pre-Inflated Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116MC	Bag Reservoir, Manometer, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116P	Bag Reservoir, PEEP Valve, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116PC	Bag Reservoir, PEEP Valve, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116PCM	Bag Reservoir, Manometer, PEEP Valve, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116PD5	Bag Reservoir, PEEP Valve, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116PM	Bag Reservoir, Manometer, PEEP Valve, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116PMT	Bag Reservoir, Manometer, PEEP Valve, Medium Adult Mask, and 14-ft. Oxygen Tubing	5
CPRM1116PW	Bag Reservoir, Pre-Attached PEEP Valve, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1116T	Bag Reservoir, Medium Adult Mask, and 14-ft. Oxygen Tubing	5
CPRM1116WM	Bag Reservoir, Pre-Attached Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126	Tubing Reservoir, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126C	Tubing Reservoir, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126F	Tubing Reservoir, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126FC	Tubing Reservoir, Filter, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126FL	Tubing Reservoir, Filter, Large Adult Mask, and 7-ft. Oxygen Tubing	6
CPRM1126FM	Tubing Reservoir, Manometer, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126FP	Tubing Reservoir, PEEP Valve, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126FPC	Tubing Reservoir, PEEP Valve, Filter, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126FPM	Tubing Reservoir, Manometer, PEEP Valve, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126FPWM	Tubing Reservoir, Manometer, Pre-Attached PEEP Valve, Filter, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126FPWT	Tubing Reservoir, Pre-Attached PEEP Valve, Filter, Medium Adult Mask, and 14-ft. Oxygen Tubing	5
CPRM1126M	Tubing Reservoir, Manometer, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126P	Tubing Reservoir, PEEP Valve, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM1126PC	Tubing Reservoir, PEEP Valve, CO <sub>2</sub> Detector, Medium Adult Mask, and 7-ft. Oxygen Tubing	5

### Small Adult Resuscitation Bags

Item No.	Description	Mask Size
CPRM4416	Bag Reservoir, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM4416M	Bag Reservoir, Manometer, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM4426	Tubing Reservoir, Medium Adult Mask, and 7-ft. Oxygen Tubing	5
CPRM4426M	Tubing Reservoir, Manometer, Medium Adult Mask, and 7-ft. Oxygen Tubing	5

If you do not see your preferred configuration, **Medline can customize a resuscitator bag to your preference.**

## Manual Resuscitator Ordering Guide

### Adult Resuscitation Bags

Item No.	Reservoir Style		Tubing Length		Pressure Relief Valve	CO <sub>2</sub> Detector	CO <sub>2</sub> FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
	Bag	40" Tubing	7-ft.	14-ft.							
CPRM1116	•		•								6/cs
CPRM1116B	•		•								6/cs
CPRM1116C	•		•			•					6/cs
CPRM1116C2	•		•			•					6/cs
CPRM1116F	•		•						•		6/cs
CPRM1116FC	•		•			•			•		6/cs
CPRM1116FM	•		•						•	•	6/cs
CPRM1116FMT	•			•					•	•	6/cs
CPRM1116FP	•		•					•	•		6/cs
CPRM1116FPC	•		•			•		•	•		6/cs
CPRM1116FPM	•		•					•	•	•	6/cs
CPRM1116L	•		•								6/cs
CPRM1116M	•		•							•	6/cs
CPRM1116MC	•		•			•				•	6/cs
CPRM1116P	•		•					•			6/cs
CPRM1116PC	•		•			•		•			6/cs
CPRM1116PCM	•		•			•		•		•	6/cs
CPRM1116PD5	•		•					•			6/cs
CPRM1116PM	•		•					•		•	6/cs
CPRM1116PMT	•			•				•		•	6/cs
CPRM1116PW	•		•					•			6/cs
CPRM1116T	•			•							6/cs
CPRM1116WM	•		•							•	6/cs
CPRM1126		•	•								6/cs
CPRM1126C		•	•			•					6/cs
CPRM1126F		•	•						•		6/cs
CPRM1126FC		•	•			•			•		6/cs
CPRM1126FL		•	•						•		6/cs
CPRM1126FM		•	•						•	•	6/cs
CPRM1126FP		•	•					•	•		6/cs
CPRM1126FPC		•	•			•		•	•		6/cs
CPRM1126FPM		•	•					•	•	•	6/cs
CPRM1126FPWM		•	•					•	•	•	6/cs
CPRM1126FPWT		•		•					•		6/cs
CPRM1126M		•	•							•	6/cs
CPRM1126P		•	•					•			6/cs
CPRM1126PC		•	•			•		•			6/cs

### Small Adult Resuscitation Bags

Item No.	Reservoir Style		Tubing Length		Pressure Relief Valve	CO <sub>2</sub> Detector	CO <sub>2</sub> FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
	Bag	40" Tubing	7-ft.	14-ft.							
CPRM4416	•		•								6/cs
CPRM4416M	•		•							•	6/cs
CPRM4426		•	•								6/cs
CPRM4426M		•	•							•	6/cs

## Resuscitation Bags

### Pediatric Resuscitation Bags

Item No.	Description	Mask Size
CPRM2216	Bag Reservoir, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216C	Bag Reservoir, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216F	Bag Reservoir, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216FC	Bag Reservoir, Filter, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216FM	Bag Reservoir, Manometer, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216FM3	Bag Reservoir, Manometer, Filter, Infant, Child and Small Adult Masks, and 7-ft. Oxygen Tubing	2, 3, 4
CPRM2216FP	Bag Reservoir, PEEP Valve, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216FPC	Bag Reservoir, PEEP Valve, Filter, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216FPM	Bag Reservoir, Manometer, PEEP Valve, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216FPM3	Bag Reservoir, Manometer, PEEP Valve, Filter, Infant, Child and Small Adult Masks, and 7-ft. Oxygen Tubing	2, 3, 4
CPRM2216M	Bag Reservoir, Manometer, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216M3	Bag Reservoir, Manometer, Infant, Child and Small Adult Masks, and 7-ft. Oxygen Tubing	3, 4
CPRM2216P	Bag Reservoir, PEEP Valve, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216P3M3	Bag Reservoir, Manometer, 30mm PEEP Valve, Infant, Pediatric, and Small Adult Masks, and 7-ft. Oxygen Tubing	2, 3, 4
CPRM2216PC	Bag Reservoir, PEEP Valve, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216PCM	Bag Reservoir, Manometer, PEEP Valve, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216PM	Bag Reservoir, Manometer, PEEP Valve, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2216PMT	Bag Reservoir, Manometer, PEEP Valve, Pediatric Mask, and 14-ft. Oxygen Tubing	3
CPRM2226	Tubing Reservoir, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226B	Tubing Reservoir, Pre-Inflated Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226C	Tubing Reservoir, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226F	Tubing Reservoir, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226F2	Tubing Reservoir, Filter, Pediatric and Small Adult Masks, and 7-ft. Oxygen Tubing	3, 4
CPRM2226FC	Tubing Reservoir, Filter, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3, 4
CPRM2226FM	Tubing Reservoir, Manometer, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226FM3	Tubing Reservoir, Manometer, Filter, Infant, Child and Small Adult Masks, and 7-ft. Oxygen Tubing	2, 3, 4
CPRM2226FP	Tubing Reservoir, PEEP Valve, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226FPC	Tubing Reservoir, PEEP Valve, Filter, CO <sub>2</sub> Detector, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226FPM	Tubing Reservoir, Manometer, PEEP Valve, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226FPWM	Tubing Reservoir, Manometer, Pre-Attached PEEP Valve, Filter, Pediatric Mask, and 7-ft. Oxygen Tubing	3
CPRM2226M	Tubing Reservoir, Manometer, Pediatric Mask, and 7 ft Oxygen Tubing	3
CPRM2226P	Tubing Reservoir, PEEP Valve, Pediatric Mask, and 7 ft Oxygen Tubing	3
CPRM2226PC	Tubing Reservoir, PEEP Valve, CO <sub>2</sub> Detector, Pediatric Mask, and 7 ft Oxygen Tubing	3
CPRM2226PW	Tubing Reservoir, Pre-Attached PEEP Valve, Pediatric Mask, and 7 ft Oxygen Tubing	3

If you do not see your preferred configuration, Medline can customize a resuscitator bag to your preference.



## Manual Resuscitator Ordering Guide

### Pediatric Resuscitation Bags

Item No.	Reservoir Style		Tubing Length		Pressure Relief Valve	CO <sub>2</sub> Detector	CO <sub>2</sub> FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
	Bag	40" Tubing	7-ft.	14-ft.							
CPRM2216	•		•		40 cm H <sub>2</sub> O						6/cs
CPRM2216C	•		•		40 cm H <sub>2</sub> O	•					6/cs
CPRM2216F	•		•		40 cm H <sub>2</sub> O				•		6/cs
CPRM2216FC	•		•		40 cm H <sub>2</sub> O	•			•		6/cs
CPRM2216FM	•		•		40 cm H <sub>2</sub> O				•	•	6/cs
CPRM2216FM3	•		•		40 cm H <sub>2</sub> O				•	•	6/cs
CPRM2216FP	•		•		40 cm H <sub>2</sub> O			•	•		6/cs
CPRM2216FPC	•		•		40 cm H <sub>2</sub> O	•		•	•		6/cs
CPRM2216FPM	•		•		40 cm H <sub>2</sub> O			•	•	•	6/cs
CPRM2216FPM3	•		•		40 cm H <sub>2</sub> O			•	•	•	6/cs
CPRM2216M	•		•		40 cm H <sub>2</sub> O					•	6/cs
CPRM2216M3	•		•		40 cm H <sub>2</sub> O					•	6/cs
CPRM2216P	•		•		40 cm H <sub>2</sub> O			•			6/cs
CPRM2216P3M3	•		•		40 cm H <sub>2</sub> O			•		•	6/cs
CPRM2216PC	•		•		40 cm H <sub>2</sub> O	•		•			6/cs
CPRM2216PCM	•		•		40 cm H <sub>2</sub> O			•		•	6/cs
CPRM2216PM	•		•		40 cm H <sub>2</sub> O			•		•	6/cs
CPRM2216PMT	•			•	40 cm H <sub>2</sub> O			•		•	6/cs
CPRM2226		•	•		40 cm H <sub>2</sub> O						6/cs
CPRM2226B		•	•		40 cm H <sub>2</sub> O						6/cs
CPRM2226C		•	•		40 cm H <sub>2</sub> O	•					6/cs
CPRM2226F		•	•		40 cm H <sub>2</sub> O				•		6/cs
CPRM2226F2		•	•		40 cm H <sub>2</sub> O				•		6/cs
CPRM2226FC		•	•		40 cm H <sub>2</sub> O	•			•		6/cs
CPRM2226FM		•	•		40 cm H <sub>2</sub> O				•	•	6/cs
CPRM2226FM3		•	•		40 cm H <sub>2</sub> O				•	•	6/cs
CPRM2226FP		•	•		40 cm H <sub>2</sub> O			•	•		6/cs
CPRM2226FPC		•	•		40 cm H <sub>2</sub> O	•		•	•	•	6/cs
CPRM2226FPM		•	•		40 cm H <sub>2</sub> O			•	•	•	6/cs
CPRM2226FPWM		•	•		40 cm H <sub>2</sub> O			•	•	•	6/cs
CPRM2226M		•	•		40 cm H <sub>2</sub> O					•	6/cs
CPRM2226P		•	•		40 cm H <sub>2</sub> O			•			6/cs
CPRM2226PC		•	•		40 cm H <sub>2</sub> O	•		•			6/cs
CPRM2226PW		•	•		40 cm H <sub>2</sub> O			•			6/cs

## Resuscitation Bags

### Infant Resuscitation Bags

Item No.	Description	Mask Size
CPRM3312	Bag Reservoir, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312C	Bag Reservoir, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312F	Bag Reservoir, Filter, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312FC	Bag Reservoir, Filter, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312FM	Bag Reservoir, Manometer, Filter, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312FM123	Bag Reservoir, Manometer, Filter, Neonatal, Infant and Pediatric Masks, and 7-ft. Oxygen Tubing	1, 2, 3
CPRM3312FP	Bag Reservoir, PEEP Valve, Filter, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312FPC	Bag Reservoir, PEEP Valve, Filter, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312FPM	Bag Reservoir, Manometer, PEEP Valve, Filter, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312FPM3	Bag Reservoir, Manometer, PEEP Valve, Filter, Small Premie, Neonatal and Infant Masks, and 7-ft. Oxygen Tubing	1, 2
CPRM3312M	Bag Reservoir, Manometer, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312M3	Bag Reservoir, Manometer, Neonatal, Infant and Pediatric Masks, and 7-ft. Oxygen Tubing	2
CPRM3312P	Bag Reservoir, PEEP Valve, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312P3M3	Bag Reservoir, Manometer, 30mm PEEP Valve, Premie Pre-Inflated, Neonatal Valved, and Infant Valved Masks, and 7-ft. Oxygen Tubing	0, 1, 2
CPRM3312PC	Bag Reservoir, PEEP Valve, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312PM	Bag Reservoir, Manometer, PEEP Valve, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312PM1	Bag Reservoir, Manometer, PEEP Valve, Neonatal Mask, and 7-ft. Oxygen Tubing	2
CPRM3312PM2	Bag Reservoir, Manometer, PEEP Valve, Pre-Inflated Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3312PM3	Bag Reservoir, Manometer, PEEP Valve, Neonatal, Infant and Pediatric Masks, and 7-ft. Oxygen Tubing	1, 2, 3
CPRM3312PMC2	Bag Reservoir, Manometer, PEEP Valve, CO <sub>2</sub> Detector, Neonatal and Infant Masks, and 7-ft. Oxygen Tubing	1, 2
CPRM3312PMTD1	Bag Reservoir, Manometer, PEEP Valve, Neonatal Mask, and 14-ft. Oxygen Tubing	2
CPRM3312PMX	Bag Reservoir, Manometer, PEEP Valve, Without Mask, and 7-ft. Oxygen Tubing	2
CPRM3312PWM	Bag Reservoir, Manometer, Pre-Attached PEEP Valve, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322	Tubing Reservoir, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322C	Tubing Reservoir, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322F	Tubing Reservoir, Filter, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322FC	Tubing Reservoir, Filter, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322FM	Tubing Reservoir, Manometer, Filter, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322FP	Tubing Reservoir, PEEP Valve, Filter, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322FPC	Tubing Reservoir, PEEP Valve, Filter, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322M	Tubing Reservoir, Manometer, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322MC	Tubing Reservoir, Manometer, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322MG	Tubing Reservoir, Manometer, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322P	Tubing Reservoir, PEEP Valve, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322PC	Tubing Reservoir, PEEP Valve, CO <sub>2</sub> Detector, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3322PM	Tubing Reservoir, Manometer, PEEP Valve, Infant Mask, and 7-ft. Oxygen Tubing	2
CPRM3412P	Bag Reservoir, PEEP Valve, Premie and Infant Masks, and 7-ft. Oxygen Tubing	0, 2

If you do not see your preferred configuration, Medline can customize a resuscitator bag to your preference.

# Manual Resuscitator Ordering Guide

## Infant Resuscitation Bags

Item No.	Reservoir Style		Tubing Length		Pressure Relief Valve	CO <sub>2</sub> Detector	CO <sub>2</sub> FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
	Bag	40" Tubing	7-ft.	14-ft.							
CPRM3312	•		•		40 cm H <sub>2</sub> O						12/cs
CPRM3312C	•		•		40 cm H <sub>2</sub> O	•					12/cs
CPRM3312F	•		•		40 cm H <sub>2</sub> O				•		12/cs
CPRM3312FC	•		•		40 cm H <sub>2</sub> O				•		12/cs
CPRM3312FM	•		•		40 cm H <sub>2</sub> O				•	•	12/cs
CPRM3312FM123	•		•		40 cm H <sub>2</sub> O				•	•	12/cs
CPRM3312FP	•		•		40 cm H <sub>2</sub> O			•	•		12/cs
CPRM3312FPC	•		•		40 cm H <sub>2</sub> O	•		•	•		12/cs
CPRM3312FPM	•		•		40 cm H <sub>2</sub> O			•	•	•	12/cs
CPRM3312FPM3	•		•		40 cm H <sub>2</sub> O			•	•	•	12/cs
CPRM3312M	•		•		40 cm H <sub>2</sub> O					•	12/cs
CPRM3312M3	•		•		40 cm H <sub>2</sub> O					•	12/cs
CPRM3312P	•		•		40 cm H <sub>2</sub> O			•			12/cs
CPRM3312P3M3	•		•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3312PC	•		•		40 cm H <sub>2</sub> O	•		•			12/cs
CPRM3312PM	•		•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3312PM1	•		•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3312PM2	•		•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3312PM3	•		•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3312PMC2	•		•		40 cm H <sub>2</sub> O	•		•		•	12/cs
CPRM3312PMTD1	•			•	40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3312PMX	•		•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3312PWM	•		•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3322		•	•		40 cm H <sub>2</sub> O						12/cs
CPRM3322C		•	•		40 cm H <sub>2</sub> O	•					12/cs
CPRM3322F		•	•		40 cm H <sub>2</sub> O				•		12/cs
CPRM3322FC		•	•		40 cm H <sub>2</sub> O	•			•	•	12/cs
CPRM3322FM		•	•		40 cm H <sub>2</sub> O				•		12/cs
CPRM3322FP		•	•		40 cm H <sub>2</sub> O			•	•		12/cs
CPRM3322FPC		•	•		40 cm H <sub>2</sub> O	•		•	•		12/cs
CPRM3322M		•	•		40 cm H <sub>2</sub> O					•	12/cs
CPRM3322MC		•	•		40 cm H <sub>2</sub> O	•				•	12/cs
CPRM3322MG		•	•		40 cm H <sub>2</sub> O					•	12/cs
CPRM3322P		•	•		40 cm H <sub>2</sub> O			•			12/cs
CPRM3322PC		•	•		40 cm H <sub>2</sub> O	•		•			12/cs
CPRM3322PM		•	•		40 cm H <sub>2</sub> O			•		•	12/cs
CPRM3412P	•		•		40 cm H <sub>2</sub> O			•			12/cs



## Colorimetric CO<sub>2</sub> Detectors

Medline's colorimetric CO<sub>2</sub> detectors offer a simple way to indicate the presence of carbon dioxide during a patient's breathing cycle.

It's quick response of changes in color between blue and yellow make it easy to visualize and determine whether exhaled CO<sub>2</sub> is present.

- Our CO<sub>2</sub> detectors are long lasting and can be used for up to 24 hours
- Both CO<sub>2</sub> detectors have a small internal dead space. Item No. CPRM2400 is 16 mL and CPRM2401 is 3 mL
- Offered in two sizes to meet your needs in supporting patients weighing 1 kg or greater

Item No.	Description	Pkg.
CPRM2400	Colorimetric CO <sub>2</sub> Detector. For patients who weigh ≥10 kg (≥22-lbs.)	10/cs
CPRM2401	Colorimetric Pediatric CO <sub>2</sub> Detector. For patients who weigh 1–15 kg (2.2–33-lbs.)	10/cs



## Resuscitation Masks and Accessories

### Resuscitation Masks

Medline's resuscitation masks provide an ergonomic cushion design for an optimal seal and clear cone to visualize patient breathing.

- Removable hook ring for hand-held or head-strap application (hook rings included on sizes 3–6)
- No valve prevents mask deflation
- Available in all sizes, including preemie

Item No.	Description	Pkg.
CPRMASK11R	Neonatal, Size 1	40/cs
CPRMASK12R	Infant, Size 2	40/cs
CPRMASK13R	Child, Size 3	40/cs
CPRMASK14R	Small Adult, Size 4	40/cs
CPRMASK15R	Medium Adult, Size 5	20/cs
CPRMASK16R	Large Adult, Size 6	20/cs

Medline offers a wide variety of masks which are not shown here. Please contact us for more information.



### Resuscitation Accessories

Medline's high quality resuscitation accessories are designed to work with all Medline resuscitation devices and allow you to customize to patient need.

Item No.	Description	Pkg.
CPRMP19	19 mm Adjustable Peep Valve	20/cs
CPRMP19MRI	19 mm Adjustable Peep Valve, MRI Compatible	100/cs
CPRMP30	19 mm Adjustable Peep Valve with 30 mm Adapter	20/cs
CPRMFILTER	Expiratory Filter	50/cs
CERFLOCAP	Maxtec FloCap® CO <sub>2</sub> Detector with Flow Indicator	24/cs
CPRMCO2	CO <sub>2</sub> Detector	10/cs



CPRMP19



CPRMFILTER



CERFLOCAP

### Mouth-To-Mask Resuscitator Mask

This ready-to-use CPR mask is a highly effective one-way filtering system that shields the rescuer and victim from cross-contamination.

- A clear plastic design allows for easy placement, and the soft inflated cuff improves the fit and patient comfort
- Standard 15 mm fittings are compatible with any resuscitator
- Standard mask size accommodates children and adults.
- Easy-to-understand illustrations printed on insert
- Includes hard plastic transport case

Item No.	Description	Pkg.
HCS64184	With One-Way Valve	12/cs

Medline offers a wide variety of masks which are not shown here. Please contact us for more information.



HCS64184






## Hyperinflation Systems

Help prevent or treat lung collapse with a hyperinflation system.

- Bag design uses a “dipped” process to enhance tactile sensitivity and improve control
- Uniform inflation and the highest level of performance
- 100% compliant with ISO 5362 standards
- Variety of sizes to more accurately fit every patient



## Hyperinflation Systems

7-ft. (2.13 m) of O<sub>2</sub> tubing included\*

Item No.	Breathing Bag	O <sub>2</sub> Connector	Mask	Manometer	Pkg.
DYNJAAHY005	0.5 L	Standard	N/A		20/cs
DYNJAAHY600	0.5 L with Tail-End Bleed Valve	Universal	N/A		50/cs
DYNJAAHY30	0.5 L	Standard	1 Neonatal		20/cs
DYNJAAHY31	0.5 L	Standard	2 Infant		20/cs
DYNJAAHY29	0.5 L	Universal	2 Infant (Pre-Inflated)		20/cs
DYNJAAHY01	1 L	Standard	N/A		20/cs
DYNJAAHY601	1 L with Tail-End Bleed Valve	Universal	N/A		50/cs
DYNJAAHY01M2	1 L	Universal	2 Infant		20/cs
DYNJAAHY02	2 L	Standard	N/A		20/cs
DYNJAAHY02M3	2 L	Universal	3 Child		20/cs
DYNJAAHY03	3 L	Standard	N/A		20/cs
DYNJAAHY03M5	3 L	Universal	5 Adult		20/cs
DYNJAAHY5005	0.5 L	Standard	N/A	●	50/cs
DYNJAAHY5053P	0.5 L	Universal	N/A	●	50/cs
DYNJAAHY32	0.5 L	Standard	1 Neonatal	●	50/cs
DYNJAAHY33	0.5 L	Standard	2 Infant	●	50/cs
DYNJAAHY35	0.5 L	Universal	2 Infant	●	50/cs
DYNJAAHY34	0.5 L	Universal	2 Infant (Pre-Inflated)	●	50/cs
DYNJAAHY501	1 L	Universal	N/A	●	50/cs
DYNJAAHY501M2	1 L	Universal	2 Infant	●	50/cs
DYNJAAHY502	2 L	Universal	N/A	●	50/cs
DYNJAAHY502M3	2 L	Universal	3 Child	●	50/cs
DYNJAAHY503	3 L	Universal	N/A	●	50/cs
DYNJAAHY503M5	3 L	Universal	5 Adult	●	50/cs

\*Custom options are available



# Suction Catheters and Trach Care

Closed Suction	72-74
Open Suction Catheters	75
Open Suction Catheters, Sleeved	75
Open Suction Kits	76
Open Suction Kits with Solutions	77
Sterile Solutions	77
Trach Care & Cleaning Trays	78
Mucus/Specimen Traps	79

NEW

# ClearPro™ Closed Suction Catheter

The power to set—and follow—best practice.

ClearPro provides the performance you expect without a time-based change-out recommendation on the label. This gives you the freedom to set your own change-out policy to meet both clinical and regulatory recommendations—and drive down costs by as much as 57%.

Plus, it includes advanced features to help provide better care:

- Isolated cleaning chamber for rapid irritation
- Pressure wiper seal to help maintain PEEP

Additional Details:


- Kitted with all necessary adapters and connectors
- Does not contain DEHP
- Available in French sizes 10–16 with catheter depth markings and color-coded by ISO standards

## Endotracheal


Item No.	Description	Pkg.	Thumb Valve Color
DYNCPDS10	Double Swivel Elbow, 10 Fr	20/cs	Black
DYNCPDS12	Double Swivel Elbow, 12 Fr	20/cs	White
DYNCPDS14	Double Swivel Elbow, 14 Fr	20/cs	Green
DYNCPDS16	Double Swivel Elbow, 16 Fr	20/cs	Orange
DYNCTP10	T-Piece Adapter, 10 Fr	20/cs	Black
DYNCTP12	T-Piece Adapter, 12 Fr	20/cs	White
DYNCTP14	T-Piece Adapter, 14 Fr	20/cs	Green
DYNCTP16	T-Piece Adapter, 16 Fr	20/cs	Orange

## Tracheostomy

Item No.	Description	Pkg.	Thumb Valve Color
DYNCPDS10T	Double Swivel Elbow, 10 Fr	20/cs	Black
DYNCPDS12T	Double Swivel Elbow, 12 Fr	20/cs	White
DYNCPDS14T	Double Swivel Elbow, 14 Fr	20/cs	Green
DYNCPDS16T	Double Swivel Elbow, 16 Fr	20/cs	Orange
DYNCTP10T	T-Piece Adapter, 10 Fr	20/cs	Black
DYNCTP12T	T-Piece Adapter, 12 Fr	20/cs	White
DYNCTP14T	T-Piece Adapter, 14 Fr	20/cs	Green
DYNCTP16T	T-Piece Adapter, 16 Fr	20/cs	Orange



## Colored Thumb Valves



## Accessories



## Medline Closed Suction Catheter Systems

### Adult

- Labeled for facilities to determine change out policy
- Soft tip helps reduce damage to mucous membranes
- Transparent connector for easy observation
- High-quality, translucent EVA sleeve helps reduce cross-contamination and facilitates operation when advancing the catheter

### Endotracheal

Item No.	Description	Pkg.	Thumb Valve Color
DYNCSDS12	Double Swivel Elbow, 12 Fr	10/cs	White
DYNCSDS14	Double Swivel Elbow, 14 Fr	10/cs	Green
DYNCSTP12	T-Piece Adapter, 12 Fr	10/cs	White
DYNCSTP14	T-Piece Adapter, 14 Fr	10/cs	Green

### Tracheostomy

Item No.	Description	Pkg.	Thumb Valve Color
DYNCSDS12T	Double Swivel Elbow, 12 Fr	10/cs	White
DYNCSDS14T	Double Swivel Elbow, 14 Fr	10/cs	Green
DYNCSTP12T	T-Piece Adapter, 12 Fr	10/cs	White
DYNCSTP14T	T-Piece Adapter, 14 Fr	10/cs	Green



DYNCSDS12T



Also available with T-Piece adapter

### Pediatric

- **NEW!** Colored depth markings for clear viewing in low-light settings
- Available with an elbow or Y-style modular patient adapter
- Three adapter sizes to ensure compatibility
- Labeled for facilities to determine change out policy
- Compatible with all modes of ventilation

### Elbow Adapter

Item No.	Description	Pkg.	Thumb Valve Color
DYNCPY5	Elbow Adapter, 5 Fr	10/cs	Black
DYNCPY6	Elbow Adapter, 6 Fr	10/cs	Light Green
DYNCPY7	Elbow Adapter, 7 Fr	10/cs	Purple
DYNCPY8	Elbow Adapter, 8 Fr	10/cs	Light Blue

### Y-Connector

Item No.	Description	Pkg.	Thumb Valve Color
DYNCPY5	Y-Connector, 5 Fr	10/cs	Black
DYNCPY6	Y-Connector, 6 Fr	10/cs	Light Green
DYNCPY7	Y-Connector, 7 Fr	10/cs	Purple
DYNCPY8	Y-Connector, 8 Fr	10/cs	Light Blue



DYNCPY6



Elbow Adapter & Y-Connector products come with three adapter sizes in each kit



## Open Suction Catheters

Medline's D419 series catheters offers a full range of straight packaging.

- Available in a complete range of continuous adult and pediatric sizes
- All Medline suction catheters feature easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while minimizing trauma potential

### Single Sterile-Straight Pack

Item No.	French Size	Pkg.
DYND41906	6 Fr	100/cs
DYND41908	8 Fr	100/cs
DYND41900	10 Fr	100/cs
DYND41901	12 Fr	100/cs
DYND41902	14 Fr	100/cs
DYND41903	16 Fr	100/cs
DYND41904	18 Fr	100/cs



DYND41902




DYND41904

### Coiled

Item No.	French Size	Pkg.
DYND41906C	6 Fr	50/cs
DYND41908C	8 Fr	50/cs
DYND41900C	10 Fr	50/cs
DYND41901C	12 Fr	50/cs
DYND41902C	14 Fr	50/cs
DYND41903C	16 Fr	50/cs
DYND41904C	18 Fr	50/cs




DYND41902C




### Poly-Cath

Item No.	French Size	Pkg.
DYND40990	10 Fr	100/cs
DYND40992	14 Fr	100/cs



DYND40992



DYND41962


### Continuous Valve

Item No.	French Size	Pkg.
DYND41960	10 Fr	100/cs
DYND41962	14 Fr	100/cs
DYND41964	18 Fr	100/cs

## Sleeved Suction Catheters

Our D407 series open-suction catheters are contained in a sleeve to help avoid cross-contamination.

Item No.	French Size	Pkg.
DYND40708F	8 Fr	50/cs
DYND40700F	10 Fr	50/cs
DYND40701F	12 Fr	50/cs
DYND40702F	14 Fr	50/cs



DYND40708F

## Kits and Trays

Medline's D409 series kits and trays provide all the necessities for airway clearance.

- Irrigation solutions can be poured into the tray or the pop-up cup
- All kits and trays contain Medline Aloetouch® 3G gloves
- Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while minimizing trauma potential

### Suction Catheter Mini Tray with Gloves

Item No.	French Size	Pkg.
DYND40986	6 Fr	100/cs
DYND40988	8 Fr	100/cs
DYND40980	10 Fr	100/cs
DYND40981	12 Fr	100/cs
DYND40982	14 Fr	100/cs
DYND40983	16 Fr	100/cs
DYND40984	18 Fr	100/cs

### Suction Catheter Pack with Gloves and Pop-up cup

Item No.	French Size	Pkg.
DYND40976	6 Fr	50/cs
DYND40978	8 Fr	50/cs
DYND40970	10 Fr	50/cs
DYND40971	12 Fr	50/cs
DYND40972	14 Fr	50/cs
DYND40973	16 Fr	50/cs
DYND40974	18 Fr	50/cs

### Suction Catheter Mini Tray with Gloves, Poly-cath

Item No.	French Size	Pkg.
DYND48980A	10 Fr	100/cs
DYND48982	14 Fr	100/cs

Medline manufactures over 125,000 sterile procedure kits every day.



DYND40982



DYND40973



DYND48982

## Suction Catheter Trays with Sterile Solutions

Our D414 series suction catheters come with everything needed for airway clearance, including a choice of irrigation fluids.

- All kits and trays contain Medline Aloetouch® 3G gloves
- Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while minimizing trauma potential

### Water

Item No.	Description	Pkg.
DYND41476	Suction Catheter Tray with Sterile Water, 100 mL; 6 Fr	36/cs
DYND41478	Suction Catheter Tray with Sterile Water, 100 mL; 8 Fr	36/cs
DYND41470	Suction Catheter Tray with Sterile Water, 100 mL; 10 Fr	36/cs
DYND41471	Suction Catheter Tray with Sterile Water, 100 mL; 12 Fr	36/cs
DYND41472	Suction Catheter Tray with Sterile Water, 100 mL; 14 Fr	36/cs
DYND41474	Suction Catheter Tray with Sterile Water, 100 mL; 18 Fr	36/cs

### Saline

Item No.	Description	Pkg.
DYND41478SA	Suction Catheter Tray with Sterile Saline, 100 mL; 8 Fr	36/cs
DYND41470SA	Suction Catheter Tray with Sterile Saline, 100 mL; 10 Fr	36/cs
DYND41442	Suction Catheter Tray with Sterile Saline, 100 mL; 14 Fr	36/cs
DYND40442	Suction Catheter Tray with Saline and Drape; 14 Fr	36/cs

## Sterile Solutions

These convenient, ready-to-use containers reduce the contamination risks and potential waste associated with bulk containers of sterile water and saline.

- Ideal for device irrigation, suction, dressing moistening and wound debridement
- Cup size provides the proper amount for suctioning

### Water

Item No.	Description	Pkg.
RDI30295	Sterile Water, 100 mL	48/cs
DYND40570	Sterile Water, 110 mL	48/cs
PCS1550	Sterile Water, 250 mL	24/cs
DYND500MLW	Sterile Water, 500 mL	18/cs
DYND1000MLW	Sterile Water, 1,000 mL	6/cs

### Saline

Item No.	Description	Pkg.
RDI30296	Sterile Saline, 0.9%, 100 mL	48/cs
DYND40540	Sterile Saline, 0.9%, 110 mL	48/cs
PCS1650	Sterile Saline, 0.9%, 250 mL	24/cs
DYND500MLS	Sterile Saline, 0.9%, 500 mL	18/cs
DYND1000MLS	Sterile Saline, 0.9%, 1,000 mL	6/cs



DYND41472



DYND40442



PCS1650

DYND1000MLS

RD130296

DYND500MLS

## Tracheostomy Care and Cleaning Trays

Medline offers a wide variety of tracheostomy care and cleaning trays in standard and customized configurations to meet diverse clinical needs.

- Available in economy, basic and full-compartment kits
- May be customized

### Did You Know?

We can build a custom kit to your exact specifications. Medline is a leading manufacturer of custom procedure kits.

	DYND40582	DYND40622	DYND40610	DYND40620	DYND40589	DYND40615	DYND40580	DYND40650	DYND40511	DYND40594	CC3T3030A	CC3T4691A	CC4682A	CC4681A
Tray Compartments	3	1	1	1	3	3	3	3	3	1	3	3	3	3
AloeTouch Gloves	2		2	2	2	2	2	2	2		2			2
Drape	1		1	1	1	1		1	1		1	1	1	1
Pipe Cleaners		2	2						3	3	3	3	3	3
Non-woven Gauze 4 x 4	4	2	2	2	2	2	2	1	3	4	1	3	3	3
Cotton Tip Applicator	2	2	5	5	2	5	2	4	2	2	2	3	3	3
Pliable Cleaning Brush	1	1	1		1	1	1	1	1	1	1	1	1	1
Twill Tape 34"	1	1	1	1	1	1	1	1		1	1	1	1	1
6-ply Trach Dressing	1		1	1	1	1	1	1	1	1	1	1	1	1
Hydrogen Peroxide					2	2	1		1					
Saline 100 mL					1			1	1					
Water 100 mL							1							
Suction Catheter 14 fr							1		1					
Forceps												1	1	
Pkg	20/cs	30/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	30/cs	20/cs	20/cs	20/cs	20/cs



## Mucus/Specimen Traps

Our mucus and specimen traps have leak-proof transport caps for more secure specimen collection.


- Sterile plastic vials, graduated in 1cc increments
- Includes threaded cap, tubing not made from natural rubber latex and patient ID label

### Sterile DeLee Mucus Trap

Item No.	Description	Pkg.
DYND44108	Trap with Catheter and Cap, Capacity 20 cc; 8 Fr	50/cs
DYND44110	Trap with Catheter and Cap, Capacity 20 cc; 10 Fr	50/cs

### Sterile Mucus Specimen Trap

Item No.	Description	Pkg.
DYND44140	Trap With Transport Cap, Capacity; 40 cc	50/cs
DYND44180	Trap With Transport Cap, Capacity; 80 cc	50/cs





# Top Respiratory Products for Seasonal Flu Demand

We've stocked up on the top respiratory products for seasonal flu demand.



## Cannulas

Item No.	Description
<b>Soft-Touch Nasal Cannula</b>	
HCS4514	Soft-Touch Nasal Cannula with 7' Tubing and Standard Connector, Adult
HCS4518	Soft-Touch Nasal Cannula with 7' Tubing and Standard Connector, Pediatric
HCSU4518	Soft-Touch Nasal Cannula with 7' Tubing and Universal Connector, Pediatric
<b>SuperSoft Nasal Cannula</b>	
HCSS4514S	SuperSoft Nasal Cannula with 7' Tubing and Standard Connector, Adult
HCSU4518S	SuperSoft Nasal Cannula with 7' Tubing and Universal Connector, Pediatric
<b>High Flow SuperSoft Nasal Cannula</b>	
HCS4514HFG	High Flow SuperSoft Green Nasal Cannula, 7' Tubing and Standard Connector, Adult



HCS4514HFG

## Masks

Item No.	Description
<b>OxyMask™</b>	
OM11258ML	OxyMask with 7' Universal Tubing, Adult
OP11258ML	OxyMask with 7' Universal Tubing, Adult Plus
OK11258ML	OxyMask with 7' Universal Tubing, Pediatric
OT10258ML	OxyMask with 7' Universal Tubing, Infant
<b>Medium Concentration Mask</b>	
HCS4600B	Medium Concentration Mask with 7' Tubing and Standard Connector, Adult
HCSU4600B	Medium Concentration Mask with 7' Tubing and Universal Connector, Adult
HCS4601B	Medium Concentration Mask with 7' Tubing and Standard Connector, Pediatric
HCSU4601B	Medium Concentration Mask with 7' Tubing and Universal Connector, Pediatric
<b>Standard Non-Rebreather Mask with Safety Vent</b>	
HCS4640B	Standard Non-Rebreather Mask with 7' Tubing and Standard Connector, Adult
HCSU4640B	Standard Non-Rebreather Mask with 7' Tubing and Universal Connector, Adult
HCS4642B	Standard Non-Rebreather Mask with 7' Tubing and Standard Connector, Pediatric
HCSU4642B	Standard Non-Rebreather Mask with 7' Tubing and Universal Connector, Pediatric
<b>Oxygen Mask 3-in-1</b>	
HCS4680B	Oxygen Mask 3-in-1 with 7' Tubing and Standard Connector, Adult
HCSU4680B	Oxygen Mask 3-in-1 with 7' Tubing and Universal Connector, Adult
HCS4685B	Oxygen Mask 3-in-1 with 7' Tubing and Standard Connector, Pediatric
HCSU4685B	Oxygen Mask 3-in-1 with 7' Tubing and Universal Connector, Pediatric
<b>Elongated Aerosol Mask</b>	
HCS4630	Elongated Aerosol Mask with Swivel Connector, Adult
HCS4631B	Elongated Aerosol Mask with Swivel Connector, Pediatric



OM11258ML

## Nebulizers

Item No.	Description
<b>Nebulizer Kit with Tee Mouthpiece</b>	
HCS4482	Nebulizer Kit with Tee Mouthpiece, 7' Tubing and Standard Connector
HCSU4482	Nebulizer Kit with Tee Mouthpiece, 7' Tubing and Universal Connector
HCS4483	Nebulizer Kit with Tee Mouthpiece, 6" Reservoir, 7' Tubing and Standard Connector
HCSU4483	Nebulizer Kit with Tee Mouthpiece, 6" Reservoir, 7' Tubing and Universal Connector
<b>Nebulizer Kit with Mask</b>	
HCS4485	Nebulizer Kit with Mask, 7' Tubing and Standard Connector, Adult
HCSU4485	Nebulizer Kit with Mask, 7' Tubing and Universal Connector, Adult
HCS4486	Nebulizer Kit with Mask, 7' Tubing and Standard Connector, Pediatric
HCSU4486	Nebulizer Kit with Mask, 7' Tubing and Universal Connector, Pediatric
<b>BAN Nebulizer Kit</b>	
HCS4433BAN	BAN Nebulizer Kit, Tee Piece with One-Way Valve, 7' Tubing and Standard Connector
<b>Nebulizer Compressor</b>	
HCS70004	Nebulizer Compressor with Disposable Nebulizer Medicine Cup, 7' Tubing and Mouthpiece



HCS4485

## Incentive Spirometers

Item No.	Description
HCS93250	Incentive Spirometer, 2,500 mL
HCS93500	Incentive Spirometer, 5,000 mL



HCS93250

## MDI Spacers

Item No.	Description
<b>OptiChamber Diamond MDI Spacers</b>	
1098213ML	OptiChamber Diamond MDI Spacer
1098214ML	OptiChamber Diamond MDI Spacer with Small Mask
1098215ML	OptiChamber Diamond MDI Spacer with Medium Mask
1098216ML	OptiChamber Diamond MDI Spacer with Large Mask
<b>MDI Spacers</b>	
HCSA100	MDI Spacer
HCSA112	MDI Spacer with Infant Mask
HCSA113	MDI Spacer with Pediatric Mask
HCSA114	MDI Spacer with Adult Mask



1098213ML

## HMEs/Filters

Item No.	Description
DYNJAABV1	Bacterial/Viral Filter
DYNJAAHEPA	Hepa Filter
DYNJAAHME5	Large HME, Filtered, No Port
DYNJAAHME12	Large HME, Non-Filtered, No Port



DYNJAAHME5

## Prefilled Bubble Humidifiers

Item No.	Description
HCS00350	Prefilled Bubble Humidifier, 350 mL
HCS00550	Prefilled Bubble Humidifier, 550 mL



HCS00350

Retail packaging available on a variety of products,  
ask your Medline Representative for more details.



**A**

Adapters  
 Aerosol, 54  
 MDI, 32  
 Valved Tee, 32  
 Aeromist BAN Nebulizer, 27  
 Aerosol  
 Masks, 33, 47, 54  
 Tubing, 52  
 Air Flowmeters, 13  
 Airway Management, 4–9  
 Aspirator, Suction, 18  
 Asthma Management, 26–35

**B**

Bacterial Filters, 20  
 Bags, Resuscitation, 62–67  
 Breathing Circuits,  
 Non-Heated, 24  
 Bronchial Hygiene, 10–11  
 Bubble Humidifiers, 55, 81  
 Bypass Heat and Moisture  
 Exchanger, 21

**C**

Canister Kits, Suction, 18  
 Cannulas, 8, 80  
 CO<sub>2</sub>, 50–51  
 Oxygen, 48–49  
 Capnoxygen Masks, 51  
 Catheters, Suction  
 ClearPro, 72–73  
 Closed, 72–74  
 Open, 75–77  
 Kits, 76–77  
 Chin Straps, CPAP, 41  
 Circuit Connector, Ventilator, 8  
 Clips, Nose, 25  
 CO<sub>2</sub>  
 Cannulas, 50–51  
 Detectors  
 Colorimetric, 68  
 FloCap, 69  
 Disposables, 46–55  
 Masks, 47, 51  
 Colorimetric CO<sub>2</sub> Detectors, 68  
 Compressors, Nebulizer, 34–35  
 Portable, 34  
 Conserver, Oxygen, 15  
 Corrugated Aerosol Tubing, 52  
 CPAP  
 Accessories, 41  
 Cleaning Supplies, 41  
 DreamPort Mask  
 Solution, 37  
 Masks, 37–40  
 Full-Face, 39  
 Nasal, 39  
 Go-PAP, O<sub>2</sub>-Max Nasal Pillow  
 Systems, 41  
 Systems, 43–44  
 CPR Mask, 69  
 Cuff Pressure Controller, 7  
 Cylinder Wrench, Oxygen, 15

**D**

Detectors, CO<sub>2</sub>  
 Colorimetric, 68

FloCap, 69  
 Disposable Canister Kits, 18  
 DreamPort CPAP Mask  
 Solution, 37

**E**

Endotracheal  
 Holder, 5  
 Suction Catheters, 72, 74  
 Tubes, Subglottic, 6  
 Exchangers, Heat and  
 Moisture, 22–23  
 Bypass, 21  
 Filtered, 23  
 Trach, 22

**F**

Filters, 81  
 Bacterial/Viral, 20  
 Heat and Moisture  
 Exchangers, 22–23  
 Pulmonary Function  
 Testing (PFT), 25  
 Filtration, 19–25  
 Finger Pulse Oximeters, 57  
 FloCap CO<sub>2</sub> Detector, 69  
 Flowmeters, Air & Oxygen 13  
 Flow meters, Peak 31  
 Foam Ear Cover Cannulas, 49

**G**

Gas  
 Delivery, 12–18  
 Hose Assemblies, 14  
 Go-PAP CPAP Systems, 43

**H**

Handheld Oximeters, 58  
 Heat and Moisture  
 Exchangers, 22–23  
 Bypass, 21  
 Filtered, 22  
 Trach, 22  
 HMEs, 81  
 High Flow Oxygen Cannulas, 49  
 Hose Assemblies, Gas, 14  
 Humidification, 19–25  
 Humidifiers, Bubble, 55  
 Prefilled, 81  
 Hyperinflation Systems, 70

**I**

Incentive Spirometers, 11, 81  
 IntelliCuff Pressure Monitor, 7  
 Inner Cannula, 8  
 Interfaces, Non-Invasive, 36–44

**M**

Manual Resuscitators, 61–67  
 Masks, 30, 39, 80  
 Aerosol, 33, 54  
 CO<sub>2</sub>, Capnoxygen, 51  
 CPAP, 38–40  
 Solution, 37  
 Gel, Full-Face, 38  
 Mouth-to-Mask  
 Resuscitator, 69

Nasal, 38–39  
 Non-Invasive Ventilation, 42  
 Oxygen, 46–47, 53  
 Resuscitation, 69  
 Southmedic OxyMask,  
 46–47  
 Face Tent, 54  
 Trach, 54  
 mBrace ET Tube Holder, 5  
 MDI  
 Adapters, 32  
 Spacers and Masks, 31, 81  
 Medication Delivery, 26–35  
 MicroFilter CO<sub>2</sub> Cannulas, 51  
 Mouthpiece, Disposable, 25  
 Mouth-to-Mask Resuscitator  
 Mask, 69  
 MRI-Conditional Transport  
 Ventilator Circuit, 24  
 Mucus Traps, 79

**N**

Nasal Pillow Systems, 41  
 Nebulizers, 27–29, 81  
 Compressors, 34–35  
 Handheld, 27–29, 81  
 Non-Heated Breathing  
 Circuits, 24  
 Non-Invasive  
 Interfaces, 36–44  
 Ventilation Masks, 42  
 Nose Clips, 25

**O**

O<sub>2</sub> Disposables, 46–55  
 O<sub>2</sub>-Max CPAP Systems, 43–44  
 Open Suction Catheters, 75  
 OptiChamber Diamond, 30  
 Oximeters  
 Finger Pulse, 57  
 Handheld, 58  
 Pulse, 56–59  
 Oxygen  
 Accessories, 55  
 Cannulas, 48–49  
 Conserver, 15  
 Cylinder Wrench, 15  
 Flowmeters, 13  
 Masks, 46–47, 53  
 Regulators, 15  
 Tubing, 52

**P**

Peak Flow Meter, 31  
 PEEP Valves, 69  
 Planogram, 82  
 Pneumatic Conserver, 15  
 Portable Nebulizer  
 Compressor, 34  
 Probes  
 Oximeter, 59  
 Replacement, 58  
 Pulmonary Function Testing  
 (PFT) Filters, 25  
 Pulse Oximeters, 57–59  
 Finger, 57  
 Probes, 59

**R**

Regulators  
 Oxygen, 15  
 Suction, 16–17  
 Respiratory Flu  
 Code List, 80–81  
 Resuscitation, 60–70  
 Accessories, 69  
 Bags, 62–67  
 Masks, 69  
 Resuscitators  
 Manual, 61–67  
 Mouth-To-Mask, 69

**S**

SideStream Nebulizer, 28  
 Single-Lumen Pneumatic  
 Conserver, 15  
 Sleeved Suction Catheters, 75  
 Solutions, Sterile, 77  
 Southmedic OxyMask, 46–47  
 Spacers, MDI, 31, 81  
 Specimen Traps, 79  
 Spirometers, Incentive, 11, 81  
 Sterile Solutions, 77  
 SubG Endotracheal Tubes, 6  
 Suction  
 Aspirator, 18  
 Canister Kits, 18  
 Catheters, 71–79  
 Devices, 12–18  
 Regulators, 16–17

**T**

Tent Masks, 54  
 Ties, Tracheostomy, 9  
 Tracheostomy, 72, 74  
 Care & Cleaning Trays, 78  
 Trach HME, 22  
 Masks, 54  
 Ties, 9  
 Traps  
 Mucus, 79  
 Specimen, 79  
 Vacuum Trap, 17  
 Tube, Endotracheal  
 Holder, 5  
 SubGlottic, 6  
 Tubing  
 Corrugated Aerosol, 52  
 CPAP, 41  
 Holder, 5  
 Oxygen, 52

**V**

Vacuum Traps, 17  
 Valved Tee Adapters, 32  
 Ventilation Masks, 42  
 Ventilator Circuit, 24  
 Connector, 8  
 Viral Filters, 20

**W**

Washer, Brass Yoke, 15  
 Wrench, Oxygen Cylinder, 15



**Medline Industries, Inc.** Three Lakes Drive, Northfield, IL 60093

**Medline United States**  
**1-800-MEDLINE (633-5463)**  
medline.com | info@medline.com

**Medline Canada**  
**1-800-396-6996**  
medline.ca | canada@medline.com

**Medline México**  
**01-800-831-0898**  
medlinemexico.com | mexico@medline.com

FOLLOW US    

REFERENCES 1. Black JM, Cuddigan JE, Walko MA, Diddier LA, Lander MJ, Kelp MR. Medical device related pressure ulcers in hospitalized patients. International Wound Journal. 2010; 7(5):358-365. 2. American Lung Association website. Measuring Your Peak Flow Rate. Available at: <http://www.lung.org/lung-health-and-diseases/lung-disease-lookup/asthma/living-with-asthma/managing-asthma/measuring-your-peak-flow-rate.html>. Accessed September 27, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2020 Medline Industries, Inc. IntelliCuff is the trademark of Hamilton Medical. Respironics and Sidestream are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. LiteTouch Respironics and LiteTouch are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. OptiChamber Respironics and OptiChamber are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. Respironics and Personal Best are trademarks of Koninklijke Philips N.V. and its affiliates, all rights reserved. Bleep, Dreamports, DreamWay, and DreamPreps are trademarks of Bleep LLC. AIRgel, Ascend, VeraSeal 2, Mojo, MiniMe, MiniMe 2, and Innova are trademarks of Sleepnet. AIRgel, Ascend, VeraSeal 2, Mojo, MiniMe, MiniMe 2, and Innova are trademarks of Sleepnet. Citrus II CPAP Mask Wipes is a trademark of Beaumont Products Inc. Go-PAP and O<sub>2</sub>-Max are trademarks of Pulmodyne. FloCap is a trademark of Sunmed Medical. The OxyMask trademark is the property of Southmedic Inc., and is used under license by Medline Industries, Inc. MKT19W734883 / LIT085 / 55