

Excel 2000 | STANDARD

- Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place
- Carbon steel frame with rust- and chip-resistant chrome plating
- Padded and upholstered armrests and calf pads are standard
- Threaded seat screw inserts reduce stripping from upholstery changes
- Dual axle hemi-height adjustable
- Chart pocket on back
- Easy-to-clean vinyl upholstery
- Weighs 38 lbs (17 kg)

ACCESSORIES

MDS85181U

O₂ Holder

MDS85196

Anti-Fold Device

MDS85183

I.V. Pole

MDS85197A

Overhead Anti-Theft Device

MDS85190

O₂ / I.V. Combo

MDS85189

Anti-Tip Device

WCA806991

Tinnerman Legrest Locks

RBY (Ruby)

NAVY (Navy)

D (Black)



MDS806300D

CHAIR	ARM STYLE	LEG STYLE	SEAT W X D INCHES (CM)	HCPCS	WT. CAP. LB (KG)
MDS806100*	Full Length Permanent	Permanent Footrest	18 x 16 (46 x 41)	K0001 or K0002	300 (136)
MDS806150*	Full Length Permanent	Detachable Footrest	18 x 16 (46 x 41)	K0001 or K0002	300 (136)
MDS806200*	Full Length Permanent	Detachable Elev. Legrest	18 x 16 (46 x 41)	K0001 or K0002 K0195 or E0990	300 (136)
MDS806250*	Desk Length Removable	Detachable Footrest	18 x 16 (46 x 41)	K0001 or K0002	300 (136)
MDS806250DFLA**	Full Length Removable	Detachable Footrest	18 x 16 (46 x 41)	K0001 or K0002	300 (136)
MDS806300*	Desk Length Removable	Detachable Elev. Legrest	18 x 16 (46 x 41)	K0001 or K0002	300 (136)
MDS806300DFLA**	Full Length Removable	Detachable Elev. Legrest	18 x 16 (46 x 41)	K0001 or K0002	300 (136)
MDS806400**	Desk Length Removable	Detachable Footrest	20 x 16 (51 x 41)	K0001 or K0002 E2201	300 (136)
MDS806450**	Desk Length Removable	Detachable Elev. Legrest	20 x 16 (51 x 41)	K0001 E2201 K0195	300 (136)

*To order, add the letter color code shown above designating color to the end of the item number. **Available only in standard black upholstery.