



DIAGNOSTIC PRODUCTS



Equipment and Testing Supplies
for Enhanced Detection of Disease

Simple Solutions. Significant Results.

We take pride in providing reliable, accurate and durable products that aid in earlier diagnosis and prevention of disease. Review the catalog to see the latest additions to our popular line of diagnostic instruments and equipment, as well as our testing and laboratory products.

 Diagnostic Instruments	3-8
 Patient Monitors	9-13
 EMS	14-15
 Blood Pressure Cuffs & Connectors	16-20
 Blood Pressure Monitors & Accessories	21-26
 Thermometers	27-30
 Stethoscopes	31-33
 Defibrillators	34
 Electrodes	35-40
 Point-of-Care Testing	41-48
 Medication Delivery & Disposal	49-55

Diagnostic Wall System


Our fully modular system can be customized to your needs and expanded or upgraded at any time.

- » Non-slip, ergonomically designed handles with patented technology for 100% light output with variable light adjustment and automatic safety shut-down after 2 minutes
- » LED or xenon illumination options for instrument heads, with a new **EliteVue macro-otoscope that exceeds the magnification of the market leader¹**
- » Wall systems available with and without premounted gray plastic mounting board that measures 31" x 16" x 1/4"

Customizable
to fit your
needs



Corrective lenses		L2
D+	1-10 12 15 20 40	
D-	1-10 15 20 25 30 35	



Model Featured:
MDS36533018B

Diagnostic Wall System

Components chart

	Wall Transformer	Otosopes				Ophthalmoscopes		Thermometer		Specula Dispenser	Blood Pressure	
Item Number	MDS3653300	MDS10565301 L2 LED	MDS10580301 L2 Xenon	MDS10512301E LED EliteVue Otoscope	MDS10513301E Xenon EliteVue Otoscope	MDS10571301 L2 Xenon	MDS10573301 L3 Xenon	MDS3656301 Tympanic Therm	MDS3658 Covidien Oral Predictive Therm	MDS3654	MDSL3655109 Aneroid	MDS1784003 Automatic
MDS36533001	•	•				•						
MDS36533002	•		•			•						
MDS36533003	•	•				•				•		
MDS36533004	•		•			•				•		
MDS36533005	•	•				•				•	•	
MDS36533006	•		•			•				•	•	
MDS36533007	•	•					•			•	•	
MDS36533008	•	•				•		•		•	•	
MDS36533009	•	•				•		•			•	
MDS36533010	•	•					•	•		•	•	
MDS36533011	•	•				•				•		•
MDS36533012	•		•			•				•		•
MDS36533013	•	•					•			•		•
MDS36533014	•	•				•			•	•	•	
MDS36533015	•	•				•			•	•		•
MDS36533016	•	•					•		•	•	•	
MDS36533017	•	•					•		•	•		•
MDS36533018	•				•	•			•	•	•	
MDS36533018B	•				•	•			•	•	•	
MDS36533019	•				•	•				•	•	
MDS36533019B	•				•	•				•	•	
MDS36533020	•			•		•			•	•	•	
MDS36533020B	•			•		•			•	•	•	
MDS36533021	•			•		•				•	•	
MDS36533021B	•			•		•				•	•	
MDS36533022	•				•	•						
MDS36533023	•				•	•				•		

■ These models come pre-mounted to a wallboard.

Get any other model above mounted to a wallboard by adding MDS12717 to your order.

Choose between LED and xenon illumination in your otoscopes

LED Illumination

- » Uniform illumination of the examination field in pure white light
- » LED illumination shows colors as they really are, which is helpful for efficient diagnosis
- » 100,000 lifetime hours per bulb, meaning less replacements and lower replacement cost over time

Xenon Illumination (XL)

- » Xenon illumination has a yellower tint to the lighting
- » Offers between 20,000 and 30,000 hours per bulb

Don't see what you're looking for?
Customize your own wall system at
<https://diagnostics.medline.com>



MDS3725

ri-scope Ophthalmoscopes

- » Available with rechargeable lithium-ion batteries, lithium batteries or C cell alkaline batteries
- » LED or Xenon (XL) lighting options
- » Filter wheel engageable for all apertures with symbol display, red-free filter, blue filter and polarization filter
- » Fully enclosed optics prevent any dust particles from entering the ophthalmoscope head

L2 Ophthalmoscope

- » Diopter wheel with 29 corrective lenses Plus 1-10, 12, 15, 20, 40; Minus 1-10, 15, 20, 25, 30, 35
- » Easy-to-operate aperture hand-wheel with semi-circle, small/medium/large circle, fixation star, and slit

Item No.	Description	Pkg.
MDS3725	3.5V XL for Lithium Batteries, C-handle	1/ea
MDS3725550	3.5V XL for Lithium-Ion Batteries, C-handle	1/ea

L3 Ophthalmoscope

- » Diopter wheel with 89 corrective lenses—Plus 1-45 in single steps, Minus 1-44 in single steps
- » Easy-to-operate aperture hand-wheel with semi-circle, small/medium/large circle, fixation star, slit and grid
- » Focusing device with zero retention force for the quick correction of ametropias

Item No.	Description	Pkg.
MDS3726	3.5V XL for Lithium Batteries, C-handle	1/ea
MDS3726550	3.5V XL for Lithium-Ion Batteries, C-handle	1/ea



MDS3703

Otoscopes

- » Available with rechargeable lithium-ion batteries, lithium batteries or C cell alkaline batteries
- » Fiber optic LED or Xenon (XL) lighting options available for optimal light transmission
- » Bi-directional swiveling, sealing optical glass with 3x magnification (L2 and L3)
- » Suitable for pneumatic tests (supplied without ball)

L2 Otoscope

- » Optional operation lens with 4x magnification available

Item No.	Description	Pkg.
MDS3716	2.5V XL, Xenon, for Alkaline Batteries, C-handle	1/ea
MDS3703	3.5V LED for Lithium Batteries, C-handle	1/ea
MDS3703550	3.5V LED for Rechargeable Li-ion Batteries	1/ea

L3 Otoscope

- » Operation lens with 4x magnification and glass mini-optic ease the insertion of instruments
- » Easy-to-operate SpecEjec specula ejection button on the rear side of instrument head

Item No.	Description	Pkg.
MDS3701	2.5V Xenon for 2 Alkaline Batteries, C-handle	1/ea
MDS3704	3.5V LED for 2 Lithium Batteries, C-handle	1/ea
MDS3704550	3.5V LED for Rechargeable Li-ion Batteries	1/ea

Otoscopes

Elite Vue Otoscope

- » More than 5-fold magnification with twice the field of view
- » Complete view of tympanic membrane, outperforming the market leader¹
- » Wheel to adjust focal depth
- » Suitable for pneumatic otoscopy (supplied without ball)
- » Comes standard with SpecEjec feature at the rear side of the otoscope to help avoid contamination
- » Must use SpecEjec speculas MDS10803533 (3mm) or MDS10803534 (4mm)

Item No.	Description	Pkg.
MDS2200202E	2.5V, XL, for 2 Alkaline Batteries, C-handle	1/ea
MDS2200204E	2.5V, LED, for 2 Alkaline Batteries, C-handle	1/ea



MDS2200202E

ri-scope Diagnostic Sets

- » Set contains one C cell handle, one otoscope head, one ophthalmoscope head, 1 sleeve of 10 specula and case

L2 Otoscope and Ophthalmoscope Diagnostic Sets

Item No.	Description	Pkg.
MDS3746004	3.5V LED/XL with C handle for Li-batteries/Xenon Opth.	1/ea
MDS3750003	3.5V LED/XL with plug-in handle with ri-accu L/Xenon Opth.	1/ea

L3 Otoscope and Ophthalmoscope Diagnostic Sets

Item No.	Description	Pkg.
MDS3746006	3.5V LED/XL with C-handle for Li-batteries/Xenon Ophthalmoscope	1/ea
MDS3750004	3.5V LED/XL with plug-in handle with ri-accu L/Xenon Ophthalmoscope	1/ea

Elite Vue Otoscope and L2 Ophthalmoscope Diagnostic Sets

Item No.	Description	Pkg.
MDS2210202E	2.5V XL EliteVue with 1 C-handle for 2 alkaline batteries/XL L2 Ophthalmoscope	1/ea
MDS2210203E	3.5V LED EliteVue with 1 C-handle for rechargeable Li-ion battery/XL L2 Ophthalmoscope	1/ea
MDS3762E	3.5V XL EliteVue otoscope with 1 C-handle for rechargeable battery/XL L2 Ophthalmoscope	1/ea



MDS3762E



MDS3763E

Desk Set Configurations

- » All desk set configurations come with desk charger, two rechargeable batteries, two handles, otoscope and ophthalmoscope heads
- » 3 hour charging time to full charge

Item No.	Description	Pkg.
MDS107032	LED L2 Otoscope/Xenon L2 Ophthalmoscope set, 3.5V, C	1/ea
MDS107042	LED L3 Otoscope/Xenon L2 Ophthalmoscope set, 3.5V, C	1/ea
MDS107052	LED L3 Otoscope/Xenon L3 Ophthalmoscope set, 3.5V, C	1/ea
MDS2212203E	LED EliteVue/LED L2 Ophthalmoscope set, 3.5V, C	1/ea
MDS3763E	Xenon EliteVue/Xenon L2 Ophthalmoscope set, 3.5V, C	1/ea

Build your diagnostic set.

Desktop. Portable. Wall systems.
Replacement parts. We have it all.

Visit today at <https://diagnostics.medline.com>



e-scope Otoscope/Ophthalmoscope Sets (Pocket Instruments)

e-scope pocket instruments are available with Xenon or LED light sources.

- » Innovative 3.7V LED illumination provides a higher contrast for a more efficient diagnosis
- » Diopter disc with 18 corrective lenses (+/- 1, 2, 3, 4, 6, 8, 10, 15, 20 diopters)
- » Easy-to-use aperture wheel with six different apertures (fixation star, large circle, small circle, red free filter, blue filter, semi circle)
- » High performance optics with aspherical condenser lens
- » Spectacle protection
- » Parallel optical path

Item No.	Description	Pkg.
----------	-------------	------

Fiber optic-2.5V Xenon, Otoscope/Ophthalmoscope

MDS2130202	White, hard case	1/ea
MDS2131202	Black, hard case	1/ea

Fiber optic-3.7V LED, Otoscope/Ophthalmoscope

MDS2130203	White, hard case	1/ea
MDS2131203	Black, hard case	1/ea

* Also available as single instruments. Contact your Medline representative for more information.



Otoscope Accessories

Item No.	Description	Pkg.
MDSK10321*	2.5 mm for L2 Otoscope and E-Scope Pocket Otoscope Only	8,500/cs
MDSK10331*	4.25 mm for L2 Otoscope and E-Scope Pocket Otoscope Only	8,500/cs
MDS10803533	3 mm SpecEjec Speculas for EliteVue and L3 Otoscopes Only	1,000/bg
MDS10803534	4 mm SpecEjec Speculas for EliteVue and L3 Otoscopes Only	1,000/bg
MDS109601	Bulb for Pneumatic Test	1/ea
MDS10961**	Connector for Pneumatic Test	1/ea
MDS10708002	Plug-In Charger with Li-ion Battery for C Handle	1/ea

* Available in child items H (34/pk) and Z (850/bg)

** Only needed with e-scope pocket otoscopes



Vital Signs Monitor with Medline CareConnection®

We offer an integrated solution that allows Long-Term Care and Skilled Nursing Facility caregivers to measure vital signs and auto-document the data in just two minutes.

- » Monitors NIBP, SpO₂ and option of oral or tympanic temperature
- » Compatible with multiple EMR systems:
 - PointClickCare
 - MatrixCare SNF
 - MatrixCare LC
 - American HealthTech
 - Cantata Health
- » Speeds clinical information to your EMR system using the wifi-enabled CareHub tablet with touchscreen[†]

Item No.	Description	Pkg.
MDSM3SNTC	M3 Vital Signs Monitor, Bp, SpO ₂ , Oral Temp	1/ea
MDSM3SNTT	M3 Vital Signs Monitor, Bp, SpO ₂ , Ear Temp	1/ea
CQCMDLN0002	CareHub Stand, Basket and Power Supply	1/ea
CQCTABLET	Touchscreen Tablet and Tablet Holder	1/ea

Additional Parts

CQCTEMPORAL	Temporal Thermometer for Constant Care Vitals Cart	1/ea
CQCKEYBOARD	Keyboard for Constant Care Vitals Cart	1/ea
MDSSPO2ADL	Adult Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2DISP**	Adult Pulse Oximeter Finger Sensor, Disposable	1/ea
MDSSPO2PED**	Pediatric Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2PEDDSP**	Pediatric Pulse Oximeter Finger Sensor, Disposable	1/ea
MDS9910HPRE*	Blood Pressure Cuff, Infant (9–14.8 cm)	1/ea
MDS9911HPRE*	Blood Pressure Cuff, Child (13.8–21.5 cm)	1/ea
MDS9913HPRE*	Blood Pressure Cuff, Adult (27.5–36.5 cm)	1/ea
MDS9914HPRE*	Blood Pressure Cuff, Large Adult (35.5–46 cm)	1/ea
MDS9915HPRE*	Blood Pressure Cuff, Thigh (45–56.5 cm)	1/ea
MDSM3AWRRNTY	Additional 1-Year Extended Warranty for Edan M3 Monitor only (not for CareHub Tablet or Stand)	1/ea
KDL500500	Filac FasTemp Thermometer	500/cs
KDL500500Z	Probe Covers	20/bx
MDSM3A7208058	Ear Probe Covers	200/bx

* Additional Medline cuff options also available. Refer to page 20 for 1-Tube HP Connector Cuffs.

** Will need to purchase MDSSPO2EXT.

† Additional software component required. For more details contact careconnection@medline.com.



CQCMDLN0002

Vital Signs Monitor

Medline's Vital Signs Monitor gives you the ability to provide spot and continuous measurement for NIBP, SpO₂ and temperature in one monitor.

- » Portable and lightweight
- » Easy-to-read and user-friendly monitor
- » Available in different configurations to meet your clinical needs

Choose Your Configuration

Item No.	BP	SpO ₂	Oral Temp	Ear Temp	Printer	Stand	Pkg.
MDSM3ASN1	•	•					1/ea
MDSM3ASNTC	•	•	•				1/ea
MDSM3ASNTCP	•	•	•		•		1/ea
MDSM3ASNTTP	•	•		•	•		1/ea
MDSM3ASNTT	•	•		•			1/ea
MDSM3ASN	•	•				•	1/ea
MDSM3ASNT	•	•	•			•	1/ea
MDSM3ASNTP	•	•	•		•	•	1/ea
MDSM3ASNTPIR	•	•		•	•	•	1/ea
MDSM3ASNTTR	•	•		•		•	1/ea

Choose Your Accessories

Item No.	Description	Pkg.
MDSM3ASTAND01	Rolling Stand w/Storage Basket	1/ea
MDSM3AWALL01	Wall Mount w/Storage Basket	1/ea
MDSSPO2ADL	Adult Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2DISP**	Adult Pulse Oximeter Finger Sensor, Disposable	1/ea
MDSSPO2PED**	Pediatric Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2PEDDSP**	Pediatric Pulse Oximeter Finger Sensor, Disposable	1/ea
MDSSPO2EAR**	Edan Pulse Oximeter Finger Sensor	1/ea
MDSM3ACUFFI*	Edan Blood Pressure Cuff, Infant (14–21 cm)	1/ea
MDSM3ACUFFC*	Edan Blood Pressure Cuff, Child (20.5–28 cm)	1/ea
MDSM3ACUFFA*	Edan Blood Pressure Cuff, Adult (25–35 cm)	1/ea
MDSM3ACUFFL*	Edan Blood Pressure Cuff, Large Adult (33–47 cm)	1/ea
MDSM3ACUFFT*	Edan Blood Pressure Cuff, Thigh (42–54 cm)	1/ea
MDSM3AWRRNTY	Additional 1-Year Extended Warranty	1/ea
KDL500500	Filac FasTemp Thermometer Probe Covers	500/cs
KDL500500Z	Filac FasTemp Thermometer Probe Covers	20/bx
MDSM3A7208058	Ear Probe Covers	200/bx
MDSM3A78035	Printer Paper	1/ea

*Medline cuffs also available. Refer to page 20 for 1-Tube HP Connector Cuffs.

** Will need to purchase MDSSPO2EXT

Additional Parts

Item No.	Description	Pkg.
MDSBPTUBE	NIBP Extension Tube, 3 meters	1/ea
MDSSPO2EXT	SpO ₂ Extension Cable, 2 meters	1/ea
MDSM3A101927	Replacement Single Caster Wheel for Vital Signs Stand	1/ea
MDSPOWERCORD	Standard Power Cord for Vital Signs Monitor	1/ea
MDSBATTERY42	Rechargeable Lithium Ion Battery, 14.8V, 2.5AH	1/ea
MDSM3ABASKET1	Primary Basket for Vital Signs Stand	1/ea
MDSM3ABASKET	Secondary Basket for Vital Signs Stand	1/ea



MDSM3ASNTP

DUS60 Ultrasonic Diagnostic Imaging System by Edan

- » Enhanced efficiency with practical tools, such as 8-segment TGC, data transfer via USB or DICOM, and multiple-pseudo-color imaging
- » Powerful and advanced for diagnostic confidence
- » Compact and light for mobility
- » Intuitive user-friendly design

Item No.	Description	Pkg.
MDSDUS60V11	Edan DUS60 System: Edan DUS60 Ultrasonic Diagnostic Imaging System	1/ea



MDSDUS60V11

U60 Ultrasound Systems by Edan

- » Compact, lightweight design and integrated battery to aid mobility, enabling operation anywhere
- » A full suite of transducers helps ensure that you can scan more patients across more applications
- » Advanced imaging technologies include speckle reduction and spatial compounding to enhance contrast resolution, as well as excellent color imaging to aid diagnostics
- » Monitor specs: 15" display on the U60

Item No.	Description	Pkg.
MDSU60	U60 Diagnostic Ultrasound System	1/ea



MDSU60

Acclarix AX7 Compact Ultrasound Systems

- » Excellent detail resolution in the infraclavicular view
- » Continuous, automatic image optimization makes seeing fluids easier
- » Highly sensitive color Doppler images for fast and accurate detection of blood flow
- » Enhanced border definition and contrast for finding potential problems

Item No.	Description	Pkg.
MDSAX7	Edan Acclarix AX7 Ultrasound System	1/ea



MDSAX7

Medline 12 Channel Interpretive ECG Machine

This 12 channel interpretive unit streamlines patient monitoring with efficient and accurate recordings that will help properly diagnose and monitor patients.

- » Easy-to-read, high resolution color display screen
- » Simplifies patient information entry with one-touch operation and shortcut keys
- » Comprehensive filters and high CMRR reduce reading contamination caused by outside noise
- » 2-year warranty on the monitor and 1-year warranty on the accessories
- » Unit moves easily between rooms with the built-in handle or on the optional mobile stand, available separately

Item No.	Description	Pkg.
MDSECGSE1201	Monitor, ECG, 12 Lead, SE 1201	1/ea
MDSECGSTAND01	Stand, ECG, F/SE 1201	1/ea
GCO32017349	Z-Fold ECG Replacement Paper	20/cs
MDS0157107048	Replacement Leadwire Cable	1/ea
MDS0157040173	Replacement Alligator Clips for Leadwire Cable (Resting)	10/pk
MDSM616101	ECG Resting Tab Electrodes	5,000/cs
MDS0157040172	Replacement Banana Clips for Leadwire Cable (Resting/Snap)	10/pk



MDSECGSE1201

Medline Handheld Fetal and Vascular Doppler

This compact doppler unit provides important, noninvasive testing to help clinicians diagnose many serious conditions, such as blood clots, heart valve defects, congenital heart disease, blocked or bulging arteries and decreased blood circulation. They are extremely popular tools found in most clinical settings.

- » Noise resistant speakers
- » Lightweight, portable design
- » LCD screen displays FHR and probe frequency
- » Interchangeable waterproof probe design
- » Headphone plug

Item No.	Description	Pkg.
MDSDOPPLERF	Doppler, Fetal, 2 MHZ	1/ea
MDSDOPPLERVF	Doppler, Fetal, 2 MHZ and Vascular 4 MHZ	1/ea
MDSDOPPLERV	Doppler, Vascular, 4 MHZ	1/ea
MDSDOPPLER8	Doppler, Vascular, 8 MHZ	1/ea
MDSPROBEF	Probe, Fetal, 2 MHZ	1/ea
MDSPROBEF3	Probe, Fetal, 3 MHZ	1/ea
MDSPROBEV	Probe, Vascular, 4 MHZ	1/ea
MDSPROBEV5	Probe, Vascular, 5 MHZ	1/ea
MDSPROBEV8	Probe, Vascular, 8 MHZ	1/ea



MDSDOPPLERF

Ri-medec Blood Pressure Machines

- » One-touch automated and manual blood pressure
- » Clinical grade technology for fast and reliable readings
- » BP independently validated to British Hypertension Society, ESH and AAMI SP10 specifications for adults and children (from 3 years onwards)
- » Manual SphygMode with auto-deflate provides consistent patient BP measurement and allows clinicians to verify automated readings without removing the cuff
- » Fast readings, typically in 30 to 40 seconds

Item No.	Description	Pkg.
MDS1780	Monitor, Digital, BP, Tabletop (no basket)	1/ea
MDS1782	Monitor, Digital, BP, Mobile (no basket)	1/ea
MDS1784	Monitor, Digital, BP, Wall (no basket)	1/ea
MDS10393	Basket, Large, For Mobile Stand	1/ea
MDS10392	Basket, Small, For Mobile Stand	1/ea

*Works with Medline cuffs. Refer to page 20 for 1-Tube MQ Connector Cuffs.



MDS1782

iM Patient Monitors

- » Edan iM50 and iM70 patient monitors provide essential measurement of patient's vital signs
- » 8.4" (iM50) or 12.1" (iM70) high resolution LCD screen
- » Portable design with bedrail on all models and optional rolling stand
- » Standard parameters are 3/5 Lead ECG, SpO₂, NIBP, Pulse Rate, Dual Temp
- » Optional Parameters include sidestream LoFlo or mainstream Respironics CO₂ modules, O₂ module (for iM70 monitors only) and printer

Item No.	Description	Pkg.
MDSIM50P	iM50 Patient Monitor w/ Printer	1/ea
MDSIM50PT	iM50 Patient Monitor w/ Printer and Touchscreen	1/ea
MDSIM50T	iM50 Patient Monitor w/ Touchscreen	1/ea
MDSIM70PTW	iM70 Patient Monitor w/ Printer, Touchscreen, and Wifi	1/ea
MDSIM70TW	iM70 Patient Monitor w/ Touchscreen and Wifi	1/ea
MDSIMRESP	Sidestream LoFlo CO ₂ Module	1/ea
MDSIMRESPMAIN	Mainstream Respironics CO ₂ Module	1/ea
MDSIMISAOR	Masimo Anesthesia Gas/O ₂ Module, F/iM70 monitors only (not compatible with iM50)	1/ea
MDSM3ASTAND01	Rolling Stand w/ Storage Basket	1/ea
MDSM3AWALL01	Wall Mount w/ Storage Basket	1/ea
MDSPLATEIM	Mounting Plate for iM70 (Needed to attach iM70 monitor to wall mount and/or mobile stand)	1/ea



MDSIM50



MDSIM70

Halogen Handheld Flashlight

- » Flashlight has a lightweight plastic body, weighing less than 0.5 lb. for easy transport
- » Easy single switch operation allows for quick use
- » Ring attached to the base for convenient carrying
- » Halogen bulb offers pure white light with 9 lumens
- » Runs on 2 D batteries (not included) for up to five hours

Item No.	Description	Pkg.
MDSFLASHLT	Halogen Flashlight	12/cs
MDSFLASHLTH		1/ea



MDSFLASHLT

DataTherm® II Continuous Temperature Monitor

- » Measures temperature in real time, updating every 4 seconds
- » Disposable, self-calibrating probes allow for multiple site application
- » Portable and water-resistant design
- » Easily programmable high and low temperature alarms
- » Kit includes DataTherm II, 2 temperature probes, a PC software disk, batteries, and Quick Start Guide

Item No.	Description	Pkg.
RGD501501	DataTherm II Temperature Monitor w/ Two 1-Meter Probes	1/ea
RGD501506	DataTherm II Temperature Monitor w/ Two 2-Meter Probes	1/ea
RGD58000D5	Replacement 1-Meter Probes	5/pk
RGD58000D15		15/pk
RGD58000D52M	Replacement 2-Meter Probes	5/pk
RGD58000D152M		15/pk



RGD501501

High-Visibility EMS Supply Bag

- » Vinyl-lined interior for fluid resistance
- » Adjustable compartments, elastic slots and large zippered interior
- » Detachable strap with rubber shoulder cushion
- » Bright, reflective tape
- » Two zippered exterior side pockets

Item No.	Description	Pkg.
MDSEMSBAG	Orange EMS Trauma Bag	8/cs
MDSEMSBAGH		1/ea



MDSEMSBAG

Geratherm Mini Rescue II Warming System

- » Provides an effective method for regulating patients' temperatures, especially during patient transport by air or ground to medical facilities
- » Safe alternative to chemical warming blankets
- » Lightweight and portable design (<7 lbs)
- » Two power source options: Portable 2590 Li-Ion battery (included in kit) and Medical grade 110-DC on-board power supply (sold separately)
- » Kit includes Mini Rescue II warming blanket (19 x 33"), 2590 Li-Ion battery w/ harness and carabiner, battery cable and charger, and blue soft carrying case

Item No.	Description	Pkg.
RGD82310KIT	Geratherm Warming Kit	1/ea
RGD82009	Optional Medical Grade 110-DC Power Supply	1/ea




RGD82310KIT

Disposable Blood Pressure Cuffs and Connectors

- » Medline's single patient use cuffs are made from a soft, comfortable woven polyester fabric
- » Each cuff has rounded corners and finished edges for optimal comfort
- » Hook-and-loop closures ensure a lasting and snug fit for every patient
- » Cuffs come in a wide range of infant to adult sizes, each color-coded for convenience. All sizes meet the American Heart Association recommendations for accuracy
- » A complete line of connectors and adapters connect cuffs to virtually every non-invasive blood pressure monitor

Choosing the Right Size Cuff

The chart below will help you determine the correct cuff size based on the patient's upper arm circumference. Please note that if you choose a cuff that is not the appropriate size, readings may not be accurate.



	Dinamap Connector (DM)		Bayonet-Style Connector (HP)	Slip Luer Connector (SL)	
Cuff Size (Limb Circumference)	1-Tube	2-Tube	1-Tube	1-Tube	2-Tube
Infant (9–14.8 cm)	MDS9910DM	MDS9920DM	MDS9910HP	MDS9910SL	MDS9920SL
Infant Long (9–14.8 cm)	Contact Representative	MDS9920DML	MDS9910HPL	Contact Representative	Contact Representative
Child (13.8–21.5 cm)	MDS9911DM	MDS9921DM	MDS9911HP	MDS9911SL	MDS9921SL
Child Long (13.8–21.5 cm)	Contact Representative	MDS9921DML	MDS9911HPL	Contact Representative	Contact Representative
Small Adult (20.5–28.5 cm)	MDS9912DM	MDS9922DM	MDS9912HP	MDS9912SL	MDS9922SL
Small Adult Long (20.5–28.5 cm)	Contact Representative	MDS9922DML	MDS9912HPL	Contact Representative	Contact Representative
Adult (27.5–36.5 cm)	MDS9913DM	MDS9923DM	MDS9913HP	MDS9913SL	MDS9923SL
Adult Long (27.5–36.5 cm)	MDS9913DML	MDS9923DML	MDS9913HPL	Contact Representative	MDS9923SLL
Large Adult (35.5–46 cm)	MDS9914DM	MDS9924DM	MDS9914HP	MDS9914SL	MDS9924SL
Large Adult Long (35.5–46 cm)	MDS9914DML	MDS9924DML	MDS9914HPL	MDS9914SLL	MDS9924SLL
Thigh (45–56.5 cm)	MDS9915DM	MDS9925DM	MDS9915HP	MDS9915SL	MDS9925SL

Cuffs constructed of vinyl are available on items listed above.



Male Luer Connector (ML)		Marquette Female Connector (MQ)		Marquette Male Connector (MQML)	Marquette Male & Female Connector (MQMF)
1-Tube	2-Tube	1-Tube	2-Tube	1-Tube	2-Tube
MDS9910ML	MDS9920ML	MDS9910MQ	MDS9920MQ	MDS9910MQML	MDS9920MQMF
Contact Representative	Contact Representative	Contact Representative	MDS9920MQL	Contact Representative	Contact Representative
MDS9911ML	MDS9921ML	MDS9911MQ	MDS9921MQ	MDS9911MQML	MDS9921MQMF
Contact Representative	Contact Representative	Contact Representative	MDS9921MQL	Contact Representative	Contact Representative
MDS9912ML	MDS9922ML	MDS9912MQ	MDS9922MQ	MDS9912MQML	MDS9922MQMF
Contact Representative	Contact Representative	Contact Representative	Contact Representative	Contact Representative	Contact Representative
MDS9913ML	MDS9923ML	MDS9913MQ	MDS9923MQ	MDS9913MQML	MDS9923MQMF
MDS9913MLL	MDS9923MLL	MDS9913MQL	MDS9923MQL	MDS9913MQMLL	MDS9923MQMFL
MDS9914ML	MDS9924ML	MDS9914MQ	MDS9924MQ	MDS9914MQML	MDS9924MQMF
MDS9914MLL	MDS9924MLL	MDS9914MQL	MDS9924MQL	MDS9914MQMLL	MDS9924MQMFL
MDS9915ML	MDS9925ML	MDS9915MQ	MDS9925MQ	MDS9915MQML	MDS9925MQMF

Greensmart Blood Pressure Cuff by Medline

- » Plastic tubing disposed with cuff to reduce cross-contamination risk
- » Adult-sized cuff reduced waste by 2,401 lb./yr at an average hospital*
- » All come with HP Bayonet-Style connectors
- » 19% less material and 20% lower weight than standard disposable cuffs*

*Data available upon request.

Item No.	Size	Pkg.
MDS9910HPNG	Infant	5/bg
MDS9910HPLNG	Infant Long	5/bg
MDS9911HPNG	Child	5/bg
MDS9911HPLNG	Child Long	5/bg
MDS9912HPNG	Small Adult	5/bg
MDS9912HPLNG	Small Adult Long	5/bg
MDS9913HPNG	Adult	5/bg
MDS9913HPLNG	Adult Long	5/bg
MDS9914HPNG	Large Adult	5/bg
MDS9914HPLNG	Large Adult Long	5/bg
MDS9915HPNG	Thigh	5/bg

greensmart®



Blood Pressure Cuff Covers and Barriers

- » Help reduce the potential of cross-contamination by using protective, disposable blood pressure cuff covers and barriers
- » Manufactured without latex
- » Made of Tyvek® material

Item No.	Description	Pkg.
Cuff Covers		
MDS9160	Child Cuff Cover	20/bx
MDS9161	Adult Cuff Cover	20/bx
MDS9162	Large Adult Cuff Cover	20/bx
Cuff Barriers		
MDS917113	Small Cuff Barrier	300/cs
MDS917114	Medium Cuff Barrier	300/cs
MDS917115	Large Cuff Barrier	150/cs



Contour Cuffs

- » Curved profile offers a proper fit, helping ensure reliable BP measurements
- » Conical shape accommodates upper arm circumference of 38-54 cm
- » Available in reusable and disposable models
- » Same helpful features as our standard cuffs

Item No.	Description	Pkg.
MDS9914BC	Disposable Contour Cuff	5/bg
MDS9914BCRE	Reusable Contour Cuff	5/bg



Reusable Blood Pressure Cuffs

- » Our reusable cuffs are made with an easy-to-clean material and include bulb, valve, and slip luer connector
- » Each cuff has rounded corners and finished edges for optimal comfort
- » Hook-and-loop closures ensure a lasting and snug fit for every patient
- » Cuffs come in a wide range of infant to adult sizes, each color-coded for convenience
- » All sizes meet the American Heart Association recommendations for accuracy



	Bayonet-Style Connector (HP)	Screw Connector (DM)	
	1-Tube	1-Tube	2-Tube
Infant (9–14.8 cm)	MDS9910HPRE	MDS9910DMRE	MDS9920DMRE
Infant Long (9–14.8 cm)	MDS9910HPLRE	MDS9910DMLRE	MDS9920DMLRE
Child (13.8–21.5 cm)	MDS9911HPRE	MDS9911DMRE	MDS9921DMRE
Child Long (13.8–21.5 cm)	MDS9911HPLRE	MDS9911DMLRE	MDS9921DMLRE
Small Adult (20.5–28.5 cm)	MDS9912HPRE	MDS9912DMRE	MDS9922DMRE
Small Adult Long (20.5–28.5 cm)	MDS9912HPLRE	Contact Representative	MDS9922DMLRE
Adult (27.5–36.5 cm)	MDS9913HPRE	MDS9913DMRE	MDS9923DMRE
Adult Long (27.5–36.5 cm)	MDS9913HPLRE	MDS9913DMLRE	MDS9923DMLRE
Large Adult (35.5–46 cm)	MDS9914HPRE	MDS9914DMRE	MDS9924DMRE
Large Adult Long (35.5–46 cm)	MDS9914HPLRE	MDS9914DMLRE	MDS9924DMLRE
Thigh (45–56.5 cm)	MDS9915HPRE	MDS9915DMRE	MDS9925DMRE

	Slip Luer Connector (SL)	Marquette Connector (MQ)		Marquette Male & Female Connector (MQMF)
	1-Tube	1-Tube	2-Tube	1-Tube
Infant (9–14.8 cm)	MDS9910SLRE	MDS9910MQRE	MDS9920MQRE	MDS9920MQMFR
Infant Long (9–14.8 cm)	MDS9910SLLRE	MDS9910MQLRE	MDS9920MQLRE	MDS9920MQMFLR
Child (13.8–21.5 cm)	MDS9911SLRE	MDS9911MQRE	MDS9921MQRE	MDS9921MQMFR
Child Long (13.8–21.5 cm)	MDS9911SLLRE	MDS9911MQLRE	MDS9921MQLRE	MDS9921MQMFLR
Small Adult (20.5–28.5 cm)	MDS9912SLRE	MDS9912MQRE	MDS9922MQRE	MDS9922MQMFR
Small Adult Long (20.5–28.5 cm)	Contact Representative	Contact Representative	MDS9922MQLRE	MDS9922MQMFLR
Adult (27.5–36.5 cm)	MDS9913SLRE	MDS9913MQRE	MDS9923MQRE	MDS9923MQMFR
Adult Long (27.5–36.5 cm)	MDS9913SLLRE	MDS9913MQLRE	MDS9923MQLRE	MDS9923MQMFLR
Large Adult (35.5–46 cm)	MDS9914SLRE	MDS9914MQRE	MDS9924MQRE	MDS9924MQMFR
Large Adult Long (35.5–46 cm)	MDS9914SLLRE	MDS9914MQLRE	MDS9924MQLRE	MDS9924MQMFLR
Thigh (45–56.5 cm)	MDS9915SLRE	MDS9915MQRE	MDS9925MQRE	MDS9925MQMFR



MDS1001



MDS3001PLUS



MDS4001



MDS5001

Medline **Elite**

Automatic Digital Blood Pressure Monitor

- » Revolutionary design and proprietary programming provide automatic inflation and error detection
- » Customized inflation levels for each individual patient
- » Color gauge for quick detection
- » Provides quick digital readings of systolic/diastolic pressure and pulse rate
- » Easy-to-use one-push button design
- » Stores the 90 most recent measurements
- » Large digital readout
- » Includes the automatic blood pressure unit, cuff and four AA batteries

Item No.	Description	Pkg.
MDS1001	Monitor with Adult Cuff (small & compact)	1/ea
MDS1001U	Monitor with Universal Cuff	1/ea
MDS3001	Monitor with Adult Cuff (large display screen)	1/ea
MDS3001PLUS	Monitor with Adult & Large Adult Cuffs	1/ea
MDS3001U	Monitor with Universal Cuff	1/ea
MDS4001	Monitor with Adult Cuff (backlight)	1/ea
MDS4001LAT	Talking Monitor with Large Adult Cuff	1/ea
MDS4001PLUS	Monitor with Adult & Large Adult Cuffs	1/ea
MDS5001	Talking Bedside Monitor with Universal Cuff, AC Adapter, Cuff Holder, Alarm and Extra Large Backlit Display	1/ea

Replacement Cuffs; Fit all Medline brand fully-automatic BP monitors

MDS9970	Small Adult Cuff (fits arms 17–22 cm)	1/ea
MDS9971	Adult Cuff (fits arms 22–30 cm)	1/ea
MDS9972	Large Adult Cuff (fits arms 30–42 cm)	1/ea
MDS9973	Extra Large Adult Cuff (fits arms 42–48 cm)	1/ea
MDS9974	Universal Sized Cuff (fits arms 22–42 cm)	1/ea

MedlinePro**Semi-Automatic Digital Blood Pressure Monitor**

- » Reliable readings from a manually inflated digital blood pressure unit
- » Provides quick digital readings of systolic/diastolic pressure and pulse rate
- » Easy-to-use one-push button design
- » Stores the 90 most recent measurements
- » Unit includes the semi-automatic blood pressure unit, cuff and four AA batteries



MDS3002

Item No.	Description	Pkg.
MDS3002	Monitor with Adult Cuff	1/ea
MDS3002LA	Monitor with Large Adult Cuff	1/ea
MDS9770	Small Adult Cuff (fits arms 17–22 cm)	1/ea
MDS9771	Adult Cuff (fits arms 22–30 cm)	1/ea
MDS9772	Large Adult Cuff (fits arms 30–42 cm)	1/ea
MDS9773	Extra Large Adult Cuff (fits arms 42–48 cm)	1/ea

MedlinePlus**Digital Wrist Blood Pressure Monitor**

- » Easy-to-use automatic wrist blood pressure monitor
- » Automatically inflates and deflates
- » Provides quick digital readings of systolic/diastolic pressure and pulse rate
- » Memory with date and time
- » Color gauge for quick detection
- » Includes two AAA batteries



MDS4003

Item No.	Description	Pkg.
MDS4003	Digital Wrist Blood Pressure Monitor	1/ea

Pedometer and Heart Rate Watch

- » Takes pulse measurements
- » Offers three pedometer display options:
 - Number of steps and current time
 - Number of steps during exercise
 - Calories burned with distance walked
- » Automatically stores last 50 days of steps
- » LCD display and four-second backlight



MDSP3044

Item No.	Description	Pkg.
MDSP3044	Fitness Watch	1/ea

Medline Elite**Medline Handheld Aneroids**

Our finest handheld aneroids.

- » Luminescent dial is easy to read, even in dim light
- » Shatter-resistant acrylic faceplate surrounds the luminescent dial
- » Exceptional accuracy from a top-quality, zinc alloy precision gauge
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included

Item No.	Description	Pkg.
MDS9411	Medline Elite Handheld Aneroid, Child	1/ea
MDS9410	Medline Elite Handheld Aneroid, Adult	1/ea
MDS9413	Medline Elite Handheld Aneroid, Large Adult	1/ea

Easy Product Selection

Differences among various aneroids can be subtle. The names of some of our popular aneroid products reflect our good, better and best in each line.



MDS9410

Medline Pro**Neoprene Handheld Aneroids**

High-quality handheld aneroids.

- » Strong neoprene tubing and bladder
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included
- » Not made with natural rubber latex

Item No.	Description	Pkg.
MDS9386LF	Medline Pro Handheld Aneroid, Infant	1/ea
MDS9387LF	Medline Pro Handheld Aneroid, Child	1/ea
MDS9380LF	Medline Pro Handheld Aneroid, Adult	1/ea
MDS9388LF	Medline Pro Handheld Aneroid, Large Adult	1/ea
MDS9389LF	Medline Pro Handheld Aneroid, Thigh	1/ea



MDS9380LF

Medline Plus**PVC Handheld Aneroids**

Great quality. Great price.

- » PVC tubing and bladder
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included
- » Not made with natural rubber latex

Item No.	Description	Pkg.
MDS9386	Medline Plus Handheld Aneroid, Infant	1/ea
MDS9387	Medline Plus Handheld Aneroid, Child	1/ea
MDS9380	Medline Plus Handheld Aneroid, Adult	1/ea
MDS9388	Medline Plus Handheld Aneroid, Large Adult	1/ea



MDS9380

MedlinePlus**Palm Aneroids**

- » Hand-held palm aneroid swivels for both left and right hand users
- » Quick connect HP Bayonet connector allows you to efficiently change out the cuff size you need depending on the patient
- » Compatible with all Medline 1 tube HP/Bayonet disposable and reusable cuffs

Item No.	Description	Pkg.
MDSANEROIDHP	Handheld Palm Aneroid with HP Connector	1/ea

Single Tube Reusable Blood Pressure Cuff w/ Bayonet Connector

MDS9911HPREH	Child	1/ea
MDS9912HPREH	Small Adult	1/ea
MDS9913HPREH	Adult	1/ea
MDS9914HPREH	Large Adult	1/ea



MDSANEROIDHP

Wall-Mount/Mobile Aneroids

- » Large, raised numbers and a luminescent coating make the wall-mount/mobile aneroid easy to read in dim light
- » Eight-foot coiled tubing and pivoting bracket extend service area
- » Inflation system and bulb assembly store easily and safely behind gauge
- » PVC and Neoprene are manufactured without latex
- » Includes Medline's nylon range finder cuff

Item No.	Description	Pkg.
----------	-------------	------

MedlinePro**Neoprene Aneroids**

MDS9400LF	Aneroid, Wall Mount, Adult Cuff	1/ea
MDS9407LF	Aneroid, Mobile, Adult Cuff with Stand	1/ea
MDS9418	Mobile Stand	1/ea

MedlinePlus**PVC Aneroids**

MDS9400	Aneroid, Wall Mount, Adult Cuff	1/ea
MDS9407	Aneroid, Mobile, Adult Cuff with Stand	1/ea
MDS9418	Mobile Stand	1/ea



MDS9400LF

MedlinePlus**Home Blood Pressure Kits**

- » Easy-to-read aneroid gauge
- » Includes complete instructions and a recording chart
- » Packed in a matching nylon carrying case

Item No.	Description	Pkg.
MDS9300	Home Blood Pressure Kit, Detached Stethoscope	1/ea
MDS9301	Home Blood Pressure Kit, Attached Stethoscope	1/ea



MDS9300

Medline **Elite** **Compli-Mates®**

Kit includes an aneroid sphygmomanometer, matching 22" (56 cm) stainless steel dual-head stethoscope and matching nylon carrying case.

Item No.	Description	Pkg.
MDS9140	Steel/Black (Stainless Steel)	1/ea



MDS9140

Medline **Deluxe** **Compli-Mates®**

Kit includes an aneroid sphygmomanometer, matching 22" (56 cm) Sprague-Rappaport stethoscope and matching nylon carrying case.

Item No.	Description	Pkg.
MDS9121	Magenta	1/ea
MDS9123	Purple	1/ea
MDS9124	Royal Blue	1/ea
MDS9125	Black	1/ea
MDS9128	Hunter Green	1/ea



MDS9124

Medline **Pro** **Compli-Mates®**

Kit includes an aneroid sphygmomanometer, matching dual-head stethoscope and matching nylon carrying case.

Item No.	Description	Pkg.
MDS9116	Magenta	1/ea
MDS9117	Teal	1/ea
MDS9118	Purple	1/ea
MDS9119	Royal Blue	1/ea
MDS9120	Black	1/ea
MDS9133	Hunter Green	1/ea
MDS9136	Red	1/ea



MDS9133

Medline **Plus** **Compli-Mates®**

Kit includes an aneroid sphygmomanometer and matching nylon carrying case.

Item No.	Description	Pkg.
MDS9104	Red	1/ea
MDS9106	Navy Blue	1/ea
MDS9108	Hunter Green	1/ea
MDS9109	Orchid	1/ea
MDS9111	Magenta	1/ea
MDS9112	Teal	1/ea
MDS9113	Purple	1/ea
MDS9114	Royal Blue	1/ea



MDS9106

Aneroid Accessories



MDS9414



MDS9630

Handheld Aneroids			Wall-Mount Aneroids	
MedlineElite	MedlinePro	MedlinePlus	MedlinePro	MedlinePlus
MDS9410, MDS9411, MDS9413	MDS9386LF, MDS9387LF, MDS9380LF, MDS9388LF	MDS9386, MDS9387, MDS9380, MDS9388	MDS9407LF, MDS9400LF	MDS9400, MDS9407

Premier Gauge and Parts

Gauge (MDS9414)	•				
Glass Face (MDS9415)	•				

Standard Gauge and Parts

Aneroid Gauge (MDS9381)		•	•		
-------------------------	--	---	---	--	--

Bulbs and Valves

PVC Bulb (MDS91481)			•		•
Neoprene Bulb (MDS91481LF)		•		•	
Premier bulb and valve assembly (MDS9630)	•				
PVC Bulb and Valve Assembly (MDS91482)		•	•	•	•

Two-Tube Inflation Bag and Nylon Range Finder Cuff

Neoprene, Infant (MDS91422LF)		•		•	
Neoprene, Child (MDS91421LF)		•		•	
Neoprene, Adult (MDS91420LF)		•		•	
Neoprene, Large Adult (MDS91423LF)		•		•	
PVC, Infant (MDS91422)			•		•
PVC, Child (MDS91421)			•		•
PVC, Adult (MDS91420)			•		•
PVC, Large Adult (MDS91423)			•		•
PVC, Thigh (MDS91424)			•		•

Two-Tube Seamless Inflation Bag

Neoprene, Adult (MDS9505LF)		•		•	
Neoprene, Large Adult (MDS9506LF)		•		•	
PVC, Adult (MDS9505)			•		•
PVC, Large Adult (MDS9506)			•		•

Nylon Range Finder Hook-and-Loop Cuffs

Child (MDS9421)		•	•	•	•
Adult (MDS9420)		•	•	•	•
Large Adult (MDS9423)		•	•	•	•

8' Coiled Tubing with Connector

PVC (MDS9475)					•
Neoprene (MDS9475LF)				•	

Miscellaneous Parts

Metal Connector/Coupler (MDS9481)				•	•
Basket for Wall-Mounted or Mobile Aneroids (MDS9472)				•	•
Mobile Stand (MDS9418)				•	•

RightTemp Tympanic Thermometer

- » This premium easy-to-use ear thermometer is perfect for practically every department and patient
- » Simple-to-read backlit LCD displays readings in Fahrenheit or Celsius
- » Accuracy: $\pm 0.4^{\circ}\text{F}$, $\pm 0.2^{\circ}\text{C}$
- » Stores 10 readings
- » Rigid disposable probe covers help reduce the risk of cross-contamination; unit is supplied with 20 covers; additional are sold separately

Item No.	Description	Pkg.
MDS8700	RightTemp Tympanic Thermometer	1/ea
MDS8701	RightTemp Tympanic Thermometer Probe Cover	5000/cs
MDS8701Z		200/bx
MDS8701H		20/pk



MDS8700

Large Display Digital Thermometer

- » Digital thermometers feature a large, clear LCD with a back light for quick, easy viewing
- » Provide highly accurate readings within 20 seconds
- » With a water-resistant probe tip for patient use: oral, axillary or rectal
- » Offer temperatures in both Fahrenheit and Celsius

Item No.	Description	Pkg.
MDS9953	Large Display F/C Thermometer	200/cs
MDS9953H		1/ea



MDS9953H

No Touch Forehead Thermometer

Innovative, quick-response forehead thermometer features:

- » One-touch operation
- » Takes measurements when holding the probe parallel to the forehead and approximately 1-inch away
- » Helps lower risk of cross contamination by eliminating patient contact
- » No probe covers needed
- » Displays green light for healthy temperature and red light for fever
- » Has human mode (head symbol) and surface mode (bottle symbol)

Item No.	Description	Pkg.
MDSNOTOUCH	Forehead Thermometer	1/ea



MDSNOTOUCH

Digital Temple Thermometer

The Medline instant read digital temple thermometer accurately measures body heat using R.A.T.E.™ technology (Rapid, Accurate, Temperature Establishment).

- » Sensors rapidly track heat flow generated from blood vessels to the skin's surface and convert measurements to body temperature
- » Easy-to-use on asleep or alert patients
- » Meets accuracy standards
- » Large digital display
- » Uses two AAA batteries
- » Memory stores previous reading

Item No.	Description	Pkg.
MDS9698L	Temple Thermometer	1/ea
MDS9699	Probe Covers	250/bx



MDS9698L

Talking Ear and Forehead Thermometer

Versatile thermometer developed for accurate, safe and fast human body temperature measurements via the ear and forehead.

- » Celsius/Fahrenheit option
- » LCD backlit display
- » Fever indicator: green (normal), slow blinking red (raised), fast blinking red (high)
- » Memory feature recalls past nine measurements
- » For home-use only

Item No.	Description	Pkg.
MDSTH1002	Ear and Forehead Thermometer	1/ea



MDSTH1002

Tympanic Thermometer

Innovative, quick-response ear thermometer with quick probe release.

- » One-touch operation
- » Mechanical probe cover quick-release allows easy probe cover change
- » Safe to use on patients of any age
- » Readings in Fahrenheit and Celsius
- » Takes an accurate, non-invasive temperature from the tympanic membrane

Item No.	Description	Pkg.
MDS9700	Thermometer	1/ea
MDS9701	Probe covers	100/bx



MDS9700

Exergen® Temporal Scanner® Thermometer

Simply scan across the temporal artery area of the forehead or touch behind the earlobe and an accurate reproducible temperature is measured in about three seconds.

- » Convenient one-handed operation
- » Completely non-invasive
- » Stainless steel probe for highest durability
- » Easy snap-on probe caps
- » Three-year warranty for EXGTAT2000; Lifetime warranty for EXGTAT5000

Item No.	Description	Pkg.
EXGTAT5000	Infrared Temporal Artery Thermometer	1/ea
EXGTAT2000	Infrared Temporal Artery Thermometer	1/ea
EXGTAT2000C*	Infrared Temporal Artery Thermometer, Consumer Model	1/ea
EXG134203	Disposable Probe Caps	1,000/cs
EXG134203H	Disposable Probe Caps for EXGTAT2000/5000	50/pk

*Disposable probe caps are not compatible with EXGTAT2000C



EXGTAT2000



EXGTAT5000

Disposable Forehead Thermometer

- » Convenient solution to monitoring patients' temperature trends, especially during minor surgical procedures.
- » 95-104° F and 35-40° C measurement range
- » Easy box dispenser with Hook and Loop closure for wall attachment
- » Unique initial temperature marking area

Item No.	Description	Pkg.
MDSTS101	Disposable Forehead Thermometer Strips	5,000/cs
MDSTS101Z		500/bx
MDSTS101H		100/pk



MDSTS101

30-Second Digital Thermometers

- » Provide accurate temperatures in about 30 seconds, depending on model
- » Large LCD screen for easy readings
- » Simple push-button start
- » Disposable sheaths ensure sanitary use

Item No.	Description	Pkg.
Oral		
MDS9950	Oral °F & °C	144/cs
MDS9950Z		12/bx
MDS9950H		1/ea
MDS9950B	Oral °F & °C, Bulk	200/cs
MDS9950BZ		20/bx
MDS9607	Oral Sheaths	100/bx
MDS9928	Oral Kit w/20 Sheaths, °F & °C	1/ea

Rectal

MDS9952	Rectal °F & °C	144/cs
MDS9952Z		12/bx
MDS9952H		1/ea
MDS9952B	Rectal °F & °C Bulk	200/cs
MDS9952BZ		20/bx
MDS9608	Pre-lubricated Rectal Sheaths	100/bx
MDS9929	Rectal Kit w/20 Sheaths, °F & °C	1/ea

Flex-tip—Used for Oral, Rectal or Axillary Measurements

MDS99901	Flex-tip, Water-Resistant, °F & °C	240/cs
MDS99901Z		12/bx
MDS99901H		1/ea
MDS99902	Flex-tip, Water-Resistant, °F & °C	240/cs
MDS99902Z		12/bx
MDS99902H		1/ea

Disposable Clinical Thermometer Probe Cover

- » Medline brand, FDA-approved disposable probe covers for use with the WelchAllyn SureTemp® and SureTemp Plus®

Item No.	Description	Pkg.
MDS7701	SureTemp Probe Covers	7,500/cs
MDS7701Z		250/bx
MDS7701H		25/pk

Disposable Penlight

- » Lightweight and disposable
- » High-intensity exam light is perfect for the physician, nurse and EMT
- » Constructed of high-impact molded styrene and can be clipped on pockets

Item No.	Description	Pkg.
MDS131040	Disposable Penlight, White	6/pk
MDS131040H		1/ea



MDS9950



MDS9952



MDS99901



MDS99902



MDS7701



MDS131040

Synergy® Dual Frequency Stethoscope

- » High end stainless steel chest piece helps trap sound
- » Strong, thick wall tubing eliminates sound interference
- » Soft, narrow ear tips for everyday comfort
- » Elevated diaphragm provides excellent auscultation and acoustics
- » Apply minimal pressure to hear low frequencies and increase pressure for higher frequencies
- » All models come with an extra set of ear tips
- » Available in adult and cardiology styles

Item No.	Description	Pkg.
MDS926501	Black	1/ea
MDS926502	Blue	1/ea
MDS926503	Hunter Green	1/ea
MDS926504	Burgundy	1/ea
MDS926601	Black, Single-head	1/ea
MDS926801	Black, Cardiology	1/ea



MDS926501

Stainless Steel Stethoscope

- » Handcrafted from high tensile strength stainless steel with durable 22" (56 cm) tubing
- » Color coordinated non-chill rings provide excellent patient comfort
- » Available in adult, pediatric, neonatal and cardiology styles
- » Cardiology stethoscope has a deep, cone-shaped bell—a unique feature that provides excellent sound transmission and acoustic sensitivity.

Item No.	Description	Pkg.
Adult		
MDS92250	Gray	1/ea
MDS92260	Black	1/ea
MDS92270	Blue	1/ea
MDS92280	Burgundy	1/ea
MDS92290	Pink Breast Cancer Awareness	1/ea
Pediatric		
MDS92200	Gray	1/ea
Neonatal		
MDS92150	Gray	1/ea
MDS92160	Pink	1/ea
Cardiology		
MDS92500	Black	1/ea



MDS92500



Sprague-Rappaport Stethoscope

- » Sprague-Rappaport style incorporates versatility into an excellent sound-transmitting device.
- » Each stethoscope comes complete with a component kit that includes two sets of ear tips, large and small diaphragm, adult bell, pediatric bell and infant bell
- » The dual 22" (56 cm) tubing design transmits even the faintest sounds

Item No.	Description	Pkg.
MDS926301	Black	1/ea
MDS926302	Gray	1/ea
MDS926303	Blue	1/ea
MDS926305	Lavender	1/ea
MDS926306	Red	1/ea
MDS926307	Hunter	1/ea
MDS926311	Yellow	1/ea
MDS926312	Orange	1/ea



MDS926303

Single- and Dual-Head Stethoscope

- » Both stethoscopes feature 22" (56 cm) tubing and lightweight design
- » Thick-walled tubing minimizes background interference
- » Dual-heads feature a non-chill ring on the bell for increased patient comfort

Item No.	Description	Pkg.
----------	-------------	------

Single-head

MDS926101	Black	1/ea
MDS926101CS	Black	100/cs
MDS926102	Gray	1/ea
MDS926103	Blue	1/ea
MDS926103CS	Blue	100/cs
MDS926105	Lavender	1/ea
MDS926106	Red	1/ea
MDS926106CS	Red	100/cs
MDS926107	Hunter Green	1/ea
MDS926108	Yellow	1/ea
MDS926108CS	Yellow	100/cs
MDS926109	Pink	1/ea
MDS926110	Orange	1/ea

Dual-head

MDS926201	Black	1/ea
MDS926201CS	Black	100/cs
MDS926202	Gray	1/ea
MDS926203	Blue	1/ea
MDS926203CS	Blue	100/cs
MDS926205	Lavender	1/ea
MDS926206	Red	1/ea
MDS926207	Hunter Green	1/ea
MDS926208	Yellow	1/ea
MDS926209	Pink	1/ea
MDS926210	Orange	1/ea



MDS926108



MDS926206



Pediatric/Neonatal Aluminum Stethoscope

- » Both stethoscopes feature thick-wall 22" (56 cm) tubing and lightweight design
- » All features of the dual-head stethoscopes but sized to fit the needs of the individual practitioner and patient
- » Thick-walled tubing minimizes background interference

Item No.	Description	Pkg.
MDS9557	Gray, Pediatric	1/ea
MDS9565	Gray, Neonatal	1/ea

Disposable Stethoscope

- » Designed for single patient use
- » Helps reduce cross-contamination in infectious areas
- » High visibility yellow

Item No.	Description	Pkg.
MDS9543	Disposable Stethoscope	1/ea
MDS9543CS	Disposable Stethoscope	100/cs

Teaching/Training Stethoscope

- » Ideal for all instructional applications
- » Dual-head design is perfect for all clinical and educational conditions
- » Includes two binaurals, two pieces of Y-tubing and an extension tube

Item No.	Description	Pkg.
MDS9750	Teaching/Training Stethoscope	1/ea

Bowles Stethoscope

- » Thick-walled 22" (56 cm) tubing
- » Chrome-plated binaural and raised-metal-stem chest piece

Item No.	Description	Pkg.
MDS9515	Black, Adult	1/ea

Glow-in-the-Dark Stethoscope

- » Add flair and fun to your work by emitting a green glow when your surroundings darken
- » Feature 22" (56 cm) thick-walled tubing to minimize background interference

Item No.	Description	Pkg.
MDS9262GLOW	Dual-Head	1/ea
MDS9263GLOW	Sprague-Rappaport	1/ea



MDS9262GLOW

Philips HeartStart

- » Designed for the ordinary person in the extraordinary moment—easy-to-use device that arrives ready-to-rescue
- » Built with on-demand CPR coaching to provide guidance during resuscitation

Item No.	Description	Pkg.
DEFM5066HOME	Philips OnSite AED with Slim Carrying Case	1/ea
DEFM5066PACK	Philips OnSite AED Package with Wall Cabinet and Carrying Case	1/ea

Philips Replacement Accessories

DEFM5071A	Adult Pads	1/ea
DEFM5072A	Pediatric Pads	1/ea
DEF5070A	OnSite Battery Pack	1/ea



DEFM5066HOME

Defibtech® Lifeline Series

- » Built to the demanding standards of professional emergency personnel, Lifeline AEDs also meet US military “drop and shock” specifications.

Item No.	Description	Pkg.
DFTDCF1310EN	Lifeline VIEW AED (Semi-Automatic)	1/ea
DFTDCFA100EN	Lifeline AED (Semi-Automatic)	1/ea
DFTDCFA120	Lifeline AUTO AED (Fully-Automatic)	1/ea

Defibtech Replacement Accessories

DFTDDP100	Lifeline AED Adult Pads for adult use only (eight years or older)	1/ea
DFTDDP200P	Lifeline AED Pediatric Pads for infants and children under eight years old	1/ea
DFTDDP2001	Lifeline VIEW AED Adult Pads for adult use only (eight years or older)	1/ea
DFTDDP2002	Lifeline VIEW AED Pediatric Pads for infants and children under eight years old	1/ea
DFTDCF200	Lifeline AED Battery Pack	1/ea
DFTDCF2003	Lifeline VIEW AED Battery Pack	1/ea



DFTDCF1310EN

Zoll AED Plus Defibrillators

Equipped with Real CPR Help, the AED Plus provides real-time feedback on the quality of CPR.

Item No.	Description	Pkg.
ZOL2110201CC	AED Plus Carrying Case Package (Semi-Automatic)	1/ea
ZOL2110201PK	AED Plus Wall-Cabinet Package (Semi-Automatic)	1/ea

Zoll AED Replacement Accessories

ZOL8900080001	Adult Pads	1/ea
ZOL8900081001	Pediatric Pads	1/ea
ZOL8000080701	Battery	1/ea



ZOL2110201CC

General Monitoring Electrodes

- » Clear tracing every time with our high performance, specially formulated gel with extended surface area
- » Electrodes stay securely on the skin, yet remove cleanly with our high-quality adhesive
- » Flexible lead wire attachment with our perforated strip backing
- » Prolonged electrode freshness with our resealable bag (50 pack only)
- » Easy removal with full-edge lift tab

MedGel® General Monitoring Foam Electrode

Recommended for: ICU/CCU, Emergency, O.R., Patient Rooms
Our best-selling electrode.

- » Flexible, foam material

Item No.	Description	Pkg.
MDSM611503	General Monitoring, Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611505		5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611550		50/pk, 1,000/ct, 6,000/cs



MDSM611550

MedGel® General Monitoring Cloth Electrode

Recommended for: ICU/CCU, Patient Rooms
Great for patients with sensitive skin.

- » Enhanced comfort by conforming to the shape of the body with our non-woven cloth material

Item No.	Description	Pkg.
MDSM611303	General Monitoring, Cloth	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611305		5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611350		50/pk, 1,000/ct, 6,000/cs



MDSM611350

MedGel® General Monitoring Tape Electrode

Recommended for: ICU/CCU, Emergency, O.R., Patient Rooms
Ideal for seeing directly under the electrode without having to remove it.

- » Provides visibility under the electrode with our transparent tape material

Item No.	Description	Pkg.
MDSM611403	General Monitoring, Tape	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611405		5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611450		50/pk, 1,000/ct, 6,000/cs



MDSM611450

Radiotranslucent Electrodes

- » Can be used during X-rays and MRIs in areas where metal is prohibited—Radiology, Cath Lab, Surgery, etc.—with our radiotranslucent snap
- » Clear tracings every time with our high performance, specially formulated gel and extended surface area
- » Stays securely on the skin yet removes cleanly with our high-quality adhesive
- » Prolonged electrode freshness and life with our resealable bag (50 pack only)
- » Flexible lead wire attachment with our perforated strip backing
- » Easy removal with full-edge lift tab

MedGel® Radiotranslucent Foam Electrode

Recommended for: ICU/CCU, Emergency, O.R., Cath Lab, Radiology, Patient Rooms
Our most popular radiotranslucent electrode.

- » Flexible, foam material

Item No.	Description	Pkg.
MDSM611803R	Radiotranslucent Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611805R		5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611850R		50/pk, 1,000/ct, 6,000/cs



MDSM611850R

MedGel® Radiotranslucent Cloth Electrode

Recommended for: ICU/CCU, Cath Lab, Radiology, Patient Rooms
Great for patients with sensitive skin.

- » Enhanced comfort by conforming to the shape of the body with our non-woven cloth material

Item No.	Description	Pkg.
MDSM611250R	Radiotranslucent Cloth	50/pk, 1,000/ct, 6,000/cs




MDSM611250R

MedGel® Radiotranslucent Tape Electrode

Recommended for: ICU/CCU, O.R., Cath Lab, Radiology, Patient Rooms
Ideal for seeing directly under the electrode without having to remove it.

- » Provides visibility under the electrode with our transparent tape material

Item No.	Description	Pkg.
MDSM611650R	Radiotranslucent Tape	50/pk, 1,000/ct, 6,000/cs



MDSM611650R

MedGel® Pediatric Foam Electrode

Recommended for: PICU

Smaller footprint fits pediatric patients better.

- » Greater percentage of gel allows for gentle use and removal
- » Clear tracings every time with our high performance, specially formulated gel
- » Prolonged electrode freshness and life with our resealable bag
- » Flexible lead wire attachment with our perforated strip backing
- » Easy removal with lift tab

Item No.	Description	Pkg.
MDSM611050	Pediatric, Foam	50/pk, 1,000/ct, 6,000/cs



MDSM611050

MedGel® Pediatric Radiotranslucent Electrode

Recommended for: Cath Lab, Radiology, PICU

Radiotranslucent foam and small size for pediatric patients.

- » Can be used during X-Rays and MRIs in areas where metal is prohibited—Radiology, Cath Lab, Surgery, etc.—with our radiotranslucent snap
- » Prolonged electrode freshness and life with our resealable bag
- » Flexible lead wire attachment with our perforated strip backing
- » Easy removal with lift tab

Item No.	Description	Pkg.
MDSM6112250R	Pediatric, Radiotranslucent Foam	50/pk, 1,000/ct, 6,000/cs



MDSM6112250R

MedGel® Pre-Wired Neonatal Electrodes

Recommended for: NICU, Pre-wired baby-sized electrodes.

- » Gentle tape material—ideal for neonatal patients
- » Ouchless removal with 100% gel coverage
- » Three electrodes come pre-wired for convenience
- » Clear tracings every time with our high performance, specially formulated gel

Item No.	Description	Pkg.
MDSM61123	Neonatal Clear Tape	3/pk, 300/ct, 3,600/cs
MDSM61143R	Neonatal, Radiotranslucent Clear Tape	3/pk, 300/ct, 3,600/cs



MDSM61143R

MedGel® EKG Resting Tab Electrodes

Recommended for: Emergency, Cardiology, Patient Rooms

- » Stable, reliable EKG signals with silver/silver chloride gel
- » Maximum contact with the skin from a large solid conductive area

Item No.	Description	Pkg.
MDSM616101	Tab Electrode, Small	100/pk, 500/bx, 5,000/cs
MDSM616201	Tab Electrode, Adult	100/pk, 500/bx, 5,000/cs



MDSM616101




MDSM616201

MedGel® Stress/Diaphoretic Foam Electrodes

Recommended for: Cardiology, Stress Test

- » Specially designed for situations where a strong adhesive is required
- » Clear tracings every time with our high performance, specially formulated gel and extended surface area

Item No.	Description	Pkg.
MDSM611903	ECG Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611905		5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611950		50/pk, 1000/ct, 6,000/cs




MDSM611950

MedGel® Radiotranslucent Stress/Diaphoretic Foam Electrodes

Recommended for: Cardiology, Cath Lab, Radiology, Stress Test

- » Radiotranslucent and designed for situations where a strong adhesive is required
- » Can be used during X-Rays and MRIs in areas where metal is prohibited—Radiology, Cath Lab, Surgery, etc.—with our radiotranslucent snap
- » Clear tracings every time with our high performance, specially formulated gel and extended surface area

Item No.	Description	Pkg.
MDSM611703R	Radiotranslucent Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611705R		5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611750R		50/pk, 1,000/ct, 6,000/cs



MDSM611750R

MedGel® Wet Gel Stress/Diaphoretic Foam Electrodes

Recommended for: Cardiology, Stress Test, EMS

- » Quick reading and stabilization for situations of particularly heavy patient sweating
- » Instant, high quality trace pickup with our liquid gel

Item No.	Description	Pkg.
MDSM618303	Wet Gel Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM618305		5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM618350		50/pk, 1,000/ct, 6,000/cs
MDSM619305R	Wet Gel Radiotranslucent Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM619350R		50/pk, 1,000/ct, 6,000/cs



MDSM618350

Medline Plus Repositionable Electrodes

Recommended for: ICU/CCU, Emergency, O.R., Cath Lab, Radiology, PICU, Patient Rooms

- » Easy to reposition at least once with a conductive hydrogel
- » Gentle on sensitive skin
- » Carbon snap for standardization throughout your facility














Item No.	Description	Pkg.
MDSM611300	Adult, Solid Gel	5/pk, 100/bx, 1,000/cs
MDSM611200	Pediatric, Solid Gel	5/pk, 100/bx, 1,000/cs



MDSM611300

MedGel® Electrodes

Recommended Electrodes for Optimal Utilization

		ICU/CCU	Emergency	O.R.	Cardiology	Cath Lab	Radiology	Stress Test	PICU	NICU	Patient Rooms	EMS
General Monitoring Foam		•	•	•							•	
General Monitoring Cloth		•									•	
General Monitoring Tape		•	•	•							•	
Radiotranslucent Foam		•	•	•		•	•				•	
Radiotranslucent Cloth		•				•	•				•	
Radiotranslucent Tape		•		•		•	•				•	
Stress/Diaphoretic Foam					•			•				
Radiotranslucent Stress/Diaphoretic Foam					•	•	•	•				
Wet Gel Stress/Diaphoretic Foam					•			•				•
Pediatric Foam									•			
Radiotranslucent Pediatric Foam						•	•		•			
Neonatal Electrodes								•		•		
Resting EKG Electrodes			•		•						•	

This list is not comprehensive and should be used as a recommended guideline only.

Defibrillator Pad Electrodes

- » Multi-function defibrillator pad electrodes are intuitive and provide reliable outcomes
- » Color coded packaging to eliminate product selection mistakes
- » 3 year shelf-life intended to reduce product waste
- » 60" lead wire for defibrillator placement flexibility
- » Available in radiotranslucent and radiolucent options
- » Available with leads-in and lead-out packaging



Phillips Medical



Physio-Control



Zoll Medical



Zoll

Item No.	Description	Pkg.	Compatibility
MDSZOL200LO	Pad, Defib, Adult, Radiotran, Leads Out	10/cs	M-Series, R-Series, E-Series, X-Series
MDSZOL200	Pad, Defib, Adult, Radiotran, Leads In	10/cs	M-Series, R-Series, E-Series, X-Series
MDSZOL100LO	Pad, Defib, Adult, Radiolucent, Leads Out	10/cs	M-Series, R-Series, E-Series, X-Series
MDSZOL100	Pad, Defib, Adult, Radiolucent, Leads In	10/cs	M-Series, R-Series, E-Series, X-Series
MDSPEDZOL	Pad, Defib, Pedi, Radiolucent, Leads In	10/cs	M-Series, R-Series, E-Series, X-Series

Medline defibrillation electrodes are NOT compatible with Zoll AED Plus and Zoll AED Pro

Physio-Control

Item No.	Description	Pkg.	Compatibility
MDSPHY200LO	Pad, Defib, Adult, Radiotran, Leads Out	10/cs	LP-9P, LP-10C, LP-12, LP-15, LP-20, LP-500, LP-500 DPS, LP-1000
MDSPHY200	Pad, Defib, Adult, Radiotran, Leads In	10/cs	LP-9P, LP-10C, LP-12, LP-15, LP-20, LP-500, LP-500 DPS, LP-1000
MDSPHY100LO	Pad, Defib, Adult, Radiolucent, Leads Out	10/cs	LP-9P, LP-10C, LP-12, LP-15, LP-20, LP-500, LP-500 DPS, LP-1000
MDSPHY100	Pad, Defib, Adult, Radiolucent, Leads In	10/cs	LP-9P, LP-10C, LP-12, LP-15, LP-20, LP-500, LP-500 DPS, LP-1000
MDSPEDPHY	Pad, Defib, Pedi, Radiolucent, Leads In	10/cs	LP-9P, LP-10C, LP-12, LP-15, LP-20

Pediatric defibrillation electrodes are NOT compatible with LP-500, and LP-500 DPS, LP-1000

Philips Medical

Item No.	Description	Pkg.	Compatibility
MDSPHL200LO	Pad, Defib, Adult, Radiotran, Leads Out	10/cs	XLT, HeartStart XL, HP Codemaster XL, MRX, MRX ALS, HeartStart FR2+, HeartStart Forerunner
MDSPHL200	Pad, Defib, Adult, Radiotran, Leads In	10/cs	XLT, HeartStart XL, HP Codemaster XL, MRX, MRX ALS, HeartStart FR2+, HeartStart Forerunner
MDSPHL100LO	Pad, Defib, Adult, Radiolucent, Leads Out	10/cs	XLT, HeartStart XL, HP Codemaster XL, MRX, MRX ALS, HeartStart FR2+, HeartStart Forerunner
MDSPHL100	Pad, Defib, Adult, Radiolucent, Leads In	10/cs	XLT, HeartStart XL, HP Codemaster XL, MRX, MRX ALS, HeartStart FR2+, HeartStart Forerunner
MDSPEDPHL	Pad, Defib, Pedi, Radiolucent, Leads In	10/cs	XLT, HeartStart XL, HP Codemaster XL, MRX, MRX ALS

Pediatric defibrillation electrodes are NOT compatible with HeartStart FR2+ and HeartStart Forerunner

Medline defibrillation electrodes are NOT compatible with HeartStart FR3 and HeartStart FRx

Anderson

Item No.	Description	Pkg.	Compatibility
MDST100LOAC	Pad, Defib, Adult, Radiolucent, Leads Out	10/cs	Compatible with defibrillators standardized to Anderson connector
MDST100AC	Pad, Defib, Adult, Radiolucent, Leads In	10/cs	Compatible with defibrillators standardized to Anderson connector
MDSPED100	Pad, Defib, Pedi, Radiolucent, Leads In	10/cs	Compatible with defibrillators standardized to Anderson connector
MDSC200LOAC	Pad, Defib, Adult, Radiotran, Leads Out	10/cs	Compatible with defibrillators standardized to Anderson connector
MDSC200AC	Pad, Defib, Adult, Radiotran, Leads In	10/cs	Compatible with defibrillators standardized to Anderson connector

Cardiac Science

Item No.	Description	Pkg.	Compatibility
MDSAED100CS	Pad, Defib, Adult, Leads In	10/cs	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro

Point-of-Care Testing

Medline Point-of-Care testing products allow for quick monitoring, diagnoses and appropriate treatments for your patients without the inconvenience of waiting for results from the lab. These tests are fast, easy-to-use, accurate and manufactured under strict quality controls, so you get the reliability and accuracy your patients expect.



MPH1540



MPH1550



MPH3540



MPH3550

EvenCare® G2® Blood Glucose Monitoring System

Advanced technology makes diabetic care easier.

- » **Easy-to-use:** Voice assistance guides you through the testing process and provides results in either English or Spanish
- » **No coding required:** Minimizes errors and inaccurate results
- » **Fast test time:** Results in 6 seconds
- » **Small blood sample size:** 1 µL
- » **Alternate site testing:** Option to test on the forearm, palm or fingertips
- » **Memory option:** Stores 300 test times and can provide 7-, 14- and 30-day averaging
- » **Alarm option:** Set reminders for up to 3 different times

Item No.	Description	Pkg.
MPH1540	Glucose Meter	1/ea
MPH1550	Test Strips	50/bx, 12 bx/cs
MPH1550Z		50/bx
MPH1560	Control Solution High/Low	1/bx, 6 bx/cs
MPH1560Z		1/bx
MPHBAAAZ	Replacement Battery, AAA	24/bx

EvenCare® G3 Professional Blood Glucose Monitoring System

The EVENCARE G3 Professional Blood Glucose Monitoring System is designed to meet the unique needs of healthcare facilities and their patients.

- » **Quality control lock mode:** Ensures quality control testing is performed
- » **Test strip ejector:** Helps reduce cross contamination by eliminating contact with used strip
- » **No coding required:** Minimizes errors and inaccurate results
- » **Fast test time:** Results in 6 seconds
- » **Small blood sample size:** 0.8µL
- » **Alternate site testing:** Option to test on the forearm, palm or fingertips

Item No.	Description	Pkg.
MPH3540	Glucose Meter	1/ea
MPH3550	Test Strips	50/bx, 12 bx/cs
MPH3550Z		50/bx
MPH3550FL	Foil Wrapped Test Strips	50/bx, 12 bx/cs
MPH3550FLZ		50/bx
MPH3560	Control Solution High/Low	1/bx, 6 bx/cs
MPH3560Z		1/bx
RAYCR2032	Replacement Battery, 3V	10/bx

EvenCare® ProView™ Blood Glucose Monitoring System

The EVENCARE ProView blood glucose meter provides simple and accurate diabetes management for resident care.

- » **Large backlit display:** Can be clearly read in low light
- » **Ergonomic design:** Has a comfortable, non-slip grip
- » **Test strip ejector:** Eliminates contact with used strip

Item No.	Description	Pkg.
MPH4540	Glucose Meter	1/ea
MPH4550	Test Strips	50/bx, 12 bx/cs
MPH4550Z		50/bx
MPH4560	Control Solution High/Low	1/bx, 6 bx/cs
MPH4560Z		1/bx
RAYCR2032	Replacement Battery, 3V	10/bx
MPHEMSGLCS	EMS Cases for Blood Glucose Monitoring Systems	240/cs
MPHEMSGLCSH		1/ea
MPH4540EMSCS	ProView EMS Protective Sleeve	1/ea



Harmony™ Blood Glucose Monitoring System

The Harmony blood glucose meter utilizes the newest technology in blood glucose testing in order to deliver more accurate results.

- » Test strips are designed to meet greater accuracy standards and limit interference with therapies such as oxygen
- » Wider hematocrit range allows for testing a wider patient population
- » **Test strip ejector:** Eliminates contact with used strip
- » **Large backlit display:** can be clearly read in low light

Item No.	Description	Pkg.
MPH6540	Glucose Meter	1/ea
MPH6550	Test Strips	50/bx, 12 bx/cs
MPH6550Z		50/bx
MPH6560	Control Solution High/Low	1/bx, 6 bx/cs
MPH6560Z		1/bx
RAYCR2032	Replacement Battery, 3V	10/bx
MPHEMSGLCS	EMS Case for Blood Glucose Monitoring Systems	240/cs
MPHEMSGLCSH		1/ea
MPH6540EMSCS	Harmony EMS Protective Sleeve	1/ea



Hemocue® Glucose 201

The Hemocue Glucose 201 system puts lab equivalent answers in health professionals' hands when they're needed most—at the point of care. Not only does the unique cuvette technology enable lab accuracy in just three simple steps, but it reduces serious contamination risks.

- » For diagnosing and monitoring diabetes mellitus
- » Reduces risk of missed hypoglycemia in newborns
- » No clinically significant lot-to-lot variation
- » Fixed analyzer calibration
- » Individually wrapped microcuvettes avoid contamination and maximize shelf life

Item No.	Description	Pkg.
HMQ120706	MG/DL, Hemocue Glucose 201 Analyzer	1/ea
HMQ121422	MG/D, Hemocue Glucose 201 DM Analyzer	1/ea
HMQ110706	Glucose 201 Vial Microcuvettes	100/bx
HMQ110723	Glucose 201 Individually Packaged Microcuvettes	50/bx
HMQ180013002	GlucoTrol-AQ Aqueous High & Low Control	1/kit
HMQ120762	Glucose DMU docking Station	1/ea



HMQ121422

Coag-Sense™ PT/INR System

The Coag-Sense PT/INR (Prothrombin Time/International Normalized Ratio) system provides a more direct approach to INR testing as the only portable monitoring system that directly detects clotting endpoint, providing greater confidence in results. Features convenient testing with no coding required, rapid and reliable results and a wide hematocrit range of 15–60%.

Item No.	Description	Pkg.
ABO03P7002	PT2 Meter, Professional User's Manual, Professional Quick Reference Guide, Auto Single Use 21G Lancets - starter bag of 8, AC wall adapter, 54 Sample Transfer Tubes with Plungers, carrying case	1/ea
ABO03PTESTKIT	Test Strip Kit—includes 54 transfer tubes (ABO03P5254Z)	50/bx
ABO03P6910Z	Control Strip Kit—includes 5 high and 5 low control strips and 1 vial of solution	10/bx
ABO03P5254	Transfer Tubes	54/vial



ABO03P7002

Adjustable-Depth Safety Lancets

- » Innovative rotating dial allows for three easy depth adjustments: 1.3mm, 1.8mm, 2.3mm
- » Pressure activated
- » Compares to the Accu-Chek Safe-T-Pro Plus

Item No.	Description	Pkg.
MPHADJDEPTH	Adjustable-Depth Safety Lancet	1,000/cs
MPHADJDEPTHZ		200/bx



MPHADJDEPTH

Safety Lancets

- » Simple and safe for both healthcare professionals and patients
- » Needle retracts automatically and cannot be re-used protecting patients and healthcare professionals from needlestick injuries
- » Different sizes allow healthcare professionals to obtain the proper amount of blood for capillary testing
- » Quick and easy-to-use
- » Double spring action and controlled needle depth provides greater patient comfort



Item No.	Description	Pkg.
Push Button		
MPHSAFETY28	28 Gauge, 1.6 mm	200/bx, 18 bx/cs
MPHSAFETY28Z		200/bx
MPHSAFETY285		50/bx, 24 bx/cs
MPHSAFETY285Z		50/bx
MPHSAFETY281		100/bx, 20 bx/cs
MPHSAFETY281Z		100/bx
MPHSAFETY23	23 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHSAFETY23Z		200/bx
MPHSAFETY21	21 Gauge, 1.8 mm	200/bx, 18 bx/cs
MPHSAFETY21Z		200/bx
MPHSAFETYBL	Blade, 1.6 mm	200/bx, 18 bx/cs
MPHSAFETYBLZ		200/bx


Pressure Activated

MPHST28	28 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHST28Z		200/bx
MPHST21	21 Gauge, 2.2 mm	200/bx, 10 bx/cs
MPHST21Z		200/bx

Item No.	Description	Pkg.
Side Button		
MPHSTIDE28	28 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHSTIDE28Z		200/bx
MPHSTIDE21	21 Gauge, 2.4 mm	200/bx, 10 bx/cs
MPHSTIDE21Z		200/bx

Contact Activated

MPHPRESS30	30 Gauge, 1.5 mm	100/bx, 20 bx/cs
MPHPRESS30Z		100/bx
MPHPRESS28	28 Gauge, 1.8 mm	100/bx, 20 bx/cs
MPHPRESS28Z		100/bx
MPHPRESS26	26 Gauge, 1.8 mm	100/bx, 20 bx/cs
MPHPRESS26Z		100/bx
MPHPRESS23	23 Gauge, 2.2 mm	100/bx, 20 bx/cs
MPHPRESS23Z		100/bx
MPHPRESS21	21 Gauge, 2.2 mm	100/bx, 20 bx/cs
MPHPRESS21Z		100/bx



MSC0002



MPH100030



MPHLANCINGPEN

Monofilament

- » Reliable, cost effective method of testing for the presence or absence of sensation
- » Helps identify patients at risk of developing neuropathy or diabetic foot ulcers
- » 5.07/10 g monofilament

Item No.	Description	Pkg.
MSC0002	Foot Monofilament,	25/bx
MSC0005Z	Diabetic Neuropathy Test	5/bx

General Purpose Lancets

- » Ultra sharp lancets have a fine-point tip for smooth and fast procedures
- » Can be used with most universal lancing devices

Item No.	Description	Pkg.
MPH100030	30 Gauge	100/bx, 25 bx/cs
MPH100030Z		100/bx
MPH1505515	23 Gauge	200/bx, 100 bx/cs
MPH1505515Z		200/bx

Lancing Device

- » Adjustable lancing device fits all general-purpose lancets
- » Non-safety

Item No.	Description	Pkg.
MPHLANCINGPEN	Medline Adjustable Lancing Device	1/ea



Heel Incision Devices

- » High-quality steel blades make small, precise incisions for reduced bruising
- » Disposable, one-time use devices
- » Easy trigger activates blade quickly, safely and precisely
- » Color-coded for easy identification

Item No.	Description	Pkg.
MPHMPREEMIE2	Micro-Preemie Heel Lancet, Light Blue	50/bx, 4 bx/cs
MPHMPREEMIE2Z		50/bx
MPHPREEMIE2	Preemie Heel Lancet, White	50/bx, 4 bx/cs
MPHPREEMIE2Z		50/bx
MPHNEWBORN2	Newborn Heel Lancet, Pink	50/bx, 4 bx/cs
MPHNEWBORN2Z		50/bx
MPHTODDLER2	Toddler Heel Lancet, Purple	50/bx, 4 bx/cs
MPHTODDLER2Z		50/bx



Pen Needles

- » Easy medication delivery
- » Universal fit
- » Thin wall technology allows for easy flow of medication during injection

Item No.	Description	Pkg.
MPHPN2912	29 Gauge, 12 mm	100/bx, 10 bx/cs
MPHPN2912Z		100/bx
MPHPN315	31 Gauge, 5 mm	100/bx, 10 bx/cs
MPHPN315Z		100/bx
MPHPN318	31 Gauge, 8 mm	100/bx, 10 bx/cs
MPHPN318Z		100/bx
MPHPN324	32 Gauge, 4 mm	100/bx, 10 bx/cs
MPHPN324Z		100/bx



MDS090735

Sterile Alcohol Prep Pads

Item No.	Description	Pkg.
MDS090735	Sterile 2-Ply Alcohol Prep Pads, Size M	3000/cs
MDS090670	Sterile 2-Ply Alcohol Prep Pads, Size L	1000/cs

Ultrasound Gel

- » Water soluble, non-greasy ultrasound gel
- » Available in 8.5-oz. squeeze bottles or 5 L jugs that include one 8.5-oz. empty bottle per 5 L jug

Item No.	Description	Pkg.
MDS092005	Ultrasound Gel, 8.5-oz. Squeeze Bottle, Blue	12/cs
MDS092005H		1/ea
MDS092010	Ultrasound Gel, 5 L Jug, Blue	4/cs
MDS092010H		1/ea
MDS092010C	Ultrasound Gel, 5 L Jug, Clear	4/cs
CTR001305	Ultrasound Gel, 8.5-oz. Squeeze Bottle, Clear	12/cs
CTR001305H		1/ea
CTR001201	Nonsterile Ultrasound Gel, 0.7-oz. Pouch, Clear	100/cs



Controlled Substances Log Book

- » Hardcover log book used for manual recordkeeping of prescription drugs
- » Nursing facilities are required by federal law to maintain records of Schedule II controlled substance inventories and administration
- » This log book will help facilities to comply with these regulations, including 21 U.S.C. 827 (b), 21 U.S.C. 827 (c)-DEA, and 42 CFR 483.60-CMS
- » Consolidate documentation in one readily accessible place for auditing-DEA requires this data be retained for at least 2 years

Item No.	Description	Pkg.
MDSMDLCSL2	Controlled Substances Log Book	1/ea



MDSMDLCSL2

BioCon 750 Bladder Scanner with BladderPoint

Real-Time Bladder Targeting

Quickly and confidently locate the bladder for accurate measurement.

You shouldn't have to do multiple scans to get an acceptable result. BioCon 750 makes it easier than ever to find the bladder with two industry-firsts:

- » Proprietary BladderPoint™ technology puts a bull's-eye in the middle of the bladder
- » Pre-scan shows a live view of the ultrasound image right from the start

Dependable Performance

Reliability is designed into every BioCon 750.

- » Standard five-year warranty includes drop protection
- » Built-in self-diagnostics alert you of any hardware or maintenance issues
- » Tough polycarbonate casing promotes robust and uninterrupted service

Versatile Care

Accommodates various patient populations.

- » One-click, 360-degree measurement requires no probe rocking
- » Probe has a rear grip for easier use on bariatric patients
- » Can be used on pediatric patients

Staff Convenience

With BioCon 750, bladder scanning is quick, simple and efficient.

- » Data and images can be easily uploaded to any EMR system
- » Integrated thermal printer provides documentation at the push of a button
- » Onboard training guide walks staff through scanning technique



Item No.	Description	Pkg.
MDSB10750	BioCon 750 Bladder Scanner w/Printer	1/ea
MDSB10700STND	Mobile Stand w/Basket	1/ea
MDSB10750CK	Calibration Kit	1/ea
MDSBSPHANTOM	Phantom Bladder	1/ea
MDSB10500PAPE	Thermal Paper	23/cs
MDS092005	Ultrasound Gel, 8.5-oz. Squeeze Bottle	12/cs

BioCon 900 Ultrasonic Bladder Scanner with Printer

Real-Time Bladder Viewing

Easily locate the bladder position for accurate measurement. You shouldn't have to do multiple scans to get an acceptable result. BioCon 900 makes it easier to find the bladder with two industry-firsts:

- » BladderPoint™ technology puts a bull's-eye in the middle of the bladder
- » Pre-scan shows a live view of the ultrasound image

Built-In Versatility

Bladder scanning is quick, simple and efficient.

- » Measure bladder volume in 3 seconds
- » Increased scan depth to industry best **230 mm** for individuals with higher BMIs
- » Data and images can be easily uploaded to your EMR system or printed from integrated printer
- » Less cleaning necessary for contact isolation patients



Dependable Performance

BioCon 900 features self-diagnosis and five-year warranty.

- » Durable probe construction helps minimize the chance for damage
- » **Five-year** warranty with drop protection coverage comes standard
- » Self-diagnosis quickly tells you if the probe is working correctly
- » No calibration required

Industry-Leading Support

We provide ongoing education, training and service throughout the device's lifespan.

- » Staff is trained on how to use the equipment through online videos, competency testing and onsite in-servicing
- » Toll-free technical service hotline: 1-866-274-2131
- » Service includes free equipment diagnosis and loaner units



Locking mechanism for safe storage and transport.




Item No.	Description	Pkg.
MDSB10900	BioCon 900 Bladder Scanner w/ Printer	1/ea
MDSMRC900	BioCon 900 Stand	1/ea
MDSBSPHANTOM	Phantom Bladder	1/ea
MDSB10500PAPE	Thermal Paper	12/ea, 12/cs
MDS092005	Ultrasound Gel, 8.5-oz. Squeeze Bottle	12/cs

Safety Syringes with Needles

- » Preassembled needle-and-syringe combinations with features to help ensure safety for nurses and patients

Item No.	Description	Pkg.
SYRS101235FZ	1 ml 23G x 1" Safety Syringe	100/bx
SYRS101255FZ	1 ml 25G x 1" Safety Syringe	100/bx
SYRS101272FZ	1 ml 27G x 1/2" Safety Syringe	100/bx
SYRS101392FZ	1 ml 30G x 1/2" Safety Syringe	100/bx
SYRS103205Z	3 ml 20G x 1" Safety Syringe	100/bx
SYRS103207Z	3 ml 20G x 1 1/2" Safety Syringe	100/bx
SYRS103215Z	3 ml 21G x 1" Safety Syringe	100/bx
SYRS103217Z	3 ml 21G x 1 1/2" Safety Syringe	100/bx
SYRS103225Z	3 ml 22G x 1" Safety Syringe	100/bx
SYRS103227Z	3 ml 22G x 1 1/2" Safety Syringe	100/bx
SYRS103235Z	3 ml 23G x 1" Safety Syringe	100/bx
SYRS103253Z	3 ml 25G x 5/8" Safety Syringe	100/bx
SYRS103255Z	3 ml 25G x 1" Safety Syringe	100/bx
SYRS105207Z	5 ml 20G x 1 1/2" Safety Syringe	100/bx
SYRS105215Z	5 ml 21G x 1" Safety Syringe	100/bx
SYRS105217Z	5 ml 21G x 1 1/2" Safety Syringe	100/bx
SYRS110207Z	10 ml 20G x 1 1/2" Safety Syringe	50/bx
SYRS110217Z	10 ml 21G x 1 1/2" Safety Syringe	50/bx
ULT101253Z	1 ml 25G x 5/8" Tuberculin Safety Syringe	100/bx



SYRS103225Z



SYRS103205Z



ULT101253Z

Insulin Safety Syringes

- » Needles can be temporarily covered by pulling the shield forward; after injection, the shield is permanently locked in place by pulling forward and twisting
- » Compliance with OSHA Bloodborne Pathogens Standard for Engineering Controls
- » Large barrel flanges for a good grip and control
- » Low dead-space design

Item No.	Description	Pkg.
SYRS101292Z	Insulin Safety Syringe with Needle, 1 ml, 29G x 0.5"	100/bx
SYRS193292Z	Insulin Safety Syringe with Needle, 0.5 ml, 29G x 0.5"	100/bx



SYRS101292Z

SYRS193292Z

Blunt Tip Needles

- » Designed for drawing up medication from vials, high flow rate, non-coring

Item No.	Description	Pkg.
SYR110022Z	18G x 1.5"	100/bx
SYR110022FZ	18G x 1.5", Filter	100/bx



SYR110022Z

Standard Hypodermic Syringes with Needles

- » Syringe with needle
- » Available in various mL and needle sizes

Item No.	Description	Pkg.
SYR101272Z	1 ml Luer Slip Disposable Syringe with 27G x ½" Needle	100/bx
SYR103205Z	3 ml Luer Lock Disposable Syringe with 20G x 1" Needle	100/bx
SYR105205Z	5 ml Luer Lock Disposable Syringe with 20G x 1" Needle	100/bx



SYR101272Z



SYR103205Z

Safety Hypodermic Needles

- » Advanced features ensure safety for nurses and patients

Item No.	Description	Pkg.
SYRS100185Z	18G x 1" Safety Needle	100/bx
SYRS100187Z	18G x 1½" Safety Needle	100/bx
SYRS100205Z	20G x 1" Safety Needle	100/bx
SYRS100207Z	20G x 1½" Safety Needle	100/bx
SYRS100225Z	22G x 1" Safety Needle	100/bx
SYRS100227Z	22G x 1½" Safety Needle	100/bx
SYRS100253Z	25G x ¾" Safety Needle	100/bx
SYRS100255Z	25G x 1" Safety Needle	100/bx
SYRS100257Z	25G x 1½" Safety Needle	100/bx
SYRS100272Z	27G x ½" Safety Needle	100/bx
SYRS100273Z	27G x ¾" Safety Needle	100/bx
SYRS100276Z	27G x 1¼" Safety Needle	100/bx
SYRS100302Z	30G x ½" Safety Needle	100/bx



SYRS100185Z



SYRS100227Z



SYRS100257Z

Luer Lock Syringes

- » Syringe only with luer lock
- » Available in various mL sizes
- » Sterile and individually packaged
- » Luer slip versions also available

Item No.	Description	Pkg.
SYR101010Z	1 ml Syringe	100/bx
SYR103010Z	3 ml Syringe	100/bx
SYR105010Z	5 ml Syringe	100/bx
SYR110010Z	10 ml Syringe	100/bx
SYR120010Z	20 ml Syringe	40/bx
SYR130010Z	30 ml Syringe	40/bx
SYR160010Z	60 ml Syringe	40/bx
SYRTCAPZ	Luer Lock Tip Caps	200/bx



SYR160010Z



SYR110010Z



SYR101010Z

SYRTCAPZ

Allergy Syringe Trays

- » Designed to deliver small doses of medication easily and accurately for efficient patient care
- » Smooth-gliding rubber stopper promotes accurate medication measurement
- » Low dead-space minimizes allergen waste
- » Easy-to-read number scale helps you prepare precise dose
- » Sterile packaged 25 per tray for easy accessibility; trays stack for ease and convenience

Item No.	Description	Pkg.
SYRA101272	27G x 0.5", 1 ml	1,000/cs
SYRA101261	26G x 3/4", 1 ml	1,000/cs



SYRA101272

ZR® Prefilled Flush Syringes

- » Help minimize the risk of catheter occlusion to enhance patient safety and help improve outcomes
- » Clearly marked saline or heparin labels help lower the risk of medication error
- » Not made with natural rubber latex and preservative free
- » Terminally sterilized for added safety
- » Barcoded to speed inventory control and minimize errors

Item No.	Description	Pkg.
EMZE010001	10 ml Saline in 10 ml Syringe	600/cs
EMZE010055	5 ml Saline in 10 ml Syringe	600/cs
EMZE010035	3 ml Saline in 10 ml Syringe	800/cs
EMZE060055	5 ml of 100U/ml Heparin in 10 ml Syringe	600/cs
EMZE050055	5 ml of 10U/ml Heparin in 10 ml Syringe	600/cs
EMZ10091280	Sterile Field 10 ml Saline in 10 ml Syringe	280/cs



EMZE010001



EMZE060055

SwabFlush Pre-filled Saline Flush Syringe

- » SwabFlush drives compliance and standardization of disinfection protocols
- » SwabCap®, conveniently integrated in the syringe, protects access valve in between line accesses
- » SwabCap disinfects valve in 5 minutes and keeps valve surface disinfected for up to 7 days if not removed

Item No.	Description	Pkg.
EMZE010301	SwabFlush 10/10 ml with SwabCap	600/cs
EMZE010335	SwabFlush 3/10 ml with SwabCap	800/cs
EMZE010355	SwabFlush 5/10 ml with SwabCap	600/cs



EMZE010301

OnGuard™ Blood Collection Holder

- » Helps your institution comply with OSHA regulations
- » Compatible with most blood collection sets, needles and tubes
- » Grooves for easy handling

Item No.	Description	Pkg.
DYNDBCHLDRZ	Blood Collection Tube Holder	250/bag



DYNDBCHLDRZ

OnGuard™ Blood Collection Sets

- » Large anti-skid wings for enhanced control
- » Sturdy non-slip grip to engage the safety feature
- » Compliant with OSHA Bloodborne Pathogens Standard

Item No.	Description	Pkg.
DYNDTNBCS21G	Blood Collection Set, 21 G x 0.75"	200/cs
DYNDTNBCS23G	Blood Collection Set, 23 G x 0.75"	200/cs
DYNDTNBCS25G	Blood Collection Set, 25 G x 0.75"	200/cs



DYNDTNBCS23G

Patient Room Sharps Containers

- » Lids designed to prohibit access to contents
- » Lid design facilitates one-handed operation
- » Easy-to-see fill levels
- » Multiple cabinet and bracket options

Item No.	Description	Pkg.
MDS705153	5-qt, Red, Wall-Mount Sharps Container	20/cs
MDS705154	5-qt, Clear, Wall-Mount Sharps Container	20/cs
MDS705253	5.4-qt, Red, Wall-Mount Sharps Container	20/cs
MDS705254	5.4-qt, Clear, Wall-Mount Sharps Container	20/cs
MDS705205	3-gal, Red, Sharps Container	12/cs
MDS705110	1-qt, Red, Phlebotomy Sharps Container	100/cs
MDS705115	1.5-qt, Red, Phlebotomy Screw Top Lid	20/cs
MDS707DART	Portable Sharps Container Snap Lid	24/cs
MDS707953	Locking Cabinet for 5-qt Sharps Container	5/cs
MDS707930	Locking Cabinet for 3-gal Sharps Container	5/cs
MDS707212	Glove Box Holder for 5-qt and 3-gal Cabinets	10/cs
MDS707315	Locking Wire Bracket for 5-qt Sharps Container	5/cs
MDS707316	Non-Locking Wire Bracket for 5-qt Sharps Container	5/cs
MDS707011	Non-Locking Wire Bracket for 1-qt Sharps Container	5/cs



Large PG-II Containers

- » Tested to comply with Department of Transportation (DOT) requirements to eliminate the need for secondary packaging upon disposal
- » Lids feature a gasket seal to resist leaking
- » Integrated side handles make containers easy to carry

Item No.	Description	Pkg.
----------	-------------	------

Chemotherapy Waste

MDS706308	8-gal., Yellow, Port Lid	9/cs
MDS706312	12-gal., Yellow, Port Lid	8/cs
MDS706218	18-gal., Yellow, Flat Lid	7/cs
MDS706318	18-gal., Yellow, Port Lid	7/cs

Pharmaceutical Waste

MDS707308	8-gal., White, Port Lid	9/cs
MDS707318	18-gal., White, Port Lid	7/cs

RCRA Hazardous Waste

MDS708308	8-gal., Black, Port Lid	9/cs
MDS708318	18-gal., Black, Port Lid	7/cs



MDS706312

MDS706308

Multi-purpose Sharps Containers

- » Variety of lids and size options
- » Nestable containers save storage space
- » Integrated handles for easy transport
- » Rolling carts available

Item No.	Description	Pkg.
MDS705201	1-gal., Red, Wall/Free, Slide Lid	32/cs
MDS701201	1-gal., Translucent Red, Wall/Free, Slide Lid	32/cs
MDS705202	2-gal., Red, Wall/Free, Slide Lid	20/cs
MDS701202	2-gal., Translucent Red, Wall/Free, Slide Lid	20/cs
MDS705208HN	8-gal., Red, Hinged Lid	10/cs
MDS705208SL	8-gal., Red, Slide Lid	10/cs
MDS705208ST	8-gal., Red, Star Lid	10/cs
MDS705210HN	10-gal., Red, Hinged Lid	6/cs
MDS705210SL	10-gal., Red, Slide Lid	6/cs
MDS705210ST	10-gal., Red, Star Lid	6/cs
MDS705212HN	12-gal., Red, Hinged Lid	10/cs
MDS705212SL	12-gal., Red, Slide Lid	10/cs
MDS705212ST	12-gal., Red, Star Lid	10/cs
MDS705218HN	18-gal., Red, Hinged Lid	5/cs
MDS705218SL	18-gal., Red, Slide Lid	5/cs
MDS705218ST	18-gal., Red, Star Lid	5/cs
MDS707CART	Rolling Wire Cart for 8, 10, 12 or 18-gal. Containers	1/ea



SureSite IV® Auto Safety IV Catheter

- » Push-button safety activation
- » Fully encapsulating needle shield
- » Notched needle for initial flash observation in sizes 20g-24g
- » Flashback chamber for clear visual confirmation
- » Manufactured with PUR (Polyurethane)

Item No.	Description	Pkg.
DYNCA14175	14G x 1.75", SureSite Auto	200/cs, 50/bx
DYNCA16177	16G x 1.77", SureSite Auto	200/cs, 50/bx
DYNCA18116	18G x 1.16", SureSite Auto	200/cs, 50/bx
DYNCA18188	18G x 1.88", SureSite Auto	200/cs, 50/bx
DYNCA20100	20G x 1", SureSite Auto	200/cs, 50/bx
DYNCA20116	20G x 1.16", SureSite Auto	200/cs, 50/bx
DYNCA20188	20G x 1.88", SureSite Auto	200/cs, 50/bx
DYNCA22100	22G x 1", SureSite Auto	200/cs, 50/bx
DYNCA24075	24G x 0.75", SureSite Auto	200/cs, 50/bx



DYNCA24075

SureSite IV® Slide Safety IV Catheter

- » Clicking lock safety mechanism
- » Self-sheathing needle shield
- » Flashback chamber for clear visual confirmation
- » Manufactured with PUR (Polyurethane)

Item No.	Description	Pkg.
DYNCS14114	14G x 1 ¼", SureSite Slide	200/cs, 50/bx
DYNCS16114	16G x 1 ¼", SureSite Slide	200/cs, 50/bx
DYNCS18114	18G x 1 ¼", SureSite Slide	200/cs, 50/bx
DYNCS20114	20G x 1 ¼", SureSite Slide	200/cs, 50/bx
DYNCS20100	20G x 1", SureSite Slide	200/cs, 50/bx
DYNCS22100	22G x 1", SureSite Slide	200/cs, 50/bx
DYNCS24340	24G x ¾", SureSite Slide	200/cs, 50/bx




DYNCS22100



DYNDTC5081



DYNDTN2522



DYNDTN2524

IV Extension Sets

- » Clear, flexible tubing
- » Male and female luer locks
- » Not made with DEHP

Item No.	Description	Pkg.
DYNDTN5077	6.5" Set with Nonbonded Needle-Free Valve, Standard Bore, 0.7 ml	100/cs
DYNDTN5081	8" Set with Nonbonded Needle-Free Valve, Small Bore, 0.3 ml	100/cs
DYNDTC5002	7" Set, Bonded Needle-Free Valve, Micro Bore, 0.1 ml	100/cs
DYNDTC5082	8" Set with Needle-Free T Connector, Small Bore, 0.4 ml	100/cs
DYNDTN2522	20" Set, Flow Controller and 1 Injection Site, Standard Bore, 3.0 ml	50/cs
DYNDTN1603	60" Set, Small Bore, 2.0 ml	100/cs
DYNDTN0007	60" Set, Mini Bore, 0.6 ml	100/cs

Gravity IV Administration Sets

- » Crystal clear tubing and drip chambers
- » Flat, clear needle-free Y-sites
- » Not compatible with infusion pumps

Item No.	Description	Pkg.
DYNDTN1545	15 drops/ml, 3 Y-site, 109"	50/cs
DYNDTN1512	15 drops/ml, 2 Y-site, 96"	50/cs
DYNDTN1529	15 drops/ml, 1 Y-site, 92"	50/cs
DYNDTN0555	60 drops/ml, 3 Y-site, 111"	50/cs
DYNDTN0512	60 drops/ml, 2 Y-site, 96"	50/cs
DYNDTN2524	Flow Controller, 1 Y-site, 92"	50/cs
DYNDTN1540A	Secondary, 40" w/Hanger	50/cs



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)

We reserve the right to correct any errors that may occur within this brochure. Some products may not be available for sale in Mexico or Canada.

©2019 Medline Industries, Inc. All rights reserved. Coag-Sense is a trademark of Coagusense, Inc. e-scope is a registered trademark of Excellus Technologies, Inc. Philips is a registered trademark of Koninklijke Philips Electronics N.V. Limited Liability Company. R.A.T.E. is a trademark of BestMed. TemporalScanner and Exergen is a registered trademark of Exergen Corporation. Tyvek is a registered trademark of E. I. Du Pont De Nemours And Company. SwabCap is a registered trademark of ICU Medical, Inc. Harmony and ProView are trademarks and Compli-mates, SureSite IV, Synergy, EVENCARE, G2, MedGel and Medline are registered trademarks of Medline Industries, Inc. CareConnection is a registered trademark of Constant Care Technology

MKT19116748 / LIT007A / 9M / QG / 23