

Personal Protective Apparel

Head-to-Toe Protection for Healthcare



**Know
you're
covered.**



Protect your people **and** **your bottom line.**

You want to know you're protected. Our protective apparel goes through rigorous testing and meets strict FDA guidelines. We follow the highest quality standards to make sure you get the protection you need.

As a leader in protective apparel, we also have the most comprehensive product line and low acquisition costs. You can find the exact products you need—and maintain cost efficiencies.

AAMI Level
isolation gowns
are color-coded
by neck binding

Recommended Gowns 4-5

AAMI Isolation Gowns 6-7

Cover Gowns 8-9

Protective Gowns 10-11

Footwear 12-13

Headwear and Beard Cover 14-15

Scrubs 16

Coveralls 17

Lab Coats 18-19

**Breast Cancer Month
Awareness PPE** 20-21

Face Masks 22

Patient Wear 23

Choose the right gown.

Our recommendations for your workplace environment.

A Housewide, Med/Surg and ICU

Our QuickComply™ AAMI Level 2 gowns feature an over-the-head style which improves compliance while reducing donning and removing time. Side ties eliminate the difficulties of tying behind the back, while thumb loop wrists allow for a good gown and glove fit.

Item No.	Description
NONLV240	Regular/Large
NONLV240XL	X-Large

Additional items on page 7

B Endoscopy and Emergency

AAMI Level 3 gowns are made from comfortable and breathable premium multi-ply material. Features a classic tie waist, hook and loop neck and knit cuffs.

Item No.	Description
NONLV325	Regular/Large
NONLV325XL	X-Large

Additional items on page 7



A Housewide, Med/Surg and ICU

B Endoscopy and Emergency



C Decontamination



D Chemotherapy

C Decontamination

Made from premium polyethylene coated materials that have been tested against ASTM F1670 and ASTM F1671. An over-the-head, open back style reduces donning and removing time, while thumb loop wrists allow for a good gown and glove fit.

Item No.	Description
NON26276	Regular/Large
NON26276XL	X-Large
NON26276XXL	XX-Large

Additional items on page 11

D Chemotherapy

Protective, polyethylene-coated materials & seams tested against top 11 chemotherapy drugs per ASTM F739-12. Features hood and loop neck, waist ties and knit cuffs. Taped seams to meet USP 800 guidelines.

Item No.	Description
NON25457M	Medium
NON25457	Regular/Large
NON25457XL	X-Large
NON25457XXL	XX-Large

Additional items on page 11

Why AAMI?

AAMI Guidelines: A Reliable Benchmark for Protection

Under AAMI Guidelines, protective apparel must meet minimum performance standards for strength, barrier protection and fluid resistance to ensure the safety of healthcare workers and patients.



Levels of Barrier Protection – AAMI PB70:20121

AAMI Level	Test	Requirements at 4% AQL	Anticipated Risk of Exposure
Level 1 💧	AATCC 42:2000	Water impact ≤ 4.5 g	Low
	AATCC 42:2000		
Level 2 💧💧	AATCC 127:1998	Spray impact ≤ 1.0 g	Moderate
	AATCC 42:2000	Hydrostatic Pressure ≥ 20 cm	
Level 3 💧💧💧	AATCC 42:2000	Spray Impact ≤ 1.0 g	Moderate to High
	ASTM F3352-19	Hydrostatic pressure ≥ 50 cm	

AAMI Level 1 Isolation Gown

Recommended for: Med/Surg, Visitors, Contact Precautions

- Made from multi-ply material
- Side ties eliminate the difficulties of tying behind the back
- Choose between over-the-head or neck tie styles

Item No.	Description	Pkg.
Medium Weight, Neck Ties, Elastic Wrists		
NONLV100	Regular/Large, Yellow	100/cs
NONLV100XL	X-Large, Yellow	100/cs
Medium Weight, Over-the-Head, Thumb Loop Wrists		
NONLV140	Regular/Large, Yellow	100/cs
NONLV140XL	X-Large, Yellow	100/cs

AAMI Level 2 Isolation Gown

Recommended for: Med/Surg, ICU, Emergency, Airborne and Contact Precautions

- Made from multi-ply material with side ties to eliminate the difficulties of tying behind the back
- Choose between over-the-head or tape tab neck styles
- Available with elastic wrists or thumb loops

Item No.	Description	Pkg.
Medium Weight, Tape Tab Neck, Elastic Wrists		
NONLV200M	Medium, Yellow	100/cs
NONLV200	Regular/Large, Yellow	100/cs
NONLV200XL	X-Large, Yellow	100/cs
NONLV200XXL	XX-Large, Yellow	100/cs
NONLV200XXXL	XXX-Large, Yellow	100/cs
Medium Weight, Tape Tab Neck, Thumb Loops, Individually Folded, Dispensing Box		
NONLV200TL	Regular/Large, Yellow	40/cs
NONLV200TLXL	X-Large, Yellow	40/cs
Medium Weight, Over-the-Head, Thumb Loops		
NONLV240	Regular/Large, Yellow	100/cs
NONLV240XL	X-Large, Yellow	100/cs
NONLV240XXL	XX-Large, Yellow	100/cs
NONLV240XXXL	XXX-Large, Yellow	100/cs
Lightweight, Over-the-Head, Thumb Loops		
NONLV242	Regular/Large, Yellow	100/cs
NONLV242XL	X-Large, Yellow	100/cs

AAMI Level 2 Isolation Gown, continued

Poly Coated recommended for: Decontam/SPD, Endoscopy, Trauma and ER

- Polyethylene coated material tested to ASTM F1670 for blood and bodily fluid penetration
- Available over-the-head
- Available with thumb loops or knit cuffs

Item No.	Description	Pkg.
Polyethylene Coated, Over-the-Head, Knit Cuff		
NON27340	Regular/Large, Yellow	100/cs
NON27340XL	X-Large, Yellow	100/cs
Polyethylene Coated, Over-the-Head, Thumb Loops		
NON27360	Regular/Large, Yellow	100/cs
NON27360XL	X-Large, Yellow	100/cs
Polyethylene Coated, Over-the-Head, Thumb Loops with Elastic Wrists		
NON27380	Regular/Large, Blue	100/cs
NON27380XL	X-Large, Blue	100/cs

AAMI Level 3 Isolation Gown

Recommended for: ICU, Emergency/Trauma, Airborne, Contact and Droplet Precautions

- Comfortable and breathable materials with classic tie waist
- Choose between tape tab, hook and loop and over the head styles
- Available with elastic wrists, knit cuffs, or thumb loops

Item No.	Description	Pkg.
Heavyweight Multi-Ply, Tape Tab Neck, Elastic Wrists		
NONLV315	Regular/Large, Yellow	100/cs
NONLV315XL	X-Large, Yellow	100/cs
Premium Multi-Ply, Hook and Loop Neck, Knit Cuffs		
NONLV325	Regular/Large, Blue	50/cs
NONLV325XL	X-Large, Blue	50/cs
Medium Weight Multi-Ply, Over-The-Head, Thumb Loop Wrists		
NONLV340	Regular/Large, Yellow	50/cs
Microporous, Tape Tab Neck, Knit Cuffs		
NONLV350	Regular/Large, White	50/cs
NONLV350XL	X-Large, White	50/cs

Gown Features



Side ties

Available on AAMI Level 1, 2 and 3 gowns



Neck tie style

Available on AAMI Level 1 gowns



Thumb loop wrists

Available on AAMI Level 1, 2 and 3 gowns



Tape tab neck style

Available on AAMI Level 2 and 3 gowns



CRI4000



Neck ties help ensure the gown remains securely in place.

Classic Cover Gown

Recommended for: Visitors, Food Services

- Made from lightweight polypropylene and multi-layer materials
- Generously sized for full coverage and flexibility
- Side ties eliminate difficulties of tying behind the back
- Neck ties and elastic wrists

Item No.	Description	Pkg.
Lightweight Polypropylene		
CRI4000	Regular/Large, Yellow	50/cs
CRI4000B	Regular/Large, Blue	50/cs
CRI4001	X-Large, Yellow	50/cs
CRI4001G	X-Large, Green	50/cs
CRI4003	XXX-Large, Yellow	50/cs
Lightweight Multi-Layer		
NON27236	Regular/Large, Yellow	50/cs
NON27236XL	X-Large, Yellow	50/cs
Medium Weight Multi-Layer		
NON27SMS2	Regular/Large, Yellow	100/cs
NON27SMS2XL	X-Large, Yellow	100/cs
NON27SMS2XXXL	XXX-Large, Yellow	100/cs
NON27SMS3	Regular/Large, Blue	100/cs
NON27SMS3XL	X-Large, Blue	100/cs



NON27SMS5



Convenient over-the-head, open back style reduces donning and removing time.

Over-the-Head, Open Back Cover Gown

Recommended for: Visitors, Food Services

- Available in lightweight polypropylene, multi-layer and medium weight multi-layer materials
- Generously sized for coverage and flexibility
- Over-the-head style reduces donning and removing time
- Thumb loop wrists

Item No.	Description	Pkg.
Lightweight Polypropylene		
CRI4200	Regular/Large	50/cs
CRI4201	X-Large	50/cs
Lightweight Multi-Layer		
NON27SMS25	Regular/Large	100/cs
NON27SMS25XL	X-Large	100/cs
Medium Weight Multi-Layer		
NON27SMS5	Regular/Large	100/cs
NON27SMS5XL	X-Large	100/cs

Polyethylene Gowns

Recommended for: Decontamination, Emergency, ICU, Food Services, Airborne and Contact Precautions and Droplet Precautions

- Material tested against ASTM F1670 for blood and bodily fluid penetration
- Over-the-head style reduces donning and removing time
- Thumb loop wrists help to prevent sleeves from sliding
- Available in bulk or individually folded in a dispenser box

Item No.	Description	Pkg.
----------	-------------	------

Polyethylene, Blue

CRI5000	Regular/Large	75/cs
CRI5010	Universal	75/cs
CRI5001	X-Large	75/cs
CRI5003	XXX-Large	75/cs

Premium Polyethylene, Blue

NONTH150	Regular/Large	75/cs
NONTH180	Universal	75/cs
NONTH200	X-Large	75/cs

Over-the-Head Microporous Gown

Recommended for: Endoscopy, GI Lab, ICU, Emergency, Maternity

- Made of microporous breathable material tested against ASTM F1670
- Over-the-head style reduces donning and removing time
- Thumb loop wrists prevent sleeve slide
- Individually folded and in a dispenser box

Item No.	Description	Pkg.
----------	-------------	------

NONTH350	Regular/Large	40/cs
NONTH400	X-Large	40/cs

Learn more
about bouffant
caps and shoe
covers.

See pages 12-14



CRI5010

NONTH350

Polyethylene-Coated Gowns

Recommended for: Decontamination, ICU, Endoscopy, GI Lab

- Material tested against ASTM F1670 for blood and bodily fluid penetration
- Neck and side ties
- Available with elastic wrists or knit cuffs

Item No.	Description	Pkg.
Polyethylene Coated, Yellow		
CRI4010	Elastic Wrists, Regular/Large	50/cs
CRI4011	Elastic Wrists, X-Large	50/cs
CRI4010CF	Knit Cuffs, Regular/Large	50/cs
CRI4011CF	Knit Cuffs, X-Large	50/cs
CRI4012CF	Knit Cuffs, XX-Large	50/cs
Premium Polyethylene Coated, Blue		
NON27114	Elastic Wrists, Regular/Large	50/cs
NON27114XL	Elastic Wrists, X-Large	50/cs
NON27116	Knit Cuffs, Regular/Large	50/cs
NON27116XL	Knit Cuffs, X-Large	50/cs
NON27116XXL	Knit Cuffs, XX-Large	50/cs

Chemotherapy Gowns

Recommended for: Oncology, Radiology, Lab, Pharmacy

- Laminate ChemoAire™ material tested against ASTM F1670 & F1671
- Materials and seams tested against top 11 chemotherapy drugs per ASTM F739-12 and ASTM D6978
- Hook and loop neck, waist ties and knit cuffs
- Taped seams to meet USP 800 guidelines

Item No.	Description	Pkg.
Protective Polyethylene Coated AAMI Level 2		
NON25457M	Medium	30/cs
NON25457	Regular/Large	30/cs
NON25457XL	X-Large	30/cs
NON25457XXL	XX-Large	30/cs
ChemoAire™ Premium Breathable Material		
NON29457	Regular/Large	30/cs
NON29457XL	X-Large	30/cs

Chemotherapy gowns tested against the following chemotherapy drugs per ASTM F739-12:

- Dacarbazine
- Carmustine (BiCNU)
- Cyclophosphamide (Cytoxan)
- Doxorubicin Hydrochloride
- 5-Fluorouracil
- Cisplatin
- Etoposide (Toposar)
- Paclitaxel (Taxol)
- Thio-Tepa
- Mitomycin C
- Vincristine Sulfate

Over-the-Head Polyethylene-Coated Gowns

Recommended for: Decontamination, Emergency, ICU, Maternity

- Polyethylene Coated materials tested against ASTM F1670
- Premium Polyethylene Coated materials tested against ASTM F1670, ASTM F1671
- Over-the-head style reduces donning and removing time
- Available with thumb loop wrists or knit cuffs

Item No.	Description	Pkg.
Polyethylene Coated		
NON272790H	Knit Cuffs, Regular/Large, Blue	50/cs
NON272790HXL	Knit Cuffs, X-Large, Blue	50/cs
NON272790HTL	Thumb Loops, Regular/Large, Blue	50/cs
NON272790HTLX	Thumb Loops, X-Large, Blue	50/cs
CRI4500	Thumb Loops, Regular/Large, Yellow	100/cs
Premium Polyethylene Coated		
NON26278	Knit Cuffs, Regular/Large, Blue	50/cs
NON26278XL	Knit Cuffs, X-Large, Blue	50/cs
NON26278XXL	Knit Cuffs, XX-Large, Blue	50/cs
NON26276	Thumb Loops, Regular/Large, Blue	50/cs
NON26276XL	Thumb Loops, X-Large, Blue	50/cs
NON26276XXL	Thumb Loops, XX-Large, Blue	50/cs

Meet
USP 800
requirements



NON29457



NON29859



CRI2010

Multi-Layer Shoe Covers

- Made of multi-layer material
- Heavyweight material features stripe non-skid pattern for increased traction
- Range of sizes to fit all shoe types

Item No.	Description	Pkg.
NON29858	Heavyweight, Regular/Large	100/bx, 300/cs
NON29859	Heavyweight, X-Large	100/bx, 200/cs
NON28858	Medium Weight, Regular/Large	100/bx, 300/cs
NON28859	Medium Weight, X-Large	100/bx, 200/cs
NON28852	Medium Weight, XX-Large	100/bx, 200/cs
CRI2006	Lightweight, Regular/Large	100/bx, 300/cs
CRI2007	Lightweight, X-Large	100/bx, 200/cs



CRI2002



CRI2000

Spunbond Polypropylene Shoe Covers

- Made of spunbond polypropylene
- Choice of non-skid pattern or smooth bottom

Item No.	Description	Pkg.
NON28758	Non-Skid, Regular/Large	100/bx, 300/cs
NON28759	Non-Skid, X-Large	100/bx, 200/cs
NON28752	Non-Skid, XX-Large	100/bx, 200/cs
CRI2002	Non-Skid, Regular/Large	100/bx, 300/cs
CRI2003	Non-Skid, X-Large	100/bx, 200/cs
CRI2000	Smooth Bottom, Regular/Large	100/bx, 300/cs
CRI2001	Smooth Bottom, X-Large	100/bx, 200/cs

Polyethylene Shoe Cover

- Smooth bottom for situations where traction is not essential
- Made from polyethylene material designed to withstand high levels of fluids

Item No.	Description	Pkg.
CRI2010	One Size Fits Most	1,000/cs

Polyethylene-Coated Shoe Covers

- Provide high-level fluid protection per ASTM F1670 for blood and bodily fluid penetration
- Medium-weight coated spunbond polypropylene
- Available with seam on bottom or seamless
- Non-skid treated for better traction

Item No.	Description	Color	Pkg.
Non-Skid/Bottom Seam			
NON30858	Regular/Large	Blue	300/cs
NON30859	X-Large	Blue	300/cs
Non-Skid/Seamless			
NON31858	Regular/Large	Blue	300/cs
NON31859	X-Large	Blue	300/cs
NON31852	XX-Large	Blue	300/cs



NON31852

Knee-High Boot Covers

- Premium breathable material tested against viral penetration per ASTM F1671 and ASTM F739-12 for top 11 chemotherapy drugs
- Blue laminate material tested against blood and bodily fluids per ASTM F1670 and viral penetration per ASTM F1671
- Poly-coated material tested against blood and bodily fluids per ASTM F1670
- Choose between non-skid bottom or foam strips for increased traction
- Available with hook and loop closure or elastic ankle

Item No.	Description	Pkg.
NON27348P	Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27348PXL	Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27349PXL	Thigh High Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27200	Blue Laminate, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27200XL	Blue Laminate, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27144	Poly-coated, Non-skid, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27144XL	Poly-coated, Non-skid, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27150	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27150XL	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27143	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, Regular/Large	50/bx, 150/cs
NON27143XL	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, X-Large	50/bx, 150/cs



NON27348P



NON27144



NON27200



NON27143




NON28229



CRI1008




NON30233B



NON28226



CRI1004



CRI1005



NON28233F



NON28230R



CRI1009



CRI1002

Multi-Layer Bouffant Caps

- Variety of prints and sizes; all come 100/bx, 500/cs
- Choice of heavy or medium weight multi-ply or heavyweight single-ply

Item No.	Description	Pkg.
NON30232B	21", Blue	500/cs
NON30233B	24", Blue	500/cs
NON30223B	28", Blue	500/cs
NON28226	24", US Flag	500/cs
NON28228	21", Brush Stroke	500/cs
NON28229	24", Brush Stroke	500/cs
NON28230R	24", Breast Cancer Ribbon	500/cs
NON28233F	24", Butterfly	500/cs
NON28533B	24", Blue, Heavyweight	500/cs

Polypropylene Bouffant Caps

- Variety of colors and sizes, all come 100/bx, 500/cs
- Choice of medium weight polypropylene or lightweight polypropylene

Item No.	Description	Pkg.
CRI1000	21", Green	500/cs
CRI1001	21", Blue	500/cs
CRI1002	21", White	500/cs
CRI1003	24", White	500/cs
CRI1004	24", Blue	500/cs
CRI1005	24", Green	500/cs
CRI1008	24", Pink	500/cs
CRI1009	24", Red	500/cs
CRI1202	24", Blue, Medium Weight	500/cs
CRI1203	24", White, Heavyweight	500/cs

**Surgical headwear
that meets
AORN Guidelines
for Surgical
Headwear**

Surgical Headwear

- Better coverage from the back of the neck to sideburn area
- One size fits most—provides all-day breathable comfort and a secure fit
- No elastic band in front means cap won't leave a mark
- Quality, longer-lasting design creates less waste than bouffant caps

Item No.	Description	Pkg.
NON61950	Premium Blue Material	600/cs
NON61990	Scrim Material	500/cs
NON61980	Multi-layer Material	500/cs



NON61980



NON61950



NON61990

Surgeon Caps

- Available in scrim, multi-layer or spunlace materials
- One size fits most
- Ties in the back
- Polypropylene crown for breathability

Item No.	Description	Pkg.
NON28625	Scrim Material for Absorbency	500/cs
NON28626	Multi-layer Material for Fluid-resistance	500/cs
NON28628	Spunlace Material for Comfort	500/cs



NON28625



NONSH100C



NONSH200



NONSH600



NONSH400

Surgeon Hood and Beard Covers

- Made of multi-layer or polypropylene material
- One size fits most

Item No.	Description	Pkg.
NONSH100C	Surgeons Hood with Multi-Layer Material, Polypropylene Crown	300/cs
NONSH200	Multi-layer Hood	300/cs
NONSH600	Polypropylene Hood	300/cs
NONSH400	Polypropylene Beard Cover	1,000/cs

Scrubwear

ComfortSure™ Scrubs

- V-neck top with 3 pockets
- Drawstring pant with pockets
- Made from premium ComfortSure™ material that provides soft, comfortable protection

Item No.	Description	Pkg.
ComfortSure™ Scrub Top, V-Neck, Green		
NON47202XS	Top, Extra Small	30/cs
NON47202S	Top, Small	30/cs
NON47202M	Top, Medium	30/cs
NON47202L	Top, Large	30/cs
NON47202XL	Top, XL	30/cs
NON47202XXL	Top, XXL	30/cs

ComfortSure™ Scrub Pant, Drawstring Waist, Green

NON47203XS	Bottom, Extra Small	30/cs
NON47203S	Bottom, Small	30/cs
NON47203M	Bottom, Medium	30/cs
NON47203L	Bottom, Large	30/cs
NON47203XL	Bottom, XL	30/cs
NON47203XXL	Bottom, XXL	30/cs



NON47203M

NON27212L

NON37202M

Traditional Scrubs

- Choice of v-neck or round neck top with chest pocket
- Choice of elastic waist or drawstring pant with two pockets
- Made from multi-ply material

Scrub Top, Round Neck, Blue

NON27212S	Small	30/cs
NON27212M	Medium	30/cs
NON27212L	Large	30/cs
NON27212XL	X-Large	30/cs
NON27212XXL	XX-Large	30/cs
NON27212XXXL	XXX-Large	30/cs

Scrub Top, V-Neck, Blue

NON27202S	Small	30/cs
NON27202M	Medium	30/cs
NON27202L	Large	30/cs
NON27202XL	X-Large	30/cs
NON27202XXL	XX-Large	30/cs

Scrub Pant, Elastic Waist, Blue

NON27213S	Small	30/cs
NON27213M	Medium	30/cs
NON27213L	Large	30/cs
NON27213XL	X-Large	30/cs
NON27213XXL	XX-Large	30/cs
NON27213XXXL	XXX-Large	30/cs

Scrub Pant, Drawstring Waist, Blue

NON27203S	Small	30/cs
NON27203M	Medium	30/cs
NON27203L	Large	30/cs
NON27203XL	X-Large	30/cs
NON27203XXL	XX-Large	30/cs

Scrub Top, Round Neck, Wine Red

NON37212S	Small	30/cs
NON37212M	Medium	30/cs
NON37212L	Large	30/cs
NON37212XL	X-Large	30/cs
NON37212XXL	XX-Large	30/cs
NON37212XXXL	XXX-Large	30/cs

Scrub Top, V-Neck, Wine Red

NON37202S	Small	30/cs
NON37202M	Medium	30/cs
NON37202L	Large	30/cs
NON37202XL	X-Large	30/cs
NON37202XXL	XX-Large	30/cs
NON37202XXXL	XXX-Large	30/cs

Scrub Pant, Elastic Waist, Wine Red

NON37213S	Small	30/cs
NON37213M	Medium	30/cs
NON37213L	Large	30/cs
NON37213XL	X-Large	30/cs
NON37213XXL	XX-Large	30/cs
NON37213XXXL	XXX-Large	30/cs

Scrub Pant, Drawstring Waist, Wine Red

NON37203S	Small	30/cs
NON37203M	Medium	30/cs
NON37203L	Large	30/cs
NON37203XL	X-Large	30/cs
NON37203XXL	XX-Large	30/cs

Coveralls

Prohibit™ Coveralls

Recommended for: Lab, Facilities Management, Maintenance and Staff

- Comfortable yet protective for high fluid or viral hazards
- Breathable tri-laminate material passes ASTM F1670 & F1671

Item No.	Description	Pkg.
NONCVP2M	Small/Medium	25/cs
NONCVP2L	Large	25/cs
NONCVP2XL	X-Large	25/cs
NONCVP2XXL	XX-Large	25/cs
NONCVP2XXXL	XXX-Large	25/cs
NONCVP2XXXXL	XXXX-Large	25/cs

Microporous Coveralls

Recommended for: Lab, Maintenance, and Visitors

- Choose between elastic wrists and elastic ankles or straight wrists and straight ankles

Item No.	Description	Pkg.
Elastic Wrists and Ankles		
NONCV700L	Large	25/cs
NONCV700XL	X-Large	25/cs
NONCV700XXL	XX-Large	25/cs
NONCV700XXXL	XXX-Large	25/cs
Straight Wrists and Ankles		
NONCV300XL	X-Large	25/cs
NONCV300XXL	XX-Large	25/cs

Multi-Layer Coveralls

Recommended for: Staff and Visitors

- Made from multi-layer polypropylene
- Elastic wrist and straight ankle

Item No.	Description	Pkg.
NONCV250XL	X-Large	25/cs
NONCV250XXL	XX-Large	25/cs



NONCVP2L

NONCV150L

Spunbond Coveralls

Recommended for: Staff and Visitors

- Available in medium weight or heavyweight polypropylene material
- Choose between elastic wrist and straight ankle or straight wrist and straight ankle

Item No.	Description	Pkg.
Medium Weight, Straight Wrists, Straight Ankles		
CRI3000	X-Large	25/cs
CRI3001	XX-Large	25/cs
Heavyweight, Elastic Wrists, Straight Ankles		
NONCV150L	Large	25/cs
NONCV150XL	X-Large	25/cs
NONCV150XXL	XX-Large	25/cs
NONCV150XXXL	XXX-Large	25/cs
NONCV150XXXXL	XXXX-Large	25/cs

ComfortSure™ Lab Jackets with Knit Collar and Cuffs

- Premium Comfortsure™ material provides soft, comfortable protection
- Two reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
NONRP602S	Small	30/cs
NONRP602M	Medium	30/cs
NONRP602L	Large	30/cs
NONRP602XL	X-Large	30/cs
NONRP602XXL	XX-Large	30/cs

Premium Lab Jackets with Knit Collar and Cuffs

- Premium multi-layer material provides comfort and protection
- Two reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
NONRP601S	Small	30/cs
NONRP601M	Medium	30/cs
NONRP601L	Large	30/cs
NONRP601XL	X-Large	30/cs
NONRP601XXL	XX-Large	30/cs



Classic Lab Coats with Knit Collar and Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
White		
NONSW500S	Small	30/cs
NONSW500M	Medium	30/cs
NONSW500L	Large	30/cs
NONSW500XL	X-Large	30/cs
NONSW500XXL	XX-Large	30/cs
NONSW500XXXL	XXX-Large	30/cs
Blue		
NONSW600S	Small	30/cs
NONSW600M	Medium	30/cs
NONSW600L	Large	30/cs
NONSW600XL	X-Large	30/cs
NONSW600XXL	XX-Large	30/cs

Classic Lab Jackets with Knit Collar and Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
White		
NONRP500S	Small	30/cs
NONRP500M	Medium	30/cs
NONRP500L	Large	30/cs
NONRP500XL	X-Large	30/cs
NONRP500XXL	XX-Large	30/cs
Blue		
NONRP600S	Small	30/cs
NONRP600M	Medium	30/cs
NONRP600L	Large	30/cs
NONRP600XL	X-Large	30/cs
NONRP600XXL	XX-Large	30/cs

Classic Lab Coats with Traditional Collar and Knit Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
White, Snap Front		
NONSW100S	Small	30/cs
NONSW100M	Medium	30/cs
NONSW100L	Large	30/cs
NONSW100XL	X-Large	30/cs
NONSW100XXL	XX-Large	30/cs
NONSW100XXXL	XXX-Large	30/cs
Blue, Snap Front		
NONSW400S	Small	30/cs
NONSW400M	Medium	30/cs
NONSW400L	Large	30/cs
NONSW400XL	X-Large	30/cs
NONSW400XXL	XX-Large	30/cs
Blue, Zipper Front		
NONSW480M	Medium	30/cs
NONSW480L	Large	30/cs
NONSW480XL	X-Large	30/cs
NONSW480XXL	XX-Large	30/cs
NONSW480XXXL	XXX-Large	30/cs

Protective Sleeve Covers

- Provides maximum amounts of protection in high fluid scenarios, or if extra protection is needed
- Available in premium breathable laminate material, which passes ASTM F1670 & F1671, or microporous which passes ASTM F1670
- Elastic bands at both ends secure the covers at the wrists and upper arms



NON18SLEEVE1

Item No.	Description	Pkg.
NON18P2SLEEVE	20" (50.8 cm) Long Blue Breathable Laminate Material	200/cs
NON18SLEEVE1	20" (50.8 cm) Long White Microporous Breathable Material	200/cs

Breast Cancer Awareness Month Pink PPE

Protection with a cause.



Staff Apparel

Knit Cuff Cover Gown

- Pink gown made from spunbond material
- Neck tie style with side waist ties
- Knit cuffs

Item No.	Description	Pkg.
NONDPL	Regular/Large	50/cs
NONDPXL	X-Large	50/cs

Lab Jacket

- Pink lab jacket made from multi-layer material
- Snap front closure
- Knit collar and cuffs

Item No.	Description	Pkg.
NONRP700S	Small	30/cs
NONRP700M	Medium	30/cs
NONRP700L	Large	30/cs
NONRP700XL	X-Large	30/cs
NONRP700XXL	XX-Large	30/cs

Classic Cover Gown

- Pink isolation gown made from multi-layer material
- Neck and waist ties
- Elastic wrists

Item No.	Description	Pkg.
NON27SMS10	Regular/Large	100/cs

Lab Coat

- Pink lab coat made from multi-layer material
- Hits right below the knee
- Snap front closure
- Knit collar and cuffs

Item No.	Description	Pkg.
NONSW700S	Small	30/cs
NONSW700M	Medium	30/cs
NONSW700L	Large	30/cs
NONSW700XL	X-Large	30/cs
NONSW700XXL	XX-Large	30/cs



NONDPL



NONRP700



NON27SM10



NONSW700

Patient Apparel

Patient Cape

- Patient cape made from multi-layer material
- Sleeveless with tape tab neck
- 31" x 21"

Item No.	Description	Pkg.
NON27038	One Size Fits Most	100/cs

Patient Gown

- Patient gown made from multi-layer material
- Sleeveless with neck and waist ties

Item No.	Description	Pkg.
NON27048	Regular/Large	50/cs

NON27038



NON27048



Bouffant Caps

Pink Ribbon Bouffant Cap

- Bouffant cap made with lightweight polypropylene material

Item No.	Description	Pkg.
NON28230R	Size 24"	500/cs



NON28230R



CRI1008

Pink Bouffant Cap

- Bouffant cap made with lightweight polypropylene material

Item No.	Description	Pkg.
CRI1008	Size 24"	500/cs

Procedure Face Masks

Face Mask with Ear Loops

- Face mask with ear loops

Item No.	Description	Pkg.
NON27712EL	Pink/White	300/cs



NON27712EL

Surgical Face Masks

Face Mask with Ties

- Face mask with ties

Item No.	Description	Pkg.
NONRIBBON	Pink Ribbon	300/cs
NON27712	Pink/White	300/cs



NONRIBBON



NON27712

Face Masks

Basic Earloop Face Masks



NON27122 Isolation Face Mask
50/bx, 6 bx/cs

- Bacterial filtration efficiency: $\geq 98\%$
- Particulate filtration efficiency (0.1 micron): $\geq 97\%$
- Spunbond polypropylene outer, thermal bond polypropylene inner



NON27378 Procedure Face Mask
50/bx, 6 bx/cs

- Bacterial filtration efficiency: $\geq 98\%$
- Particulate filtration efficiency (0.1 micron): $\geq 97\%$
- Spunbond polypropylene outer, thermal bond polypropylene inner



NON27375 Procedure Face Mask
50/bx, 6 bx/cs

- Bacterial filtration efficiency: $\geq 98\%$
- Particulate filtration efficiency (0.1 micron): $\geq 97\%$
- Cellulose inner, spunbond polypropylene outer facings

Basic Surgical Face Masks



NON27402 Surgical Face Mask
50/bx, 6 bx/cs

- Bacterial filtration efficiency: $\geq 98\%$
- Particulate filtration efficiency (0.1 micron): $\geq 97\%$
- Cellulose inner, spunbond polypropylene outer facings



NON27385 Hypoallergenic Surgical Face Mask
50/bx, 6 bx/cs

- Bacterial filtration efficiency: $\geq 98\%$
- Particulate filtration efficiency (0.1 micron): $\geq 97\%$
- Cellulose inner and outer facings



NON27386 Sensitive Skin Surgical Face Mask
50/bx, 6 bx/cs

- Bacterial filtration efficiency: $\geq 99\%$
- Particulate filtration efficiency (0.1 micron): $\geq 97\%$
- Spunbond polypropylene outer, thermalbond polypropylene inner

Patient Wear

Thoughtfully designed in a variety of sizes to fit every shape.

Multi-Layer Patient Gowns

- Disposable material has a fabric-like feel for comfort
- Opening in front or back for easy access
- Attached ties for a custom fit

Item No.	Description	Pkg.
NON27046	Gown, Sleeveless, Waist Tie, Regular/Large, Blue	50/cs
NON27048	Gown, Sleeveless, Waist Tie, Regular/Large, Pink	50/cs
NON27146	Gown, Sleeveless, Neck and Waist Tie, Regular/Large, Blue	50/cs
NON27146SL	Gown, Short Sleeves, Neck and Waist Tie, Regular/Large, Blue	50/cs
NON27146SLXL	Gown, Short Sleeves, Neck and Waist Tie, XL, Blue	50/cs
NON27346SL	3 Arm Hole Gown, Sleeveless, Regular/Large, Blue	50/cs



NON27048

NON27146SL

Reinforced Tissue Patient Gowns

- 3-ply tissue-poly-tissue (TPT) material aids in fluid absorption and protection
- Scrim material provides added strength and absorption
- Opening in front or back for easy access

Adjustable ties for a custom fit

Item No.	Description	Material	Pkg.
NON24356	Gown, Waist Tie, 30" x 42", Blue	TPT	50/cs
NON24355	Gown, Waist Tie, 30" x 42", White	TPT	50/cs
NON24252	Gown, Waist Tie, 30" x 42", Mauve	TPT	50/cs
NON24354	X-Ray Gown, Waist Tie, 30" x 42", Blue	TP	50/cs
NON254572	Bariatric Gown, Waist Tie, 72" x 45", Blue	Scrim	25/cs



NON24356



NON24355



NON24252

Medline PPE Companion™

Protection that's always by your side.

Streamline caregiving. Move PPE closer to the point of care. This unique isolation cart has all you need to protect yourself—and your patients.

Isolation Cart

A Precaution Identification

Multilingual instructions show the clinically appropriate level of personal protection.

B Organized Storage

Multilingual labels on each drawer.

C Hand Protection

Three glove sizes. More storage behind.

D Face Protection

At-a-glance supply of face masks with shields and earloop masks, ample storage for N95 respirators and more.

E Body Protection

Gowns dispense one at a time with a touch and a pull.

F Gown Storage

Additional storage for a minimum of 50 gowns.

G Quick Disinfection

Sanitize hands. Disinfect surfaces.



NONISOORGCART

Disposable Isolation Station

- Keeps hallways free from isolation cart congestion
- Organizer stores flat-packed isolation gowns, three gloves boxes, and one box of face masks
- 3 rows of pockets with eyelets for hanging
- Made of SMS material



NONISOCADDY

Know you're protected. Contact your Medline representative for more information.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093
Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996 | medline.ca | canada@medline.com
Medline México
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    

REFERENCE:

1. Considerations for selecting protective clothing used in healthcare for protection against microorganisms in blood and body fluids. Centers for Disease Control and Prevention website. Table 3: ANSI/AAMI PB 70:12 Classification of barrier performance of surgical gowns, other protective apparel, surgical drapes and drape accessories. Available at: <http://www.cdc.gov/niosh/npptl/topics/protectiveclothing/>. Accessed May 3, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

© 2019 Medline Industries, Inc. Companion and QuickComply are trademarks and Medline is a registered trademark of Medline Industries, Inc. MKT19W719051 / LIT171 / 1.5M / TC / 57