



Invacare® IVC™ Tracer® SX5 Wheelchair

The Invacare IVC Tracer SX5 lightweight frame weighs less than 36 pounds, which makes this chair perfect for rental and long- or short-term use. The Tracer SX5 now combines the design and technology of the 9000 series wheelchairs, making it possible to interchange components. The Tracer SX5 offers multiple seat widths, depths, armstyles, sectional back and seat-to-floor heights by utilizing the dual axle position frame base. The 20" and 22" widths include a semi-pneumatic caster and chrome handrim to accommodate a weight capacity of 300 lb.

FEATURES

- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Dual axle positions allow variation of seat-to-floor heights
- 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Sectional back available in 16" and 18" back heights
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Meets RESNA standards for manual wheelchairs in product performance

Specifications

Dimensions	
Seat sizes	
Widths	14", 16", 18", 20", 22"
Depths (optional)	16" standard, 18" by seat extension
Seat-to-floor heights	17.5", 19.5" 15.5", 16.5", 18.5" with different wheel and caster packages
Back height	16", 18" sectional
Front riggings	
Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
Articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles	
Full-length fixed height, padded	
Desk-length fixed height, padded	
Full-length adjustable height, padded	
Desk-length adjustable height, padded	
(All arms are flip back non-removable)	
Weight capacity	250 lb. (300 lb. 20" and 22" width)
Product weight (TRSX5/WD86/28FB)	34 lb.
Warranty	lifetime on frame and crossbraces