


Graphite



Maple

Dione PR446



Comfort Zone
by GOLDEN

Shown in Maple

Inspired by the moon, the DeLuna™ Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep. The Dione is offered in both medium and large sizes and features a stylish, overstuffed seamed back and plush chaise pad. It comes standard with adjustable lumbar and headrest, programmable AutoDrive™ hand control, and USB charging port for your phone or tablet. Also standard is our patented four-zone Comfort Zone positioning technology, which provides a stress-free and supportive reclining experience.

Standard Features:

- Adjustable headrest
- Adjustable lumbar
- MaxiComfort® Technology
- Programmable AutoDrive™ hand control
- USB charging port



PR446 Dione



Model Number	PR446M	PR446L
Size	Medium	Large
Weight Capacity	375 LBS	375 LBS
Positioning	Infinite	Infinite
Comfort Zone(s)	4	4
Twilight Technology	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	YES	YES
Back Type	Seam	Seam
Removable Back Type	Bracket	Bracket
Head Rest Cover	NO	NO
Arm Rest Covers	NO	NO
Overall Width	34"	35"
Overall Height	43.5"	46"
Width Between Arms	22"	22.5"
Floor to Top of Seat	20.5"	22"
Seat Depth	21"	21.5"
Seat to Top of Back	28"	30"
Seat to Top of Arm	7"	7.5"
Distance Required from Wall Reclined	24"	26"
Extension Reclined from Back of Seat to End of Foot Rest	42.5"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.
Hand Control Routing	Inside Arm	Inside Arm
Extra Pocket Available	Standard	Standard
Heat & Massage Available	NO	NO
Articulating Headrest Available	Standard	Standard
Adjustable Lumbar Available	Standard	Standard
AutoFlex Available	NO	NO
Foot Rest Extension Available	NO	NO
Straight Lift Option	YES	YES
USB Charging Port	Outside Arm	Outside Arm
Place of Origin	China	China

DeLuna™
SERIES



Shown in Graphite



Convenient built-in USB charging port for your phone or tablet located on the outside arm.

Battery Backup & SmartTek Diagnostic System STANDARD on all Golden models.


Designed and engineered in Old Forge, PA. Assembled in China using Golden components. Crossover fabrics not available on this model.



Comfort Zone
by GOLDEN

FOUR ZONE POSITIONING SYSTEM

Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



Take a picture with your smartphone to learn more.



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

SHT-PR446-121520