

Portable Bath Step



FIGURE A



FIGURE B



RETAIL BOX

- Makes getting in and out of the bath tub easier and safer
- Unique modular design allows the multiple steps to be stacked (Figure A)
- Step height is 4"
- When stacking, each additional step adds 2"
- Bath step comes with four removable covers on the surface
- When stacking, remove covers to allow steps to lock together (Figure B)
- Attractive, textured surface is slip resistant and comfortable
- Retail Packaged

PRODUCTS

Description

RTL12068	Contains 1 Bath Step per box, 1/cs
----------	------------------------------------

SPECIFICATIONS

Description

Dimensions	19.5" (W) x 4" (H) x 16" (D)
Weight Limit	450 lbs.
Warranty	Limited One Year