

# Jazzy<sup>®</sup>

## ELITE HD

The Jazzy<sup>®</sup> Elite HD features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.



### LARGE 14" DRIVE WHEELS



### PERFORMANCE-MATCHED ELECTRONICS



### COMPACT CASTERS FOR SUPERIOR INDOOR MANEUVERABILITY

## AMERICA'S FAVORITE POWER CHAIR

Learn more online.



[gojazzy.com](http://gojazzy.com)

# Specifications

<b>Weight Capacity</b>	450 lbs.
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	3"
<b>Turning Radius</b>	24"
<b>Overall Length</b>	35.5"
<b>Overall Width</b>	24.25"
<b>Drive Wheels</b>	14" solid
<b>Front Wheels</b>	3" solid anti-tips
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, front-wheel drive
<b>Braking System</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard Electronics</b>	60A, Dynamic Shark controller
<b>Battery Charger</b>	Off-board, 5 amp
<b>Per-Charge Range</b>	Up to 9.8 miles
<b>Battery Requirements</b>	(2) 12 volt, NF-22
<b>Battery Weight</b>	37 lbs. each
<b>Base Weight</b>	99 lbs.
<b>Standard Seat Weight</b>	52 lbs. (high-back)
<b>Standard Seat</b>	18" x 18"-20", 20" x 18"-20", 22" x 20"-22"
<b>Max. Seat Dimensions</b>	24" x 24"
<b>Warranty</b>	Lifetime limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

## Colors



Jazzy Red



Jazzy Blue

## Models

**Jazzy Elite HD-2HD-SS** | Code: K0824

**Jazzy Elite HD-2HD-C** | Code: K0825

\*Additional charges may apply.

# Features and Benefits

- Depth adjustable, limited recline comfort high back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 70" lap belt

## Seating Options\*

- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 18"x18" or 18"x20"
- Synergy® Seat

## Options

- Bell shaped footrest
- Elevating leg rests
- Heel loops
- Height adjustable foot platform bracket
- Residual limb support
- Swing- away leg rests
- Swing-away joystick

## Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

## Recommended Lifts

- Backpacker® Series
- Outlander Full Platform
- Commander 450



Find a compatible lift for your vehicle online:  
[PRIDEMOBILITY.COM/LIFT-FINDER](http://PRIDEMOBILITY.COM/LIFT-FINDER)



(US) 800-800-8586 • (Canada) 888-570-1113 • [pridemobility.com](http://pridemobility.com) •



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.