

NEW!



BARIATRIC
700 LB
WEIGHT CAPACITY

BARIATRIC **Steel Folding Walker, Two Button** 10220HD Heavy Duty, Adult, 2/cs

FEATURES:

- For bariatric use (700 lb. weight capacity).
- Sturdy, steel construction.
- Vinyl contoured hand grips.
- Easy push-button mechanisms may be operated by fingers, palm or side of hand.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10220HD	19.5"	22.5"	25.5"	22"	39.5"	32.5"	7.8 lbs.	18.2 lbs.	700 lbs.

bariatric steel folding walker

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM



**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**

Side Walker

10240-1 Adult, (Hemi-Walker), 1/cs

10240-2* Adult, (Hemi-Walker), 2/cs

FEATURES:

- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Height adjustment 29.5" - 37".
- Folds easily with one hand.
- Ideal for users who do not require a walker but need a wider base for support.
- Limited Lifetime Warranty.



Bariatric Side Walker

10242 Adult, (Hemi-Walker), 1/cs

FEATURES:

- **Manufactured with reinforced tubing.**
- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Height adjustment 30" - 36".
- Folds easily with one hand.
- Ideal for users who do not require a walker but need a wider base for support.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Overall Width	Overall Depth	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10240-1/2	17.5"	15.5"	37"	29.5"	3.5 lbs.	5.5/9 lbs.	300 lbs.
10242	16.5"	16"	36"	30"	10 lbs.	13 lbs.	500 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Deluxe Folding Walker, Two Button with 5" Wheels & Tennis Ball Glides

10210TB-2 Adult, 2/cs

FEATURES:

- Comes with 5" wheels and Tennis Ball glides. (item #10119 with additional glides, 10107 and 10107C).
- Comes with wheels. For replacement wheels, refer to item # 10109.
- Easy push-button mechanisms may be operated by fingers, palm or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Unique U-shape frame design creates greater clearance.
- The patient is able to bring the walker closer to assist in standing.
- Comes with rear tennis ball glides.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Vinyl contoured handgrip for added comfort.
- Limited Lifetime Warranty.



New & Improved Bariatric (Wider & Deeper Frame Design) Aluminum Folding Walker, Two Button

10220-1 Adult, 1/cs
10220-2 Adult, 2/cs
10223-2 Junior, 2/cs
10220-2WW Adult, with Heavy Duty Wheels, 2/cs

BARIATRIC
500 LB
WEIGHT CAPACITY



THIS ITEM
IS AVAILABLE
IN A 1-PACK!

FEATURES:

- For Bariatric use.
- Steel legs and side braces.
- Wider and deeper frame design accommodates individuals up to 500 lbs.
- Vinyl contoured hand grip.
- Item # 10220-2WW comes with wheels. For replacement wheels, refer to item # 10118SV.
- Precision design provides additional strength while adding minimal weight.
- Limited Lifetime Warranty.



TWO BUTTON BARIATRIC WALKER WITH WHEELS FEATURE MATCHING SILVER VEIN "U" BRACE AND SIDE FRAMES & LEGS!

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10210TB-2	17"	19"	25.5"	18"	39"	32"	6.9 lbs.	15.8 lbs.	350 lbs.
10220-1	20"	23"	24"	21"	39"	32"	7.5 lbs.	8.8 lbs.	500 lbs.
10220-2	20"	23"	24"	21"	39"	32"	7.5 lbs.	17 lbs.	500 lbs.
10223-2	20"	23"	24"	20"	32"	25"	7.3 lbs.	16.6 lbs.	500 lbs.
10220-2WW	20"	23"	27"	24"	39"	32"	11 lbs.	24.2 lbs.	500 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

walkers

drive

To order call toll free: 877.224.0946

WWW.DRIVEMEDICAL.COM

BARIATRIC
500 LB
WEIGHT CAPACITY



**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**

**EASY HEIGHT
ADJUSTMENT
WITH SELF
THREADING
KNOB!**




10215RD-1
SEE PAGE 52 & 53 FOR
REPLACEMENT SEATS.

Go-Lite Bariatric Steel Rollator, Padded Seat, 8" Casters with Loop Locks

10215BL-1 Blue, (Has replaced 685 series and 466), 1/cs
10215RD-1 Red, (Has replaced 685 series), 1/cs

FEATURES:

- Large 8" casters are ideal for indoor and outdoor use.
- Comes with large 21" (W) x 9" (D) x 6" (H) basket that can be mounted under seat (item # 10215BSK).
- Tool Free removable padded back rest for comfort.
- Strong Steel reinforced frame supports 500 lbs.
- Soft padded oversized seat (18" (W) x 13" (D)).
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Handles are adjustable in height.
- Limited Lifetime Warranty.
Lifetime Warranty on brake cable.



**PARTS
DIAGRAMS
ON PAGE 67**

Adjustable Width Winnie Walker, Bariatric EXP Steel Rollator, 8" Casters with Loop Locks

466* Light Blue & Chrome, 1/cs

FEATURES:

- Large 8" casters are ideal for indoor and outdoor use.
- 4" adjustable width allows entry into narrow doorways.
- Strong Steel reinforced frame supports 400 lbs.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Soft padded oversized seat.
- Handles are adjustable in height.
- Comes with oversized 14" x 12" x 7" wire basket (item # 828).
- Limited Lifetime Warranty.
Lifetime Warranty on brake cable.

BARIATRIC
400 LB
WEIGHT CAPACITY

**EASY TO
REACH KNOB
ADJUSTS
ROLLATOR
WIDTH 4" TO
ALLOW ENTRY
INTO NARROW
DOORWAYS.**



466*
SEE PAGE 52 & 53 FOR
REPLACEMENT SEATS.

ADJUST FRAME WIDTH UP TO 4" WHEN USING LEVER



**SEE ROLLATOR ACCESSORIES ON PAGES 50-57.
ROLLATOR PARTS DIAGRAMS ON PAGE 58-67.
ROLLATOR COLOR PAINT SWATCHES ON PAGE 68.**

SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
GO-LITE BARIATRIC	39"	34"	30"	29"	13"	18"	22"	26 lbs.	30.5 lbs.	500 lbs.
EXP/466	38"	34"	25" - 29"	30"	10"	20"	25"	28 lbs.	32.5 lbs.	400 lbs.

rollators

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM