



## Free-standing Toilet Safety Rail

## Portable Bath Step

Item #	UOM
RTL12079	2 retail boxes/cs

- Free standing around toilet.
- Easy, tool-free assembly using wing nuts.
- Fits standard and elongated toilets.
- Comes with padded armrests for added comfort.
- Comes with a magazine rack for added convenience.

specifications	
Width Between Arms	21"
Overall Width	26.5"
Weight	9.0 lbs.
Carton Shipping Weight	28 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Item #	Product Description	UOM
RTL12068	Contains 1 Bath Step per box	1/cs

- Makes getting in and out of the bath tub easier and safer.
- Unique modular design allows the multiple steps to be stacked. (Figure A)
- Step height is 4". When stacking, each additional step adds 2".
- Bath step comes with four removable covers on the surface. When stacking, remove covers to allow steps to lock together. (Figure B)
- Attractive, textured surface is slip resistant and comfortable.
- Retail Packaged.

specifications	
Dimensions	19.5" (W) x 4" (H) x 16" (D)
Weight Limit	450 lbs.
Warranty	Limited One Year