

Innerspring Mattress

HCPCS: E0271



- High quality innerspring design
- The unit has a 1" high-density, firm, combustion-modified fiber topper on each side
- Multi-layered with two 1" high-density, firm, combustion modified fiber toppers on the top and one 1" high density, firm, combustion-modified fiber topper on the bottom (15006EF only)
- Vinyl cover is waterproof and anti-static for easy use and care

PRODUCTS	Description
15006	Mattress 80" Innerspring, 1/EA
15006EF	Mattress,Inner Spring,ExtraFir,1/ea
15014	Mattress,Inner Spring36"x84"XF,1/ea

SPECIFICATIONS	Description
Warranty	1 Year
Weight Capacity	350 lbs