



DRIVE

Universal Walker Wheels with Rear Glide Caps

ITEM #	SIZE	TUBING	UNITS	MFG #
1081A	5"	Grey	Pair	10109
1081B	3"	Grey	Pair	10111
1081D	5"	Blue	Pair	10109-BL

- Universal walker wheels with "2-Sided" positioning holes fit most leading manufacturers walkers, including Invacare®, Guardian®, Carex® & Lumex®
- Converts folding walker into wheeled walker
- Allows for 8" height adjustments
- Rubber wheel allows walker to glide easily and smoothly over irregular surfaces
- Rear glide caps (included) provide easy mobility over most surfaces
- Weight Capacity: 350 lbs



DRIVE

5" Walker Wheels and Extension Legs Combo Pack

ITEM #	DESCRIPTION	MFG #
1086D	Combo Pack	10108WC
1086A	Extension Legs Only	10108

- Allows for height adjustments from 36" to 43" (walker handgrip to floor), (Adds 4" to height of walker).
- 5" Rubber wheel allows walker to roll easily and smoothly over irregular surfaces.
- For use with all Drive and most leading manufactures walkers
- Includes rear glide caps, (item #10107) and glide covers (item #10107C) allowing use on all surfaces
- Weight Capacity: 350 lbs



DRIVE

5" Swivel Wheels with Rear Glide Caps

ITEM #	MFG #
1102A	10103

- Includes rear glide caps (1085A) and glide covers (1085B)
- Weight Capacity: 350 lbs




DRIVE

Rear Glide Brakes with 3" Wheels

ITEM #	UNITS	MFG #
1084A	Pair	10106

- For use on back of wheeled walker
- 3" wheel glides along surface until weight is applied to hand grip, which activates rubber tip brake
- Allows for 8" height adjustments
- Fits 1" diameter walker tubing
- Weight Capacity: 350 lbs



DRIVE

Tall Extension Legs

ITEM #	UNITS	MFG #
1086A	Set/4	10108

- Allows for height adjustments from 36" to 43" (walker handgrip to floor)
- 5" Rubber wheel allows walker to roll easily and smoothly over irregular surfaces
- For use with all Drive and most leading manufacturers walkers
- Includes rear glide caps and glide covers allowing use on all surfaces
- Weight Capacity: 350 lbs



DRIVE

Rear Glide Brake

ITEM #	UNITS	MFG #
1083A	Pair	10104

- For use on back of wheeled walker
- Nylon tip glides along surface until weight is applied to hand grip, which activates rubber tip brake
- Allows for 8" height adjustments
- Fits 1" diameter walker tubing
- Weight Capacity: 350 lbs



DRIVE

Glide Cap

ITEM #	UNITS	MFG #
1085A	Pair	10107

- Allows walker to glide easily and smoothly over irregular surfaces
- Replaces rubber tip



DRIVE Glide Cover

ITEM #	UNITS	MFG #
1085B	Pair	10107C

- Tip style glide cover for use over existing rubber tip
- Allows walker to "glide" across most surfaces
- Easily added to walker
- Increases life of walker tips



DRIVE Walker Shoe Glides

ITEM #	MFG #
1085K	RTL100014

- Stylish tennis shoe glides personalize the walker
- Allows walker to glide easily and smoothly over most surfaces
- Easy, tool free installation
- Replaces rubber tip



CAREX Walker Glides with Raised Points

ITEM #	SIZE	UNITS
RA835	1"	Pair
RA836	1-1/8"	Pair

- These glides have three points of contact for better stability and longer wear
- For use on the back of a wheeled walker



COMPLETE MEDICAL SUPPLIES, INC

Walker Glide Skis for 1" Walker Tubing

ITEM #	COLOR	UNITS
1085C	White	1 Pair
1085D	Black	1 Pair

- For use with all 1" folding walkers
- Easily added to walker
- Durable composite construction ensures long life



CAREX Walker Glide Skis for 1-1/8" Walker Tubing

ITEM #	UNITS
A80900	Pair

- Glides smoothly over most surfaces, including grass and cement
- The glides insert into extension legs on the back of wheeled walkers
- Fits most walkers
- Fits 1-1/8" walker tubing



DRIVE Tennis Ball Glides with 1 Pair Replaceable Glide Pads

ITEM #	DESCRIPTION	MFG #
1088A	Tennis Ball Glides	10119
1088B	Repl. Glide Pads	10123

- Provides a quiet, smooth and durable glide experience when used with a walker
- Easy and safe to install
- Comes with an additional pair of glide pads
- Last longer than plastic glide cap
- Warranty 1 year



Replacement Walker Tips

ITEM #	WALKER DIAM.	COLOR	UNITS
6102G	1 1/8"	Gray	1 Pair
1552	1"	Gray	1 Pair

- Measurement must be taken at base of walker
- Provides stability to walker
- Retail boxed



Use with Walkers, Rollators
and Wheelchairs!



Universal Cup Holder

ITEM #
10679

- New clamp allows cup holder to swivel to keep contents level.
- Easily attaches to wheelchairs, rollators or walkers with tubing diameter from 5/8" to 1"
- Special clamp insert can be added or removed as necessary to accommodate different size tubing
- Made with durable, lightweight plastic
- Dimensions: 3" (L) x 3" (W)
- Limited 1 Year Warranty



Blue Print **Burgundy Print** **Solid Gray** **Solid Burgundy**

MADDAK

Double-Sided Quilted Walker Bag

ITEM #	COLOR
1947A	Print, Blue
1947B	Print, Burgundy
1947D	Solid, Gray
1947C	Solid, Burgundy

- Designer inspired, decorative cotton walker bag has roomy pockets in the front and back so you can take your "stuff" wherever you go
- Bag conveniently hangs over the walker crossbar and is held in place by hook-and-loop closures
- Pockets are 13" x 9" (33 x 23 cm)
- Machine wash and tumble dry



DRIVE

Carry Pouch for Walker

ITEM #	SIZE
10258-1	Small, 8.5" x 9" x 1"
10258	Large, 10.5" x 14" x 2"

- Easily attaches to walker
- For use with most folding walkers
- Large zipped pouch
- 3 Mesh pouches
- Made of durable, easy-to-clean nylon
- 90 Day Warranty



DRIVE

Side Walker Carry Pouch

ITEM #
10255

- Easily attaches to side frame of walker
- For use with most folding walkers
- 2 large pockets with velcro type closure
- Made of durable, easy-to-clean nylon
- Detachable mesh pocket
- Size: 14.5" x 15" x 2.5"
- 90 Day Warranty

Bariatric
500 lbs



DRIVE

Platform Walker/Crutch Attachment

ITEM #	WT. CAPACITY	MFG #
1087A	300 lbs	10105-2
1087C	500 lbs	10105HD-2

ACCESSORIES

ITEM #	DESCRIPTION
1087D	Replacement Clamps for 1087A
1087B	Bracket for Platform Crutch Attachment


- Designed for those who are unable to grip the walker or crutch
- Cradles forearm with soft vinyl padding universally adjustable
- Adjustable strap securely holds arm in place
- Fits most manufacturers' walkers and aluminum crutches
- In order for 1087A to be used with crutches, the bracket (1087B) needs to be purchased as well
- **HCPSCS Suggested Code: E0154**



TO₂TE Oxygen Tank Carrier™ For Wheeled Walkers

ITEM #	DESCRIPTION
8037A	Carrier for - E Cylinder
8037B	Carrier for - D Cylinder
8037C	Carrier for - M6 Cylinder
8037D	Carrier for - M9 Cylinder

- Oxygen Cylinder Holder increased mobility and freedom for individuals who are required to use oxygen and a wheeled walker
- Secures to the walker with 2 hook and loop fasteners and a safety buckle strap at the top, and a centering strap near the bottom of the walker
- Convenient pocket for the wrench
- Manufactured and sold under U.S. Patent #6, 672, 321



CAREX

Universal Attachable Flashlight

ITEM #
A850

- Easily attaches to most walkers, canes, wheelchairs and rollators
- Adjustable swivel attachment allows you to direct light where you need it most
- Lights up to 17 feet away
- Batteries included (3 AAA)
- Extra bright light from 14 LED bulbs
- Easy to install on most mobility products or can be used as a portable flashlight



COMPLETE MEDICAL SUPPLIES, INC.
Walker Rest Seat

ITEM #
1090

- The Walker Rest Seat provides a comfortable, fold-down padded seat when you become tired and need to rest
- Two inch foam padded vinyl cover can be easily cleaned with damp cloth
- Plywood base with wood frame supports
- Two hook and loop fastener straps secure the seat to the uprights or horizontal rails
- Manufactured and sold under U.S. Patent #D457,111



COMPLETE MEDICAL SUPPLIES, INC.
Walker Fleece Covers

ITEM #
1547A

- Protects against soreness and irritation
- Machine washable/polyester
- Easy on and off
- Fastens securely with hook and loop
- Two-ply cushioning for added comfort



DRIVE
Walker Carry Basket

ITEM #	DESCRIPTION
A830	Walker Basket with Insert
A8301	Additional Insert Only

- Walker basket with insert for use with all 1" folding walkers
- Comes with plastic insert tray and cup holder
- Dimensions: 16" L x 5" W x 7" D



No-Wire™ Walker Basket

ITEM #
10614

- Firmly attaches to the horizontal and vertical bars of the walker
- Flexible, retractable stabilizing arms lock onto the legs of the walker and eliminate excess swinging
- The front of the basket has a slot for an identification card and the inside of the basket has a built-in glass, can or cup holder, a small compartment for coins or medications and a trough for utensils, pens, straws or eyeglasses
- No additional tray is needed
- Assembled basket measures 16.25" x 8.75" x 8.25" (41.3 cm x 22.2 cm x 21 cm)
- Made of polypropylene



CAREX
Collapsible Canvas Walker Basket

ITEM #
RA829

- Attach this basket to your walker for carrying personal items
- Elastic straps secure items
- Collapsible hinges allow flat storage



COMPLETE MEDICAL SUPPLIES, INC.
Universal Walker Tray

ITEM #
1094

- Allows you to carry personal items from room to room easily
- Made of easy to clean, durable plastic
- Two cup holders
- Fits most manufacturers walkers
- Easy to install
- Size: 21"W x 16"D
- Limited 1 Year Warranty



DRIVE
Walker Folding Flip Tray with Cup Holders

ITEM #	SIZE	MFG #
1092	16" W x 12" D	10125

- Fits most manufacturers 2 button walkers
- Easily folds down when not in use
- Made of easy to clean, durable plastic
- Two cup holders
- Easy to install
- Limited 1 Year Warranty

Equipment Bags



Plastic Bags for Walkers

ITEM #	BAG SIZE	UNITS
EBW3871	30" x 12" x 45"	100/roll

- Protects your walkers

Able TRIPOD Cane Tip

ITEM # CANE DIAM.

10349	3/4"
10349A	5/8"
10349B	7/8"
10349C	1"



- Replaces existing cane tip
- Provides greater support
- For right or left handed base
- Fits most manufacturers canes
- To determine which size tip you need, remove your cane or forearm crutch tip and accurately measure the diameter of the cane shaft
- 10349 fits #1624, 1626 & 1628 amongst other canes



1753

1754



1742-1752A

Cane Tips

RETAIL BOX

ITEM #	SHAFT DIAM.	COLOR	UNITS
1753	3/4"	Black	Pk/2
1754	5/8"	Black	Pk/4
1755A	1"	Black	Pk/2

UNBOXED - BULK PLASTIC BAG

ITEM #	SIZE	SHAFT DIAM.	UNITS
1742	16	1/2"	Bg/12
1742A	16	1/2"	Bg/4
1744	17	5/8"	Bg/12
1744A	17	5/8"	Bg/4
1746	18	3/4"	Bg/12
1746A	18	3/4"	Bg/4
1748	19	7/8"	Bg/12
1748A	19	7/8"	Bg/4
1750	20	1"	Bg/12
1750A	20	1"	Bg/4
1752	21	1 1/8"	Bg/12
1752A	21	1 1/8"	Bg/4

- Black
- Ridged bottom to prevent slipping
- Hard rubber for long use
- Assorted sizes for any cane or walker




Ice Cane Attachment

ITEM #

1756

- For safer walking in ice and snow
- Fits on all canes
- Strong steel five-prong grip



Cane Clip

ITEM #

A702

- This double clip device can be used on wheelchairs
- One side can be locked onto a wheelchair via our rubber ring system, the other clip is plastisol coated to grip canes, reachers, etc



COMPLETE MEDICAL SUPPLIES, INC.

Cane Holder

ITEM #

A703

- Allows cane to lean on the table
- Slip-resistant foam disc helps hold cane

Crutches



DRIVE

Push Button Aluminum Adjustable Crutches

ITEM #	PATIENT HEIGHT	MFG #
1455D	4'0" - 4'6"	10416-1
1455A	4'6" - 5'2"	10401-8
1455B	5'2" - 5'10"	RTL10400
1455C	5'10" - 6'6"	RTL10402

- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Push-pin adjustment makes it easy to adjust crutch height in 1" increments
- Non-skid, jumbo size rubber tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable rubber underarm pad and hand grip
- Limited Lifetime Warranty
- Sold as a pair
- Weight capacity: 300 lbs
- HCPCS Suggested Code: E0114

Crutch Rack



SEE PAGE 81