



bellavita

Auto Bath Lifter

477200252 1/cs

The Bellavita is our new ultra modern automatic bath lift. The Bellavita is easy to transport and stow away. It is ideally suited for use in the home or away. The new stylish design was created with the assistance of healthcare professionals and includes the following features...



FEATURES:

- The lightest bathlift on the market at only 20 1/2 lbs.
- At 2.3" the Bellavita is the lowest bath lift available.
- Seat height reaches a record 18.8". Excellent for deeper model tubs
- The easy to use quick charge lithium-ion battery is accomodated in the Hand Control.
- The backrest reclines to a market leading 50 degrees at its lowest position.
- Completely tool free set up makes it easy to install and operate.
- The completely water tight hand control floats, making it easy to access the control buttons while bathing.
- Two piece constrution makes the lift easy to assemble/disassemble for transport.
- When not in use, back easily folds down for space-saving storage.
- Padded seat and backrest provides added comfort.
- Washable hygienic covers available in white, blue or grey.
- Hygienic cutaway and slimline profile for easy handling.
- Built in safety controls.
- 4 quick release suction cups with option for an additional 2 for anti-slip bath base.
- Headrest available as an optional extra.
- 308 lbs. weight capacity.
- 2 Year Warranty on battery/handset.
- Lifetime Warranty on frame & motor.



The battery is located in the hand control and can be easily unplugged from the lift for charging in an outlet.

bath safety

To order call toll free: 877.224.0946

drive
MEDICAL DESIGN & MANUFACTURING

REPLACEMENT ACCESSORIES

SEAT AND BACKREST COVERS

- 460900252 Cover-Set Comfort, White
- 460900312 Cover-Set Comfort, Blue
- 460900432 Cover-Set Comfort, Grey

SPARE PARTS

- 460900403 Charger for Bellavita
- 415252 Plug + Washer for Cover (1 ea.), White
- 415312 Plug + Washer for Cover (1 ea.), Blue
- 415432 Plug + Washer for Cover (1 ea.), Grey
- 460900213 Suction Cup 1/ea 2.95 " (75 mm) with clip, Grey
- 460900220 Set of Suction Cups 4/pk (4 x 2.95 " 75 mm) with Clips, Grey)
- 460900600 Hand Control Including Storage Battery

ACCESSORIES

- 101385 Rotating & Transfer Aid, White 1/bg
- 102385 Rotating Aid, White 1/bg
- 410200252 Headrest, Padded, White
- 410200312 Headrest, Padded, Blue
- 410200432 Headrest, Padded, Grey
- 460900205 Side Flap (1 ea.)



bellavita



HIGHEST SEAT POSITION



LIGHTEST BATHLIFT AVAILABLE



RECLINING BACKREST LOWEST BATHLIFT AVAILABLE

SPECIFICATIONS:

Overall length (with backrest in a reclined position)	32.6"
Overall width (seat incl. side flaps)	27.5"
Width of seat	13.7"
Depth of seat	16.5"
Overall length of seat	22.4"
Sitting height min.	2.3"
Sitting height max.	18.8"
Height of backrest	25.1"
Width of backrest	12.5" - 14.1"
Backrest reclines	50°
Width of baseplate	12.2"
Length of baseplate	22"
Load capacity max.	308 lbs.
Total weight (excl. hand control)	20.5 lbs.