

Gel "E" 2" General Use Gel/Foam Wheelchair Cushion

HCPCS: E2601




FIGURE A



FIGURE B

- The Medicare approved General Use 2" gel-foam wheelchair cushion will provide comfort and support for the wheelchair bound patient
- High-Density foam encapsulates a viscous gel bladder that provides proper immersion and envelopment for pressure redistribution
- The cover is fluid resistant low-shear nylon zippered and removable with a durable vinyl bottom
- Carry handle included for convenient transport
- Product Dimensions: 18" (W) X 16" (D) X 2" (H)
- Meets flammability standard TB117

PRODUCTS

Description

14888	W/C Cushion,2in Gel/Foam,18x16x2in,1/cs
-------	-----------------------------------------

SPECIFICATIONS

Description

Cover Material	Fluid-Resistant Nylon Waterproof Vinyl Bottom
Warranty	1 Year Limited
Weight Capacity	250 lbs
Weight Capacity (13046 13046-HD)	250 lbs.