

## Deluxe Steel Drop-Arm Commode with Padded Seat

HCPCS: E0165



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Easy-to-assemble without tools in 4 easy steps (Figure A - D)
- Padded open front vinyl toilet seat
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode
- Easy-to-clean gray powder-coated steel finish
- Padded armrests and back provide extra comfort and support
- Comes with 12 qt commode bucket with carry handle and splash shield
- Made using recycled materials

PRODUCTS	Description
11125PSKD-1	Commode,DropArm,Padded Seat,KD,1/cs

  

SPECIFICATIONS	Description
Carton Shipping Weight	23 lbs.
Height	19" - 23"
Seat Dimensions	15.25" (W) x 16.75" (D)
Warranty	Limited Lifetime
Weight	21 lbs.
Weight Capacity	300 lbs
Width	19" between arms   23.75" Outside

  

PARTS	Description
11125PSKD-AL	LeftArmWithReleaseLever11125PSKD,1ea, EA
11125L	Replacement leg, 11125,1/ea, EA
12002PB	Push Button for Bath Bench,Walker,1/ea, EA
11125PSKD-B	Back for 11125PSKD, 1/ea, EA