

INVACARE BED OVERVIEW

Homecare and LTC Bed Systems

MARC GELLIDA | PRODUCT MANAGER



Yes, you can.®

ICC/LTC Beds

This document is intended for the internal use only and may not be distributed externally or reproduced for external distribution in any form



ICC/LTC Bed Portfolio



Yes, you can.®




ICC/LTC BEDS

820 DLX



SC900 DLX



Carroll CS3



Carroll CS5



Carroll CS7



Carroll CS9 FX600





Product Flow | CS Series



Universal CS Features

CS3



CS5



CS7



CS9 FX600



Cover general nursing, memory care, short term rehab and long term acute care.



Pure vertical
rise




Side-to-side
pendant
adaptability



Slat deck



Tool-free
assembly



Tubular
frame



Wood or
laminate panels
in variety of
styles/colors



76" to 80" to
84" sleep
surface



INVACARE CS3




Yes, you can.®

Ergonomic
Hand Pendant



Optional fixed 36", 39" or
42" deck widths



450lb total
weight capacity



Optional battery
back-up



Analog
electronics

27" to 8.25"



Mobile only at
low height

Universal CS
Features



Pure vertical
rise




Side-to-side
pendant adaptability



Slat deck



Tool-free
assembly



Tubular
frame



Wood or laminate
panels in variety
of styles/colors



76" to 80" to
84" sleep
surface



INVACARE CS5



Yes, you can.®

Ergonomic
Hand Pendant



Optional fixed 36", 39" or
42" deck widths



500lb total
weight capacity



360° Mobility
at any height



Analog
electronics

30" to 9"

Optional One-Touch
Floor Lock



Universal CS
Features



Pure vertical
rise




Side-to-side
pendant adaptability



Slat deck



Tool-free
assembly



Tubular
frame



Wood or laminate
panels in variety
of styles/colors



76" to 80" to
84" sleep
surface



INVACARE CS7



Yes, you can.®


Auto Contour




Optional fixed 36",
39" or 42" deck widths



500lb total
weight capacity



Lowest bed: 30" to 7"



Optional attendant
control with
Trendelenburg



Advanced electronics
with auto-leveling
motors



Open Bus
Electronic
(AC/TR/AL)



Mobile at any height
with floor lock

Universal CS
Features



Pure vertical
rise




Side-to-side
pendant adaptability



Slat deck



Tool-free
assembly



Tubular
frame



Wood or laminate
panels in variety
of styles/colors



76" to 80" to
84" sleep
surface



INVACARE CS9 FX600





Yes, you can.®

Built in adjustable width
from 36", 39" and 42"

Bariatric 600lb total
weight capacity




6- Function
Ergonomic
Hand Pendant



Adjustable width
Frame/Bed ends

30" to 9"



Optional attendant
control with
Trendelenburg



Advanced electronics
with auto-leveling
motors



Open Bus
Electronic
(AC/TR/AL)



Mobile at any height
with floor lock

Universal CS
Features



Pure vertical
rise




Side-to-side
pendant adaptability



Slat deck



Tool-free
assembly



Tubular
frame



Wood or laminate
panels in variety
of styles/colors



76" to 80" to
84" sleep
surface



INVACARE CS line comparison



Yes, you can.®



Key Specifications	CS3	CS5	CS7	CS9
Travel range	8.25" – 27"	8.75" – 30"	7" – 30"	9" – 30"
Weight Capacity	450 Lbs	500 Lbs		600 Lbs
Bed Width	Fixed 36"/39"/42"			Adjustable
Bed Length	76"/80"			
Bed Length Extender (80"- 84")	Optional			
Bed Ends	Wood or Laminate			
Mobile Height	Low	Any		
360° Maneuverability	Yes			
Locking Mechanism	Locking Casters (One-Touch optional)		One-Touch Floor Lock	
Vertical Hi-Lo Travel	Pure Vertical			
Adjustable Bed Ends	No			Yes
Deck Design	Slat			



INVACARE CS line comparison



Yes, you can.®



Electronics	CS3	CS5	CS7	CS9
Actuators/Motors	4 DC			
Pendant Style	Rubberized			
Pendant buttons	6-button		8-button	
Side to Side Adaptability	Yes			
Auto Contour	No		Yes	
Trendelenbrg w/ Reverse Trend	No		Optional (in ACP)	
Auto Leveling	No		Yes	

Others	CS3	CS5	CS7	CS9
Warranty	2 years on electronics 10 years on frame	3 years on electronics 10 years on frame	2 years on electronics 10 years on frame	2 years on electronics 10 years on frame



INVACARE CS Series




Yes, you can.®



CS Rail Options

Assist Rail



- ✓ IHCSRLAS
IHCSRLAS-39
IHCSRLAS-42
- ✓ 3 Positions:
 - 10 Guard
 - 10 Assist for Transfers
 - 10 Clear for Transfers
- ✓ Available for 39" or 42" decks.
- ✓ Available without "guard" position.

Thinksoft Positioning Device



- ✓ IHCSTSPD
- ✓ Blow molded device provides positioning support at the head end.
- ✓ Mounts to head deck, moves with articulation of head.
- ✓ Fits all widths.

Assist Bar



- ✓ IHCSBAAS
IHCSBAAS-39
IHCSBAAS-42
- ✓ Folds flat for transfers, provides handle for transfers.
- ✓ Available for 39" or 42" decks.

1¾ or Partial Rail



- ✓ IHCSRLPA
- ✓ 68" rail leaves area at the foot end for access to care or transfers.



INVACARE CS and DLX Series



CS and DLX Bed End Options



Shapes



Kensington Maple




Solar Oak



Amber Cherry




Tuscan Walnut




Williamsburg Cherry



Biltmore Cherry



African Walnut



Cafelle

Finishes