



CS7 bed featured with Invacare® ThinkSoft® Positioning Device, Drake Bed Ends and Invacare® Solace® Mattress

## CARROLL CS SERIES™ CS7™ BED



**Yes, you can.®**

The innovative **Carroll CS Series™ CS7™ Bed** is designed to provide the ultimate versatility for any facility or resident. The ultra low height works well for Memory Care Wings, while the ratcheting foot lift and optional Trendelenburg offer solutions for the Rehab Wing. 30" Working height gives the caregiver comfort while performing routine tasks for every resident.

The Auto Contour feature is comfort for the resident by simultaneously raising both the head and knees, which is designed to help keep residents comfortably in place while being better positioned for reading, watching TV or talking to family members.

Invacare® beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the recommended FDA Guidelines for bed rail entrapment. Visit [www.invacare.com/bedsafety](http://www.invacare.com/bedsafety) for additional safety and entrapment information.

## Features

**Advanced Electronics** – Designed with synchronized motors to keep the bed level at all times, providing stability throughout the travel range.

**One-Touch Floor Lock** – Patented one-touch floor lock secures the bed at any height.

**Pure Vertical Rise** – Raises straight up and down throughout the travel range, with zero horizontal movement of the casters, minimizing potential damage to the wall or flooring.

**Exclusive Slat Deck Design** – Designed to provide a comfortable sleep surface, allows for improved infection control and helps reduce wear and tear on the mattress.

**Ergonomic Hand Pendant** – Water-resistant<sup>1</sup>, rugged hand pendant designed to fit comfortably in the resident's hand.

### Customizable Options –

- 39" or 42" Wide Deck
- 4" Bed Extender (84" Length)
- Attendant Control Panel
- Trendelenburg / Reverse Trendelenburg / Comfort Chair positioning
- Battery Backup on board or portable
- Hi-Lo Lockout
- 4 Rail options



1. Tool-less Assembly  
2. One-Touch Floor Lock  
3. Ergonomic Hand Pendant  
4. Slat Deck Design

**Invacare Corporation**  
[www.invacare.com](http://www.invacare.com)

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2020

Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

Form No. 20-032  
200204

Model #	Height Range	Sleep Surface	Weight Capacity*	Regulation Compliance	Limited Warranty
IHCS7	7" - 30" H	76" or 80" L 36", 39", or 42" W	500 lb.	†Classified to: UL 60601-1 IEC 60601-2-38 CAN/CSA C22.2 No. 601.1	10 Years Frame 5 Years Mechanical 3 Years Electronics

NOTE: Whenever bed rails are used as part of any hospital or homecare bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

\* Maximum safe working load.

1. Pendant is water-resistant to IPX6 per IEC 60529.

†Testing conducted to select portions of UL/IEC Standards.

