
REVOLUTION SPINEBOARD ATTACHMENT



THE REVOLUTION LIFT CAN BE PURCHASED WITH AN OPTIONAL SPINEBOARD ATTACHMENT (WITH SPINEBOARD) FOR ADDITIONAL FUNCTIONALITY. THE SPINEBOARD ATTACHMENT WILL WORK IN BOTH ABOVE-GROUND AND IN-GROUND CONFIGURATIONS, AND WILL ALSO CLEAR A 34" [864mm] SPA WALL WITH THE STANDARD SEAT-ARM. WITH THE OPTIONAL SPA SEAT-ARM THE SPINEBOARD ATTACHMENT WILL CLEAR A 48" [1219mm] WALL. THE SPINEBOARD MEASURES 73" [1854mm] LONG BY 18-1/4" [464mm] WIDE BY 2-1/2" [64mm] THICK, AND IS MADE OF ROTATIONALLY MOLDED LLDPE. THE ATTACHMENT IS MADE OF TYPE 304 STAINLESS STEEL, AND IS POLYESTER POWDER COATED.

*NOTE: WHEN USING THIS ATTACHMENT WITH A REVOLUTION LIFT, THE LIFT IS DE-RATED TO 350 POUNDS FOR THE STANDARD REVOLUTION AND 300 POUNDS FOR THE DEEP DRAFT REVOLUTION!



9889 GARRY MORE LANE MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602

NOTE:

*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE *ALL DIMENSIONS ARE IN INCHES

PRODUCT:	F-734RSA	
DWG#: F-734RSA-1		
SCALE:	1:30	DATE: 7/8/2020
DRAWN BY:	KR	REVISION: