

MILFORD / SMART TRANSFER FEATURE SHEET



THREE WAYS TO USE THE MILFORD



IN YOUR VEHICLE

- Compatible with nearly every vehicle on the road
- A single Milford Lift can be used with multiple vehicles by installing Mounting Posts in each vehicle
- Allows customer to keep their existing vehicle, OEM vehicle seat and wheelchair
- Eliminates virtually all effort from the operator
- The only permanent fixture is the Mounting Post, which was designed to be as compact and unobtrusive as possible
- 330lb Maximum Lifting Capacity
- Weight of Milford Arm is about 25lbs

AROUND THE HOME

- The Smartbase accessory allows the Milford to go mobile
- Includes a rechargeable battery for long-lasting power
- Wheels lock for safety and stability, and the leg width can be adjusted to fit the space
- Folds flat for easy storage and great portability
- Ideal for use in larger rooms in the home or for use during travel
- 330lb Maximum Lifting Capacity
- Smartbase w/Battery weight is about 60lbs





IN SMALL SPACES

- The Wall Bracket accessory is ideal for transfer in tight spaces such as the bathroom, staircase, bedroom or favorite chair
- With extensions, the Milford on a Wall Bracket can reach from 39" to 70" (each extension will reduce lifting weight capacity)
- The Milford can be transferred from vehicle to the home, no need to purchase more than one Milford
- Wide variety of slings available to suit the customers needs