

There is Power Mobility. And then there's

POWERFUL MOBILITY



CityCruzer



There is Power Mobility. And then there's

POWERFUL MOBILITY

CityCruzer

With a 285 lb load capacity and 4 wheels, the CityCruzer can transport higher weighted individuals with stability not usually found in portable scooters. The folding mobility seat has deep cushions and armrests that flip up and/or can be adjusted laterally.



MAX SPEED
4.0
MPH



LOAD CAPACITY
285
LBS



BATTERY
12 Ah
SLA



RANGE
8
MILES

Overall Length	42.5"
Overall Width	19.5"
Overall Height (at Maximum Tiller Height)	37"
Legroom - Floorboard Length x Width	11.5" x 16"
Deck Height - Floorboard to Floor	4.7"
Ground Clearance	2"
Turning Radius	53"
Wheelbase	32"
Front Section	45 lbs
Rear Section	26 lbs
Battery Pack	30 lbs
Seat	25 lbs
Weight with Batteries	127 lbs
Weight without Batteries	97 lbs

Seat Type	Fish on seat
Seat - Width & Depth & Height	16" x 14" x 14"
Seat - Cushion Thickness	1.5" - 2"
Seat Height - Cushion to Deck	17" - 19"
Seat Height - Cushion to Floor	21" - 24"
Tire Type	PU
Rim Type Rear	Plastic
Rim Type Front	Plastic
Outside Diameter Rear	8"
Tire Width Rear	2"
Outside Diameter Front	8"
Tire Width Front	2"
Battery Type & Power	12AH (standard) 20AH (optional)

Brake System	Electromagnetic brake
Climbing Angle - Maximum Recommended	8°
Charger Type	Input AC 100-240V, Output 24V, 2 Amp
Charging Time	Up to 8 hours
Controller	Dynamic 50A
*Driving Range	Up to 8 miles (12AH) Up to 13 miles (20AH)
Load Capacity	285 lbs
Motor Size	270 Watts
*Speed - Maximum	4 mph
Colors	Pearl Green, Blue, Silver and Black
Tiller Type	Delta
Throttle Control Type	speed control lever
Anti-Tipper Wheels	Rear

(*) SPECIFICATION PERFORMANCE IS DETERMINED IN THE LABORATORY SETTING BASED ON IDEAL CONDITIONS (MINIMAL FRICTION, NO WIND AND CONTROLLED TEMPERATURE AND PRESSURE WITH AN AVERAGE STANDARD WEIGHT). ACTUAL PERFORMANCE WILL BE DIFFERENT THAN RESULTS OBTAINED DURING LABORATORY TESTING.

Call us at 1-800-944-1997 | sales@evrider.com | www.evrider.com