



Uplift[®]
Safe Mobility Solutions



Safe Patient Handling and Falls Management Catalog

Stay in sync with upcoming ANA recommendations.

The power of Uplift

Uplift is our breakthrough program designed to help protect patients and clinicians. It mirrors the American Nurses Association (ANA) safe patient handling and mobility (SPHM) guidelines to integrate safe patient handling, falls management and early mobility into one united initiative. Inside this catalog, you will find a wide array of products that support the Uplift program and your patient mobility goals.




Table of Contents

Products that count as qualifying spend for the Uplift Program

Uplift® Program

Uplift Program	4–5
----------------	-----

Repositioning and Lateral Transfer

Comfort Glide	
Repositioning Sheets	6
Air-Assisted Repositioning	7
Fitted Sheets	8
Offloading Foam Wedges	8
Drypads	8
Inflatable Lateral Transfer	10–11
Disposable Slide Sheets	9

Patient Lifts

Overhead Lift System	12
Powered Total Lift	13
Manual Patient Lift Cradle	13
Hoyer® Presence® Patient Lift	14
Hoyer Stature® Patient Lift	14
Electric Stand Assist Lift	15
Hoyer Stand Assist Lift	15
Raizer® Mobile Lifting Chair	16–17
Reflex Shower Chair	18

Patient Slings

U-Shaped Patient Slings	19–22
U-Shaped Clip Patient Slings	23
Repositioning Patient Slings	24–25
Full Body Patient Slings	26
Stand Assist Patient Slings	27–28
Specialty Patient Slings	29–30
Sling Materials Chart	31

Patient Alarms

Fall Alarm Monitors	32–33
Crest Call Cords	34
Nurse Call Cables	34
Alarm Accessories	34
Reusable Bed/Chair Sensor Pads	35
Premium Alarm Sensor Pads	35
Economy/Single Patient Use Sensor Pads	36
Advantage® Magnetic Patient Alarm	37
Anti-Wander Alarm Kits and Accessories	37
Smart Caregiver® Wireless Monitoring System	38
Sensible® Care Proactive Patient Monitoring System	39

Gait Belts and Safety Essentials

Gait Belts	40–41
Hand Protectors	42
Multi-Purpose Arm and Leg Splints	43
Activity Apron	43
Activity Overlay	43
Limb Holders	44–47
Repositioning Belts	48
Breezeline Pelvic Holders	49
Plaid Safety Vests and Jackets	49–50
Pulse Oximeter Probe Wraps	50
Infant ID Bracelets	50
Circumstraint Board Straps	50

Products that do not qualify for the Uplift Program

Patient Alarm

Infrared Patient Bed Alarm	38
----------------------------	----

Identification and Protective Products

Falls Prevention Slippers	51
PerforMAX™ Falls Prevention IV Gowns	52
Falls Prevention Kits	53
Fall Mats	54
Hip Protectors	54
Men's Premium Hip Protector Briefs	55
Hipsters Men's Briefs	55–56
Hipsters Unisex Sweatpants	56
Unisex Protective Caps	56
Knockdown Transfer Bench with Microban®	57

Knockdown Bath Bench with Microban	57
Commodes	57–58
Grab Bars	59
Toilet Safety Rail	59
Raised Toilet Seats	60
Toilet Sensor Pads	60
Urinals	61
Spil-Pruf Male Urinals	61
UltraFlex® Silicone Male External Catheters	61
Seizure Side Rail Pads	62
Torso Wedge	62
Bed Assist Bar	63
Side Rail Protector	63
Advantage Mattress Bolster	63



Uplift®

Safe Mobility Solutions

Make a move to help protect both patients and caregivers with safe mobility.

The challenge

Until now, there's been no single program to address the integrated facets of safe patient handling and falls management to promote patient and caregiver safety through early mobility.

Mobility-related injuries drive up costs on both ends of the spectrum.



Patients

1M

An estimated 700,000 to 1 million hospital patients fall each year.¹

\$57K

The average cost of treating a patient with Ventilator Acquired Pneumonia.²

\$40K

Added cost for a pressure injury.³

\$9,600

Added cost for readmission.⁴

Caregivers

#1 HIGHEST INJURY RATES

Caregivers suffer the highest musculoskeletal injury rates, resulting in the most lost work days of any profession.⁵

\$37,000

The average back injury claim.⁶

20%

Nurses who leave direct patient care positions due to risks associated with work.⁷

The solution

Uplift your team with safe mobility made simple.

Introducing Uplift, the first integrated solution designed to keep caregivers and patients safe through every stage of progressive mobility in one easy-to-implement program.

Safe Patient
Handling

Falls
Management

Early Mobility

Uplift Program

Every aspect of safe mobility is interconnected, so we treat them as a single holistic solution.

Uplift covers all three stages of progressive mobility:



1. Discovery Assessment

We understand that every facility is different. That's why we conduct discovery assessments to understand your unique challenges and opportunities from every angle.




Boosting and repositioning



2. Custom Recommendations

Based on those findings, we develop a solution that works for you, from policy reviews to environmental changes to product recommendations. No two plans are alike.



Sitting to standing



3. Implementation

We make the transition to Uplift simple and seamless with practical tools that get everybody on board. A comprehensive onboarding kit, onsite clinical training and ongoing education help promote a culture of safety.



Walking

Ask your Medline Representative how you can qualify for the Uplift program.

The Comfort Glide® Solution

Take the pain out of repositioning.

Comfort Glide repositioning products work together to ease the stress of turning on both caregivers and patients. These include repositioning sheets, wedges, drypads, and fitted sheets.

Comfort Glide Repositioning Sheets, Single-Patient Use

- Low-friction nylon bottom surface requires 35% less effort than a draw sheet to initiate and maintain sliding⁸
- Multiple handles help optimize weight distribution and leverage
- Air-permeable materials allow the sheet to remain under the patient for an extended period of time and are always available to aid repositioning⁸
- Soft, breathable microfiber covers the entire top surface for patient comfort and is compatible with low air-loss mattresses⁸
- Optional retention straps anchor the sheet to the bed
- Safety rated for up to 575-lbs.; bariatric sheet rated to 1,000-lbs.



Item No.	Description	Pkg.
CGLIDESHEETNS	Glide Sheet 35 x 55"	12/cs
CGLIDESHEET	Glide Sheet 35 x 55", with Straps	12/cs
CGSHEETXLD	Glide Sheet, X-Large, 42 x 78", Two Optional Straps	12/cs
MSC50013	Glide Sheet, 35 x 55", Tacky Top	12/cs
CGLIDESHEETM	Glide Sheet, Mini, 30 x 35"	12/cs

Comfort Glide Repositioning Sheet and Sling in One (see pg. 24)




MDSC40018

Item No.	Description	Size	Pkg.
MDSC40018	Patent-Pending 2-in-1 Repositioning Sheet and Sling, Disposable, 1,000-lbs. Weight Capacity	42 x 78"	12/cs
MDSC40017	Patent-Pending 2-in-1 Repositioning Sheet and Sling, Disposable, 1,000-lbs Weight Capacity	42 x 65"	12/cs

Two products in one



Premium glide sheet




Sling for patient lifts

Comfort Glide Air

Inflatable repositioning and transfer sheet gives caregivers a safe and physically easy way to move patients.

- Cushion of air to help transfer and reposition patients without lifting
- Easy access handles allow for ergonomic grip and repositioning
- Wider size helps optimize leverage
- Air-permeable so it can remain under patients to aid repositioning
- Requires 80% less effort than a standard draw sheet to move a patient⁸
- Safety rated for patients up to 1,000-lbs.⁸

Item No.	Description	Pkg.
MSC60011	Comfort Glide Air Inflatable Sheet	12/cs
MSC60090	Comfort Glide Air Blower	1/ea
MSC60099	Comfort Glide Air Cart	1/ea





CGLIDEBASE36

Comfort Glide Fitted Sheets

Low-friction, breathable polyester fabric to help reduce the effort to slide a patient.

- Contoured elastic perimeter hugs the mattress for a close fit without awkward buckles
- Compatible with low air-loss mattresses
- Can be laundered up to 50 times
- Available in four different sizes to fit a wide range of mattresses

Item No.	Description	Pkg.
CGLIDEBASE36	Fitted Base Sheet, 36 x 80", Blue	12/cs
CGLIDEBASE42	Fitted Base Sheet, 42 x 80", Purple	12/cs
CGLIDEBASE50	Fitted Base Sheet, 50 x 80", Black	12/cs
CGLIDEBASE54	Fitted Base Sheet, 54 x 80", Black	12/cs

Comfort Glide Foam Wedges

High-density foam helps keep patients in the recommended 30° tilt position, reducing pressure on bony prominences.

- Non-slip base grips the surface of any mattress or linen
- "Down" arrows ensure proper placement
- Soft, composite fabric covering can be wiped clean
- Holds patients in a 30-degree tilt position
- High density foam doesn't bottom out



MSC50023

Item No.	Description	Pkg.
MSC50023	30° Foam Wedges, Non-Slip Base: 11 x 5 x 16"	2/pk, 6/pk/cs
MSC50024	30° Foam Wedges, Non-Slip Base: 15 x 6½ x 16"	2/pk, 2/pk/cs
MSC50025A	30° Foam Wedges, Non-Slip Base: 24 x 12 x 7"	2/cs
MSC50070	15° Foam Wedges, Non-Slip Base: 16 x 11 x 3"	2/pk, 6 pk/cs

Comfort Glide Drypads

These ultra-absorbent level drypads are designed to manage moisture over an extended period of time and keep patients 26 times drier than the leading competitor.⁸

- Draws in moisture, locks it away from skin and feels dry to the touch within minutes⁸
- Air-permeable for optimal patient comfort and zero strikethrough⁸
- Compatible with low air-loss mattresses⁸
- Sizes designed to completely cover the repositioning sheets

57" Standard Length
CGLIDEPAD

Item No.	Description	Pkg.
CGLIDEPAD	Drypads, 36 x 57"	5/bg, 60/cs
CGLIDEPADXL	X-Large Drypads, 40 x 80"	5/bg, 40/cs

Disposable Slide Sheets

Makes lateral transfers and patient repositioning easier.

- Variety of styles to support your facility's repositioning and lateral transfer protocols
- Single use designs are available when you need them and avoid infection risk
- Convenient storage options available, including wall mounted holder

Item No.	Description	Pkg.
MDSSLIDE5848	Disposable Slide Sheet, Tube, 58 x 48"	50/cs
MDSSLIDEWALL	Optional Wall Holder for MDSSLIDE5848	1/ea
MDSMEDSLIDE	Disposable Slide Sheet, Tube, 33 x 39"	60/cs
MDS30180	Disposable Slide Sheet, Tube, 28 x 28"	30/ea




MDSSLIDE5848



MDSSLIDEWALL




MDSMEDSLIDE



MDS30180

Comfort Glide Lateral Transfer

Comfort Glide LT provides caregivers with a safe and easy way to laterally transfer patients of all sizes from one surface to the next by raising the patient in one smooth motion on a cushion of air, without the need for lifting. The reduced exertion helps lower the risk of injury and requires less manpower to complete a move. The pre-printed surface instructions makes the process easy and the breathable fabric helps prevent shear and friction forces on the patient's skin, for a simple and comfortable transfer experience.



56%

Reduction in injury risk following safe patient handling and mobilization program implementation in 44 intervention sites.⁹

75%

Reduction in effort compared to standard of care.⁸

Standard

MSC60042

Standard transfer of any patient.



Bariatric

MSC600LTXL

Transfer larger patient to a specialty bed.



Narrow

MSC600LT34

Transfer any patient to a more narrow surface.



Half Length

MSC634LTHF

For transfers wherever the bed/table can drop off.



- Comfort Glide LT supports the body for comfort during transfer
- Sheet can remain under the patient throughout their hospital stay
- Perimeter and extended handles allow multiple caregivers to assist with the move
- Helps reduce effort by 75% compared to standard of care⁸

**Less effort
means less risk**

Contact your Medline Representative today.

Comfort Glide LT

Item No.	Description	Size	Pkg
MSC60042	Comfort Glide LT Inflatable Sheet, With Drypad	40 x 80"	12/cs
MSC60041	Comfort Glide LT Inflatable Sheet	40 x 80"	12/cs
MSC600LTXL	Comfort Glide LT Inflatable Sheet, Bariatric	52 x 80"	6/cs
MSC600LT34	Comfort Glide LT Inflatable Sheet, Narrow	34 x 80"	12/cs
MSC634LTHF	Comfort Glide LT Inflatable Sheet, Half Length	34 x 46"	12/cs

Accessories

Item No.	Description	Size	Pkg
CGLIDEPAD	Drypad, Super Absorbent Polymer	36 x 57"	60/cs
MSC60090	Comfort Glide Air Blower	—	1/ea
MSC60099	Comfort Glide Air Blower Mobile Cart with Basket	—	1/ea
DYNJE4250	Hose Protection Sleeve	5 x 96"	25/cs
MSC600925	MRI-Safe Hose	25'	1/ea

Patient Lifts

Elevate patient and caregiver safety.

The Association of Safe Patient Handling Professionals (ASPHP) recommends one mobile lift for every eight patients who need one. Our lifts can help you reduce risks of caregiver injuries and lift patients with exceptional safety, convenience and dependability.

Overhead Lift System

This system makes it simple to own and install overhead lifts in almost any facility.

- Single-motor has a weight capacity of up to 600-lbs., or two motors can be combined with reinforced track to handle up to 1,000-lbs.
- Auto-locking cradle design guarantees that slings stay secured to the cradle during lifting
- Track lengths may be trimmed to any room size that is 7 meters or less
- Industry-leading 3-year warranty on the lift motor



Motor Components

Item No.	Description	Pkg.
ERGLUNA	600-lbs. Capacity Overhead Lift Motor	1/ea
ERGHANDCTRL	Hand Control Pendant for Overhead Lift	1/ea
ERGSPREADER	Spreader Bar for Overhead Lift	1/ea
ERGCHARGER	Charging Station for Overhead Lift	1/ea

Contact
your Medline
representative
to quote overhead
patient lifts



ERGLUNA



ERGHANDCTRL



ERGSPREADER

Mobile Lifts

Easy-to-use and reliable mobile lifts designed for all essential clinical tasks including patient transfers and repositioning.

- Emergency release provides for smooth lowering in case the battery dies unexpectedly.
- State-of-the-art control panel displays battery charge level and more
- Removable battery pack reduces downtime with off-board charging
- Can be lowered to floor height for lifting a patient who has fallen
- Compatible with lift scales for convenient weighing of patients
- Six-point cradle for maximum flexibility and safety
- Warranty: three years on all mechanical and electrical parts, one year on the battery



MDS700EL
(Shown with optional
gait training arms)

Medline Mobile Patient Lifts

Item No.	Description	Pkg.
MDS450EL	450-lbs. Capacity Lift with Powered Base	1/ea
MDS700EL	700-lbs. Capacity Lift with Powered Base	1/ea
MDSMELSCALE	700-lbs. Capacity Scale Accessory for Mobile Lift	1/ea
MDS700ARMS	Gait Training Arms Accessory for MDS700EL	1/ea

Medline Manual Patient Lift Cradle

Enables caregivers to position patients with one hand, providing better control and helping prevent injuries.

- Helps to eliminate grabbing and pulling
- Improves patient confidence
- Makes the caregiver's job easier
- Compatible with Medline mobile patient lifts MDS450EL and MDS700EL



Item No.	Description	Pkg.
MDS700MCRD	Patient Lift Cradle (Use with Clip Slings)	1/ea

Hoyer® Lifts

Feature-rich total lifts with lightweight aluminum frames from the trusted Hoyer® brand of lifting products.



Hoyer® Presence® Patient Lift

- Intuitive push pad assists with initiating movement
- Smart monitor technology drives user compliance and lift longevity
- Lightweight aluminum construction
- Low leg option kit for the Hoyer Presence lift is designed to fit under beds, stretchers, and other furnishings with low clearance

Item No.	Description	Weight Capacity	Pkg.
JRN500PELLB	Lift, Presence, Low Base, Electric	500-lbs.	1/ea
JRN500PELSLB	Lift, Presence, Low Base, Electric with Scale	500-lbs.	1/ea




Hoyer® Stature® Patient Lift

- Intuitive Push Pad assists with initiating movement
- Swan Neck style legs allow Stature to get closer to the resident
- True Vertical Lifting
- Lightweight aluminum construction
- Powered clip style positioning cradle (optional)

Item No.	Description	Weight Capacity	Pkg.
JRN500SEL	Lift, Stature, Electric	500-lbs.	1/ea
JRN500SELS	Lift, Stature, Electric with Scale	500-lbs.	1/ea
JRN500CR	Powered Clip Cradle for JRN500SEL/SELS	500-lbs.	1/ea

Stand Assist Lifts

For transferring activities with patients who have some ability to bear weight. Or to provide support for patients to stand and practice bearing weight.



Medline Electric Stand Assist Lift

- Adjustable knee pad for finding the perfect fit for anyone
- State-of-the-art control panel displays battery charge level and more
- Removable battery pack reduces downtime with off-board charging
- Warranty: three years on all mechanical and electrical parts, one year on the battery

Item No.	Description	Pkg.
MDS500SA	500-lbs. Capacity Stand Assist with Powered Base	1/ea



Hoyer® Stand Assist Lift

Stand assist with integrated scale.

- Built-in scale to seamlessly weigh patients during transport
- Intuitive Push Pad assists with initiating movement
- Smart Monitor technology drives user compliance and lift longevity
- One hand adjustable kneepad with infinite settings
- Angled, removable foot tray

Item No.	Description	Weight Capacity	Pkg.
JRN440SAS	With Scale	440-lbs.	1/ea



Service plans for extra confidence

Keep patients and caregivers safe while protecting your facility from costly last-minute repairs with your choice of two service plans:

- **Preventive Maintenance Plan**—full inspection and lubrication of moving parts to ensure proper lift operation
- **Full-Service Plan**—all the advantages of the Preventive Maintenance Plan, plus unlimited troubleshooting and repair hours—ideal for equipment that's past warranty

Raizer® Mobile Lifting Chair

Now there's a better way to lift someone after a fall.

This sleek, innovative, battery-operated chair enables one assistant to safely raise a large adult to a sitting or standing position in just a few minutes.

Raizer Mobile Lifting Chair for Safe Patient Handling by LiftUp

- Easy to carry—Two handy bags hold seven lightweight components
- Simple setup—Color-coding guides the placement of components, and a loud click confirms each connection
- Gentle on caregiver—Requires virtually no physical effort besides a supportive hand
- Fast—In minutes, one assistant can set up and raise a large adult from the floor to a standing position
- Easy to clean—Cleaning's a snap with the smooth surfaces on all Raizer components

Item No.	Description	Pkg.
MDSRAIZER	Raizer Mobile Lifting Chair	1/ea



Mechanically lift using the console buttons or remote control

Powerful rechargeable battery provides up to 40 lifts on one charge

Small footprint fits in tight spaces

Patients are lifted safely and comfortably from the floor in minutes



Add this innovative chair to your safe patient handling program and prevent patient and staff injuries.

All in 4 easy steps:



1 Place the seat motor under the person's knees



2 Attach components to the seat motor



3 Mechanically lift the patient using the console buttons or remote control



4 Transfer the patient to a walker or wheelchair



Color-coding makes setup simple



Two handy carrying bags hold seven lightweight components

Technical Specifications

Weight Capacity	330-lb.
Lifting Time	20–30 secs.
Product Weight:	28-lb. total
Charger	Two pieces: Soniel 12v, 100–240V 12V car adapter
Number of lifts on a full charge	40 at maximum weight capacity
Charging Time	6 hours

Lifting Shower Chair

Make bathroom tasks easier, more efficient and safer for the caregiver.



Use with a bedpan or over a toilet



Easily tilt backward for feet access



Adjust height and tilt with simple controls

Reflex Shower Chair

This lifting shower chair is designed with OSHA requirements in mind to help you reduce staff risks. Plus, working at the optimal height makes bathing tasks easier, faster and more efficient.

- Smooth, energy efficient electric height and tilt adjustments
- Intuitive, easy-to-use control console
- Multifunctional for showering and toileting
- Comfortable foam padding and adjustable headrest
- Arms lock down for safety or fold away for easier transfers

Item No.	Description	Pkg.
MDSREFLEX	300-lbs. Capacity Shower Chair Lift	1/ea



Patient Slings

Choose the right slings for the right task in the right fabric

U-Shaped Patient Slings

Medline's U-shaped patient lift slings are designed for comfort, versatility, and ease of application. The underbody crossing design helps distribute weight, keeping the patient comfortable. This style of sling is easier to get underneath patients, particularly if they are in a seated position.

Single-Patient Use U-Shaped Slings

Standard Sling, Polypropylene Fabric

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDS2PTDHBM	Medium	38 x 48"	700-lbs.	85–220 lbs.	10/cs
MDS2PTDHBL	Large	39½ x 54"	700-lbs.	170–400 lbs.	10/cs
MDS2PTDHBXL	X-Large	42 x 56"	700-lbs.	320–700 lbs.	10/cs



MDS2PTDHBL



MDS2PTDLBL

Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDS2PTDLBM	Medium	38 x 39"	700-lbs.	85–220 lbs.	10/cs
MDS2PTDLBL	Large	40 x 44"	700-lbs.	170–400 lbs.	10/cs
MDS2PTDLBXL	X-Large	41 x 47"	700-lbs.	320–700 lbs.	10/cs

Comfort Sling, Polypropylene Fabric

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMDHS1	Small	33 x 49"	700-lbs.	100–250 lbs.	10/cs
MDSMDHS2	Medium	38½ x 54"	700-lbs.	200–450 lbs.	10/cs
MDSMDHS3	Large	42½ x 59"	700-lbs.	300–600 lbs.	10/cs
MDSMDHS4	X-Large	47 x 73"	700-lbs.	500–700 lbs.	10/cs



MDSMDHS1



MDSMD3

Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMD2	Medium	39 x 47"	700-lbs.	200–450 lbs.	10/cs
MDSMD3	Large	43 x 48"	700-lbs.	300–600 lbs.	10/cs



Reusable U-Shaped Slings

Standard Sling, Polyester Fabric

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDS2PTHBM	Medium	38 x 48"	700-lbs.	85–220 lbs.	1/ea
MDS2PTHBL	Large	39½ x 54"	700-lbs.	170–400 lbs.	1/ea
MDS2PTHBXL	X-Large	42 x 56"	700-lbs.	320–700 lbs.	1/ea

Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDS2PTLBM	Medium	38 x 39"	700-lbs.	85–220 lbs.	1/ea
MDS2PTLBL	Large	40 x 44"	700-lbs.	170–400 lbs.	1/ea
MDS2PTLBXL	X-Large	41 x 47"	700-lbs.	320–700 lbs.	1/ea



Comfort Sling, Polyester Mesh Fabric

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMMHS1	Small	33 x 49"	450-lbs.	85–150 lbs.	1/ea
MDSMMHS2	Medium	38½ x 54"	450-lbs.	100–250 lbs.	1/ea
MDSMMHS3	Large	42½ x 59"	450-lbs.	200–450 lbs.	1/ea

Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMM2	Medium	39 x 47"	450-lbs.	100–250 lbs.	1/ea
MDSMM3	Large	43 x 48"	450-lbs.	200–450 lbs.	1/ea

Comfort Sling, Padded Knitted Fabric

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMPHS1	Small	31½ x 50"	450-lbs.	85–150 lbs.	1/ea
MDSMPHS2	Medium	38 x 54"	450-lbs.	100–250 lbs.	1/ea
MDSMPHS3	Large	44 x 64"	450-lbs.	200–450 lbs.	1/ea
MDSMPHS4	Bariatric	48 x 73"	700-lbs.	400–700 lbs.	1/ea



Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMP2	Medium	38 x 50"	450-lbs.	100–250 lbs.	1/ea
MDSMP3	Large	44½ x 51"	450-lbs.	200–450 lbs.	1/ea
MDSMP4	Bariatric	46 x 59"	700-lbs.	400–700 lbs.	1/ea



Wide Leg Sling, Padded Knitted Fabric

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMHC70012	Medium	37 x 47"	500-lbs.	100–210 lbs.	1/ea
MDSMHC70011	Large	41 x 57"	500-lbs.	200–350 lbs.	1/ea
MDSMHC70010	Bariatric	46 x 60"	700-lbs.	290–700 lbs.	1/ea



Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMHC70002	Medium	34 x 36"	500-lbs.	100–210 lbs.	1/ea
MDSMHC70001	Large	41 x 57"	500-lbs.	200–350 lbs.	1/ea
MDSMHC70000	Bariatric	46 x 46"	700-lbs.	290–700 lbs.	1/ea





MDSPUNVM

Comfort Sling, Spacer Fabric

Recommended use: Toileting and routine lifting, especially for patients with pressure ulcers or limited upper-body function and tone.

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSPUNVS	Small	33½ x 41"	450-lbs.	85–195 lbs.	1/ea
MDSPUNVM	Medium	40 x 48"	450-lbs.	150–300 lbs.	1/ea
MDSPUNVL	Large	44 x 53"	450-lbs.	200–350 lbs.	1/ea
MDSPUNVXL	X-Large	47 x 60"	450-lbs.	300–450 lbs.	1/ea



MDSPHAML

Hammock Sling, Spacer Fabric

Recommended use: Toileting and routine lifting, especially for patients with pressure ulcers or limited upper-body function and tone.

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSPHAMS	Small	33½ x 37½"	450-lbs.	85–195 lbs.	1/ea
MDSPHAMM	Medium	38 x 42"	450-lbs.	150–300 lbs.	1/ea
MDSPHAML	Large	41 x 43"	450-lbs.	200–350 lbs.	1/ea
MDSPHAMXL	X-Large	46½ x 46½"	450-lbs.	300–450 lbs.	1/ea



Pediatric Sling, Polyester Mesh Fabric

Recommended use: Bathing and toileting of pediatric patients over 35-lb. (16 kg).

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSPM111	One Size Fits Most	21 x 33"	100-lbs.	1/ea



Pediatric Sling, Padded Knitted Fabric

Recommended use: Transfers of pediatric patients over 35-lb. (16 kg).

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSPPP121	One Size Fits Most	21 x 32"	100-lbs.	1/ea

U-Shaped Clip Patient Slings

Used with clip-style cradles, Medline's Clip Slings are designed for ease of positioning. Head support and padded leg supports offer patients comfort and protection.

Single-Patient Use Clip Sling

Comfort Clip Sling, Polypropylene Fabric



Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSCLIPSLNGD2	Medium	35 x 47"	600-lbs.	85–264 lbs.	10/cs
MDSCLIPSLNGD3	Large	38 x 51"	600-lbs.	220–440 lbs.	10/cs
MDSCLIPSLNGD4	X-Large	42 x 55"	600-lbs.	350–600 lbs.	10/cs



MDSCLIPSLNGD2

Reusable Clip Sling

Comfort Clip Sling, Polyester Fabric



Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSCLIPSLING2	Medium	35 x 47"	600-lbs.	85–264 lbs.	1/ea
MDSCLIPSLING3	Large	38 x 51"	600-lbs.	220–440 lbs.	1/ea
MDSCLIPSLING4	X-Large	42 x 55"	600-lbs.	350–600 lbs.	1/ea



MDSCLIPSLING2

Repositioning Patient Slings

Designed with breathability so they may be left under patients, Medline's Repositioning Slings are readily available for all lifting activities including turning, repositioning and transfers.

Single-Patient Use Repositioning Slings



Standard Repositioning Sling, Polyester Cotton Fabric

Item No.	Dimensions	Wt. Cap.	Pkg.
MDSMDREPO	42 x 78"	1,000-lbs.	10/cs



MDSC40018

2-in-1 Comfort Glide Sling and Repositioning Sheet, Microfiber/Nylon Fabric

2-in-1 product is both a sling that can be used with patient lifts, and a repositioning sheet for in-bed positioning or lateral transfers.

Item No.	Style	Dimensions	Wt. Cap.	Pkg.
MDSC40018	2-in-1 Repositioning Sling and Sheet	42 x 78"	1,000-lbs.	12/cs
MDSC40017	2-in-1 Repositioning Sling and Sheet	42 x 65"	1,000-lbs.	12/cs

Reusable Repositioning Slings

Standard Repositioning Sling, Polyester Cotton Fabric

Item No.	Dimensions	Wt. Cap.	Pkg.
MDSTRIRR01	42 x 78"	1,100-lbs.	1/ea



2-in-1 ComfortRise Sling and Repositioning Sheet, Spacer Fabric

Item No.	Style	Dimensions	Wt. Cap.	Pkg.
CLIFT	2-in-1 Sling and Repositioning Sheet	42 x 78"	1,000-lbs.	1/ea



Repositioning Sling, Spacer Fabric

Recommended use: Turning or repositioning patients in bed, especially for those with pressure ulcers.

Item No.	Style	Dimensions	Wt. Cap.	Pkg.
MDSPPOS	Flat	42 x 78"	440-lbs.	1/ea
MDSPPOSE	Fitted, with Elastic	42 x 78"	440-lbs.	1/ea



MDSPPOS

Full Body Patient Slings

With their comfortable and secure design, Medline's Fully Body Slings are available in various reusable materials to match all patient lifting needs.



MDSMR110

Reusable Full Body Slings

Full Body Sling, Polyester Mesh Fabric

Without Commode Opening

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMR110	Medium	30 x 47"	450-lbs.	100–195-lbs.	1/ea
MDSMR111	Large	34 x 52"	450-lbs.	150–300-lbs.	1/ea
MDSMR140	Bariatric	42 x 52"	700-lbs.	290–700-lbs.	1/ea



MDSMR114

With Commode Opening

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMR114	Medium	32 x 48"	450-lbs.	100–195-lbs.	1/ea
MDSMR115	Large	35 x 50"	450-lbs.	150–300-lbs.	1/ea
MDSMR116	X-Large	37 x 56"	450-lbs.	200–450-lbs.	1/ea
MDSMR141	Bariatric	42 x 50"	700-lbs.	290–700-lbs.	1/ea



MDSMR112

Full Body Sling, Padded Woven Flannel/Cordura Fabric

Item No.	Size	Dimensions	Wt. Cap.	Best Fit	Pkg.
MDSMR112	Medium	30 x 47"	450-lbs.	100–195-lbs.	1/ea
MDSMR113	Large	34 x 52"	450-lbs.	150–300-lbs.	1/ea
MDSMR117	X-Large	36 x 57"	450-lbs.	200–450-lbs.	1/ea

Stand Assist Patient Slings

Medline's Stand Assist slings provide support across the chest and under the arms for safe and comfortable lifting. All Medline Stand Assist slings are padded at the back and under the arms for maximum patient comfort.

Sling-Patient Use Stand-Assist Sling

Stand-Assist Sling, Polypropylene Fabric

Recommended use: Patient toileting and transfers.

Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSMDSA662	Medium	30–48"	500-lbs.	10/cs
MDSMDSA663	Large	36–60"	500-lbs.	10/cs



MDSMDSA663

Reusable Stand Assist Slings

Stand Assist Sling, Padded Knitted/Flannel Fabric

Recommended use: Patient toileting and transfers.

Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSMSA1	Small	28–42"	500-lbs.	1/ea
MDSMSA2	Medium	30–48"	500-lbs.	1/ea
MDSMSA3	Large	36–60"	500-lbs.	1/ea
MDSMSA4	X-Large	50–79"	600-lbs.	1/ea



MDSMSA1



MDSGRIPSLING3

Stand Assist Sling, Grip Neoprene Fabric

The Grip Sling adheres to clothing and is less likely to slide up under the arms and create pressure during a stand-assist transfer.

- Uniquely comfortable for patient satisfaction
- Specially textured neoprene material
- Soft padding in the underarm and lumbar areas

Recommended use: Patient toileting and transfers with stand-assist lifts for extra support for standing.

Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSGRIPSLING2	Medium	30 x 48"	500-lbs.	1/ea
MDSGRIPSLING3	Large	36 x 60"	500-lbs.	1/ea



MDSMBSB

Stand Assist Buttock Strap, Knitted Fabric

Recommended use: An accessory for extra support in addition to padded stand-assist slings.

USE WITH
STAND
ASSIST
SLINGS

Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSMBSS	Standard	34 x 7"	400-lbs.	1/ea
MDSMBSB	X-Large	58 x 11"	700-lbs.	1/ea

Specialty Patient Slings

Medline's specialty slings are designed to assist with unique patient lifting and clinical activities.

Single-Patient Use Specialty Sling

Limb Lift Sling, Polypropylene Fabric

Designed specifically for lifting heavy extremities. It supports and offers access to arms and legs for various care or therapeutic activities.



Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSLIMBLIFT	One Size Fits Most	40 x 6½"	450-lbs.	10/cs

Reusable Specialty Slings

Turning Patient Sling, Polyester Fabric

Medline's Turning Sling is designed for ease of placement under patients as well as its assistive capabilities to turn and position patients in bed with the use of a total patient lift.



Item No.	Style	Dimensions	Wt. Cap.	Pkg.
MDSTURN	Bariatric Turning Sling	48 x 24"	1,000-lbs.	1/ea



Toileting Sling, Padded Knitted/Flannel Fabric

Patient toileting activities. Support strap around chest and padding under arms offers patient comfort while the long leg straps fully support the patient's legs.

USE WITH
TOTAL
OR
STAND
ASSIST
LIFTS

Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSMT01	Small	26–42"	450-lbs.	1/ea
MDSMT02	Medium	30–48"	450-lbs.	1/ea
MDSMT03	Large	38–60"	450-lbs.	1/ea
MDSMT04	Bariatric	58–73"	700-lbs.	1/ea



Gait Training Sling, Polyester Fabric


To be used with total lifts, Medline's Gait Training Sling provides balance and support for standing, rising and gait training activities.

USE WITH
2-PT.
CRADLES

Item No.	Chest Circ.	Wt. Cap.	Pkg.
MDSMGAIT	40–54"	300-lbs.	1/ea

Sling Materials

Choose the right slings for the right task in the right fabric.




Polypropylene

Non-launderable polypropylene slings are designed for single-patient use. The right choice for facilities concerned about slings not returning from the laundry, it has the reliable strength for high weight capacities.



Hybrid Microfiber Nylon

Soft, breathable microfiber top and smooth nylon bottom. The skin-friendly, breathable material can remain in bed with the patient throughout their length of stay (single-patient use).




Knitted

This fabric has just the right amount of "give" for a secure, comfortable and easy fit. But don't let it fool you. This heavy-duty choice is designed to withstand high weight capacities and rigorous healthcare use.




Spacer

Warp-knitted polyester "spacer" fabric redistributes pressure and reduces shear to protect skin, while it is also breathable and comfortable. This material is reusable and can be left under patients.




Textured Grip

Strong and flexible, soft to the touch, neoprene adds advanced functionality to a standard sling.



Padded Flannel

Soft flannel provides a secure and easy fit as comfortable as your favorite shirt. Yet strong enough to support patients up to 450-lb.



Polyester Mesh

Ultra-strong and soft to the touch. The synthetic material retains its brilliant color, wash after wash. Skillful workmanship and reinforced stitching support high weight capacities.



Synthetic Polyester

Strong and durable, it supports high weight capacities and retains its color longer than traditional reusable fabrics.

Patient Alarms

Handle high-risk patients with care.

Support patients safely and securely with high-quality products at a great value.



MDT8510T

Tamper-Resistant Nurse Mode Alarms

- Alarm alerts staff when patient gets out of bed or chair
- Affordable monitor for use with most bed or chair sensor pads and floor mats
- Alarm and low battery lights communicate functionality, so you're always powered up
- Auto-sensing technology monitors patient movement and triggers the alarm

Item No.	Description	Pkg.
MDT8510T	Tamper-Resistant Nurse Mode Alarm with Single Port	1/ea
MDT8520T	Tamper-Resistant Nurse Mode Alarm with Dual Port	1/ea



MDT8450Y

Tamper-Resistant Personal Safety Alarm

- Auto-sensing technology monitors patient movement and sounds an alarm if patient gets up from a bed or wheelchair
- Prevent tampering and patient interference by pressing two buttons simultaneously to disengage the alarm
- Low battery indicator light and volume control
- Extra-durable wheelchair attachment clip included

Item No.	Description	Pkg.
MDT8450Y	Tamper-Resistant Alarm Monitor**	1/ea

Miscellaneous Alarm Accessories

Item No.	Description	Pkg.
MDT84YBRKT	Optional Patient Alarm Bracket for MDT84 Series	1/ea

Recordable Alarm Monitor

- Customizable with your desired message or alarm tone/sound
- When patient sits up in bed or stands up from a chair, the monitor will alert you and get the attention of patient with your message
- Protective boot for maximum protection
- Auto reset functionality so alarm is always ready to use

Item No.	Description	Pkg.
MDT800400	Recordable Alarm Monitor**	1/ea

*All alarms are universally compatible with pads listed on pages 50–51, excluding MDT85-Series alarms.

**Contact your Medline rep for cord options.



Quick Alert Nurse Call Alarms

- For use with bed, chair, seatbelt, toilet and floor sensor mats
- Distinguish between patients with 8 different tones and tunes
- Includes 3 volume levels with low battery indicator, yellow protective cover, 2 AAA batteries
- Basic alarm contains all the same features without delay option, tones, and tunes

Item No.	Description	Pkg.
MDT9100	Universal Alarm, Multi-tone**	1/ea
MDT9105	Universal Alarm**	1/ea



MDT9100



MDT85NCS



MDT85NCC



MDT851ADAPT

Crest Call Cords

Call cords specific for your custom call light system.

- Universal options in 7-ft. and 12-ft.
- Specialty Cords available upon request.
Please inquire with your Medline Representative.

Item No.	Description	Pkg.
CRE101968	Nurse Call Cord, 7-ft., 1/4" Both Ends	1/ea
CRE10196812	Nurse Call Cord, 12-ft., 1/4" Both Ends	1/ea
CRE113708	Pendant Only, Red, Clean Cord	15/pk
CRE113743	Cord, Pad Call, 1/4", Plug, 10-ft., Non-Locking	1/ea
CRE114192	Cord, Clean Cord, White, 100-ft.	1/ea
CRE114775	Cord, Pull, 6-ft. Glow, Dark White/Red	1/ea
CRE7700DB12	Cord, Call, 12-ft., Beige, 1/4" Both Ends, Phone Plug	1/ea
CRE7700DB12	Cord, Call, 12-ft. Double, Beige, 1/4" Both Ends, Phone Plug	1/ea
CRE7700DW10	Cord, Call, Double, 10-ft., White	1/ea
CREB112C	Call Station Cord, Bedside, Cornell, Double	1/ea
CREE1143C	Station, Emergency Pull Cord	1/ea

Universal Nurse Call Cables

Item No.	Description	Pkg.
MDT85NCC	Coiled Universal Nurse Call Cable	1/ea
MDTUNURSE	Coiled Universal Nurse Call Cable	1/ea
MDT85NCS	Straight Universal Nurse Call Cable	1/ea

Miscellaneous Alarm Accessories

Item No.	Description	Pkg.
MDT91NURSE	Quick Alert Nurse Call Cord, 1/4" Both Ends	1/ea
MDT9100AC	Quick-Alert Series Alarm Adapter, AC	1/ea
MDT91YELBT	Yellow Cover, Quick-Alert Series	1/ea
MDT81STRAP	Nylon Strap with Quick-Release Buckle	1/ea
MDTWSTRAP	Quick-Release Wipeable Alarm Strap	1/ea
MDT851ADAPT	Wall Mount Adapter for MDT85 Alarms	1/ea
MDTWMT	Wall Mount for MDT85 Alarms	1/ea

Reusable Bed/Chair Sensor Pads

- When attached to an alarm, the sensor pad alerts staff that the patient has taken weight off of the pad and may be at risk of falling
- Foldable, reusable and easy-to-store
- High-grade, water-resistant vinyl material provides protection against incontinence

Reusable Chair Pad

Item No.	Description	H x W	Pkg.
MDTRUC1	1-Year Warranty	12 x 15"	1/ea

Reusable Bed Pad

Item No.	Description	H x W	Pkg.
MDTRUB1	1-Year Warranty	15 x 30"	1/ea



MDTRUC1

Premium Alarm Sensor Pads

- Weight-sensitive pad activates alarm if patient leaves a chair or bed
- Water-resistant, high-grade vinyl
- Use on a chair or bed (underneath the fitted sheet)

Chair

Item No.	Description	H x W	Pkg.
MDT85130C6*	Chair Sensor Pads for MDT85 Alarms	13 x 13"	6/cs
MDT85130C1*	Chair Sensor Pad for MDT85 series alarm	12 x 15"	1/ea
MDT9130C5	Chair Sensor Pads for MDT91 Alarms	12 x 15"	5/cs
MDT9190C5	Chair Sensor Pads for MDT91 Alarms	12 x 15"	5/cs
MDT91180C5	Chair Sensor Pads for MDT91 Alarms	12 x 15"	5/cs



MDT85130C6

Bed

Item No.	Description	H x W	Pkg.
MDT85130B6*	Bed Sensor Pads for MDT85 Alarms	13 x 30"	6/cs
MDT9130B5	Bed Sensor Pads for MDT91 Alarms	15 x 30"	5/cs
MDT9190B5	Bed Sensor Pads for MDT91 Alarms	15 x 30"	5/cs
MDT91180B5	Bed Sensor Pads for MDT91 Alarms	15 x 30"	5/cs



MDT85130B6

*Pads are only compatible with MDT85 alarm series.



MDT821180C6



MDT82145B6

Economy/Single Patient Use Sensor Pads

- When pad is connected to the patient alarm, it will trigger the alarm when weight is removed from the pad
- Durable, high-grade vinyl pad is easy to store
- For single patient use on a chair or bed (underneath the fitted sheet)
- MDT82 series includes a full-size, non-slip grip sheet for maximum stabilization

Chair Sensor Pad

Item No.	Description	H x W	Pkg.
MDT83130C6	30-Day Warranty	7 x 15"	6/cs
MDT82145C6	45-Day Warranty	7 x 15"	6/cs
MDT831180C6	180-Day Warranty	7 x 15"	6/cs
MDT821180C6	180-Day Warranty	7 x 15"	6/cs
MDT831365C6	1-Year Warranty	7 x 15"	6/cs
MDT821365C6	1-Year Warranty	7 x 15"	6/cs

Bed Sensor Pad

Item No.	Description	H x W	Pkg.
MDT83130B6	30-Day Warranty	10 x 30"	6/cs
MDT82145B6	45-Day Warranty	10 x 30"	6/cs
MDT831180B6	180-Day Warranty	10 x 30"	6/cs
MDT821180B6	180-Day Warranty	10 x 30"	6/cs
MDT831365B6	1-Year Warranty	10 x 30"	6/cs
MDT821365B6	1-Year Warranty	10 x 30"	6/cs

Advantage® Magnetic Patient Alarm

- Control patient wandering with an alarm that attaches to a bed, chair or patient clothing
- Magnetic pull switch releases from any angle, helping to eliminate false alarms and clothing tears
- Recessed on/off switch limits patient interference
- Battery tester button ensures that alarm is always powered up
- 2 AAA batteries included

Item No.	Description	Pkg.
MDT5000	Alarm, Advantage, Magnetic Tether	5/cs
MDT5000Z	Alarm, Advantage, Magnetic Tether	1/ea
MDT5000XCD	Advantage Magnetic Alarm Replacement Pull Cords	3/bg



MDT5000Z


Anti-Wander Alarm Kits and Accessories

- A great solution for patients with Alzheimers or dementia who may be prone to wandering
- Kit includes three easy-to-use components: a door strip monitor, magnetic contact switch, and resident wristband
- Easy installation with no hard wiring

Item No.	Description	Pkg.
MDT86200	Anti-Wander Kit, One Door	1/ea
MDT86200H	Anti-Wander Kit, Home Setting	1/ea
MDT86200TR	Resident Wristband	1/ea
MDT86100W	Reset Button	1/ea
MDT862CMU	Central Monitoring Unit	1/ea



MDT862CMU



MDT86200



MDTIRM1

Infrared Patient Bed Alarm*

- Alarm affixes to the headboard or bed frame, so patient can move freely within bed
- If patient attempts to leave bed, the infrared beam detects motion and sends out your choice of alarm or chime tone
- Alarm head swivels for precise beam alignment
- Attach using the bracket (optional) or hook and loop strips (included)
- Requires 9V battery (not included)

Alarm

Item No.	Description	Pkg.
MDTIRM1	Infrared Alarm	1/ea

Metal Mounting Bracket for Infrared Alarm

Item No.	Description	Pkg.
MDTIRM1BRKT	Metal Mounting Bracket	1/ea

*This product does not qualify towards Uplift Program spend.

Smart Caregiver® Wireless Monitoring Systems

- Wireless monitor helps reduce tripping hazards and eliminates broken or tangled pad cords
- Cordless pads communicate with the monitor using wireless technology with both an audible and visual alert
- Mount the monitor 15–20-ft. away from the sensor pads to eliminate bedside alarm noise in the resident's room

Monitors

Item No.	Description	Pkg.
MDT8800	Wireless Monitor (Up to 2 Pads)	1/ea
MDT433CMU	Central Monitoring Unit (Up to 40 Pads)	1/ea



MDT8800

Wireless Alarm Sensor Pads and Accessories

Item No.	Description	Pkg.
MDT8830C1	Cordless Chair Alarm Pad for MDT8800	1/ea
MDT8830B1	Cordless Bed Alarm Pad for MDT8800	1/ea
MDT88365B1	Bed Sensor Mat, 1-year warranty	1/ea
MDT88365C1	Chair Sensor Mat, 1-year warranty	1/ea
SMRTIR1	Infrared Patient Bed Alarm for MDT8800	1/ea
MDT8836XWB1	Bed Sensor Mat, 1-Year Warranty, X-Wide	1/ea
MDT8800FM	Cordless Floor Mat	1/ea
MDT433NC	Nurse Call Button for MDT433CMU	1/ea
MDT433RB	Reset Call Button for MDT433CMU	1/ea
MDT433PGD	Pager with LCD for MDT433CMU	1/ea
MDT433PAGER	Pager without LCD for MDT433CMU	1/ea

Sensable® Care Proactive Patient Monitoring System

The only clinically proven technology available that can reduce bed falls and alarm fatigue.

A. Sensor Pad

Lightweight, thin, reusable pad is inserted into the fitted sheet to monitor patient movement from the shoulders to the thighs, detect potential bed exits, and reduce bed falls. It also issues repositioning alerts.

B. Alarm Monitor

When there's an unscheduled bed exit attempt, this bedside device flashes red, issues an audible message to remind the patient to wait for caregiver, and alerts staff.

C-D. Dashboard and Computer Server

Placed at the nursing station, the nursing dashboard receives real time data from the alarm monitor and can simultaneously monitor up to 35 beds in real time. The computer server collects data from sensors and devices to generate unit-specific reports.

E. Customizable Fitted Sheet

Poly sheet accommodates the sensor pad and is customized to your bed size.

Free Mobile Device App

Nurses receive potential bed exit alerts and notifications to reposition patients via the free Sensable Care App available for Apple iOS® and Android™ devices.



Sensable Care Monitoring System

- Enhance patient safety and security with Sensable Care, a state-of-the-art monitoring system that detects early patient movement patterns and alerts staff prior to a potential bed exit

Item No.	Description	Pkg.
MDT9590B1	Sensor Pad	1/ea
MDT9500	Alarm Monitor	1/ea
MDT95DASH	Dashboard	1/ea
MDT95SERVER	Computer Server	1/ea
MDT9590SHEET	Fitted Sheet	1/ea



Gait Belts and Safety Essentials



MDT821203YS



MDT821203

Cotton Gait Belts

- Belt fits securely around the patient's waist and provides caregivers with a surface to grab when helping patients ambulate or transfer
- Aid with transferring patients
- Additional support helps reduce the risk of caregiver back injuries
- Available in pediatric, standard and bariatric sizes

Item No.	Design	Length	Buckle	Pkg.
MDT821203YS	Solid Yellow Fall Management	54"	Metal	1/ea
MDT821203Y	Solid Yellow Fall Management	60"	Metal	1/ea
MDT821203YL	Solid Yellow Fall Management	72"	Metal	1/ea
MDT821203	White with Blue and Red Pin Stripes	60"	Metal	1/ea
MDT828203	White with Blue and Red Pin Stripes	60"	Metal	6/cs
MDT821203L	White with Blue and Red Pin Stripes	72"	Metal	1/ea
MDT821204	Pastel Stripes	60"	Metal	1/ea
MDT828204C	Pastel Stripes	60"	Metal	6/cs
MDT828208PASL	Pastel	72"	Metal	1/ea
MDT821203BS	Solid Black	54"	Metal	1/ea
MDT821203B	Solid Black	60"	Metal	1/ea
MDT821203BL	Solid Black	72"	Metal	1/ea
MDT828203NAT	Natural	60"	Metal	1/ea
MDT82803BLU	Navy Blue	60"	Metal	1/ea
MDT821208NATM	Natural	60"	Plastic	1/ea
MDT821208PAS	Pastel Stripes	60"	Plastic	1/ea
MDT828209PASL	Pastel	72"	Plastic	1/ea

Wipeable Gait Belts

- Versatile belts provide help with mobilizing and transferring patients
- Straps are durable and easy to keep clean
- Save laundry time—wipe clean with standard PDI wipes/hospital-grade disinfectant
- Buckle provides easy adjustment

Polypropylene—Flexible, Webbed Material

Item No.	Design	Length	Buckle	Pkg.
MDT821206BLU	Solid Blue	60"	Plastic	1/ea
MDT821206RAIN	Rainbow	60"	Plastic	1/ea
DII27051BL	Solid Blue	60"	Plastic	1/ea
DII27051YL	Solid Yellow	60"	Plastic	1/ea
DII27051BK	Solid Black	60"	Plastic	1/ea
DII24051BL	Solid Blue	60"	Metal	1/ea
DII24061BL	Solid Blue	70"	Metal	1/ea
DII27061YL	Solid Yellow	72"	Plastic	1/ea
DII27061BK	Solid Black	72"	Metal	1/ea



MDT821207BLUL

Polyurethane—Smooth, Durable Material

Item No.	Design	Length	Buckle	
MDT828203PWP	Solid Pink	60"	Plastic	1/ea
MDT828203BWP	Solid Black	60"	Plastic	1/ea
MDT821203WP	Solid Blue	60"	Plastic	1/ea
MDT828203YWP	Solid Yellow	60"	Plastic	1/ea
MDT821203WPL	Solid Blue	70"	Plastic	1/ea
MDT828203OWPL	Solid Orange	70"	Plastic	1/ea



MDT828203BWP

MDT828203PWP



MDT823256

MDT823266

Hand Protectors

- Provide a way for patients to exercise fingers during rehabilitation and to prevent scratching or picking and limit contracture in stroke patients
- Secure to bed frame or wheelchair with included arm restraint
- Machine washable polyester and cotton fabric is filled with resilient polyester pile
- One size, for left or right hands
- Safety check design has top flap for checking patients' hands
- Finger-control shoelace style mitt separates fingers and features a breathable mesh backing

Adult Hand Protectors

Item No.	Description	Straps	Closure	Pkg.
MDT823256	Flexible Palm	Double	Metal Buckle with Teeth Slider	1/pr
MDT823266	Rigid Palm	Double	Metal Buckle with Teeth Slider	1/pr
MDT823243	Safety Check	Single	Hook-and-Loop Material	1/pr
MDT823242	Finger-Control	Single	Shoelaces	1/pr



MDT823242

Toddler and Youth Mitts

Item No.	Description	L x W	Wrist Circumference	Pkg.
PSY2911M	Toddler	4 x 1¼"	5"	1/pr
PSY2911L	Youth	6½ x 1¼"	5½"	1/pr
PSY2911S	Pediatric	3½ x 1¼"	4¾"	1/pr



MDT823243

Multi-Purpose Arm and Leg Splints

- Keep patient from pulling at tubes, drains and lines with an arm splint
- Easily applied with hook and loop closure
- Vinyl mesh outer covering has a 100% cotton lining
- Zippered opening enables beads to be added or removed for patient comfort

Item No.	Style	Limb Length	Limb Circumference	Pkg.
MDT828816P	Pediatric	7"	1½–4"	1/ea
MDT828816XS	Pediatric 2–7 Years	10"	1½–4"	1/ea
MDT828816S	Small	1"	3–14"	1/ea
MDT828816M	Medium	13"	5–17"	1/ea
MDT828816L	Large	14"	8–21"	1/ea
PSY8169	Adult Leg	18"	14–27"	1/ea



MDT828816M

Activity Overlay

- Provides motor/sensory stimulation for patients with diminished cognitive function
- Bright colors and texture provide visual and tactile stimulation
- Attaches to lap or chair with sticky tape
- Machine washable

Item No.	Description	Pkg.
MDTOVERLAY	Activity Overlay	1/ea



Activity Apron

- Bright colors and texture variation provide visual and tactile stimulation for patients with decreased cognitive function
- Includes clear picture frame for photos
- Secures like an apron around waist
- Machine washable

Item No.	Description	Pkg.
MDTACTIVITY	Activity Apron	1/ea





MDT823100

Quick Release Limb Holders

- Secure wrists or ankles to bed frame with hook and loop closure
- Will not overtighten with patient movement
- Non-detachable double straps offer extra range of motion restriction
- 7 x 2½" quilted cuff with double straps for range of motion

Item No.	Description	Cuff L x W	Pkg
MDT823100	Quick-Release Ties	7 x 2½"	1/pr
MDT821301Q	Black, Quick Release	7 x 3"	1/pr
MDT821301QL	Black, Long Strap, Quick Release	7 x 4"	1/pr
MDT821302QL	Black, 60" Strap, Quick-Release Veriflex Buckle	7 x 4"	1/pr



MDT823040

Quick Control Limb Holders

- Help keep patients safe who may be a danger to themselves and staff
- 36" non-detachable straps offer additional range of motion restriction
- 9 x 3" size features a slide harness strap and synthetic fur cuff for patient comfort
- Machine washable

Item No.	Description	Cuff L x W	Pkg	Pkg
MDT823040	Orthopedic Felt, Cotton	Metal Buckle and Hook-and-Loop	10 x 3"	1/pr
MDT823070*	Synthetic Fur, Polyester	Metal Buckle with Teeth Slider	9 x 3"	1/pr

*For a case of 6 pairs, order item MDT829070C.



Disposable Limb Holder

- Single-patient use for economy and efficiency
- Secure wrists or ankles to a bed frame
- Metal buckle with teeth slider cinch closure
- 10 x 3" foam lined cuff with double 30" long non-stretch, non-detachable anchor straps

Item No.	Description	Cuff L x W	Pkg
MDT829030	Metal Buckle with Teeth Slider	10 x 3"	6 pr/cs



Reusable Padded Limb Holder

- Secure wrists or ankles of at-risk patients to a bed frame
- Will not tighten with patient movement
- Slide harness minimizes chafing during wrist rotation
- 10 x 4" fur lined cuff with 36" long single strap

Item No.	Description	Cuff L x W	Pkg
MDT829082C	Metal Buckle with Teeth Slider	10 x 4"	6 pr/cs

Foam Limb Holders

- Prevent patients from pulling at medical devices, interrupting treatments and disrupting staff members
- Padded foam cuff secures wrists or ankles to a bed frame
- Will not tighten with patient movement
- Adjustable $\frac{3}{8}$ " D-ring closure
- $13\frac{1}{2} \times 3$ " foam lined cuff with 55" long double straps

Adult

Item No.	Description	Cuff L x W	Pkg
MDT829520	Double Strap with Metal D-ring Closure	$13\frac{1}{2} \times 3$ "	1/pr
MDT829528	Single Strap with Metal D-ring Closure	$13\frac{1}{2} \times 3$ "	1/pr



MDT829520



MDT829473

Infant

Item No.	Description	Cuff L x W	Pkg
MDT829473	Hook-and-Loop	$5\frac{1}{2} \times 1$ "	24/bx

Disposable Personal Safety Limb Holder

- Thick, breathable foam patient safety devices feature adjustable hook-and-loop closures to ensure a comfortable fit
- Limits movement of wrists or ankles
- Nylon D-ring provides added security
- Anchor straps will not stretch
- $8\frac{3}{4} \times 3\frac{3}{4}$ " foam lined cuff with 45" long double strap

Item No.	Description	Cuff L x W	Pkg
NON243074	Limb Holder	$8\frac{3}{4} \times 3\frac{3}{4}$ "	48/cs
NON243074HH	Limb Holder	$8\frac{3}{4} \times 3\frac{3}{4}$ "	1/ea



NON243074



Personal Safety Limb Holder

- Cotton personal safety limb holders prevent patients from pulling out an IV line or catheter
- For use on wrists or ankles, they won't over tighten with movement
- Synthetic fur lined cotton cuff with adjustable hook-and-loop fastener and straps that tie to bed frame
- 10 x 4" fur lined cuff with 44" long double straps

Item No.	Description	Cuff L x W	Pkg
MDT829052	Synthetic Fur, Cotton	10 x 4"	6 pr/cs



Cinch-Type Patient Safety Limb Holder

- Foam limb holder limits movement of wrists or ankles
- Hook-and-loop tab adheres anywhere on the foam backing
- D-ring stress point prevents cuff from tightening when patient pulls on it
- Adjusts for adult or pediatric use
- 10 x 3" foam lined cuff with 40" long strap

Item No.	Description	Cuff L x W	Pkg
MDT829020C	Cinch-Type with D-Ring Closure	10 x 3"	6 pr/cs

Wrist and Ankle Limb Holders

- For patients who are a threat to themselves or others in locking and non-locking styles
- Neoprene cuffs are gentle on skin and adjust to fit most limb sizes
- Two universal keys included with locking cuffs
- Hook and loop closure for easy adjustment
- Machine washable

Keylock Limb Holders

Item No.	Description	Cuff L x W	Strap Length	Pkg.
MDT829098W	Wrist	12 x 2½"	46"	1/pr
MDT829099A	Ankle	14½ x 2½"	44½"	1/pr

Hook & Loop Ring Limb Holders

Item No.	Description	Cuff L x W	Strap Length	Pkg.
MDT829090W	Wrist	12 x 2½"	48½"	1/pr
MDT829091A	Ankle	14½ x 2½"	50"	1/pr

Quick-Release Limb Holders

Item No.	Description	Cuff L x W	Strap Length	Pkg.
MDT829090WQ	Wrist	12 x 2½"	45"	1/pr
MDT829091AQ	Ankle	12 x 2½"	46½"	1/pr

Stretcher Limb Holders

Item No.	Description	Cuff L x W	Strap Length	Pkg.
MDT8292750	Wrist	12 x 2½"	9 x 2"	1/pr
MDT8292755	Ankle	14½ x 2½"	9 x 2"	1/pr



MDT829098W



MDT829090W



MDT829090WQ



MDT8292755

Repositioning Belts

- Ideal for patients with minimal upper trunk control who require positioning support and are at risk of falling, or those who require a restraint for medical treatment
- For use on a bed, stretcher or chair, wheelchair
- For combative patients, utilize keylock restraints for maximum security
- Flame-retardant cotton wraparound belt has a secure hook-and-loop closure
- Roll belts are less restrictive in that they allow patients to sit up or roll from side to side in bed



MDT3658ML



MDT821231



MDT822125

Locking Belts with Keylock

(For key replacements use item no.: MDT829074)

Item No.	Description	Pkg.
MDT829334C	Safety Belt with Keylock, Cotton	1/ea
MDT829337C	Three-Piece Swedish Belt	1/ea

Wrap-around Belts with Hook and-Loop Closure

Item No.	Description	Strap Length	Pkg.
MDT3658ML	Medium/Large, Quilted Fabric	77"	1/ea
MDT3658XL	X-Large, Quilted Fabric	96"	1/ea
MDT36582XL	X-Large, Quilted Fabric	121"	1/ea
MDT36592XL	2X-Large, Foam Fabric	121"	1/ea

Roll Belts with Tie Closure

Item No.	Description	Pkg.
MDT821231	Single Strap	1/ea
MDT821135	Roll Belt, Double Strap	1/ea

Roll Belts with Quick Release Buckles

Item No.	Description	Pkg.
MDT821231Q	Single Strap	1/ea
MDT821135Q	Dual Strap	1/ea

Safety/Soft Belts

Item No.	Description	Pkg.
MDT8214125	Soft Belt with Ties	1/ea
MDT822125	Poly, 1-pc.	1/ea
MDT822126	Poly, 2-pc.	1/ea
MDT822127	Quick-Release, 2-pc.	1/ea

Breezeline Pelvic Holders

- Mesh holder helps prevent forward sliding in wheelchairs or geri-chairs
- Straps have quick-release buckles
- Machine washable

Item No.	Description	L x W	Upper Waist Strap	Lower Waist Strap	Pkg.
PSY4427S	Small	16½ x 17"	52½"	50"	1/ea
PSY4427M	Medium	19 x 17"	56½"	65"	1/ea
PSY4427L	Large	23 x 17"	60½"	56"	1/ea



PSY4427L

Quick Release Blended Safety Vests and Jackets

- Keep patient securely in bed or in an upright position in a wheelchair
- Vests are available with quick-release straps or ties and a zip back
- Long sleeve jacket has an off-center zipper back to avoid spinal irritation
- Machine washable cotton/polyester

Vests—with Quick-Release Straps, Zip Back

Item No.	Size	Pkg.
MDT828306M	Medium, 110–155 lbs.	1/ea
MDT828306L	Large, 145–190 lbs.	1/ea



MDT828306L

Vests—with Tie Straps, Zip Back

Item No.	Size	Pkg.
MDT828305S	Small, 80–120 lbs.	1/ea
MDT828305M	Medium, 110–155 lbs.	1/ea
MDT828305L	Large, 145–190 lbs.	1/ea
MDT828305XL	X-Large, 180–230 lbs.	1/ea



MDT8283360M

Jackets

Item No.	Size	Pkg.
MDT8283360S	Small, 80–120 lbs.	1/ea
MDT8283360M	Medium, 110–155 lbs.	1/ea
MDT8283360L	Large, 145–190 lbs.	1/ea



MDT828002

Koolnit Safety Vests

- Vests prevent patient from sitting up or sliding down in bed or wheelchair
- Ties in back out of patient's reach with connecting straps that attach to bed frame or wheelchair tilt bars
- Koolnit polyester mesh is comfortable and washable

Item No.	Size	Pkg.
MDT828000	Small, 80–120 lbs.	6/cs
MDT828002	Medium, 110–155 lbs.	6/cs
MDT828004	Large, 145–190 lbs.	6/cs

Pulse Oximeter Probe Wraps

- Help secure optical transducers and probe sensors without adhesive tape
- Hold optical components firmly to help prevent misalignment
- Universal design may be trimmed to fit neonatal foot, infant toe or pediatric finger



Item No.	Size	Pkg.
MDT829554	5¼ x 1¼"	12/bx



Infant ID Bracelets

- Foam ID Bracelet for infant identification
- Panel is blank for writing or printed label placement
- Easily repositioned to accommodate IV Site changes

Item No.	Size	Pkg.
MDT829464P	5¼ x 1¼"	12/dz



Circumstraint Board Straps

Infant

Item No.	Description	Cuff L x W	Pkg
MDT829474P	Hook-and-Loop	10 x 2"	24/bx

*The products in this section do not count as qualifying spend towards the Uplift Program.

Identification and Protective Products*

Best practices lead to better outcomes.

Make it easy to do the right thing. Arm your team with the tools they need to quickly identify at-risk patients.

Fall Prevention Slippers*

- Bright colored slippers help identify at-risk patients
- Double-sided, non-slip tread help patients stay on their feet
- Soft, breathable terrycloth interior keeps feet comfortable, dry and warm

Yellow—Moderate Risk

Item No.	Size	Tread Color	Pkg.
MDTPS4B01FPY	Child	White	48 pr/cs
MDTPS4B01FPYH	Child	White	1/ea
MDTPS4B03FPY	Medium	White	48 pr/cs
MDTPS4B03FPYH	Medium	White	1/ea
MDTPS4B04FPY	Large	White	48 pr/cs
MDTPS4B04FPYH	Large	White	1/ea
MDTPS4B05FPY	X-Large	White	48 pr/cs
MDTPS4B05FPYH	X-Large	White	1/ea
MDTPS4B06FPY	2X-Large	White	48 pr/cs
MDTPS4B06FPYH	2X-Large	White	1/ea
MDTPS4B07FPY	XXX-Large	White	48 pr/cs
MDTPS4B07FPYH	XXX-Large	White	1/ea

Red—High Risk

Item No.	Size	Tread Color	Pkg.
MDT211218R	One size fits most	White	48 pr/cs
MDT211218RH	One size fits most	White	1/ea
MDTPS4B01FPR	Child	Gray	48 pr/cs
MDTPS4B03FPR	Medium	Gray	48 pr/cs
MDTPS4B03FPRH	Medium	Gray	1/ea
MDTPS4B04FPR	Large	Gray	48 pr/cs
MDTPS4B05FPR	X-Large	Gray	48 pr/cs
MDTPS4B06FPR	XX-Large	Gray	48 pr/cs
MDTPS4B07FPR	XXX-Large	Gray	48 pr/cs

Identify The Risk

- Yellow—Moderate Risk
- Red—High Risk

Double-sided,
non-slip tread
helps keep patients
on their feet.



PerforMAX™ Fall Prevention IV Gowns

- Quickly identify high-risk patients with this bright yellow gown
- Waist ties are sewn into side seam and bar tacked for extra strength
- Side tie eliminates the need for double gowning
- Durable PerforMAX polyester fabric has superior color retention wash after wash

Item No.	Size	Pkg.
➤ MDTPG5ITSHRG	One Size	6 dz/cs
➤ MDTPG7ITSHRG	One Size	6 dz/cs
➤ MDTOG7ITSHRG	3X-Large	1 dz/cs
➤ MDTBG7ITSHRG	10X-Large	1 dz/cs

Identify The Risk

- Yellow—Moderate Risk
- Red—High Risk



MDTPG7ITSHRG



Falls Prevention Kits

- Choose from kits with yellow components to identify patients at moderate risk for falls or kits with red components to identify patients at high risk
- Each kit contains a color-coded 3 x 3" magnet to place on the patient's door frame, one pair of color-coded non-slip slipper socks and a falls risk ID band that easily secures around the wrist

Item No.	Description	Pkg.
➤ DYK100MFALL	Moderate Falls Risk Prevention Kit	60/cs
➤ DYK100HFALL	High Falls Risk Prevention Kit	60/cs



DYK100MFALL



DYK100HFALL

Custom Kit Component Checklist

Most Common Components

- ☐ Yellow Fall Precaution Label
- ☐ Yellow ID Band
- ☐ Yellow XL Slippers
- ☐ Yellow 3 x 3" Magnet
- ☐ Yellow Blanket
- ☐ Yellow ID Band
- ☐ Red 4 x 3" High Fall Risk Label
- ☐ Red XL Slippers
- ☐ Red 3 x 3" Magnet
- ☐ Window Cling

Additional Items Prevention

- ☐ Sensor Pad
- ☐ Gait Belt
- ☐ Urinal
- ☐ Education Materials for LTC
- ☐ Signage

Packaging Options

- ☐ KareBag*
- ☐ Zip-Lock Bag
- ☐ Drawstring Bag
- ☐ Perforated Bag



*Customization available. Lead time for custom printing on bags is 120 days.

Fall Mats

- Our most economical fall mats help reduce injury due to a fall from bed
- Absorbent, high-compression foam in bi-fold or non-folding style
- Durable, easy-to-clean vinyl cover

Bi-Fold

Item No.	Description	Pkg.
MDTFM24662B1	24 W x 66 D x 2" H	1/ea
MDTFM24663B1	24 W x 66 D x 3" H	1/ea
MDTFM36662B1	36 W x 66 D x 2" H	1/ea
MDTFM48662B1	48 W x 66 D x 2" H	1/ea

Non-Folding

Item No.	Description	Pkg.
MDTFM24662	24 W x 66 D x 2" H	1/ea
MDTFM24663	24 W x 66 D x 3" H	1/ea
MDTFM36662	36 W x 66 D x 2" H	1/ea
MDTFM48662	48 W x 66 D x 2" H	1/ea

Beveled Edge Fall Mats

- Shock-absorbent mats help protect patients from fall-related injuries
- Low profile, beveled edges help reduce trip hazard
- Non-skid material helps prevent mat from sliding on floor

Item No.	Description	Pkg.
MDTFALLMATB	Brown, 23 W x 68 D x 1" H	1/ea
MDTFALLMATG	Gray, 23 W x 68 D x 1" H	1/ea



MDTFALLMATB



MDTFALLMATG

Hip Protector

- Reduce the risk of hip fractures associated with a fall
- Breathable, comfortable and form-fitting
- Durable and machine washable
- Holds incontinence liners

Item No.	Description	Pkg.
HIPPROTSM	Small, 30–38" (76–97 cm) Waist	10/cs
HIPPROTMD	Medium, 36–44" (91–112 cm) Waist	10/cs
HIPPROTLG	Large, 40–48" (102–122 cm) Waist	10/cs
HIPPROTXL	X-Large, 44–56" (112–142 cm) Waist	10/cs
HIPPROTXXL	XX-Large, 48–60" (122–152 cm) Waist	10/cs
HIPPROTSMZ	Small, 30–38" (76–97 cm) Waist	1/ea
HIPPROTMDZ	Medium, 36–44" (91–112 cm) Waist	1/ea
HIPPROTLGZ	Large, 40–48" (102–122 cm) Waist	1/ea
HIPPROTXLZ	X-Large, 44–56" (112–142 cm) Waist	1/ea
HIPPROTXXLZ	XX-Large, 48–60" (122–152 cm) Waist	1/ea



HIPPROTM

Open Hip Protector

- Reduce the risk of hip fractures associated with a fall
- Breathable, comfortable and form-fitting
- Durable and machine washable
- Open design—perfect for frequent restroom users

Item No.	Description	Pkg.
HIPPROTOPSM	Small, 30–38" (76–97 cm) Waist	10/cs
HIPPROTOPMD	Medium, 36–44" (91–112 cm) Waist	10/cs
HIPPROTOPLG	Large, 40–48" (102–122 cm) Waist	10/cs
HIPPROTOPXL	X-Large, 44–56" (112–142 cm) Waist	10/cs
HIPPROTOPXXL	XX-Large, 48–60" (122–152 cm) Waist	10/cs
HIPPROTOPSMZ	Small, 30–38" (76–97 cm) Waist	1/ea
HIPPROTOPMDZ	Medium, 36–44" (91–112 cm) Waist	1/ea
HIPPROTOPLGZ	Large, 40–48" (102–122 cm) Waist	1/ea
HIPPROTOPXLZ	X-Large, 44–56" (112–142 cm) Waist	1/ea
HIPPROTOPXXLZ	XX-Large, 48–60" (122–152 cm) Waist	1/ea



HIPPROTOPSM

Hipsters Men's Standard Briefs

- Body-contouring, removable hip pads within the brief help provide protection throughout the critical hip fracture area
- Hip pads are impact-absorbing, perforated EVA foam that allows skin to breathe, while reducing perspiration and resulting skin breakdown
- Slip on easily over undergarments or wear as underwear
- Virtually undetectable under clothing
- Comfortable and washable polyester and cotton blend

Item No.	Size	Waist	Hips	Pkg.
PSY6031S	Small	28–30"	35–37"	1/ea
PSY6031M	Medium	30–34"	37–41"	1/ea
PSY6031L	Large	34–38"	41–45"	1/ea



PSY6031M



PSY6017SEA

Hipsters Men's Incontinent Briefs

- High impact-absorbing pads help protect hip bones against injury from falls
- Discreet, low-profile standard foam pads are virtually undetectable under clothing to help improve compliance versus hard-shell products
- Comfortable and washable polyester and cotton blend

Item No.	Size	Waist	Hips	Pkg.
PSY6017SEA	Small	28–30"	35–37"	1/ea
PSY6017MEA	Medium	30–34"	37–41"	1/ea
PSY6017LEA	Large	34–38"	41–45"	1/ea



PSY6009NM

Hipsters Unisex Sweatpants

- Removable hip pads provide protection throughout the critical fracture area to help reduce the risk of injury from a fall
- Long pants provide additional comfort and warmth
- Neutral navy color for women and men
- Comfortable and washable polyester and cotton blend

Item No.	Size	Waist	Hips	Pkg.
PSY6009NM	Medium	28–30"	35–37"	1/ea
PSY6009NL	Large	30–34"	37–41"	1/ea



PSY6040M

Unisex Protective Caps

- Absorbing, built-in soft foam pads help reduce the force of impact to the head from a fall
- Discreet, low profile baseball style cap
- Fleece interior offers comfort and reduced heat buildup
- Comfortable and washable lightweight cotton blend

Item No.	Size	Circumference	Pkg.
PSY6040M	Medium	22½"	1/ea
PSY6040L	Large	24"	1/ea

Knockdown Transfer Bench with Microban®*

- Microban antimicrobial* product protection is built-in to this transfer bench to help prevent the growth of stains and odors caused by bacteria, mold and mildew. Always clean this product thoroughly after each use.
- Suction cup foot piece secures bench to tub floor
- Lightweight, easy to assemble
- Side arm provides additional stability

* These antimicrobial properties are built-in to the transfer bench. The transfer bench does not protect users or others against bacteria, viruses, germs or other disease organisms.



MDS86960KDMB

Item No.	Description	Wt. Cap.	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
MDS86960KDMB	Transfer Bench	350-lbs.	19½–23½"	16"	N/A	26"	2/cs
MDS86960KDMBH	Transfer Bench	350-lbs.	19½–23½"	16"	N/A	26"	1/ea

Knockdown Bath Bench with Microban®*

- Microban antimicrobial* product protection is built-in to this bath bench to help prevent the growth of stains and odors caused by bacteria, mold and mildew. Always clean this product thoroughly after each use.
- Bath benches help patients who have difficulty sitting in a standard bath tub or standing in the shower
- Benches with a back give extra support and comfort; Benches without a back provide easy access and maneuverability when bathing
- Non-slip, rubber feet resist slipping on the bathtub or shower floor

* These antimicrobial properties are built-in to the bath bench. The bath bench does not protect users or others against bacteria, viruses, germs or other disease organisms.



MDS89745KDMB

Item No.	Description	Wt. Cap.	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
MDS89745KDMB	With Back	300-lbs.	13¾–20¾"	12¼"	N/A	17"	2/cs
MDS89745KDMBH	With Back	300-lbs.	13¾–20¾"	12¼"	N/A	17"	1/ea
MDS89740KDMB	No Back	300-lbs.	13¾–20¾"	12"	N/A	17"	2/cs
MDS89740KDMBH	No Back	300-lbs.	13¾–20¾"	12"	N/A	17"	1/ea

Commode with Microban®*

- Microban antimicrobial* product protection is built-in to this commode to help prevent the growth of stains and odors caused by bacteria, mold and mildew. Always clean this product thoroughly after each use.
- Tool-free assembly
- Powder-coated for easy cleaning
- Splash guard included for over-the-toilet use
- Features a convenient toilet paper holder
- Hospital quality construction

* These antimicrobial properties are built-in to the commode. The commode does not protect users or others against bacteria, viruses, germs or other disease organisms.



Item No.	Description	Wt. Cap.	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
MDS89664KDMB	Steel	350-lbs.	15½–21¼"	19"	N/A	25½"	1/cs



G30213-4F

Folding 3-in-1 Commodes

- Folding frame design easily opens and folds, and will fold flat for convenient storage and transportation
- Clip-on seat can be removed for easy cleaning
- Includes pail, lid and splashguard
- Elongated seat commode has a deeper seat to allow more room for cleaning

Item No.	Description	Wt. Cap.	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
G30213-4F	Basic, Steel	350-lbs.	16½–22½"	16"	19"	21¼"	4/cs
G30210-4F	Premium, Aluminum	250-lbs.	17½–21½"	16"	19"	21¼"	4/cs
MDS89664EL	Elongated Seat, Steel	400-lbs.	17–23"	16"	19"	22"	4/cs



MDS89664

Non-Folding Commodes

- Use in the bathroom as an elevated toilet seat or as toilet safety rails when combined with most standard toilets
- Includes pail, lid and splashguard
- Clip-on seat can be removed easily for cleaning
- Gray powder-coated steel commode is rust-resistant and durable
- Bariatric commodes are extra-wide for larger patients

Item No.	Description	Wt. Cap.	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
MDS89664	Basic, Steel	350-lbs.	15–21"	15"	19"	21¼"	4/cs
G30216B	Bariatric, Steel	500-lbs.	16½–23½"	18½"	23½"	27"	2/cs
MDS89664XW	Extra-Wide Bariatric, Steel	650-lbs.	18–23½"	15"	24"	28"	1/ea



G98202

Drop-Arm Commodes

- Arms swing out of the way for safe lateral transfers on and off the commode
- Durable, rust-resistant, powder coated reinforced frame
- Includes pail, lids and splash guard

Item No.	Description	Wt. Cap.	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
MDS89668	Basic, Steel	250-lbs.	17–21"	15"	17.5"	25½"	1/ea
G98202	Standard Non-Padded, Steel	300-lbs.	18¼–23¼"	18–20"	18½"	22"	1/ea
G98204	Premium Padded, Steel	350-lbs.	20–25"	18–20"	18½"	22"	1/ea
MDS89668XW	Bariatric, Steel	850-lbs.	20–25"	19"	28"	29"	1/ea

Grab Bars

- Help make bathrooms and stairways safer
- Chrome bars have a non-slip knurled surface
- Hanging hardware and instructions included

Chrome

Item No.	Length	Pkg.
MDS86012CHR	12"	3/cs
MDS86016CHR	16"	3/cs
MDS86018CHR	18"	3/cs
MDS86024CHR	24"	3/cs
MDS86032CHR	32"	3/cs

White

Item No.	Length	Pkg.
G87520M	12"	3/cs
G87530M	16"	3/cs
G87540M	18"	3/cs
G87550M	24"	3/cs
G87560M	32"	3/cs
G87570M	36"	3/cs



Toilet Safety Rail

- Assist patients with joint pain or at risk of falling
- Handles adjust and rotate back for a wide range of positions
- Foam armrests provide a comfortable, sure grip
- Easy-to-clean aluminum frame mounts securely onto bowl
- Adjustable bracket range of 18–24"

Item No.	Description	Wt. Cap.	Height Range	Width Range	Overall Depth Arms	Overall Depth Legs	Pkg.
G30300	Aluminum	300-lbs.	26 x 31"	18 x 24"	16¾"	12"	2/cs



G30300



G30270A

Locking Raised Toilet Seats

- Front clamping mechanism locks the seat onto the toilet, so it won't come loose or shift during use or transfer
- Wide, contoured surface raises seat level by 5"
- Fits most toilet bowls from 11 x 14"
- Optional padded armrests assist patient on and off the toilet

Item No.	Armrests	Wt. Cap.	Seat Elevation	Overall W x L	Opening W x L	Pkg.
G30260	No	350-lbs.	5"	20½ x 18"	8 x 9½"	3/cs
G30270A	Yes	350-lbs.	5"	20½ x 18"	8 x 9½"	3/cs



Economy Raised Toilet Seat

- Makes using the toilet easier for those with trouble bending
- Simple design made of polypropylene resin is easy to clean
- Non-slip pads and lip inside toilet offer additional stability

Item No.	Weight Capacity	Overall W x L	Opening W x L	Pkg.
G30250	250-lbs.	14½ x 15½"	8 x 9½"	3/cs



Compatible with
Quick-Alert
Alarms

Toilet Sensor Pads

- Compact toilet sensor pad allows caregivers to monitor patients while maintaining their dignity and privacy
- Position under the toilet seat
- 30-day warranty
- Alarm sold separately

Item No.	Pad Size	Pkg.
MDT9130TL1	6 W x 1 D x 3" H	1/ea
MDT9130TL5	6 W x 1 D x 3" H	5/cs

Urinals

- Patented splashproof lid and side storage help prevent accidental tipping and spilling for easier bedside care
- Graduated in ounces and milliliters
- 32-oz./1,000 mL clear plastic container

Item No.	Description	Pkg.
DYND80234	Mini Pack, Individually Shrink-Wrapped	12/cs
DYND80235S	Supreme Male Urinal	48/cs
DYND80235SD	Glow-in-the-dark Lid	48/cs
DYND80235SNL	Supreme, No lid	48/cs
PGA101	EZP Urinal Holder with Hook and Loop Strap	50/cs
PGA201	EZP Urinal Holder, Vinyl Strap	50/cs



DYND80235DH

Spil-Pruf Male Urinals

- Discreet and easy-to-use urinal helps increase patient confidence due to its spill and leak proof, odor-reducing features
- Offer peace of mind to bed and wheelchair bound patients
- Lightweight and durable polypropylene urinal holds up to 500 cc of liquid

Item No.	Description	Pkg.
EKJ149B	Black	12/cs
EKJ149CL	Clear	12/cs



EKJ149B

UltraFlex® Silicone Male External Catheters

- Built-in adhesive band offers ease of use and patient comfort
- Avoid allergic reactions to latex proteins with 100% silicone
- Non-sterile

Item No.	Description	Pkg.
RCH33301	Small, 25 mm	30/bx
RCH33302	Medium, 29 mm	30/bx
RCH33303	Large, 32 mm	30/bx
RCH33304	Large, 36 mm	30/bx
RCH33305	Extra-Large, 41 mm	30/bx



RCH33301



PSY5725

Seizure Side Rail Pads

- Form fitted pads fit snugly over existing bedside rails
- 1" layer of medium density foam helps safeguard against patient injury
- Sheer, lightweight nylon outer shell helps prevent patient "rub" burns and allows unobstructed movement of the side rails to a lowered position
- Open windows in upper section allow for easy access to controls on most bed models
- Includes 2 upper and 2 lower pads

Item No.	Description	Pkg.
PSY5742	Fits CareAssist and Advanta 2 Beds	1/ea
PSY5746	Fits TotalCare Med/Surg Beds	1/ea
PSY5755	Fits Secure I and II Beds with Rounded Side Rails	1/ea
PSY5763	Fits GO Bed II with Rectangular Side Rails	1/ea
PSY5725	Fits Advanta Bed, Waterproof	4/cs



NON081244

Disposable Foam Positioning Wedges

- Helps provide stability and elevation during surgeries, imaging and rehab
- Radiolucent so as not to interfere with crucial testing

Item No.	Size	Pkg.
NON081230	8 x 11 x 18"	6/cs
NON081233	16 x 11 x 5"	8/cs
NON081242	18 x 18 x 26"	4/cs
NON081244	7½ x 7½ x 17"	8/cs
NON081244H	7½ x 7½ x 17"	1/ea
NON081245	8 x 11 x 22"	8/cs
NON081245H	8 x 11 x 22"	1/ea
NON081254	10 x 7 x 2.5"	16/cs
NON081255	21¼ x 7 x 7"	16/cs
NON081265	23 x 11 x 11"	12/cs
NON081265H	23 x 11 x 11"	1/ea

Bed Assist Bar

- Provides help getting in and out of bed
- Features a soft grip for comfort, pocket for handy storage and straps for added safety
- Built to last with quick and easy tool-free assembly
- Includes a nylon bag for storage

Item No.	Description	Pkg.
MDS6800BA	Bed Assist Bar	3/cs



Side Rail Protector

- Help reduce entrapment by closing up gaps within and between the bed side rails
- Wipeable vinyl composition
- Pads secure easily with multiple hook and loop attachments

Item No.	Description	Pkg.
MDT821570	Side Rail Protector	1/ea



Advantage Mattress Bolster

- Transform any mattress to a raised edge mattress
- 8" bolsters help prevent patient falls and lower the risk of entrapment
- Quick installation and removal
- Easy-to-clean waterproof nylon cover reduces shear and friction

Item No.	Description	Pkg.
MDTMATBOLST	Mattress Bolster	1/ea



Create a culture of safety with Uplift.

Why Uplift?

- Improves caregiver satisfaction and well-being
- Ensures patient comfort and safety
- Helps retain staff and promote clinician wellness

Learn more at www.medline.com/go/uplift

Contact your Medline Representative
to schedule a free discovery assessment.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    

Reference: 1. Falls, Agency for Healthcare Research and Quality. January 2019: <https://psnet.ahrq.gov/primers/primer/40/Falls> 2. Coanour, CS et al. Cost of a Ventilator-Associated Pneumonia in a Shock Trauma Intensive Care Unit: <https://www.ncbi.nlm.nih.gov/pubmed/15865552> 3. Preventing Falls in Hospitals, Agency for Healthcare Research and Quality: <https://www.ahrq.gov/professionals/systems/hospital/pressureulcertoolkit/putool1.html> 4. Six Stats on the Cost of Readmission for CMS-Tracked Conditions. Becker's Hospital Review: <https://www.beckershospitalreview.com/quality/6-stats-on-the-cost-of-readmission-for-cms-tracked-conditions.html> 5. Safe Patient Handling—Preventing Musculoskeletal Disorders in Nursing Homes. 2012. OSHA(3108): www.osha.gov/Publications/OSHA3708.pdf 6. Worker Safety In Your Hospital, Occupational Safety and Health Administration. Available at: https://www.osha.gov/dsg/hospitals/documents/1.1_Data_highlights_508.pdf. Accessed July 18, 2019. 7. Safe Patient Handling. U.S. Department of Labor, Occupational Safety and Health Administration: <https://www.osha.gov/SLTC/healthcarefacilities/safepatienthandling.html> 8. Data on file. Available upon request. 9. NCBI, U.S. National Library of Medicine National Institute of Health. Outcomes of Safe Patient Handling and Mobilization Programs, Results. Teeple, Erin MD, MPH, Collins, Jamie E., PhD, Shrestha, Swastina, BS, Dennerlein, Jack T., PhD., Losina, Elena, PhD., and Katz, Jeffrey N., MD, MSc. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6138450/> Accessed June 20, 2019.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

© 2020 Medline Industries, Inc. All rights reserved. All product names, regardless of whether or not they appear with a trademark symbol, are the property of their respective companies unless otherwise noted. Trademark designations refer to trademarks registered or pending in the United States only. Medline and Uplift are registered trademarks of Medline Industries, Inc. MKT19W307853 / LIT053R / 2.5M / PP / 23