



# zip'r™ Traveler

- Easily disassembles into 5 pieces
- Maximum speed of up to 4.25 mph.
- Wireless, auto-latching assembly can be performed in less than 30 seconds.
- Headlight, basket and non-marking, flat-free tires are standard on all scooters.
- Upgrade to the Xtra for 4 inches of additional legroom, a glossy finish, and an upgraded seat.

Traveler disassembled



## Specifications Zip'r Traveler 3 & 4

Weight capacity	250 lbs.
Maximum speed*	4.25 mph
Per-charge range*	Up to 10 miles
Turning radius	25.4" / 37"
Ground clearance	2.5"
Headlight	24V, 5 W
Overall length	36.5"
Overall width	19.5"
Overall height	31.5"
Weight w/ batteries	99 lbs. / 104.8 lbs.
Heaviest piece	31.7 lbs. / 34.5 lbs.
Weight: Rear section	31.7 lbs.
Weight: Front section	28.7 lbs. / 34.5 lbs.
Weight: Seat	17.6 lbs.
Weight: Batteries	21 lbs. w/case
Battery	(2) 12 V, 12 Ah
Charger	Off-board, 2A
Drive	Rear
Motor specification	24 V, 155 W (.21 hp)
Brakes	Electromagnetic, regenerative
Seat type	Padded, foldable
Seat width	14"
Seat depth	12.5"
Seat back height	13"
Front wheel(s)	7.5" x 2" flat-free
Rear wheels	7.5" x 2.5" flat-free
Body colors	Red/Blue
HCPCS code	K0800

### Limited Warranty\*\*

Frame	3 years
Drivetrain	1 year
Electronic Components	1 year
Batteries	6 months

All specifications are subject to change without notification.

\*Specification may vary due to user weight, speed, battery condition and terrain.

\*\*Exclusions apply. See owner's manual for more details.