

SPECIFICATIONS

	Zip'r3	Zip'r4	Zip'r3 Xtra	Zip'r4 Xtra
Overall Dimensions: L x W x H	36.4" x 19.5" x 27"	36.4" x 19.5" x 27"	40.5" x 19" x 28"	40.5" x 19" x 28"
Front Tires	1 pcs 2.0" x 7.5"	2 pcs 2.0" x 7.5"	2 pcs 2.0" x 7.5"	2 pcs 2.0" x 7.5"
Rear Tires	2 pcs 2.5" x 7.5"	2 pcs 2.5" x 7.5"	2 pcs 2.5" x 7.5"	2 pcs 2.5" x 7.5"
Weight - Total	102 lbs	106 lbs	107 lbs	111 lbs
Weight - Heaviest Piece	29 lbs	29 lbs	33 lbs	33 lbs
Battery	12 V, 12 AH, 2pcs	12 V, 12 AH, 2pcs	12 V, 12 AH, 2pcs	12 V, 12 AH, 2pcs
Battery Charger	Off board, 2A current	Off board, 2A current	Off board, 2A current	Off board, 2A current
Power Rating of Motor	155 W, .21 hp	155 W, .21 hp	155 W, .21 hp	155 W, .21 hp
Ground Clearance	2.1"	2.1"	3"	3"
Turning Radius	30"	38"	40"	43"
Maximum Weight Capacity	250 lbs	250 lbs	250 lbs	250 lbs
Body Colors	Red/Blue	Red/Blue	Red/Blue Gloss Paint	Red/Blue Gloss Paint
Seat Type	Foldable foam filled	Foldable foam filled	Upgraded Foldable foam filled	Upgraded Foldable foam filled
Per-Charge Range (Max)	10 miles	10 miles	10 miles	10 miles
Maximum Speed	4.25 mph	4.25 mph	4.25 mph	4.25 mph