

**Medline Basic Steel Transport Chairs**

- Basic Steel Transport Chair, Hammertone, 19"**
- Back folds down for easy storage and transport
 - Features include comfortable nylon upholstery, detachable footrests and a seat belt for safety
 - Restaurant-style permanent armrests with durable plastic arm pads allow the user to pull up to most table heights
 - Optional accessories: IV pole (MDS85183FT*), O2 holder (MDS85181FT*), O2/IV combo (MDS85190FT*) antitip device (MDS85189FT), antifold antitheft device (MDS85196), Tinnerman leg rest locks (WCA806991), carrying case with strap (MDSCHAIRCASE)
- *Recommended for use with Antitip Device

Packaging	1 Case / Each
Manufacturer	Medline
Manufacturer #	MDS808200B

SKU MDS808200B

Specifications

ARM LENGTH	Full Length
BARIATRIC	No
COLOR	Black
FABRIC TYPE	Nylon
FOLDABLE	Yes
HANDBRAKES	No
LATEX FREE	No
LEG REST STYLE	Swing Away
MATERIAL	Steel
OVERALL WIDTH	23 "
QUICK RELEASE AXLES	No
REAR WHEEL SIZE	8 "
REMOVABLE ARM	No
SEAT DEPTH	16 "
SEAT WIDTH	19 "
SEAT TO FLOOR HEIGHT	18"
UNSPSC	42192210
WEIGHT CAPACITY	250 lb.

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Case (CS)	1.0 Each	0.0 / 33.05	5.167	24.0 x 12.0 x 31.0	193489002624
Each (EA)	1.0 Each	0.0 / 33.05	5.167	24.0 x 12.0 x 31.0	193489002624

California Customers Only... WARNING

HCPCS Information*HCPCS Code: NO GOVERNMENT OR SUPPLIER HCPCS CODE VERIFICATION**HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are***Documents/SDS**

No SDS content is available for this product.

VIDEOS

intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.