



FR 168-4S Series Spec

Model	FR-168-4S	FR-168-4S II
Top Speed:	4 mph	4 mph
Weight Capacity	350 lb	460 lb
Turning Radius	53.25"	
Angle Adjustable Tiller	Yes	
Angle Adjustable Seat	No	
Basket Included	Yes	
Battery Pack Weight	U1:22lb/ 50ah: 33 lb	
Battery Type	AGM	
Battery Charger Type	Off board	
Charging Port Location	Tiller	
Drive Range ¹	15 mile @ U1/ 20 miles @ 50ah	
Disassembles	No	
Delta Tiller	Yes	
Drive Type	Rear Wheel Drive 24V DC	
Flip-Back Armrests	Yes, full-long and padded	
Ground Clearance	3.5"	
Handbrake	No	
Headlight	Yes	
Height Adjustable Seat	No	
Height Adjustable Tiller	No	
Larger Battery Option	50ah AGM	
Max Incline Rating	6 degree	
Max Seat to Ground Height	22.83"	
Mirror Included	Yes, 1 unit	
Multiple Seating Options	No	
Overall Width	28"	
Overall Length	48"	
One-Hand Operation	Yes	
Front Wheel Size	10"	
Rear Wheel Size	10"	
Seat Belt Option	Optional	
Seat Type	Middle-Back	
Seat Width	20"	
Suspension	Yes (front and rear)	
Wheel Type	10" Solid PU	
Weight Without Batteries	110 lb	

Remark

1. Range will vary due to rider weight, drive surface, terrain, and battery type