




Alterra Long-Term Care Beds
Designed for safety, built to last



Find the bed that's just right.

Offering a wide range of choices, Alterra gives you the perfect balance of features and prices to meet your needs. The premium Maxx model has an expandable frame and offers 23 inches of height adjustment. Looking for something more basic? The 1100 is our most economical full-electric bed.

Make the most of your capital budget. Ask your Medline Representative about all our Alterra models, plus our individually tailored payment plans.

A Better Bedside Manner.

Developed from real-world observations of resident and caregiver needs, Alterra long-term care beds are a giant step forward in safety and performance.

RESIDENT AND CAREGIVER SAFETY

We saw residents struggling to get in and out of bed and staff straining their backs while providing care. So we engineered Alterra with essential and advanced safety features, and to meet FDA entrapment guidelines.

- » Exceptional high-low adjustment range helps reduce fall injury risks and ensure an ergonomic working position
- » Mattress retainers reduce slippage
- » Stability knobs keep head/footboards properly secured

INDUSTRY-LEADING PERFORMANCE

In our conversations with you, we learned that enhancing the comfort and well-being of both residents and staff is a top concern. Alterra is designed to optimize care and perform for years.

- » Top-of-the-line parallel electronics system delivers smooth and quiet operation
- » Synchronized motors mean no horizontal movement when raising and lowering (except model 1100)
- » Custom head/footboards, hardwood accents and beautiful finishes create a home-like feel
- » Exceptionally long warranty ensures peace of mind





MAX
height of
30.25"

LOW
height of
7.25"

Headboards and footboards
sold separately

Alterra Maxx

Expandable High-Low Bed



Width and length extenders are standard.



Five-function hand pendant comes standard with the Alterra Maxx bed.



OneTouchLock makes it easy to secure the casters for stationary positioning.



Built-in prop stand holds up bed deck head section during motor maintenance.

UNIVERSAL PATIENT FIT

Built-in extensions expand the length and width for patients of all sizes. Five-function hand pendant comes standard for exceptional ease of use. At the head and foot sections, our exclusive OneTouchLock™ secures the height with a simple tap of the foot. The smooth deck surface makes cleaning a snap.

SAFETY

A convenient LED nightlight helps patients see the floor when getting out of bed at night. In case of a power outage, the built-in battery backup keeps all functions available.

To avoid strain on caregiver backs, the deck raises to 30.25"—the highest in our line. The low height of 7.25" reduces risks of patient injury.

STRENGTH

The full 600-lb. weight capacity is the highest we offer.

APPEARANCE

Premium Alterra Maxx beds make a great impression on visitors and families.

Item number	Description
FCEALMAXX	Alterra Maxx Bed (Comes standard with integrated width and length extensions, 5-function hand pendant, mounting hardware, wall guard, LED underbed lights, head section kickstand, mattress retainers)
FCEMAXX37	Alterra Maxx Staff Control Kit

IDEAL FOR:



Skilled nursing facilities



Assisted living facilities



Nursing homes



Hospices


PRODUCT SPECIFICATIONS

Height range:	7.25" to 30.25"
Bed deck dimensions:	36"/39"/42" W x 80"/84"/88" L
Overall length:	91"
Weight capacity:	600-lb.
Movement styles:	Head, knee hinge, high-low, Trendelenburg, reverse Trendelenburg, chair
Deck style:	Rib
Rolling position:	Move at any height
Motors:	4 DC
Casters:	Locking, tracking
Horizontal movement:	No horizontal movement when raising or lowering
Battery backup:	Standard

Warranty information:

» Electrical components:	5 years
» Bed frame:	15 years
» Welds:	Lifetime

Bed Deck minimum and maximum size: ————— 42" —————





MAX
height of
29.88"

LOW
height of
7.75"

Assist bars, headboards and footboards
sold separately

Alterra 1385

High-Low Full-Electric Bed



Premium motors for smooth vertical movement.



OneTouchLock makes it easy to secure the casters for stationary positioning.



Raises to 29.88" to help reduce risks of caregiver back injury.

RESIDENT SAFETY

OneTouchLock feature secures the bed and prevents rolling. A low height of 7.75" helps prevent injuries with high fall-risk residents. All Alterra beds meet the entrapment guidelines of the FDA.

SAFER PATIENT HANDLING

Caregivers can reduce back injury risks by raising the Alterra 1385 to a full 29.88".

PERFORMANCE

Top-of-the-line parallel electronics system is exceptionally smooth and quiet. Perfectly synchronized high-low motors keep the bed deck level even after frequent, long-term use.

Battery backup helps during power outages, and the 500-lb. weight capacity accommodates heavier patients. Wall guard, mattress retainers and stability knobs for the headboard and footboard are standard.

Item number	Description
FCE1385B	Alterra 1385 Bed only
FCE1385NRQO	Alterra 1385 Economy Bed Package (Includes Alterra 1385 bed, head/footboards in Traditional style, Solar Oak finish, mounting hardware)
FCE1385SRQO	Alterra 1385 Bed Package (Includes Alterra 1385 bed, soft touch side rails, head/footboards in Traditional style, Solar Oak finish, mounting hardware)
FCE1385038	Alterra 1385 Staff Control Kit

IDEAL FOR:



Skilled nursing facilities



Assisted living facilities



Nursing homes



Hospices

PRODUCT SPECIFICATIONS

Height range:	7.75" to 29.88"
Bed deck dimensions:	35" W x 80" L
Overall length:	87.5"
Weight capacity:	500-lb.
Movement styles:	Head, knee hinge, high-low
Deck style:	Combination grid/rib
Rolling position:	Move at any height
Motors:	4 DC
Casters:	8 total, 4 locking
Horizontal movement:	No horizontal movement when raising or lowering
Battery backup:	Standard

Warranty information:

» Electrical components:	5 years
» Bed frame:	15 years
» Welds:	Lifetime



MAX
height of
26"

LOW
height of
7.25"


Side rails, headboards and footboards
sold separately

Alterra 1232

High-Low Full-Electric Bed



Lowers to just 7.25" for high fall-risk residents.



Battery backup helps during power outages.

RESIDENT SAFETY

A low height of only 7.25" helps reduce injuries from falling with high-risk patients. All Alterra beds meet the entrapment guidelines of the FDA.

PERFORMANCE

Battery backup helps during power outages. Wall guard, mattress retainers and stability knob for headboard and footboard are standard.

Item number	Description
FCE1232B	Alterra 1232 high-low bed only, with battery backup
FCE1232NRQ0	Alterra 1232 Economy Bed Package (includes Alterra 1232 bed, head/footboards in Solar Oak finish, mounting hardware)
FCE1232SRQ0	Alterra 1232 Bed Package (includes Alterra 1232 bed, head/footboards in Solar Oak finish, mounting hardware, soft touch side rails)
FCE1232Q0	Alterra 1232 Bed Package (includes Alterra 1232 bed, head/footboards in Solar Oak finish, mounting hardware, head section assist rails, 4 mattress retainers)
FCE1232038	Alterra 1232 Staff Control Kit

IDEAL FOR:



Skilled nursing facilities



Assisted living facilities



Nursing homes



Hospices

PRODUCT SPECIFICATIONS

Height range:	7.25" to 26"
Bed deck dimensions:	35" W x 80" L
Overall length:	87.5"
Weight capacity:	450-lb.
Movement styles:	Head, knee hinge, high-low
Deck style:	Combination grid/rib
Rolling position:	Move at any height
Motors:	4 DC
Casters:	8 total, 4 locking
Horizontal movement:	No horizontal movement when raising or lowering
Battery backup:	Standard

Warranty information:

» Electrical components:	4 years
» Bed frame:	15 years
» Welds:	Lifetime



FCE1100Q0 Bed Package shown.

Alterra 1100

Full-Electric Bed

IDEAL FOR:



Wall guard near the head section helps prevent damage.



Stability knobs ensure that the head and footboards are properly secured.



Mattress retainers keep the mattress in place for safer transfer in and out of bed.

COST-EFFECTIVE

Provides positioning comfort with independent and simultaneous operation of the head, knee and hi-low functions. Strong, well-built and low-maintenance electric bed is suitable for most residents in long-term care facilities.

RESIDENT SAFETY

A low height of 12" helps reduce injuries from falling with high-risk patients. All Alterra beds meet the entrapment guidelines of the FDA.

PERFORMANCE

The DC motors provide smooth, whisper-quiet operation.

PRODUCT SPECIFICATIONS

Height range:	12" to 23.5"
Bed deck dimensions:	35" W x 80" L
Overall length:	87.5"
Weight capacity:	500-lb.
Movement styles:	Head, knee hinge, high-low
Deck style:	Combination grid/rib
Rolling position:	Low position
Motors:	3 DC
Casters:	4 total, 3" (8 cm) locking
Horizontal movement:	Minimal horizontal movement when raising or lowering

Warranty information:

» Electrical components:	3 years
» Bed frame:	15 years
» Welds:	Lifetime

Item number	Description
FCE1100B	Alterra Full-Electric 1100 Bed only
FCE1100QO	Alterra Full Electric 1100 Bed Package (includes Medline Full Electric 1100 bed, head section assist rails, head/footboards in Solar Oak finish, mounting hardware, 4 mattress retainers)

At Home™ Headboards and Footboards



Help residents feel right at home. Solid wood headboards and footboards make the perfect finishing touch to a resident room. With our wide selection, you're sure to find ideal finishes to complement your décor.

- » Solid wood mounts directly to any Medline Alterra long-term care bed
- » Available in 36", 39" and 42" wide and 14 finish options
- » Available with a recessed opening to flush mount a staff control keypad (sold separately)

STANDARD FINISHES



Item number	Description
FCEHB36	36" Headboard Only
FCEFB36	36" Footboard Only
FCEFB36D	36" Footboard with Staff Control Cutout (Staff Control sold separately)
FCEFB36DN	36" Footboard with Staff Control Cutout for Alterra Maxx (FCEMAXX37 sold separately)
FCEHB39	39" Headboard Only
FCEFB39	39" Footboard Only
FCEFB39D	39" Footboard with Staff Control Cutout (Staff Control sold separately)
FCEFB39DN	39" Footboard with Staff Control Cutout for Alterra Maxx (FCEMAXX37 sold separately)
FCEHB42	42" Headboard Only
FCEFB42	42" Footboard Only
FCEFB42D	42" Footboard with Staff Control Cutout (Staff Control sold separately)
FCEFB42DN	42" Footboard with Staff Control Cutout for Alterra Maxx (FCEMAXX37 sold separately)
FCE1232038	Alterra 1232 Staff Control Kit
FCE1385038	Alterra 1385 Staff Control Kit
FCEMAXX37	Alterra Maxx Staff Control Kit

Made to order. Ask your Medline Representative about lead times.

Colors are reproduced as faithfully as printing processes allow. To ensure 100% satisfaction, please call your Medline Representative for swatch samples.



Ask your Medline Representative about our matching dressers, overbed tables, nightstands and more.

Custom Headboards and Footboards

Create a polished look and differentiate your facility with custom headboards and footboards. We offer a variety of finishes and styles, from traditional to modern, to match your room furnishings.



Avalon

Classic

Modern

Traditional

Transitional

Style	Headboard and Footboard Set			Headboard and Footboard Set with Staff Control Cutout		
	36" Wide	39" Wide	42" Wide	36" Wide	39" Wide	42" Wide
Avalon	FCET01RHFAV	FCET01RHFAVW	FCETRHFVAVXW	FCE36AVS	FCE39AVS	FCE42AVS
Classic	FCET01RHFCL	FCET01RHFCLW	FCERHFCLXW	FCE36CLS	FCE39CLS	FCE42CLS
Modern	MDRT01HFB	MDRT01HFBW	MDRT01HFBXW	FCE36MDS	FCE39MDS	FCE42MDS
Traditional	FCET01RHF	FCET01RHFW	FCET01RHFXW	FCE36TRS	FCE39TRS	FCE42TRS
				FCE36TRSMAX	FCE39TRSMAX	FCE42TRSMAX
Transitional	FCET01RHFTR	FCET01RHFTRW	FCERHFTRXW	FCE36TNS	FCE39TNS	FCE42TNS

Made to order. Ask your Medline Representative about lead times.

STANDARD FINISHES



Colors are reproduced as faithfully as printing processes allow. To ensure 100% satisfaction, please call your Medline Representative for swatch samples.

QUICK SHIP HEADBOARDS AND FOOTBOARDS

Delivered in 5-7 business days; staff control not available

Item	Style	Finish	Size
FCET01RHFAVO	Avalon	Solar Oak	36"
MDRT01HFBEM	Modern	Empire Mahogany	36"
FCET01RHFCHR	Traditional	American Cherry	36"
FCE1232RHFQ02	Traditional	Solar Oak	36"

Mounting Hardware (FCE1232HDW) included with Alterra Maxx (FCEALTMAXX), but must be ordered separately with all other headboard and footboard combinations.

Bed Accessories

Our extensive range of bed accessories complement our long-term care beds. Your Medline representative can help you determine the best options to meet your specific needs.



STAFF CONTROL
(FCE1232038/FCE1385038/
FCEMAXX37)*

This convenient control is mounted into the footboard of Alterra beds and allows staff members to lock any or all bed functions to make the bed safer for select residents.

*Staff Control is ordered as item FCEMAXX37 for Alterra Maxx bed (see pg. 5), item FCE1385038 for Alterra 1385 bed (see pg. 7), or item FCE1232038 for Alterra 1232 bed (see pg. 9). Staff control is not available for the Alterra 1100 bed.



ASSIST BAR
(FCE1232ABN)

This positioning aid has a cane shape that is easy to grip for elderly patients, when rolling over or getting in and out of bed. Also, swings away for easier patient transfer.



SOFT-TOUCH ASSIST RAILS
(FCE1232RSRN)

These steel rails have a coating that is soft to the touch and offer a less institutional look. They are comfortable to grip when rolling over or getting in and out of bed.

These assist rails can be mounted either on the head section or the foot section of the bed deck.



STANDARD ASSIST RAILS
(FCE1232RSR)

These 29" assist rails are long enough to cover the entire head section of the bed and offer additional protection for patients at a high risk for falls.

They can only be mounted to the head section of the bed deck

Ask your Medline Representative about our wide variety of mattresses, matching casegoods and linens.



QUICK-RELEASE SIDE RAILS
(FCEQRSIDE)

Quick-release side rails for Alterra beds are comfortable to grip when rolling over or getting in and out of bed.

These soft-touch, blow-molded side rails can easily be attached and detached for staff convenience.



EXTRA-TALL ASSIST RAILS
(FCE1232RSRXT)

Low air loss mattresses are generally between 10" and 12" thick, about the same height as standard bed rails. Our extra-tall rails are the perfect complement for low air loss mattresses.

These assist rails can be mounted either on the head section or on the foot section of the bed deck



TRAPEZE MOUNTING KIT
(FCE1232TPZ)

Our trapeze adapter kit is mounted on the bed's frame so that it can move up and down with the bed, allowing the trapeze to be used at any height.

Note: Trapeze (MDS80615T) sold separately.



BED EXTENDERS
(Length: FCE1232EXT/
Width: FCE42KIT)

Length extenders modify the bed from 80" to 84" to accommodate taller residents.

Width extenders modify the bed deck from 36" to 42" wide.

Note: Available for FCE1232B, FCE1100B and FCE1385B. Alterra Maxx features built-in width and length extenders.

Service Plans for Long-Term Care Beds

Pulling each bed from a resident room.
Performing maintenance and servicing.
Then putting it back. It all adds up to days
of staff time.

There's a better way.

Our bed service plans keep minor
issues from becoming larger
problems. Expert technicians from
Certified Centers of Excellence
(ISO 13485, 900) guarantee your
beds operate properly for
12 months after your service,
regardless of bed age. So you
can make the most of your
bed investment.



Choose the Plan That's Right for You

ANNUAL PREVENTIVE MAINTENANCE PLANS

- » Full inspection and lubrication of moving parts to ensure proper operation
- » Complete functional check of all electronics, mechanical components and rails

Ordering Information

Item No.	No. of Beds
FCESRVPLP1	5-20
FCESRVPLP2	21-50
FCESRVPLP3	51-100
FCESRVPLP4	100+

Hourly rates for troubleshooting or repairing a bed apply. Mileage fee applies if service is greater than 120 miles from our nearest service center. Contact your Medline Representative for more information.

ANNUAL FULL-SERVICE PLAN

All the advantages of the Preventive Maintenance Plan, plus unlimited troubleshooting and repair hours—ideal for equipment that's past warranty.

Ordering Information

Item No.	No. of Beds
FCESRVPLFS1	5-20
FCESRVPLFS2	21-50
FCESRVPLFS3	51-100
FCESRVPLFS4	100+

Mileage fee applies if service is greater than 120 miles from our nearest service center.

Ask your Medline Representative about pricing.



Ordering Guide

Choose Your Bed

- Alterra Maxx expandable high-low bed
- Alterra 1385 high-low bed
- Alterra 1232 high-low bed
- Alterra 1100 full-electric bed

Choose Your Headboard and Footboard

AT-HOME

- Headboard only
- Footboard only
- Footboard with Staff Control cutout

Width

- 36"
- 39"
- 42"

Finish

- Natural Maple
- Millwork Cherry
- Spiced Cherry
- Cordovan
- Warm Walnut
- Bayou Walnut
- Sand Dune
- Deep Forest
- Riverwood
- Asian Walnut
- Chestnut
- American Cherry
- Winchester Walnut
- Wild Cherry

CUSTOM

- Headboard and footboard set
- Headboard and footboard set with Staff Control cutout

Width

- 36"
- 39"
- 42"

Style

- Avalon
- Classic
- Modern
- Traditional
- Transitional

Finish

- Solar Oak
- Kensington Maple
- Empire Mahogany
- American Cherry
- Walnut

QUICK SHIP HEADBOARD AND FOOTBOARD SET

Width

- 36"

Style and Finish

- Avalon/Solar Oak
- Modern/Empire Mahogany
- Traditional/American Cherry
- Traditional/Solar Oak

Choose Your Hardware

Mounting Hardware for Headboards and Footboards

Must order for Alterra 1385, 1232 and 1100 headboards and footboards.

Comes standard with the Alterra Maxx.

Mounting hardware kit includes mattress retainers, "L" brackets and mounting screws for headboard and footboards.

Choose Your Accessories

	Alterra Maxx	Alterra 1385	Alterra 1232	Alterra 1100
Staff Control Kits <small>Note: Quick ship option not compatible with Staff Control kit</small>				
<input type="checkbox"/> For Alterra Maxx expandable high-low bed	X FCEMAXX37			
<input type="checkbox"/> For Alterra 1385 high-low bed		X FCE1385038		
<input type="checkbox"/> For Alterra 1232 high-low bed			X FCE1232038	
Rails and Assist				
<input type="checkbox"/> Assist bar	X	X	X	X
<input type="checkbox"/> Soft-touch assist rails	X	X	X	X
<input type="checkbox"/> Standard assist rails	X	X	X	X
<input type="checkbox"/> Quick-release side rails	X	X	X	X
<input type="checkbox"/> Extra-tall assist rails	X	X	X	X
<input type="checkbox"/> Trapeze adapter kit	X	X	X	X
<input type="checkbox"/> Trapeze	X	X	X	X
Bed Extenders				
<input type="checkbox"/> Length extender (to 84")	Extenders built in	X	X	X
<input type="checkbox"/> Width extender (to 42")		X	X	X
5-Function Hand Pendant				
<input type="checkbox"/> For Alterra 1385 high-low bed	Comes standard with bed	X FCE1385034		
<input type="checkbox"/> For Alterra 1232 high-low bed			X FCE1232034	

Alterra Maxx bed, comes standard with integrated width and length extensions, 5-function hand pendant, mounting hardware, wall guard, LED underbed lights, head section kickstand, and mattress retainers.

Choose Your Service Plan

ANNUAL PREVENTIVE MAINTENANCE PLANS

Number of Beds

- 5-20
- 21-50
- 51-100
- 100+

ANNUAL FULL-SERVICE PLAN

Number of Beds

- 5-20
- 21-50
- 51-100
- 100+

World Class Service, Delivered.

You found the beds you need. Now how do you get them delivered quickly and installed with a minimum of disruption? We can help.

From schedule coordination to old equipment removal, our team of experts will handle everything. Contact your Medline Representative for a complete list of our service offerings.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    