

**ComfortEZ 3-Position Recliners****CAL 117 Compliant ComfortEZ 3-Position Recliner, No Tray, Blueridge, Extra-Wide** • The ComfortEZ 3-position recliner is

- designed for comfort, versatility, and ease of use.
- Made with vinyl upholstery.
- 3-position headrest provides the comfort needed.
- Activity tray included on FURN10100BLT and FURN10100BET models

Packaging	1 Each / Each
Manufacturer	Medline
Manufacturer #	FURN10300BLC

SKU FURN10300BLC**Specifications**

ADJUSTMENT STYLE	Manual
ARM STYLE	Fixed
BARIATRIC	Yes
BASE STYLE	Mobile
CASTER SIZE	5 "
CASTER STYLE	All Locking/2 Swivel
CASTERS	Yes
COLOR	Blueridge
FABRIC TYPE	Vinyl
LATEX FREE	Yes
LENGTH	69 in
LENGTH INCHES	69 "
MRI SAFE	No
MATERIAL	Steel
NUMBER OF POSITIONS	3
PATIENT TRANSFER SYSTEM	No
RECLINER ACCESSORY OPTIONS	Activity Tray
RECLINER TYPE	Geri
SEAT DEPTH	20 in
SEAT DEPTH INCHES	20"
SEAT WIDTH INCHES	25"
TB-133	Yes
TRENDELENBURG	No
UNSPSC	42192102

Documents/SDS

No SDS content is available for this product.

WEIGHT CAPACITY
POUNDS

500 lb

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Each (EA)	1.0 Each	0.0 / 138.0	30.0	40.0 x 36.0 x 36.0	10888277143067

HCPCS Information

HCPCS Code: NO CODE - PRODUCT DOES NOT MEET HCPCS SPECIFICATIONS

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.