

Model	Go-Go® Endurance Li ¹
Part Number	S54XLIT
Weight Capacity	325 lbs.
Maximum Speed ²	Up to 4.2 mph
Ground Clearance ⁴	3.25" at side deck
Turning Radius ⁴	48.25"
Overall Length ⁴	42"
Overall Width ⁴	20.5"
Seat-to-Deck Height ⁴	16.25" - 17.75"
Seat-to-Ground Height ⁴	21.25" - 22.75"
Front Tires ⁴	2" x 8" solid
Rear Tires ⁴	2.5" x 8" solid
Range Per Charge ^{2,3}	Up to 7.2 miles (8AH)** Up to 14.4 miles (16AH MAXX) ⁵
Total Weight Without Batteries ⁴	93.5 lbs. w/o battery pack; with 8AH battery pack**: 101.5 lbs.; with 16AH MAXX battery pack ⁵ : 105.5 lbs.
Heaviest Piece When Disassembled ⁴	40.75 lbs. (front section)
Standard Seating	Type: Compact, foldable Material: Black vinyl Dimensions: 17" x 17"; 20" x 18" Weight: 23 lbs. (17" x 17"); 27.75 lbs. (20" x 18")
Motor Output	24-volt DC motor, sealed transaxle
Brakes	Regenerative and electromechanical
Battery Requirements ⁴	Size: 8 AH lithium ion pack** or 16 AH MAXX lithium ion pack ⁵ Weight: 8 lbs. each (8AH); 12 lbs. each (16AH MAXX)
Battery Charger	Off-board, 3.5 amp (lithium ion)
Warranty	Frame: Lifetime limited Electronics: 2-year limited Drive motors: 2-year limited Batteries: 2-year limited

HCPCS Code: K0899
POWER MOBILITY DEVICE, NOT CODED
BY DME PDAC OR DOES NOT MEET CRITERIA

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Results derived from theoretical calculations based on lithium battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Does not meet airline requirements.

PRIDE®
MOBILITY
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    



** Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go® Endurance Li and its lithium ion battery are permitted.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GOGOENDURANCELi_082621

GOGO®¹
ENDURANCE Li

DISCOVER THE
PORTABLE POWER
OF LITHIUM ION



FDA Class II Medical Device¹

LIVE YOUR BEST®

pridemobility.com |    

LITHIUM ION MEETS PORTABILITY



**LONGER
LIFE**

LITHIUM ION BATTERY


- Faster charging time
- Sustained power performance
- Lighter weight
- More useable energy
- Longer life
- Smart battery protection



CTS Suspension for greater comfort



Charger port located conveniently in tiller.



Feather-touch disassembly for quick transport or storage.



Comfort-Trac front and rear suspension for added comfort.

With the **Go-Go® Endurance Li**, traveling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With the standard 8AH airline-compliant lithium-ion battery pack**, the Go-Go Endurance is perfect for air travel. Feather-touch disassembly provides the portability you need. When you arrive at your destination, enjoy a smooth, comfortable ride with Comfort-Trac Suspension.

Additional Features

- Per charge range up to 14.4 miles w/ 16AH MAXX battery pack⁵; up to 7.2 miles w/ 8AH battery pack**
- Maximum speed up to 4.2 mph
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrists
- Convenient dual charging ports can charge battery pack on-board or off-board
- Includes two sets of easily changeable, colored panels
- Standard 8AH airline-compliant lithium-ion battery pack** or upgrade to 16AH MAXX battery pack⁵
- Two seat options: 17" x 17" and 20" x 18"

****Airline-compliant. Please see note on back page for details.**

Colors



Garnet Red (matte)

Euro Grey

Accessories*

- Cane clip attachment
- Cane/crutch holders
- Cell phone holder
- Cup holder
- Lap belt (50", 60", 70")
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Under seat storage
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Additional spare 8AH lithium-ion battery pack**
- Additional spare 16AH MAXX lithium-ion battery pack⁵

Military Patches



***Only one rear accessory can be used at a time.**