

Folding Steel Commode

HCPCS: E0163



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Easy, tool-free assembly
- Easily opens and folds for storage and transportation
- Durable plastic snap-on seat and lid
- Powder-coated steel-welded construction
- Plastic arms for added comfort
- Includes 7.5 qt commode bucket with carry handle, cover and splash shield
- Retail Packaged
- Made using recycled materials

PRODUCTS

Description

RTL11158KDR	Commode, Folding Frame, RTL,1 EA
-------------	----------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	17.5 lbs
Seat Dimensions	13.75" (W) x 16.25" (D)
Seat Height	17.0" X 23.0"
Warranty	Limited Lifetime
Weight	13 lbs
Weight Capacity	350 lbs
Width	18.125" between arms 22.25" Outside