



FUTURE MOBILITY PRODUCTS INC.

DIVISION OF ABF MFG GROUP INC.
WWW.FUTUREMOBILITY.COM



CAPELLA

TEL: 716-783-9130
FAX: 716-783-9236
yellahi@futuremobility.com

Account Information

Account #: _____
Business Name: _____
Ship to Address: _____
Contact: _____
Tag For: _____

Order Date: _____
P.O.# _____
Phone: _____
City: _____
State: _____ Zip Code: _____

The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy. Prices are subject to change without notice

1. FRAME TYPE

	Parts #	E-Code	Retail Price
<input type="checkbox"/> Capella 20° Dynamic Tilt (Max Weight 350 lbs)	109	E1161	\$2,800.00
<input type="checkbox"/> Capella 45° Tilt (Max Weight 250 lbs)	109T45	E1161	\$3,400.00
<input type="checkbox"/> Capella 40° Power Tilt (Capella 45 only)	109PT45	-	\$1,600.00

2. SEAT DUMP ANGLE

<input type="checkbox"/> 0°	109SDA1	-	Standard
<input type="checkbox"/> 2°	109SDA2	-	No Charge
<input type="checkbox"/> 4°	109SDA3	-	No Charge

3. FRAME FINISH

<input type="checkbox"/> Wet Black		109C01	-	Standard
<input type="checkbox"/> Silver		109C02	-	No Charge
<input type="checkbox"/> Sky Blue		109C22	-	No Charge
<input type="checkbox"/> Cobalt Blue		109C03	-	No Charge
<input type="checkbox"/> Green Apple		109C04	-	No Charge
<input type="checkbox"/> Burnt Orange		109C05	-	No Charge
<input type="checkbox"/> Grape Purple		109C06	-	No Charge
<input type="checkbox"/> Fire Red		109C07	-	No Charge
<input type="checkbox"/> Yellow Ribbon		109C08	-	No Charge
<input type="checkbox"/> Custom Colours		109C17	-	\$300.00

Allow 7-10 days for powder coating time for custom colors.

4. SEAT WIDTH

<input type="checkbox"/> 16"	109SW16	-	No Charge
<input type="checkbox"/> 17"	109SW17	-	No Charge
<input type="checkbox"/> 18"	109SW18	-	Standard
<input type="checkbox"/> 19"	109SW19	-	No Charge
<input type="checkbox"/> 20"	109SW20	E2201	\$395.00
<input type="checkbox"/> 20"	109SW20	E2201	\$395.00
<input type="checkbox"/> 21" (Up to 350 lbs)	109SW21	E2201	\$395.00
<input type="checkbox"/> 22" (Up to 350 lbs)	109SW22	E2201	\$395.00
<input type="checkbox"/> 23" (Up to 350 lbs)	109SW23	E2201	\$395.00
<input type="checkbox"/> 24" (Up to 350 lbs)	109SW24	E2202	\$600.00

5. SEAT DEPTH

<input type="checkbox"/> 15"	109SD15	-	No Charge
<input type="checkbox"/> 16"	109SD16	-	Standard
<input type="checkbox"/> 17"	109SD17	-	No Charge
<input type="checkbox"/> 18"	109SD18	-	No Charge
<input type="checkbox"/> 19"	109SD19	-	No Charge
<input type="checkbox"/> 20"	109SD20	E2203	\$511.00

6. SEAT PAN OPTIONS

<input type="checkbox"/> No Seat Pan	109SE2	-	No Charge
<input type="checkbox"/> Solid Seat Pan	109SE1	E2231	\$99.00
<input type="checkbox"/> Adjustable Seat Pan	109AS1	E2231	\$99.00



CAPELLA FRONT SEAT TO FLOOR HEIGHT MATRIX

Using the tables find desired seat to floor height. Then select front caster size and rear wheel size combination. All Dimensions are approximate

Capella Front STF Height Availability										
4" Caster										
Rear Wheel	12"	11.75"	12.25"							
16"	11.75"									
5" Caster										
Rear Wheel	16"	13"	13.5"	14"	14.5"					
20"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"		
22"			14"	14.5"	15"	15.5"	16"	16.5"		
24"					15"	15.5"	16"	16.5"		
6" Caster										
Rear Wheel	16"		13.5"	14"	14.5"					
20"		13.5"	14"	14.5"	15"	15.5"	16"	16.5"		
22"			14"	14.5"	15"	15.5"	16"	16.5"	17"	
24"					15"	15.5"	16"	16.5"	17"	
7" Caster										
Rear Wheel	16"									
20"					15"		16"	16.5"		
22"					15"		16"	16.5"	17"	17.5"
24"						15"	16"	16.5"	17"	17.5"
8" Caster										
Rear Wheel	22"								17"	
24"									17"	18" 19"

Front STF	Rear Wheel	Axel Pos.
11.75	16"	2
12.25	12"	4
13.0	16"	3
13.0	20"	1
13.5	16"	3
13.5	20"	2
14.0	16"	4
14.0	20"	2
14.0	22"	1
14.5	16"	4
14.5	20"	2
14.5	22"	1
15.0	20"	3
15.0	22"	2
15.0	24"	1
15.5	20"	3
15.5	22"	2
15.5	24"	1
16.0	20"	4
16.0	22"	3
16.0	24"	2
16.5	20"	4
16.5	22"	3
16.5	24"	2
17.0	22"	4
17.0	24"	3
17.5	22"	4
17.5	24"	4
18.0	24"	4
19.0	24"	4

* 1" RISER BLOCK REQUIRED

7. SEAT TO FLOOR HEIGHT

Parts #	E-Code	Retail Price
<input type="checkbox"/> 11.75"	109SF28	- No Charge
<input type="checkbox"/> 12.25"	109SF29	- No Charge
<input type="checkbox"/> 13"	109SF02	- No Charge
<input type="checkbox"/> 13.5"	109SF03	- No Charge
<input type="checkbox"/> 14"	109SF04	- No Charge
<input type="checkbox"/> 14.5"	109SF05	- No Charge
<input type="checkbox"/> 15"	109SF06	- No Charge
<input type="checkbox"/> 15.5"	109SF07	- No Charge
<input type="checkbox"/> 16"	109SF08	- Standard
<input type="checkbox"/> 16.5"	109SF09	- No Charge
<input type="checkbox"/> 17"	109SF10	- No Charge
<input type="checkbox"/> 17.5"	109SF11	- No Charge
<input type="checkbox"/> 18"	109SF12	- No Charge
<input type="checkbox"/> 19" (Uses 1" Riser Block Kit)	109SF14	- \$175.00

8. CASTER

<input type="checkbox"/> 4" x 1"	109CW13	- No Charge
<input type="checkbox"/> 5" x 1.25"	109CW02	- No Charge
<input type="checkbox"/> 6" x 1.25"	109CW04	- Standard
<input type="checkbox"/> 7" x 1"	109CW06	- No Charge
<input type="checkbox"/> 8" x 1.25"	109CW07	- No Charge

1 Buffalo River Place, Buffalo, NY 14210, United States Phone: (716) 783-9130 Fax: (716) 783-9236
E-mail: yellahi@futuremobility.com Website: www.FutureMobility.com

Prices and specs are subject to change without notice

9. REAR WHEEL

	Parts #	E-Code	Retail Price
<input type="checkbox"/> 12" Composite	109RW01	-	No Charge
<input type="checkbox"/> 16" Composite	109RW02	-	No Charge
<input type="checkbox"/> 20" Composite	109RW03	-	No Charge
<input type="checkbox"/> 22" Composite	109RW04	-	Standard
<input type="checkbox"/> 24" Composite	109RW05	-	No Charge
<input type="checkbox"/> 20" Spoke	109RW06	-	\$255.00
<input type="checkbox"/> 22" Spoke	109RW07	-	\$255.00
<input type="checkbox"/> 24" Spoke	109RW08	-	\$255.00

10. REAR WHEEL TIRE

<input type="checkbox"/> Airless (Solid Rubber)	109RT1	-	Standard
<input type="checkbox"/> Pneumatic (Only Available for 20", 22" or 24")	109RT2	E2211	\$60.00

11. AXLE TYPE

<input type="checkbox"/> Fixed Rear Axle	109AX3	-	Standard
<input type="checkbox"/> Quick Release Axle	109AX1	-	No Charge

12. HANDRIM

Avail. with 20"/22"/24" Rear Wheels

<input type="checkbox"/> Chrome	109HR1	-	Standard
<input type="checkbox"/> Plastic Coated	109HR2	-	\$90.00

13. TILT TRIGGER LOCATION


<input type="checkbox"/> Push Cane	109TH1	-	Standard
<input type="checkbox"/> Stroller Bar	109TH2	-	No Charge
<input type="checkbox"/> Foot Tilt Release - Left	109TH3	-	\$95.00
<input type="checkbox"/> Foot Tilt Release - Right	109TH4	-	\$95.00

14. ARM TYPE

<input type="checkbox"/> T-Style Fixed Height	109AR02	E0973	Standard
<input type="checkbox"/> Flip Back Fixed Height	109AR10	E0973	\$210.00
<input type="checkbox"/> Flip Back Adj. Height	109AR09	E0973	\$210.00
<input type="checkbox"/> T-Style Adj. Height Pin Release - Youth	109AR05	E0973	\$210.00
<input type="checkbox"/> T-Style Adj. Height Pin Release - Adult	109AR06	E0973	\$210.00
<input type="checkbox"/> T-Style Adj. Height Knob Release	109AR30	E0973	\$210.00

Armrest Matrix for T - Style Pin Release Arm

Armrest Style	Transfer Tube Length D1	Arm Guard Length D2	H Range
Adult	6.75"	7.50"	9.50" - 14.50"
Youth	4.50"	5.00"	8.75" - 13.75"



Arm Pad Length L		
Desk	Full	Curved
10.25"	14.0"	17.0"

15. ARM PAD TYPE

<input type="checkbox"/> Desk length Vinyl 10" (pair)	109AP1	-	No Charge
<input type="checkbox"/> Full length Vinyl 14" (pair)	109AP2	-	No Charge
<input type="checkbox"/> Contour Gel 17" (pair)	109AP3	-	Standard
<input type="checkbox"/> Contour Gel Wide 16" (pair)	109AP4	-	\$150.00

16. SEAT TO BACK ANGLE

<input type="checkbox"/> 83°	109SBA1	-	No Charge
<input type="checkbox"/> 90°	109SBA2	-	Standard
<input type="checkbox"/> 97°	109SBA3	-	No Charge
<input type="checkbox"/> 104°	109SBA4	-	No Charge
<input type="checkbox"/> 111°	109SBA5	-	No Charge
<input type="checkbox"/> 118°	109SBA6	-	No Charge

17. BACK POST

<input type="checkbox"/> Fold Down Canes 20" Height	109PC3	-	No Charge
<input type="checkbox"/> Fold Down Canes 25" Height	109PC4	-	No Charge
<input type="checkbox"/> Stroller Bars Fixed 20" Height	109SH9	-	No Charge
<input type="checkbox"/> Stroller Bars Fixed 25" Height	109SH10	-	No Charge
<input type="checkbox"/> Stroller Bars Fold Down 20" Height	109SP13	-	No Charge
<input type="checkbox"/> Stroller Bars Fold Down 25" Height	109SP14	-	No Charge
<input type="checkbox"/> Stroller Bars Adj. Hinged 20" Height	109SH5	-	No Charge
<input type="checkbox"/> Stroller Bars Adj. Hinged 25" Height	109SH7	-	No Charge
<input type="checkbox"/> Stroller Bars Adj. Hinged Fold Down Low 20" Height	109SH11	-	No Charge
<input type="checkbox"/> Stroller Bars Adj. Hinged Fold Down High 25" Height	109SH13	-	No Charge

Recline Options

<input type="checkbox"/> Stroller Bars Recline 20" Height	109SP11	E1225	\$650.00
<input type="checkbox"/> Stroller Bars Recline 25" Height	109SP12	E1225	\$650.00

A. Recline options cannot be used with the Capella 40 power tilt.

18. BACK OPTIONS

	Parts #	E-Code	Retail Price
<input type="checkbox"/> No Back	101BK49	-	Standard
<input type="checkbox"/> Prism Basic Back 14"-20" wide	BB14-20	E2620	\$615.00
<input type="checkbox"/> Prism Basic Back 22"-26" wide	BB22-26	E2621	\$750.00
<input type="checkbox"/> Prism Basic Back 27"-32" wide	BBCUS	E2621	\$775.00
<input type="checkbox"/> Prism Ultra Back 14"-20" wide	UB14-20	E2620	\$615.00
<input type="checkbox"/> Prism Ultra Back 22"-26" wide	UB22-26	E2621	\$750.00
<input type="checkbox"/> Prism Ultra Back 27"-32" wide	UBBCUS	E2621	\$775.00
<input type="checkbox"/> Prism Air Back 14"-20" wide	AB14-20	E2615	\$655.00
<input type="checkbox"/> Prism Air Back 22"-26" wide	AB22-26	E2615	\$755.00
<input type="checkbox"/> Prism Air Back 27"-32" wide	ABCUS	E2616	\$815.00
<input type="checkbox"/> Prism TrueFitt Back 14"-20" wide	TFB14-20	E2615	\$625.00
<input type="checkbox"/> Prism TrueFitt Back 22"-26" wide	TFB22-26	E2616	\$725.00
<input type="checkbox"/> Prism TrueFitt Back 27"-32" wide	TFBCUS	E2616	\$775.00
<input type="checkbox"/> Max Ultra Back 14"-20" wide	MUB14-20	E2620	\$645.00
<input type="checkbox"/> Max Ultra Back 22"-26" wide	MUB22-26	E2620	\$695.00
<input type="checkbox"/> Max Ultra Back 27"-32" wide	MUBCUS	E2620/21	\$725.00
<input type="checkbox"/> Max Air Back 14"-20" wide	MAB14-20	E2620	\$725.00
<input type="checkbox"/> Max Air Back 22"-26" wide	MAB22-26	E2620	\$755.00
<input type="checkbox"/> Max Air Back 27"-32" wide	MABCUS	E2620/21	\$775.00
<input type="checkbox"/> Max TrueFitt Back 14"-20" wide	MTFB14-20	E2620	\$675.00
<input type="checkbox"/> Max TrueFitt Back 22"-26" wide	MTFB22-26	E2620	\$725.00
<input type="checkbox"/> Max TrueFitt Back 27"-32" wide	MTFBCUS	E2620	\$775.00
<input type="checkbox"/> Lunar Back Low No Lateral 14"-20"	LBSNL14-20	E2620	\$595.00
<input type="checkbox"/> Lunar Back Low Basic Lateral 14"-20"	LBSBL14-20	E2620	\$670.00
<input type="checkbox"/> Lunar Back Mid No Lateral 14"-20"	MBSNL14-20	E2620	\$645.00
<input type="checkbox"/> Lunar Back Mid Basic Lateral 14"-20"	MBSBL14-20	E2620	\$720.00
<input type="checkbox"/> Lunar Back Mid Ultra Lateral 14"-20"	MBSUL14-20	E2620	\$795.00
<input type="checkbox"/> Lunar Back Full No Lateral 14"-20"	FBSNL16-22	E2620/21	\$695.00
<input type="checkbox"/> Lunar Back Full Basic Lateral 14"-20"	FBSUL16-22	E2620/21	\$770.00
<input type="checkbox"/> Lunar Back Full Ultra Lateral 14"-20"	FBSUL16-22	E2620/21	\$845.00
<input type="checkbox"/> Lunar Back Custom 14"-20"	LUNARCUS	E2620/21	\$950.00

W () x H ()

19. BACK PINS

<input type="checkbox"/> Longer Pins for Back Hardware 1" Growth	109BP1	E2620	\$50.00
<input type="checkbox"/> Longer Pins for Back Hardware 2" Growth	109BP2	E2620	\$50.00

20. SEAT CUSHION

<input type="checkbox"/> No Cushion	101SC20	-	Standard
<input type="checkbox"/> Serene Prime 16"-20"	PR16-24	E2601	\$100.00
<input type="checkbox"/> Serene Paramount 16"-20"	PA16-20	E2605	\$300.00
<input type="checkbox"/> Prism Ideal 16"-20"	IC16-20	E2601	\$100.00
<input type="checkbox"/> Prism Ideal 22"-26"	IC22-26	E2601	\$175.00
<input type="checkbox"/> Prism Ideal 27"-32"	ICCUS	E2601/02	\$225.00
<input type="checkbox"/> Prism Elite 16"-20"	EC16-20	E2603	\$225.00
<input type="checkbox"/> Prism Elite 22"-26"	EC22-26	E2603	\$300.00
<input type="checkbox"/> Prism Elite 27"-32"	ECCUS	E2603/4	\$375.00
<input type="checkbox"/> Prism Supreme 16"-20"	SC16-20	E2607	\$365.00
<input type="checkbox"/> Prism Supreme 22"-26"	SC22-26	E2607	\$495.00
<input type="checkbox"/> Prism Supreme 27"-32"	SCCUS	E2607/08	\$595.00
<input type="checkbox"/> Prism Supreme II Plus 16"-20"	SCP16-20	E2607	\$385.00
<input type="checkbox"/> Prism Supreme II Plus 22"-26"	SCP22-26	E2607	\$425.00
<input type="checkbox"/> Prism Supreme II Plus 27"-32"	SCPCUS	E2607/08	\$475.00
<input type="checkbox"/> Prism Supreme Air 16"-20"	SAC16-20	E2624	\$475.00
<input type="checkbox"/> Prism Supreme Air 22"-26"	SAC22-26	E2625	\$550.00
<input type="checkbox"/> Prism Supreme Air 27"-32"	SACCUS	E2624/25	\$595.00
<input type="checkbox"/> Prism Supreme Gel 16"-20"	SGC16-20	E2607	\$395.00
<input type="checkbox"/> Prism Supreme Gel 22"-26"	SGC22-26	E2607	\$425.00
<input type="checkbox"/> Prism Supreme Gel 27"-32"	SGCUS	E2607/08	\$475.00

W () x H ()



1 Buffalo River Place, Buffalo, NY 14210, United States Phone: (716) 783-9130 Fax: (716) 783-9236
 E-mail: yellahi@futuremobility.com Website: www.FutureMobility.com

Prices and specs are subject to change without notice

21. FRONT RIGGING

	Parts #	E-Code	Retail Price
<input type="checkbox"/> 70° Pin Style	109FR01	-	No Charge
<input type="checkbox"/> 80° Pin Style	109FR02	-	No Charge
<input type="checkbox"/> 60° Swingaway	109FR03	-	No Charge
<input type="checkbox"/> 70° Swingaway	109FR04	-	No Charge
<input type="checkbox"/> 90° Swingaway (Not available w/ 8" casters)	109FR05	-	No Charge
<input type="checkbox"/> 60° Dual Rotation	109FR06	-	\$230.00
<input type="checkbox"/> 70° Dual Rotation	109FR07	-	\$230.00
<input type="checkbox"/> 80° Dual Rotation	109FR08	-	\$230.00
<input type="checkbox"/> 90° Dual Rotation	109FR09	-	\$230.00
<input type="checkbox"/> Swingaway Elevating Legrest (pair)	109FR10	E0990	\$275.00
<input type="checkbox"/> Dual Swingaway Elevating Legrest (pair)	109FR25	E0990	\$290.00
<input type="checkbox"/> No Front Rigging	109FR11	-	No Charge
<input type="checkbox"/> No Footrest Hanger	109FR12	-	No Charge

22. FOOTPLATE

<input type="checkbox"/> Composite Flip Up Small	109FR01	-	No Charge
<input type="checkbox"/> Composite Flip Up Medium	109FR02	-	No Charge
<input type="checkbox"/> Aluminum Flip Up Small	109FR05	-	No Charge
<input type="checkbox"/> Aluminum Flip Up Medium	109FR06	-	No Charge
<input type="checkbox"/> Aluminum Angle Adj. Flip Up Small	109FR10	K0040	\$120.00
<input type="checkbox"/> Aluminum Angle Adj. Flip Up Medium	109FR11	K0040	\$120.00
<input type="checkbox"/> Aluminum Angle Adj. Flip Up Large	109FR12	K0040	\$120.00
<input type="checkbox"/> One Piece Footboard 15"-20"	109FR19	-	\$160.00
<input type="checkbox"/> One Piece Footboard 21"-24" Custom	109FR20	-	\$190.00

23. FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> Heel Loops	109FA01	E0951	\$40.00
<input type="checkbox"/> No Heel Loops	109FA02	-	No Charge
<input type="checkbox"/> Padded Calf Strap	109FA03	K0038	\$45.00
<input type="checkbox"/> Swingway Calf Pad (pair)	109FA04	-	\$185.00
<input type="checkbox"/> Foam Guard on Elevating Legrest (pair)	109FA05	-	\$45.00
<input type="checkbox"/> Foam Pad One Piece Footboard	109FA06	-	\$50.00
<input type="checkbox"/> Foam Pad on Footplate (pair)	109FA07	-	\$45.00

24. WHEEL LOCKS

<input type="checkbox"/> Push to Lock	109WL1	-	No Charge
<input type="checkbox"/> Push to Lock + Wheel Lock Extension	109WL3	E0961	\$60.00
<input type="checkbox"/> Pull to Lock	109WL2	-	No Charge
<input type="checkbox"/> Pull to Lock + Wheel Lock Extension	109WL3	E0961	\$60.00
<input type="checkbox"/> Attendant Wheel Lock (Foot Activated) RT	109WL8	E2206	\$195.00
<input type="checkbox"/> Attendant Wheel Lock (Foot Activated) LT	109WL9	E2206	\$195.00
<input type="checkbox"/> Unilateral Push Wheel Lock RT Activated	109WL10	-	\$250.00
<input type="checkbox"/> Unilateral Push Wheel Lock LT Activated	109WL11	-	\$250.00
<input type="checkbox"/> Unilateral Pull Wheel Lock RT Activated	109WL12	-	\$250.00
<input type="checkbox"/> Unilateral Pull Wheel Lock LT Activated	109WL13	-	\$250.00
<input type="checkbox"/> Berner Push Wheel Lock	109WL14	-	\$95.00
<input type="checkbox"/> Berner Pull Wheel Lock	109WL15	-	\$95.00

25. HEADREST HARDWARE

<input type="checkbox"/> Omit Headrest Hardware	109HE99	-	Standard
<input type="checkbox"/> Fixed Hardware	109HE01	E1028	\$195.00
<input type="checkbox"/> Fixed Headrest Curved	109HE15	E1028	\$195.00
<input type="checkbox"/> Fixed Hardware (+ 4")	109HE02	E1028	\$200.00
<input type="checkbox"/> Fixed Headrest Curved (+ 4")	109HE16	E1028	\$200.00
<input type="checkbox"/> Single Multi Axial	109HE03	E1028	\$217.00
<input type="checkbox"/> Single Multi Axial Curved	109HE17	E1028	\$175.00
<input type="checkbox"/> Single Multi Axial (+ 4")	109HE04	E1028	\$222.00
<input type="checkbox"/> Single Multi Axial Curved (+ 4")	109HE18	E1028	\$200.00
<input type="checkbox"/> Dual Multi Axial	109HE05	E1028	\$219.00
<input type="checkbox"/> Dual Multi Axial Curved	109HE19	E1028	\$225.00
<input type="checkbox"/> Dual Multi Axial (+ 4")	109HE06	E1028	\$220.00
<input type="checkbox"/> Dual Multi Axial Curved (+ 4")	109HE20	E1028	\$250.00
<input type="checkbox"/> Infinite Multi Axial	109HE21	E1028	\$225.00
<input type="checkbox"/> Infinite Multi Axial Curved	109HE22	E1028	\$250.00
<input type="checkbox"/> Infinite Multi Axial (+ 4")	109HE23	E1028	\$270.00
<input type="checkbox"/> Infinite Multi Axial Curved (+ 4")	109HE24	E1028	\$270.00
<input type="checkbox"/> Links Multi Adj. Hardware (Comes with 2 links)	109HE10	E1028	\$200.00
<input type="checkbox"/> Links Multi Adj. Curved Hardware (Comes with 2 links)	109HE11	E1028	\$240.00
<input type="checkbox"/> Additional Link	109HE12	E1028	\$50.00
<input type="checkbox"/> 2 Additional Links	109HE13	E1028	\$100.00
<input type="checkbox"/> 3 Additional Links	109HE14	E1028	\$150.00

26. HEADREST ADAPTOR BRACKETS

<input type="checkbox"/> Universal Headrest Bracket Assembly	109HEBA	E1028	\$86.00
--	---------	-------	---------

27. HEADREST FOAM PAD

<input type="checkbox"/> Omit Headrest Foam Pad	109HE09	-	Standard
<input type="checkbox"/> Steady - Small	109HF1	E0955	\$214.00
<input type="checkbox"/> Steady - Medium	109HF7	E0955	\$230.00
<input type="checkbox"/> Steady - Large	109HF3	E0955	\$250.00
<input type="checkbox"/> Curvus - Small	109HF2	E0955	\$225.00
<input type="checkbox"/> Curvus - Medium	109HF8	E0955	\$260.00
<input type="checkbox"/> Curvus - Large	109HF4	E0955	\$275.00
<input type="checkbox"/> Pro-Fit - Small	109HF12	-	\$225.00
<input type="checkbox"/> Pro-Fit - Medium	109HF13	-	\$260.00
<input type="checkbox"/> Pro-Fit - Large	109HF14	-	\$275.00

28. ANTI TIP

<input type="checkbox"/> Quick Flip Anti Tipper	109MA23	E0971	\$120.00
---	---------	-------	----------

29. POSITIONING STRAP

<input type="checkbox"/> Auto	109MA03	E0978	\$52.00
<input type="checkbox"/> Airline	109MA04	E0978	\$76.00
<input type="checkbox"/> Velcro	109MA05	E0978	\$48.00



30. ATTACHMENT TIE DOWN

* <input type="checkbox"/> Wheelchair Attachment Tie Down	109MA58	-	\$150.00
---	---------	---	----------

*Non occupied

31. MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> Padded Lap Tray 15"-16"	109MA09	E0950	\$230.00
<input type="checkbox"/> Padded Lap Tray 17"-18"	109MA10	E0950	\$230.00
<input type="checkbox"/> Padded Lap Tray 19"-20"	109MA11	E0950	\$230.00
<input type="checkbox"/> Padded Lap Tray 21"-22"	109MA12	E0950	\$385.00
<input type="checkbox"/> Padded Lap Tray 23"-24"	109MA13	E0950	\$385.00
<input type="checkbox"/> Half Padded Lap Tray, left side	109MA14	E0950	\$350.00
<input type="checkbox"/> Half Padded Lap Tray, right side	109MA15	E0950	\$350.00
<input type="checkbox"/> Elbow Stop, left	109MA16	K0108	\$94.00
<input type="checkbox"/> Elbow Stop, right	109MA17	K0108	\$94.00
<input type="checkbox"/> Shoulder Lateral Pad (each)	109MA18	-	\$90.00
<input type="checkbox"/> Adjustable Side Pad	109MA19	-	\$157.00
<input type="checkbox"/> IV Pole, Telescopic	109MA20	K0105	\$175.00
<input type="checkbox"/> Oxygen Tank Holder	109MA21	E2208	\$175.00



WARRANTY:

Future Mobility Healthcare warrants the Capella/Capella 45 to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user.

NOTES: