

## Cirrus Plus LT

HCPCS: K0822



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Foldable, lightweight frame for easy portability
- Attractive steel cross-brace frame
- Armrests are padded, adjustable and removable
- Adjustable-length controller mount
- Adjustable-tension back upholstery
- Swing-away footrests with heel loops
- Positioning belt included
- Padded back upholstery with storage pouch and cushioned seat
- Programmable controller
- Flat-free, non-marking tires
- Easy free-wheel operation
- Manual wheel locks
- General foam seat cushion included
- Anti-tip wheels

| PRODUCTS | Description                           |
|----------|---------------------------------------|
| CP18FBAN | CirrusPlusPowerChair,18 FlpBkArm,1/cs |

| SPECIFICATIONS        | Description               |
|-----------------------|---------------------------|
| Anti-Tip Wheels       | Yes                       |
| Base Weight           | 90 lbs                    |
| Batteries (pair)      | 12V x 33AH                |
| Battery Weight (pair) | 47 lbs                    |
| Brakes                | Electromagnetic           |
| Casters               | 8" x 2" solid             |
| Charger               | 4A Offboard               |
| Climbing Angle        | 6 degrees                 |
| Color                 | Black                     |
| Controller            | PG Vsi 50A                |
| Drive Wheels          | 12.5 x 2.25 solid PU tire |
| Freewheel Mode        | Yes                       |
| Ground Clearance      | 4"                        |
| Heaviest Piece        | Frame: 90 lbs             |
| Material              | Steel                     |
| Max Speed             | 5 mph                     |
| Maximum Range         | 12.5 miles                |

| SPECIFICATIONS       | Description   |
|----------------------|---|
| Motor                | 2 x (24V 320W x 4600 rpm)   |
| Number of Wheels     | 4   |
| Product Weight       | Without Battery: 90 lbs   |
| Seat Depth           | 16.5"   |
| Seat To Floor Height | 21.5"   |
| Seat Width           | 18"   |
| Turning Radius       | 30"   |
| Warranty             | Batteries: 12 Months   Drive Train Components: 24 Months   Electronics: 24 Months   Frame: Limited Lifetime |
| Weight Capacity      | 300 lbs   |
| ACCESSORIES          | Description   |
| CPELR                | Elevating Legrests, CPN, 1/pr, EA   |