

Active Care™

Fixed Height Bed



**Safely
enter and exit the
bed from a full chair position with the
push of ONE button.**

Since the inception of the hospital bed - one of the biggest challenges for caregivers and patients is **safely** getting in and out of the bed. Today, Med-Mizer has innovated the **ActiveCare** Bed to solve this challenge.

- **OneButtonSafeTurnExit** – exclusive feature that fully automates the bed to turn into a chair and allow the patient to safely enter and exit the bed.
- **FullChair Position** – full chair position to allow the patient the ability to sit in the bed for added comfort. Watch TV, eat a meal, and socialize with family and friends while remaining comfortably in bed.
- **RollBack inPlace™** – head of bed can articulate to 80 degrees with the assistance of rollback in place to alleviate shearing/friction in the lower back.
- **FixedHeight** – comes with a low height of around 19" (standard height of a wheelchair) and has the ability to have casters added (3" and 5") to raise the bed height to allow the end user an appropriate height to stand up and sit down.

ActiveCare® Fixed Height Bed SPECIFICATIONS


Deck Height Low	12.5"
Deck Height High	15.5" or 18.5"
Sleeping Surface	35"w, 42"w x 80"l
Maximum Safe Working Capacity	600 lbs

LIMITED WARRANTY


15 years on bed frame mechanics	1 year on wood, casters & handsets
2 years on drive system	Lifetime warranty on all welds



ENDBOARD OPTIONS



MATTRESS OPTIONS



Deluxe - 7" Multi Layered foam with middle die cut for superior pressure redistribution



Elite - 7" Copper gel infused memory foam for superior immersion and heat dissipation

SIDERAIL OPTIONS

Mobility Bar

Up/Down Rail



ACTIVECARE PORTFOLIO



Deluxe - Premium Aesthetics, Night Light, USB Charging, Hidden Casters, Battery Backup



Standard - Powered Height with OneButtonSafeTurn

PROUDLY MADE IN THE U.S.A.

(877) 867-7365
MED-MIZER.COM



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS PREVENTION



BARIATRIC CARE



MOBILITY



REHABILITATION



WOUND CARE