

RECONDITIONING PROCESS

It is our goal is to produce a product that meets or exceeds it's originally manufactured condition and per-formance. Our hospital beds undergo an in-depth refurbishment process to thoroughly and professionally restore them back to manufacturer's specifications.

DISASSEMBLE

All components are removed and inspected.

CLEAN

Equipment is thoroughly washed removing any dirt, grime, and tape residue to sanitize, prep for painting and repair.

REPAIR

- Damaged and/or missing parts are repaired and/or replaced with OEM replacement parts.
- Moving parts are lubricated and/ or greased.
- Casters are removed and replaced to look like new.

REPAINT

- All components are removed and disassembled for sanding, deep cleaning and repainted.
- Frame is painted for a "like new" appearance.
- Headboards and footboards are repainted if necessary.
- Metal components masked and painted if necessary.

REASSEMBLE

Equipment is reassembled.

RELABEL

- Equipment is relabeled with OEM equipment labels and decals
- OEM serial # will remain intact.

FINAL INSPECTION

- All motors and functions are biomedically checked, tuned and electrically inspected.
- All side rail functions are checked and tuned.
- All circuit boards are biomedically checked.
- All beds receive a final cleaning.