

[Documents/SDS](#)

No SDS content is available for this product.

#### Bed Extension Kits

- 4" L Extension Kit for Homecare Beds**
- 4" Bed extension kit increases bed length from 80" to 84".
  - Easy assembly, no tools required

- Bed retains weight capacity

Packaging	1 Each / Each	<a href="#">Report Incorrect Product Information</a>
Manufacturer	Medline	
Manufacturer #	MDR630EXTN	

SKU **MDR630EXTN**

#### Specifications

BARIATRIC	No
BED MODEL COMPATIBILITY	Medline Hc Beds
BED SPECIFIC	Yes
LATEX FREE	Yes
RETAIL PACKAGING	No
REUSABLE OR DISPOSABLE	Reusable
UNSPSC	42191807
UNIT COMPATIBILITY	MDR107003L

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Each (EA)	1.0 Each	0.0 / 12.787	1.79	37.795 x 6.299 x 12.992	40884389126357

California Customers Only... WARNING

#### HCPCS Information

**HCPCS Code: NO CODE - PRODUCT DOES NOT MEET HCPCS SPECIFICATIONS**

**HCPCS Disclaimer:** HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.