

**Momentum Shower Chairs**

- Momentum Shower Chair, Grey** • Our Momentum shower chair for adults features an advanced design with easy assembly, pushbutton height adjustment, and slip-resistant feet for grip on wet surfaces
- Strong aluminum frame resists rust and scratches and supports up to 250 lb.
 - Armrests and curved backrest provide greater comfort and security than a typical shower bench
 - Cushioned seat can be removed from the shower seat for easy cleaning
 - Seat cushion, armrests and backrest are infused with Microban antimicrobial* product protection to help prevent stains and odors caused by bacteria, mold and mildew
 - *These antimicrobial properties are built in to protect the product. The product does not protect users or others against bacteria, viruses, germs or other disease-causing organisms.

Packaging	1 Piece / Each	Report Incorrect Product Information
Manufacturer	Medline	
Manufacturer #	MDSMOMCHAIRG	

SKU MDSMOMCHAIRGH**Specifications**

COLOR	Gray
LATEX FREE	No
MRI SAFE	No
MATERIAL	Aluminum W/ Foam Padding
SEPARATE BACK AVAILABLE	No
UNSPSC	42211610
WEIGHT CAPACITY POUNDS	250 lb

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Each (EA)	1.0 Each	0.0 / 14.95	2.243	18.15 x 9.125 x 23.4	00888277720841
Piece (PC)	1.0 Each	0.0 / 14.95	2.243	18.15 x 9.125 x 23.4	00888277720841
Case (CS)	2.0 Each	0.0 / 32.5	4.529	18.25 x 18.25 x 23.5	20888277720845

California Customers Only... WARNING

HCPCS Information**HCPCS Code: E0240**

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

For more information on HCPCS codes, [Click here](#)

[Documents/SDS](#)

No SDS content is available for this product.

VIDEOS