

[Documents/SDS](#)

No SDS content is available for this product.

**Economy Overbed Table****Economy Overbed Table, Walnut** Ideal for activities involving hemi, super-hemi, pediatric height wheelchairs and low beds

- Contact your Medline sales representative for more economy options

Packaging	<b>1 Case / Each</b>	<a href="#">Report Incorrect Product Information</a>
Manufacturer	<b>Medline</b>	
Manufacturer #	<b>MDS104015</b>	

SKU **MDS104015****Specifications**

BASE STYLE	H-Base
CASTER STYLE	Threaded Twin Wheel
LATEX FREE	Yes
LENGTH INCHES	28 "
MAX HEIGHT INCHES	42 "
MIN HEIGHT INCHES	30 "
TABLE STYLE	Flush Edge
UNSPSC	42191801
VANITY	No
WEIGHT CAPACITY POUNDS	50 lb
WIDTH INCHES	15 "

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Case (CS)	1.0 Each	0.0 / 22.0	1.566	30.8 x 18.3 x 4.8	10080196768530
Each (EA)	1.0 Each	0.0 / 22.0	1.566	17.93 x 3.94 x 30.54	00080196768533

California Customers Only... WARNING

**HCPCS Information****HCPCS Code: E0274**

**HCPCS Disclaimer:** HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

For more information on HCPCS codes, [Click here](#)