



**Composite H-Base Overbed Table**

**Overbed Table with 30" x 15" Composite Top** Table features a 30" x 15" (76 x 38 cm) molded composite top that is easy to keep clean

- It is impervious to moisture and will not crack or peel like wood tables
- 1-hand adjustable from 30" to 42" (76 to 107 cm)
- Durable hammertone base provides stability
- Holds up to 50 lb. (23 kg)

Packaging	<b>1 Case / Each</b>	<a href="#">Report Incorrect Product Information</a>
Manufacturer	<b>Medline</b>	
Manufacturer #	<b>MDS107015</b>	

SKU **MDS107015**

**Specifications**

BASE STYLE	H-Base
CASTER STYLE	Threaded Twin Wheel
LATEX FREE	Yes
LENGTH INCHES	30 "
MAX HEIGHT INCHES	41"
MIN HEIGHT INCHES	30 "
TABLE STYLE	Composite
UNSPSC	42191801
VANITY	No
WEIGHT CAPACITY POUNDS	50 lb
WIDTH INCHES	14.75"

[Documents/SDS](#)

No SDS content is available for this product.

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Case (CS)	1.0 Each	0.001 / 20.6	1.249	30.75 x 17.55 x 4.0	00080196746333
Each (EA)	1.0 Each	0.001 / 20.6	1.249	30.75 x 17.55 x 4.0	00080196746333

California Customers Only... WARNING

**HCPCS Information**

HCPCS Code: E0274

*HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.*

*For more information on HCPCS codes, Click here*