

**Tilt Top H-Base Overbed Tables**

- H-Base Overbed Table with 30" x 15" Walnut Tilt Top**
- Tabletop measures 30" x 15" (76 x 38 cm)
 - Three-position tilt with squeeze of the handle

- Height-adjustable from 30" to 42" (76 to 107 cm)
- Separate tilt and raise/lower handles allow for safe and independent operation of table functions

Packaging	1 Case / Each	Report Incorrect Product Information
Manufacturer	Medline	
Manufacturer #	MDS104950	

SKU MDS104950

Specifications

BASE STYLE	H-Base
CASTER STYLE	Threaded Twin Wheel
LATEX FREE	No
LENGTH INCHES	30 "
MAX HEIGHT INCHES	42 "
MIN HEIGHT INCHES	30 "
TABLE STYLE	Tilt Top
UNSPSC	42191801
VANITY	No
WEIGHT CAPACITY POUNDS	50 lb
WIDTH INCHES	15 "

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Case (CS)	1.0 Each	0.0 / 27.6	1.826	21.0 x 4.54 x 33.09	10080196768547
Each (EA)	1.0 Each	0.0 / 23.81	1.515	20.472 x 3.937 x 32.48	00080196768540

California Customers Only... WARNING

HCPCS Information*HCPCS Code: E0274*

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

For more information on HCPCS codes, Click here

[Documents/SDS](#)

No SDS content is available for this product.