



Fold-Flat Overbed Table

Fold-Flat Overbed Table with Basket, Walnut An overbed table with basket accessory that folds flat for convenient, compact storage

- Durable, water-resistant composite material with spill edge and hammertone finish
- Weighs only 18.8 lbs. (8.5 kg)
- Height adjustments: 29.5 to 41" (75 to 104 cm)
- Weight capacity: 50 lbs. (23 kg)

[Documents/SDS](#)

No SDS content is available for this product.

Packaging	1 Each / Each	Report Incorrect Product Information
Manufacturer	Medline	
Manufacturer #	MDS108015	

VIDEOS

SKU **MDS108015**

Specifications

BASE STYLE	H-Base
CASTER STYLE	Dual Wheel Casters
LATEX FREE	No
LENGTH INCHES	30 "
MAX HEIGHT INCHES	41"
MIN HEIGHT INCHES	29.5"
TABLE STYLE	Composite
UNSPSC	42191801
VANITY	No
WEIGHT CAPACITY POUNDS	50 lb
WIDTH INCHES	15 "

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Each (EA)	1.0 Each	0.0 / 23.733	1.819	15.748 x 6.299 x 31.693	10888277417656

California Customers Only... WARNING

HCPCS Information

HCPCS Code: NO GOVERNMENT OR SUPPLIER HCPCS CODE VERIFICATION

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.