



**Freedom Plus Lightweight Bariatric Transport Chairs**

- **22"W Freedom Plus Lightweight Bariatric Transport Chair**
  - Back folds down for easy storage and transport
  - Handbrakes and a seat belt provide better

- control and safety
- Easy-to-clean vinyl upholstery, push bar with foam grips, cup holder and chart pocket on the back offer convenience
- Full-length permanent armrests and elevating leg rests offer comfort
- Recommended for use with antitip device

Packaging	<b>1 Case / Each</b>	<a href="#">Report Incorrect Product Information</a>
Manufacturer	<b>Medline</b>	
Manufacturer #	<b>MDS808200BAR</b>	

[Documents/SDS](#)

No SDS content is available for this product.

SKU **MDS808200BAR**

**Specifications**

ARM LENGTH	Full Length
ASSEMBLY REQUIRED	No
AXLE TYPE	Single
BACK HEIGHT	22 "
BARIATRIC	Yes
CASTER SIZE	8 "
COLOR	Red
FABRIC TYPE	Vinyl
FOLDABLE	Yes
HANDBRAKES	Yes
LATEX FREE	Yes
LEG REST STYLE	Elevating
MRI SAFE	No
MATERIAL	Aluminum
OVERALL HEIGHT	37.5 "
OVERALL WIDTH	29.5 "
OVERALL WIDTH INCHES	29.5"
PRODUCT TYPE	Transport Wheelchair
PRODUCT WEIGHT (LB)	42 lb
QUICK RELEASE AXLES	No
REAR WHEEL SIZE	12 "
REMOVABLE ARM	No
SEAT DEPTH	19 "
SEAT WIDTH	22 "
SEAT TO FLOOR HEIGHT	19.5"

UNSPSC 42192210  
WARRANTY Yes  
WEIGHT CAPACITY 400 lb.

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Case (CS)	1.0 Each	0.0 / 45.6	7.254	27.25 x 11.5 x 40.0	10080196285976
Each (EA)	1.0 Each	0.0 / 47.069	5.747	27.165 x 10.433 x 35.039	00080196285979

California Customers Only... WARNING

## HCPCS Information

*HCPCS Code: E1039*

*HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.*

*For more information on HCPCS codes, [Click here](#)*