



ALL NEW!
**Great for Indoor
& Outdoor Use!**

**BUZZ
AROUNDLT™**
3 wheel GB121

The all-new Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! This sporty design features a plush, stadium style seat, a weight capacity of 300 lbs., and increased foot room. The new Buzz LT also has a charger port in the delta tiller, a new hard plastic rugged basket, and painted color shrouds. Enjoy its tight indoor maneuverability and smooth outdoor use with 3" of ground clearance!

 **GOLDEN**



3 wheel Specifications

Model	GB121
Medicare Code	K0800
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range ¹	8 miles
Ground Clearance Under Deck	3"
Turning Radius	38"
Type Batteries	2-12AH
Colors	Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
Dimensions:	
Length	41"
Width	21.5"
Height (ground to top of seat back)	32" - 34"
Ground to Top of Seat	18" - 20"
Ground to Top of Deck	5"
Deck to Top of Seat	13.5" - 15.5"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	31.5"
Weight of Unit: (assembled)	98.5 lbs
Front Half	30.5 lbs
Rear Half	24 lbs
Seat w/Arms	23 lbs
Batteries (pack with 2 batteries)	21 lbs
Tires:	Solid, Low Profile
Front	7"
Rear	7.5"
Rear Anti-Tippers	2"
Standard Seat:	Stadium Style
Back Height - no headrest	14.375"
Width x Depth	16.5" x 15.5"
Color	Black Vinyl
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Battery	1 Year

¹. Battery range at maximum weight capacity and will vary due to rider weight, drive surface, terrain, and battery type.

². Electronics warranty excludes batteries.

Features

- Designed for indoor and outdoor use
- 300 lb. weight capacity
- Offset Delta tiller provides added clearance at the knees
- Plush stadium style seat standard
- Increased foot room on a non-slip deck
- Painted color shrouds
- Charging port on the tiller
- Disassembles quickly and easily for transport



Delta tiller folds easily.



Adjustable, flip up arm rests for customized comfort.



Conveniently located charging port on the delta tiller.



Quickly disassembles into five, lightweight pieces for easy transport!

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

