

## Bariatric Aluminum Rehab Shower Commode Chair

HCPCS:



- Versatile 3-in-1 design enables this model to easily be positioned over a standard toilet, used as a portable self-contained commode or shower chair.
- Simply and quickly adjust the height to achieve the ideal seated position.
- Padded flip-back armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing.
- Easy to empty commode bucket to reduce splashes and spills. Use commode liners for even quicker clean up and less mess!
- 5" locking casters allow for extra security when the unit is not in motion.
- Durable and easy-to-clean padded seat (22" x 17") and backrest.
- Seat opening allows for unobstructed access to the perineal area for easy and effective cleansing.
- Swing-away detachable footrests with height and angle adjustments.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 500 lbs.
- Simple, tool-free assembly.

PRODUCTS	Description
NRS185008	Bari Shower&Commode Chair,1/ea

SPECIFICATIONS	Description
Gross Weight	43.3 lbs
Net Weight	36.8 lbs
Overall Height	40" - 42"
Overall Length (with riggings)	35.04"
Overall Width	26"
Seat Dimensions	22" (W) x 17.1" (H)
Seat To Floor Height	20.5" - 22.4"
Warranty	Limited Lifetime
Weight Capacity	500 lbs

PARTS	Description
NRS185008-SEAT	Seat, NRS185008, 1/ea, EA
STDS851-SA	Velcro Seat Belt Self Attaching,1ea, EA
NRS185006	Self Propel Shower Commode Chair,1/ea, EA
NRS185008	Bari Shower&Commode Chair,1/ea, EA
NRS185006-02	SeatClipW/ScrewsForNRS185006/7, 1/ea, EA
NRS185006-05	ArmrestBracketLeftForNRS185006/7, 1/ea, EA
NRS185006-04	ArmpadForNRS185006/7, 1/ea, PR

PARTS	Description
NRS185006-06	ArmrestBracketRightForNRS185006/7, 1/ea, EA
NRS185006-08	CasterWithLegForNRS185006/7, 1/ea, EA
NRS185006-09	FrameClipForNRS185006/7, 4/pkg, PKG
NRS185006-11	SwingawayFootrestForNRS185006/7, 1/pr, PR
NRS185006-12	FootplateForNRS185006/7, 1/pr, PR
NRS185006-14	DowntubePlasticInsrtsForNRS185006/7,1/ea, EA
NRS185006-13	AngleAdjFtplateMntW/HdwrNRS185006/7,1/ea, EA
NRS185006-15	FootrstHgtAdjKnobW/BoltNRS185006/7,1/ea, EA
NRS185006-16	FootrstPivotSaddleW/BoltNRS185006/7,1/ea, EA
NRS185006-17	FootrestFrameLatchForNRS185006/7, 1/ea, EA
NRS185006-21	BackLockingKnobForNRS185006/7, 1/pr, PR
NRS185006-22	AntiRattleInsertForNRS185006/7, 1/ea, EA
11106-1	Commode Bucket,7.5qt.1/ea, EA