

ROADSTER

SCOOTERS



ROADSTER 3

HCPCS Code

K0800

Weight Capacity: 300 lbs

Range: Up to 6 miles

Heaviest Piece: 27 lbs.

FEATURES

- Extra charging port underneath the dashboard
- Detachable battery pack allows you to charge on the go
- 3 wheel provides a tighter turning radius for better indoor maneuverability
- Designed to easily disassemble into five light components for transport and storage.
- Adjustable tiller angle for ergonomic driving and easy transferring of the unit.
- 18" width deluxe stadium seat



ROADSTER 4

HCPCS Code

K0800

Weight Capacity: 300 lbs

Range: Up to 6 miles

Heaviest Piece: 32 lbs.

FEATURES

- Extra charging port underneath the dashboard
- Detachable battery pack allows you to charge on the go
- 4 wheel provides more stability
- Designed to easily disassemble into five light components for transport and storage.
- Adjustable tiller angle for ergonomic driving and easy transferring of the unit.
- 18" width deluxe stadium seat



ROADSTER DELUXE

HCPCS Code

Retail

Weight Capacity: 260 lbs

Range: Up to 8 miles

Heaviest Piece: 25 lbs.

FEATURES

- Hydraulic-adjustable tiller angle for ergonomic driving and easy transferring
- Extra charging port underneath the dashboard
- Detachable battery pack allows you to charge on the go
- Designed to easily disassemble into five light components for transport and storage.
- Convenient LED headlight and large plastic basket.
- 18" width deluxe stadium seat

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-3252

sales@meritsusa.com

merits

	Roadster 3	Roadster 4	Roadster Deluxe
Model Number	S731A	S741A	S731
Weight Capacity	300 lbs	300 lbs	260 lbs
Seat	18"x17"x15"	18"x17"x15"	18"x17"x15"
Seat-to-Deck Height	16"-18"	16"-18"	16"-18"
Overall Length	40"	38"	40"
Overall Width	21"	21"	21"
Overall Height	36"	36"	35"
Turning Radius	35.4"	43.3"	35.4"
Ground Clearance	2"	2"	2"
Weight of Heaviest Piece	35 lbs	32 lbs	25 lbs
Total Weight	99 lbs	106 lbs	101 lbs
Max Speed	4 mph	4mph	4mph
Per Charge Range	up to 6 miles	up to 6 miles	up to 8 miles
Batteries	12V/12Ah	12V/12Ah	12V/15Ah
Gradient	6°	6°	6°
HCPCS Code	K0800	K0800	Retail
Warranty	Motors and electronics - 18 months, batteries - 1 year and frame - 5 years		

OPTIONS

Oxygen Tank Holder
 Oxygen Tank Bag Holder
 Cane & Crutch Holder
 Walker Holder
 Foldable Rear Basket (Med & Large)
 Canopy
 Cover
 Mounting Bracket for Cup Holder
 Side Bag
 Safety Flag
 USB Charger
 Dock-N-Lock System with Fabric Basket
 Gray- Powerchair Cover
 Brindle- Powerchair Cover



merits

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-3252

sales@meritsusa.com