



PRIUS STRIVING TO BE THE DIFFERENCE



VESOFLOW IPCS Prophylaxis System

- **The whisper** quiet pump guarantees a good nights sleep and promotes patient comfort.
- **Three Systems In One**
The VESOFLOW IPCS system is a 3-in-1 design which provides choice of foot, calf or thigh compression, ensuring "the Right Therapy at the Right Time."
- **The VESOFLOW IPCS** pump is a completely sealed unit preventing the ingress of fluids, therefore reducing infection control risks.

VESOFLOW IPCS

Prophylaxis System

Ease Of Use

Easy to apply garments with hook and loop closures ensure snug fit and simple adjustment to different leg shapes.

Simple and secure snap-lock connection of garments to the pump prevents accidental disconnection yet promotes easy set up.

Snaplock air hoses make for easy storage and simplicity when cleaning.

Timed Therapy


Prius encourages compliance, therefore the visual timed therapy indicator allows the caregiver to monitor therapy time.

Three Systems In One

The system is a 3-in-1 design which provides choice of foot, calf or thigh compression, ensuring "the Right Therapy at the Right Time."




- **Our 100% LATEX free** garment range ensures higher patient compliance. Patients limbs remain cool and moisture free reducing all risk of skin maceration.
- **Timed therapy** enables the caregiver to monitor therapy times and therefore compliance.
- **The Safety lockout feature** guarantees correct compression pressures at all times thereby ensuring patient comfort and safety.
- **Large LED Screen** allows for a clear view of actual pressure being delivered to the patient's limb ensuring patient safety and caregiver peace of mind.



VFI - 004 (L)
LATEX Free Thigh Garment (pair), Large.
For thigh, size up to 35" circumference

VFI - 003 (M)
LATEX Free Thigh Garment (pair), Medium.
For thigh, size up to 28" circumference




VFI - 1 (XS)
LATEX Free Calf Garment (pair), Extra Small.
For calf, size up to 13.8" circumference


VFI - 001 (S)
LATEX Free Calf Garment (pair), Small.
For calf, size up to 16.9" circumference

VFI - 002 (M)
LATEX Free Calf Garment (pair), Medium.
For calf, size up to 18.9" circumference

VFI - 005(B)
LATEX Free Calf Garment (pair), Bariatric.
For calf, size up to 31.9" circumference



VFI-FOOTL (L)
LATEX Free VFI Foot garment (pair) large
Size 9 - 11



VFI-UNIF (U)
LATEX Free VFI Foot garment (pair)
Regular Size 7 - 9

VFAH 1.5m
Air hose extension of 60 inches (pair).

VFAH 3m
Air hose extension of 118 inches (pair).

VFAH 4.5m
Air hose extension of 177 inches (pair).



PUMP SPECIFICATIONS:

| | |
|-----------------------|----------------------|
| Size | 12" x 4.3" x 7.7" |
| Weight | 5.1 lbs |
| Pressure Range | 40 and 120 mmHg |
| Input Rating | 100-240V 60/50Hz |
| Fuse Rating | 1A 250V |
| Classification | Class I, Type BF |
| Operation Humidity | 30-75% |
| Operation Temperature | 59° - 95F |
| Applied Part | Garment and Air Tube |

